SR Triplex[®] Cold

Instructions for Use

Gebrauchsinformation

Mode d'emploi

Istruzioni d'uso

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Bruksanvisning

Brugsanvisning

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Bruksanvisning

Productinformatie

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Date information prepared: 2018-06-14, Rev. 1
563814/WE3

The logo graphic consists of a series of black dots of varying sizes arranged in a semi-circular pattern above the text.
ivoclar
vivadent[®]
technical

Description

SR Triplex® Cold is a cold-curing, PMMA-based denture base material for the pouring technique. It offers you ease of handling. Additionally, a veined version of the material is available. The shades are coordinated with those of the heat-curing SR Triplex Hot denture base material. To achieve optimum results, read these Instructions carefully.

Composition

Powder

Polymethyl methacrylate, catalyst, pigments

Liquid

Methyl methacrylate stab., dimethacrylate, catalyst, stabilizer

Indication

- Partial dentures
- Combination dentures
- Relining
- Repairs

Contraindication

- Direct intraoral contact of unpolymerized material.
- If the patient is known to be allergic to any of the ingredients in SR Triplex Cold

Side effects

In individual cases, local allergic reactions to polymethyl methacrylate materials have been reported.

Application procedure

Preparation

Isolate the boiled out, well-wetted, hand-warm plaster surfaces with two layers of Ivoclar Vivadent Separating Fluid and allow it to dry. To ensure an adequate bond with the denture base, well roughen the teeth, provide mechanical retention and wet with monomer.

Dosage

- **Ideal mixing ratio for one partial denture**
13 g polymer (powder) : 10 ml monomer (liquid)
- **With dosage system**
Polymer: 1st graduation on measuring cylinder
Monomer: 10 ml

Too high a monomer content alters the shade and handling properties of the material. Additionally, polymerization is impaired. This may result in irritation of the mucous membrane.

Mixing

Thoroughly mix polymer and monomer in the given ratio with a spatula. Allow the material to rest for 15 seconds to permit any trapped air to rise.

Flow phase

At room temperature (23 °C / 73 °F), the flow phase is 2.5 to 3 minutes. Pour the material into the saddle within this time span.

A higher room temperature shortens the working and the setting time.

Modelling phase

After a transition period of approx. 5 to 6 minutes, the material has set. It can be modelled during an additional 3 minutes.

Maturing time and working time depend on the room temperature.

Polymerization

Polymerization is carried out in a pressure device (eg Ivomat) at 40 °C / 104 °F and at 2 to 6 bar pressure for 15 minutes.

Finishing

Remove precast and finish in the usual manner. Subsequently, polish the surface to a high gloss.

Repair and correction possibilities of SR Triplex Cold

Corrections and repairs can be carried out with SR Triplex Cold or ProBase® Cold. Thoroughly roughen the corresponding surfaces with a bur and wet with monomer. Subsequently, apply the mixed material and polymerize according to the corresponding Instructions. The residual monomer content after polymerizing the material according to the method described is <4.5%.

Warnings

- The monomer contains methyl methacrylate (MMA).
- MMA is easily flammable and irritating (flash point +10 °C / 50 °F).
- Irritating to eyes, skin, and respiratory system.
- May cause sensitization by skin contact.
- Avoid contact of the skin with monomer or uncured material. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.
- Do not breathe vapour.
- Keep away from sources of ignition – no smoking.
- Do not empty into drains.
- Take precautionary measures against static discharges.

Storage

- Store material in a cool, dark, well-ventilated place.
Storage temperature: 2–28 °C (36–82 °F).
- Do not use materials after the indicated date of expiration.
- Keep out of the reach of children.

Classification according to EN ISO 20795-1, Type 2, Class 1

The product meets the requirements of EN ISO 20795-1:2013.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Polymethylmethacrylat, Katalysator, Pigmente

Flüssigkeit:

Methylmethacrylat stab., Dimethacrylat, Katalysator, Stabilisator

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.