# SR Nexco Paste

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
Inhibition Layer /  
Inhibitionsschicht

SR Nexco Stains +  
SR Nexco Dentin Shades

SR Nexco Incisal +  
SR Nexco Effect Shades

SR Gel

## Finish

See instructions

# SR Nexco Paste

Anterior crown framework-free / Gerüstfreie Frontzahnkrone

2x SR Nexco Liner

Inhibition Layer /  
Inhibitionsschicht

SR Nexco Stains +  
SR Nexco Dentin Shades

SR Nexco Incisal +  
SR Nexco Effect Shades

SR Gel

Finish

See instructions

# SR Nexco Paste

## Metal-supported restorations / Metallgestützte Restaurationen

Al<sub>2</sub>O<sub>3</sub>, 80-100 µm, 2-3 bar  
SR Link

2x SR Nexco Opaquer

Inhibition Layer /  
Inhibitionsschicht

SR Nexco Stains +  
SR Nexco Dentin Shades

SR Nexco Incisal +  
SR Nexco Effect Shades

SR Gel

### Finish

See instructions

# SR Nexco Paste

Zirconium oxide restorations / Zirkongestützte Restaurationen

$\text{Al}_2\text{O}_3$ , 80-100  $\mu\text{m}$ , max. 1 bar  
SR Link

2x SR Nexco Opaquer

Inhibition Layer /  
Inhibitionsschicht


SR Nexco Stains +  
SR Nexco Dentin Shades

SR Nexco Incisal +  
SR Nexco Effect Shades

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
### Important information:

Regular maintenance and functional control of the polymerization devices is required.

### Notes:

Units suitable for intermediate curing: Quick (Ivoclar Vivadent AG), Hilite pre (Heraeus), Visio Alfa (3M), Sublite V (Shofu), Steplight SL-I (GC).

\* Thinly apply a first layer of opaque and precure it. Then apply the second layer of opaque so that it covers the entire surface and cure it as indicated in the polymerization table.

\*\* Observe the maximum layer thickness! Precure, if required.

\*\*\* Stains should only be applied in extremely thin layers – a very low curing depth is achieved with dark shades.

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

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CE 0123

Rx ONLY

For dental use only!  
Made in Liechtenstein

Manufacturer:

Ivoclar Vivadent AG  
9494 Schaan/Liechtenstein  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)