Ivomat IP3

Operating Instructions
Bedienungsanleitung
Mode d'emploi
Istruzioni d'uso
Instrucciones de uso
Instruções de Uso

ivoclar
vivadent®
technical

**KONFORMITÄTSERKLÄRUNG
 DECLARATION OF CONFORMITY
 CERTIFICAT DE CONFORMITÈ
 DICHIARAZIONE DI CONFORMITÀ
 DECLARACIÓN DE CONFORMIDAD
 DECLARAÇÃO DE CONFORMIDADE**

**ivoclar
 vivadent**
 BENDERERSTR. 2
 FL-9494 LICHTENSTEIN
 TEL ++423 / 235 35 35
 FAX ++423 / 235 33 60

Produkt / Product / Produit / Prodotto / Producto / Produto

Ivomat IP3

- DE** Hiermit erklären wir in alleiniger Verantwortung, dass das oben aufgeführte Produkt den erwähnten Normen entspricht.
 Gemäss den Bestimmungen der EU-Richtlinie(n):
- GB** We herewith declare that the product listed above complies with the mentioned standards.
 Following the provisions of Directive(s):
- FR** Par la présente, nous déclarons que le produit ci-dessus indiqué est conforme aux normes énoncées.
 Conformément aux dispositions de la (des) Directive(s) CE:
- IT** Con la presente dichiariamo sotto la nostra responsabilità, che il prodotto sopra menzionato corrisponde alle norme citate.
 Secondo le disposizioni della/e Direttiva/e CEE:
- ES** Por la presente declaramos que el producto arriba indicado cumple con las normas citadas.
 Siguiendo las indicaciones de la Directiva:
- PT** Declaramos que o produto citado cumpre as normas mencionadas.
 De acordo com as especificações da(s) Diretriz(es):

87/404/EWG 73/23/EWG 89/336/EWG	EN 50081-1	1992
	EN 50082	1992
	EN 55011	1991
	EN 55022	1994
	EN 60335-1	1988
	EN 61000-3-2	1995
	EN 61000-3-3	1995

Schaan, 20.09.2002

Bürs, 20.09.2002

Dipl. Ing. Wolfgang Vogrin

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Ivomat IP3

Table of Contents

- 1. Introduction / Signs and Symbols**
 - 1.1 Preface
 - 1.2 Introduction
 - 1.3 Signs and Symbols

- 2. Safety First**
 - 2.1 Indications
 - 2.2. Health and safety instructions

- 3. Product Description**
 - 3.1 Components
 - 3.2 Functional description
 - 3.3 Hazardous areas and safety equipment

- 4. Installation**
 - 4.1 Unpacking and checking the contents
 - 4.2 Selecting the location
 - 4.3 Connections
 - 4.4 Connecting the return water canister

- 5. Start-Up**
 - 5.1 Switching on/off
 - 5.2 Filling in water

- 6. Operation**
 - 6.1 Placing the objects and starting the program
 - 6.2 End of the program and removing the objects



- 7. Maintenance, Cleaning, and Diagnosis**
 - 7.1 Monitoring and maintenance
 - 7.2 Cleaning
 - 7.3 Maintenance

- 8. What If ...**
 - 8.1 Power failure
 - 8.2 Repair
 - 8.3 Technical malfunctions

- 9. Product Specifications**
 - 9.1 Delivery form
 - 9.2 Technical data
 - 9.3 Acceptable operating conditions
 - 9.4 Acceptable transportation conditions

List of Parts

- 1) Air filter control
- 1a) Sight glass
- 1b) Release button
- 2) Safety lever
- 3) Bold knob
- 4) Pressure chamber lid
- 5) Pressure chamber
- 6) Thermostat
- 7) Manometer
- 8) Timer
- 9) Mains switch
- 10) Pilot lamp
- 11) Non-return valve
- 12) Compressed air hose
- 13) Power cord
- 14) Fuse cap
- 15) Fuse
- 16) Safety valve
- 17) Reset button
- 18) Water drain hose
- 19) Canister cap
- 20) Return water canister
- 21) Canister hose
- 22) Sealing washer for sound absorber
- 23) Sound absorber
- 24) Collar
- 25) Air hole
- 26) Collar screw
- 27) Polymerization container lid
- 28) Alligator clips
- 29) Polymerization container
- 30) Replenishing cup

1. Introduction / Signs and Symbols

1.1 Preface

Dear Customer,

Thank you for having purchased the IVOMAT IP3. This apparatus is a technically advanced product. Please read these Operating Instructions carefully and use the IVOMAT IP3 according to the instructions. Should you have any further questions, please contact your dealer or Ivoclar Vivadent directly.

1.2 Introduction

Apparatus: IVOMAT IP3

Target group: Dental professionals

The IVOMAT IP3 is suitable for the polymerization of dental resins. The Operating Instructions ensure safe, correct, and economic use of the IVOMAT IP3. They are divided into several, clearly structured chapters. This should enable you to locate specific topics quickly and easily.

To inform you about risks, dangers, important information, and contraindications, these Instructions contain corresponding signs/symbols to mark important paragraphs.

We recommend keeping the Instructions in a safe place near the apparatus to have immediate access to the information if necessary.

Should you lose the Operating Instructions, extra copies can be ordered at a normal fee from your local Ivoclar Vivadent Service Center.

1.3 Signs and symbols

The signs and symbols in these Operating Instructions facilitate the finding of important points and have the following meanings:

In the Operating Instructions

Risks and dangers

Information

Contraindication

On the apparatus

Alternating current

On

Off

2. Safety First

This chapter is especially important for individuals who work with the IVOMAT IP3 or who have to carry out maintenance or repair work. This chapter must be read and the corresponding instructions followed.

2.1 Indications

The IVOMAT IP3 has been especially developed for the polymerization of dental resins and materials and should be used for these purposes only.

Uses other than the ones stipulated are contraindicated. The manufacturer does not assume any liability for damage resulting from misuse. The user is solely responsible for any risk resulting from failure to observe these instructions.


Further instructions to assure proper use of the furnace:

- The instructions, regulations, and notes in these Operating Instructions must be observed.
- The apparatus must be operated under the indicated environmental and operating conditions (see Chapter 9)
- The IVOMAT IP3 must be properly maintained.

2.2.1

Risks and dangers


Make sure that no liquids or other foreign objects enter the air vents, since this may result in an electrical shock.

2.2.2

Risks and dangers

The housing must not be opened while it is connected to the power supply, since there is a risk of electrical shock. The connection cover may only be opened by a qualified service center.

3. Product Description

3.1 Components

The IVOMAT IP3 consists of the following components:

- Apparatus with air filter control
- Non-return valve
- 1.5 m compressed air hose
- 1.2 m water drain hose with canister cap
- Canister hose with sound absorber; polymerization container in pressure chamber, and 3 object clips secured to removable lid
- 1 return water container
- 1 replenishing cup

Special accessories

- IVOMAT refill container type IN 1
- Wallfixator

3.2 Functional description

The pressure chamber is equipped with a heater. In this way, the water can be heated to the temperature required for the corresponding material. The pressure chamber is impinged with compressed air to achieve optimum homogeneity of the material. The polymerization time can be set by means of the timer.

3.3 Hazardous areas and safety equipment

Description of the hazardous areas:

Hazardous area	Type of risk
Safety holder	Electrical shock
Air vents	Electrical shock

Description of the safety equipment

Safety equipment	Protective effect
Protective conductor	Protection from electrical shock
IP44 (system of protection)	Protection from electrical shock
Protective thermal switch	Protects the IVOMAT IP3 from overheating

4. Installation

4.1 Unpacking and checking the contents

Remove the apparatus from its packaging and check for possible transportation damage. We recommend keeping the original packaging for future transportation purposes. Use only the original packaging for transportation.

Do not carry the IVOMAT IP3 by means of the power cord. Support the bottom of the IVOMAT IP3 with both hands to carry it.

4.2 Selecting the location

Place the apparatus on a flat table using the rubber feet. Make sure that the apparatus is not placed in the immediate vicinity of heaters or other sources of heat. Furthermore, protect the apparatus from direct sunlight.

4.3 Connections

Make sure the voltage indicated on the rating plate complies with the local power supply. If this is not the case, the apparatus must not be connected.

Connecting the compressed air hose

Make sure that the hose is properly connected.

Power connection

Make sure that the power cord does not touch any parts of the apparatus that become hot during use.

The electrical installations of the room where the apparatus is located must comply with the national regulations and IEC standards.

5. Start-up

4.4 Connecting the return water canister

- Fill return water canister (20) to a depth of 3–5 cm.
- Insert the canister hose into the canister and screw cap on (19).
- Put the canister on the floor in a suitable place.
- Lift the safety lever (2) and engage. Push the bolt to the right using the knob (3) and lift the lid (4) of the vessel.
- Remove protective insert used during transportation.

5.1 Switching on / off

Take the polymerization container (29) with lid (27) out of the pressure chamber. Switch on main (9) (0 to I). The pilot lamp (10) will light up when the apparatus is ready for use.

The apparatus may be kept ready for use throughout the day (mains switch on I), but the pressure chamber lid (4) should remain either unbolted or open after each use.

5.2 Filling in water

Fill the pressure chamber with water.

- Maximum filling height: up to the market groove in the chamber
- Minimum filling height: objects to be polymerized must be completely under water

6. Operation

6.1 Placing the objects and starting the program

- Remove the lid (27) from the polymerization contains. Secure the objects to be polymerized on the alligator clips (28).
- Slowly insert the polymerization container with lid and objects, or models held by the Wallfixator, into the pressure chamber.
- Close chamber lid (4). Bolt the chamber by sliding the knob (3) to the left. Press safety lever down until it audibly clicks into place.

Chamber lid

- Set polymerization thermostat (6). Move the red indicator to the desired temperature.

Thermostat (°C)

The program automatically begins, once the timer is set. Note: For times under 5 minutes, first turn the timer to approximately 10 minutes, then immediately turn it back to the desired time. (If required, the time may be returned to zero during the automatic program in progress. The time may also be prolonged or shortened.)

Timer

- Automatic program:
 - The manometer (7) indicates the rise in pressure (check to make sure that the required operating pressure has been reached by the time the selected temperature has been achieved).

Manometer

- The water in the pressure chamber is heated. The black indicator on the thermostat registers the current temperature.
- Once the desired temperature is reached, the set polymerization time begins (timer countdown).
- After the polymerization time, the water is ejected from the pressure chamber and flows into the canister. At the same time, the pressure is released (indicated on the manometer).

7. Maintenance, Cleaning, and Diagnosis

6.2 End of the program and removing the objects

- Opening the pressure chamber lid:
Do not lift the safety lever before the manometer needle is at zero and the black thermostat indicator has fallen below 110 °C (230 °F). (Any residual pressure is blown off against the chamber lid.) The lid can then be safely unbolted by pushing the knob to the right and the chamber lid opened.
- Removing the polymerized objects:
Remove the polymerization container including the lid from the pressure chamber and allow it to cool for a few minutes before removing the polymerized objects. Cooling too quickly may result in stress within the resin materials.

- **Use softened water (not distilled) or water 4–9° hard (70–160 ppm/°US = 5–12° Brit.) (Ph 7–8/20 °C).** If tap water of more than 9° hard (16° Fr) is used, special care must be given to maintenance and cleaning or decalcification (see decalcification, page 11).
- **Make sure that no wax enters the pressure chamber.**
- **Set the polymerization time only after bolting the pressure chamber lid.**

- **Caution: Do not touch hot parts of the apparatus during use.**

For safety reasons, disconnect power cord for all cleaning and maintenance procedures.

7.1 Monitoring and maintenance

Keep pressure chamber sealing rim (5a) and sealing ring clean and undamaged. If using water as prescribed in point IV, remove the perforated plate and clean the pressure chamber, heater, outflow filter, and outflow system at least every three months. (Remove any wax residue using a suitable wax solvent.)

7.2 Cleaning

Dust the IVOMAT IP3 from time to time. Use, for example, a vacuum cleaner with a cleaning brush for that purpose.

Parts	Cleaning material
Housing	Cloth
Pressure chamber	Deliming agent

7.3 Maintenance

Please observe the relevant accident control regulations and other accepted rules regarding safety and industrial medicine. The connection cover, which is marked with the rating plate below may only be opened by a qualified service technician and with the power disconnected.

Decalcification

If hard tap water of more than 9° (160 ppm/°US = 5–12° Brit.) is used, the apparatus must be regularly checked, according to use about once a week, for possible calcareous deposits on the heating element. For that purpose, the perforated plate must be removed from the pressure chamber.

If there is a visible chalky deposit on the heating element, the device must be cleaned using a deliming agent as follows: (customary deliming agents used for coffee makers, boilers, etc. can be used).

8. What If ...

- Fill the pressure chamber with approx. 0.5 l water and dissolve the amount of deliming agent recommended in the corresponding instructions for use.
- Close and bolt the pressure chamber lid.
- Set the thermostat to 80 °C and the timer to 5 minutes.
- Once the program is completed, open the lid and check the cleaning results.
- Remove any large, loose particles from the chamber.
- If necessary, repeat the procedure until the heating element and pressure chamber are clean.
- After decalcifying, repeat the procedure with fresh water to rinse the pressure chamber and the outflow system. (Once the operating pressure is reached, the program can be interrupted by setting the timer to 0)
- Empty the return water canister and rinse with fresh water.

Deliming agents are caustic and toxic.
Avoid contact with skin and eyes. In case of accidental contact, rinse with copious amounts of fresh water.

Condensation accumulated in the sight glass (1a) can be drained off by pushing the release button (1b).

At three-months intervals, unscrew the canister cap (19) and clean (clear) the small air hole (25) of the canister hose (21). Access to the hole (25) is gained by loosening the screw (26) and removing the collar (24).

8.1 Power failure or blown fuse

- In the event of a power failure, the automatic program is interrupted. Once the power supply returns, the program will automatically continue.
- To take the objects out of the pressure chamber during prolonged power failures, proceed as described under 3a.
- If the fuse (15) of the apparatus is defective, open the fuse cap and replace the fuse.

If problems arise, which are not mentioned in the table below, please contact our Customer Service Department.

8.2 Repair

Repairs may only be carried out by a certified Ivoclac Vivadent Service Center. Please contact our Customer Service Department.

8.3 Technical malfunctions

Error or malfunction	Possible causes	Corrective action
1. Objects not fully cured, even though the program ran normally (Timer has run to "0", manometer shows "0").	<ul style="list-style-type: none"> – Too little water in the pressure chamber, or – temperature too low, or – time too short, or – air pressure too low 	1a) Repeat program (polymerization) with sufficient water, correct temperature and time, and stipulated air pressure.
2. Program fails to continue after polymerization time has been set.	Chamber lid not properly closed and bolted.	2a) Close and bolt chamber lid according to the instructions.
3. Automatic program has stopped too soon (Timer has not run to "0", manometer indicates pressure).	No or far to little water in the pressure chamber (excess heat protection of the heater engaged).	3a) Disconnect power. Release pressure by means of the safety valve. Turn air knob (16a) at the end of the safety valve (16) counter-clockwise until the manometer indicates "0" and the black needle on the thermostat registers less than 90 °C/194 °F. Open the chamber lid and allow chamber to cool for approx. 6 min. Then, turn air knob clockwise until it stops. Press red reset button (17) on the back of the apparatus.
	Heater and overheating sensor in the chamber very dirty or calcareous.	3b) Proceed as described in 3a). Remove perforated plate from the pressure chamber and clean it. Clean heater with overheating sensor and clean or delime the chamber according to the service recommendations.
	Power failure or defective fuse.	3c) See point 7.
4. Pressure falls abnormally slowly after the end of the program. (Manometer shows "0", timer has run to "0". Polymerization temperature has been reached.	Outflow filter in the pressure chamber is partly blocked.	4a) Remove perforated plate from the pressure chamber. Clean outflow filter, chamber, and heater, or delime according to the service recommendations.
	Outflow system partly blocked.	4b) Clean sound absorber (23). Unscrew water drain hose (18) from the back of the apparatus and blow compressed air through the hoses. Delime if necessary.

9. Product Specifications

9.1 Delivery form

The IVOMAT IP3 consists of the following components:

- Apparatus with air filter control
- Non-return valve
- 1.5 m compressed air hose
- 1.2 m water drain hose with canister cap
- Canister hose with sound absorber; polymerization container in pressure chamber, and 3 object clips secured to removable lid
- 1 return water container
- 1 replenishing cup

9.2 Technical data

Power supply: Single-phase alternating current

Standard version:

- 220 V / 50 Hz
- 240 V / 50 Hz

Acceptable voltage fluctuations: +10 % to –15 %

Power consumption: 1030 W

Electrical fuses:

100–118 V: 12.5 A, slow. Diameter: 6.3 x 32 mm

220 V: 6.3 A, slow. Diameter: 5 x 20 mm

240 V: 5 A, slow. Diameter 5 x 20 mm

Compressed air:

Operating pressure (already set): approx. 6 bar = 6.10 Pa
(Tolerance: green area of the scale)

System pressure: 6 to max. 12 bar = 6 to max. 12.10 Pa

Hose with an inside diameter of 6 mm

Dimensions of the closed apparatus:

Width: 312 mm, Depth: 302 mm

Height: 220 mm

Useful dimensions of pressure chamber:

Diameter: 112 mm, height: 98 mm

Max. water capacity 0.78 dm³ (Liter)

Weight: 9.75 kg

Safety information:

EN 61010, Part 1, EMC tested

9.3 Acceptable operating conditions

Acceptable temperature range

+5 °C to +40 °C (+41 °F to +104 °F)

Acceptable altitude:

The apparatus has been tested for use at altitudes of up to 2000 m

Atmospheric pressure of 500 mbar to 1060 mbar

9.4 Acceptable transportation and storage conditions

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Ivoclar Vivadent S.A.
Ivoclar Vivadent AB
Ivoclar Vivadent, Inc.