Manufacturer's notes

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
**ivoclar**
vivadent
clinical

General notes:

Please also observe the country-specific regulations and guidelines regarding the hygiene standards and reprocessing of medical devices in dental practices.

New non-sterile instruments must be subjected to a complete processing cycle before they are used for the first time. Clean and disinfect all instruments prior to each use. It is strongly recommended to additionally sterilize the instruments in a steam sterilizer!

Single-use items and single-use components of reusable items must not be cleaned, disinfected, sterilized and used again!

Users are strongly advised not to use damaged instruments. Maintenance and repair work should only be carried out by specialist staff.

Do not place instruments in NaCl solution (to avoid the risk of pitting and/or stress corrosion cracking).

Protective gloves and goggles complying with the requirements of directive 89/686/EEC must be worn to handle any used and contaminated instruments. In line with the recommendations of the Robert Koch Institute and the guideline compiled by the pertinent bodies (DGKH - German Society for Hospital Hygiene, DGSV - German Society for Sterile Supply and AKI - Working Group Instrument Preparation), automated processing in a washer-disinfector (WD) is to be generally preferred over manual processing. Manual processing should only be used if automated processing is not possible.

Only a WD (compliant with DIN EN ISO 15883) with proven effectiveness should be used to ensure reliable automated cleaning and disinfection. The operator is responsible for ensuring that validation, performance re-qualification and periodic routine checks of the processing procedures in the WD are defined, documented and carried out appropriately.

Restrictions on reprocessing:

The life span of the instruments is mainly determined by wear and damage through use. Frequent reprocessing has little impact on the performance of these instruments.

Processing at the patient treatment area:

Carefully remove large, unhardened residues from the instruments with a lint-free cellulose tissue immediately after use on the patient. Before you continue with the processing procedure, remove the modelling tip from the handheld instrument and discard it. The tips are designed for single use only. They must not be cleaned, disinfected, sterilized and re-used. The handheld instruments may either be placed in a removal container (dry removal) or directly immersed in a combination of detergent and disinfectant (wet removal), see below, for transportation to the processing and sterilization area. In any case, it is advisable to transport soiled instruments from the patient treatment area to the processing area as soon as possible. Waiting periods of several hours between contamination and processing should be absolutely avoided.

Preparation for cleaning/disinfection:

Preliminary treatment of the instruments is strongly recommended to protect personnel against infection and to prepare the instruments for the automated or manual processing procedure!

To prevent residues from drying on the surfaces and to ensure preliminary disinfection, place the instruments in a disinfection bath (with basket insert and lid) filled with aldehyde-free alkaline disinfectant (e.g. ID 212 forte, Dürr Dental, for 15 min). Make sure that the instruments are completely covered in disinfectant and that the disinfectant is free of bubbles. The disinfectant used should be suitable for the instruments in question and tested for effectiveness (e.g. DGHM or FDA certification or CE mark). Please adhere to the disinfectant manufacturer's instructions (concentration rate, temperature and contact time) to ensure that the disinfectant is used correctly.

If you use an automated processing method, carefully rinse the instruments under running tap water after preliminary disinfection to prevent residues of the cleaning and disinfecting agent from reaching the WD.

Automated cleaning and disinfection:

Requirements to carry out the validated procedure:

- WD: G 7882, Miele (A_0 -value > 3000)
- Program: Vario TD
- Cleaning and disinfecting agent: Neodisher FA, Dr. Weigert

Place the instruments in an instrument container that is suitable for your particular type of WD, according to the WD manufacturer's directions. Fill the WD with cleaning and disinfecting agent (e.g. Neodisher FA, Dr. Weigert) as appropriate for the intended process. Observe the directions of the disinfectant and equipment manufacturer. Start the cleaning and disinfection program (e.g. Vario TD) validated for the particular type of unit.

Drying:

After the cleaning and disinfection program has been completed, remove the instruments from the WD and completely dry any remaining moisture with compressed air or a clean, lint-free single-use towel.

Maintenance, checking and inspection:

Visually inspect the instruments to see if they are clean, undamaged and function properly. If macroscopically visible residual debris is present on the instruments, repeat the cleaning and disinfection process. If the instruments have only been cleaned in the WD (without automated disinfection), thermal disinfection of the unwrapped instruments in a steam sterilizer is absolutely mandatory!

Manual cleaning and disinfection:

Requirements to carry out the validated procedure:

- Soft brush: nylon instrument cleaning brush, Miltex
- Cleaning and disinfecting agent: ID 212 forte, Dürr Dental

Only use soft brushes for the manual removal for visible residues. Never use metal brushes or steel wool.

Brush the instruments with a suitable brush (e.g. nylon instrument cleaning brush, Miltex) in the disinfection bath of the preliminary treatment until you can no longer see any contamination on the surface. Remove the instruments

from the disinfection bath, leave them in the basket insert and carefully rinse them under running tap water (< 40 °C) until all residues are rinsed off. Do not use chemicals for this step. Allow remaining water to sufficiently drain off.

Next, place the instruments again in a disinfection bath (with basket insert and lid) filled with aldehyde-free alkaline disinfectant (e.g. ID 212 forte, Dürr Dental, for 15 min.) and allow the disinfectant to take action, similar to the preliminary treatment step. Make sure that the instruments are completely covered in disinfectant and that the disinfectant is free of bubbles. The disinfectant that you use should be suitable for the instruments in question and tested for effectiveness (e.g. DGHM or FDA certification or CE mark). Please adhere to the manufacturer's instructions (concentration rate, temperature and contact time) to ensure correct use of the disinfectant. Subsequently, carefully rinse the instruments in the basket under running tap water. Allow remaining water to sufficiently drain off.

Drying:

After cleaning and disinfecting, completely dry the instruments with compressed air or a clean, lint-free single-use towel.

Maintenance, checking and inspection:

Visually inspect the instruments to see if they are clean, undamaged and function properly. If macroscopically visible residual debris is present on the instruments, repeat the cleaning and disinfection process.

Thermal disinfection of the unwrapped instruments in a steam sterilizer is recommended after manual cleaning and disinfection.

Sterilization:

Requirements to carry out the validated procedure:

- Steam sterilizer: Systec HX-320
- Aluminium sterilization container
- Program: 3 x fractionated pre-vacuum, 134 °C (273 °F), 4 min.

It is strongly recommended to sterilize semi-critical instruments of Classes A & B!

Packaging:

The instruments should be autoclaved in an appropriate packaging (e.g. aluminium sterilization container).

Sterilize the items in a steam sterilizer (the unit should be compliant with EN 285 and validated according to DIN EN ISO 17665) using a fractionated pre-vacuum at 134 °C (273 °F) for at least 4 minutes. Please also observe the directions of the unit manufacturer.

If you sterilize several instruments at the same time in the steam sterilizer, make sure that you observe the unit manufacturer's stipulations on the maximum permissible load to avoid overloading.

The instruments are designed for at least 200 sterilization cycles in an autoclave.

The stipulated sterilization procedure has been validated by the manufacturer as appropriate for the purpose in question according to EN ISO 17664.

The manufacturer does not assume responsibility for other sterilization procedures. If other sterilization processes are used, the operator is responsible for validating them to ensure that they are effective and fit for the purpose.

Storage:

Store the instruments in dry conditions to avoid condensation build-up. It is recommended to store the modelling tips in their original packaging at room temperature, protected from dust, humidity, compression and contamination. Generally, storage at temperatures between 12 - 28 °C (54 - 82 °F) is appropriate. Please observe the date of expiry on the label!

Keep material out of children's reach.

For use in dentistry only.

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damage resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding. For USA: Rx only.

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www.ivoclarvivadent.pl**Ivoclar Vivadent Marketing Ltd.**
www.ivoclarvivadent.ru**Ivoclar Vivadent Marketing Ltd.**
www.ivoclarvivadent.es**Ivoclar Vivadent AB**
www.ivoclarvivadent.com**Ivoclar Vivadent Limited**
www.ivoclarvivadent.co.uk**Ivoclar Vivadent, Inc.**