Heliobond

English

EN Instructions for Use

- Light-curing bonding resin

DE Gebrauchsinformation

- Lichthärtender Haftvermittler

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ivoclar
vivadent®
clinical

Description

Light-curing, single-component bonding resin for optimizing the enamel-etch technique in combination with all light-curing restorative materials (Heliomolar®, Helio Progress®, Tetric®, Tetric EvoCeram®, Tetric EvoCeram® Bulk Fill).

Composition

Heliobond contains:

Bis-GMA

60 %wt.

Triethylene glycol dimethacrylate

40 %wt.

Indication

- Bonding agent for adhesive restorations
- Transparent sealant for fissures and restorations
- Bonding layer for the repair of resin-, crown-, or bridgework

Contraindication

Heliobond should not be used if patients are known to be allergic to any of the ingredients.

Interactions

To achieve an optimum bond to the enamel, the enamel-etch technique (e.g. Total Etch) must be applied. Contamination of the etched and dried enamel surface with saliva or blood will hamper this effect. Phenolic substances (e.g. eugenol) inhibit polymerization. The use of base materials containing these substances must be avoided.

Application

Preparation of enamel

- Depending on the case, bevel enamel surfaces to be treated using a diamond or clean them carefully with an abrasive (e.g. pumice powder and water). Do not use any oily or greasy paste. Possible residue from base materials on the enamel should be removed in this way. Then rinse and dry.

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- Apply Total Etch on the enamel surfaces to be treated and allow it to react for 15-30 s. Afterwards, carefully rinse with water spray and dry the etched surfaces with a stream of air (use only water- and oil-free air).
- A dry operating field must be ensured. Ideally, a rubber dam should be placed (e.g. OptraDam® Plus).
- The etched and dried enamel surface must not be contaminated prior to the application of Heliobond (if the surface comes into contact with saliva or blood, the etching and drying procedures must be repeated).

Dosage

- Dosage depends on the indication. When used as a bonding agent, a very thin layer of Heliobond should be applied. If necessary, use a stream of air to achieve an optimal thin layer.
- When used as a transparent fissure sealant, Heliobond must be applied using a thicker layer.

Method of application

1. Bonding agent for adhesive restorations

- Etch desired enamel surfaces (see preparation of the enamel). If necessary, cover exposures, dentin surfaces with Syntac® Adhesive (see Syntac instructions for use) or appropriate base material.
- Apply a thin layer of Heliobond onto the etched enamel surface using a brush or spherical instrument.
- An optimal, thin layer can be achieved using a stream of air.
- If used in combination with light-curing luting composites, Heliobond does not require separate polymerization. For direct restorations, separate polymerization of Heliobond is required: Light-cure for 10 s using a light intensity of 500 to 1100 mW/cm² (e.g. Bluephase®).
- Application of the composite.
- Polymerization of the composite.
- Finishing of the restoration.

2. Transparent sealant for fissures and pits

- Etch desired enamel surfaces (see preparation of the enamel).

- Apply Heliobond into the fissures using a suitable instrument or brush. Avoid entrapment of air and wait for approx. 15 s to allow penetration.
- Polymerize with light for 20 s using a light intensity of 500 to 1100 mW/cm² (e.g. Bluephase).
- After curing, remove inhibited layer. Check occlusion and grind off any excess material.

3. Preparation of resin repair work

- Roughen the smooth resin surfaces.
- If necessary, apply Monobond® Plus (see Monobond Plus instructions for use).
- Apply a thin layer of Heliobond (use a stream of air, if necessary).
- Polymerize with light for 10 s using a light intensity of 500 to 1100 mW/cm² (e.g. Bluephase).

Polymerization

- A polymerization time of 10 s is sufficient for thin layers using a light intensity of 500 to 1100 mW/cm² (e.g. Bluephase).
- Thicker layers (e.g. sealings) require 20 s curing time using a light intensity of 500 to 1100 mW/cm² (e.g. Bluephase).
- The light probe should be held less than 5 mm above the Heliobond surface and should not touch unpolymerized material

Note

- Do not expose Heliobond to intensive light during application since this will significantly shorten the working time.
- If Heliobond is used as a cover layer (e.g. sealant), a thin layer of unpolymerized, sticky material will remain on the surface after polymerization due to the oxygen inhibition. It is recommended to remove this layer with cotton rolls, pellets or by polishing.

Storage

- Close bottle immediately after use. Exposure to light may lead to premature polymerization.
- Storage: 2–28 °C / 36–82 °F.

-
-

We
Av
me
cal
ral
me
th
Ke
Fo

The
be c
accr
stip
for i
tion
binc

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- Do not use Heliobond after the date of expiration.
- Date of expiration: See note on the bottle/packaging.

Warning

Avoid contact of Heliobond with the skin/mucous membrane or eyes. Unpolymerized Heliobond may cause slight, reversible irritation, or lead to a general sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

**Keep material out of the reach of children!
For use in dentistry only!**

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.