IPS e.max[®] CAD

Abutment Solutions

EN Instructions for Use

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1 Intended use

Intended purpose

Implant-supported hybrid restorations for the replacement of single teeth

Patient target group

Adult patients with dental implants

Intended users / Special training

- Dentists (chairside fabrication of restorations; clinical workflow)
 - Dental laboratory technicians (fabrication of restorations in the dental laboratory)
- No special training required.

Use

For dental use only.

Description

IPS e.max® CAD Abutment Solutions are CAD/CAM-fabricated, implant-supported hybrid restorations (hybrid abutments and hybrid abutment crowns) for the replacement of single teeth. These hybrid restorations are individually fabricated from a lithium disilicate glass-ceramic (LS₂) block and cemented onto a titanium bonding base.

Performance characteristics

Property	Specification	Typical mean value
CTE (25–500°C) [10 ⁻⁶ /K ⁻¹]	10.1 ± 0.5	–
Flexural strength (biaxial) [MPa]	≥ 360	530*
Chemical solubility [µg/cm ²]	< 100	–
Type/Class	Type II / Class 3	–

According to ISO 6872:2015

* mean biaxial flexural strength evaluated over 10 years of quality measurements

Indications

- Partial edentulism in the anterior and posterior region

Contraindications

- Bruxism
- If the patient is known to be allergic to any of the material's constituents

Limitations of use

- Failure to observe the requirements stipulated by the implant manufacturer for the selected implant type (diameter and length of the implant must be approved for the respective position in the jaw by the implant manufacturer)
- Exceeding or falling short of the acceptable ceramic layer thicknesses
- Use of a luting composite other than Multilink Hybrid Abutment to adhesively cement IPS e.max CAD to the titanium bonding base
- Intraoral cementation of the ceramic structures to the titanium bonding base
- Temporary cementation of the crown on the hybrid abutment
-  Do not re-use

Processing restrictions

Failure to observe the following information will compromise the results achieved with IPS e.max CAD:

- Milling of the blocks in a non-compatible CAD/CAM system
- When producing a hybrid abutment crown, the opening of the screw channel must not be located in the area of the contact points. If this is not possible, it is better to produce a hybrid abutment with a separate crown.
- Cut-back technique in the fabrication of hybrid abutment crowns
- The use of IPS e.max CAD Crystall./Glaze Spray
- Combination with materials other than IPS e.max Ceram, IPS Ivocolor® or IPS e.max CAD Crystall./.
- Crystallization in a non-approved and non-calibrated ceramic furnace
- Crystallization in a high-temperature furnace
- Crystallization using deviating firing parameters
- Failure to observe the manufacturer's directions regarding the processing of the titanium bonding base.

System requirements

IPS e.max CAD Abutment Solutions must be processed with an authorized CAD/CAM system.¹ Select a suitable titanium bonding base in accordance with the implant inserted and the CAD/CAM system used. Please observe the instructions for use and the processing instructions of the respective manufacturer.

Side effects

There are no known side effects to date.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

Lithium disilicate glass-ceramic

After the manufacturing process of the glass-ceramic, a stable and inert network is formed in which the different elements are incorporated via oxygen bridges. The composition is determined as oxides.

Oxide	in wt.%
SiO ₂	57.0 – 80.0
Li ₂ O	11.0 – 19.0
K ₂ O	0.0 – 13.0
P ₂ O ₅	0.0 – 11.0
ZrO ₂	0.0 – 8.0
ZnO	0.0 – 8.0
Al ₂ O ₃	0.0 – 5.0
MgO	0.0 – 5.0
Coloring oxides (CeO ₂ , MnO ₂ , V ₂ O ₅ , Tb ₄ O ₇ , Er ₂ O ₃)	0.0 – 8.0

2 Application

Shade selection

Clean the teeth prior to shade determination. The shade is determined by the shade of the neighbouring teeth.


Minimum layer thicknesses of the ceramic structures

– Hybrid abutment:

- The wall thickness W_{HA} must be at least 0.5 mm.
- The hybrid abutment should be designed in a similar way as a prepared natural tooth:
 - Circular epi-/supragingival shoulder with rounded internal angles or a chamfer
 - In order for the crown to be cemented to the hybrid abutment using conventional or self-adhesive cementation methods, retentive surfaces and an adequate "preparation height" must be created.
 - Create an emergence profile with a right angle at the transition to the crown (see picture).
- The crown width B_{Crown} is limited to 6.0 mm from the axial height of contour to the screw channel of the hybrid abutment.
- The implant manufacturer's instructions regarding the maximum height of the hybrid abutment and separate crown must be observed.


– Hybrid abutment crown:

- The wall thickness W_{HAC} must be larger than 1.5 mm for the entire equatorial circumference.
- The opening of the screw channel must not be located in the area of contact points. If this is not possible, it is better to produce a hybrid abutment with a separate crown.
- The width of the hybrid abutment crown B_{HAC} is limited to 6.0 mm from the axial height of contour to the screw channel.
- The implant manufacturer's instructions regarding the maximum height of the hybrid abutment crown must be observed.

Block selection

The block is selected in accordance with the desired tooth shade and the selected titanium bonding base. An IPS e.max CAD MO or LT block is selected depending on the application. The range of blocks available may vary depending on the CAD/CAM machine used.

	 IPS e.max CAD MO A14 (Medium Opacity)	 IPS e.max CAD LT A14 (Low Translucency)	 IPS e.max CAD LT A16 (Low Translucency)
 IPS e.max CAD Hybrid Abutment	✓	✓	-
 IPS e.max CAD Hybrid Abutment Crown	-	✓	✓

Finishing

Observe the grinding instrument recommendations² and the minimum layer thicknesses when finishing and contouring the ceramic structures. If possible, adjustments by grinding should be carried out while the restoration is still in the pre-crystallized (blue) state, at low speed and exerting only light pressure to prevent delamination and chipping at the margins. Overheating of the ceramic must be avoided. Cut the ceramic structure from the block using a diamond separating disc. Hybrid abutment: using a separating disc, cut slightly into the attachment area at the incisal side of the abutment, then cut through the attachment point fully from the basal aspect. Carefully place the ceramic structure on the titanium bonding base and check the fit. Observe the position of the anti-rotation lock.

– Finishing the outer surface of the ceramic structure

Do not grind the shoulder of the ceramic structure so as not to affect the accurate fit on the titanium bonding base. If necessary, finish the emergence profile taking the fit to the gingiva and the minimum thickness (0.5 mm) into account.

– Hybrid abutment

Smooth out the attachment point to the block with fine diamond grinding instruments taking the shape of the emergence profile and the crown margin into account. Do not perform any individual shape adjustments, as this will negatively affect the fit of the crown on the hybrid abutment. Information regarding the crown: If the fit on the hybrid abutment is inaccurate, carry out adjustments on the crown.

– Hybrid abutment crown

Smooth out the attachment point to the block with fine diamond grinding instruments taking the shape of the emergence profile and the proximal contacts into account. Gently grind the entire occlusal surface with a fine diamond to smooth out the surface structure created by the CAD/CAM process. Check the proximal and occlusal contacts. Create the surface textures.

Always clean the ceramic structures with ultrasound in a water bath or with the steam jet before further processing. Make sure to thoroughly remove any milling additive residue of the CAD/CAM milling unit. Residue of the milling additive remaining on the surface may result in bonding problems and discolouration. Do not blast the ceramic structure with Al₂O₃ or glass polishing beads.

Optional: Clinical try-in of the blue restorations

A clinical try-in to check the accuracy of fit can be conducted prior to further processing. Clinical try-in may also take place at a later stage, i.e. with the crystallized, tooth-coloured IPS e.max CAD ceramic structure.

– Provisional fixation of the ceramic structure on the titanium bonding base

To facilitate the intraoral handling and prevent damage to the pre-crystallized ceramic structure, the components must be temporarily attached to one another using silicone impression material, e.g. Virtual[®] Extra Light Body Fast Set.

The untreated titanium bonding base and the ceramic structure are steam-cleaned and then dried with compressed air. The ceramic structure is placed on the titanium bonding base (which is screwed onto the model analog) and the relative position of the components is marked with a waterproof pen. This step makes it easier to attain the correct position when the parts are temporarily assembled. The screw channel of the untreated titanium bonding base is sealed with a foam pellet. Virtual Extra Light Body Fast Set is applied to the titanium bonding base and directly into the ceramic structure. The titanium bonding base is inserted into the ceramic structure. The alignment of the two components must be checked (anti-rotation lock/markings). The components must be held firmly in the correct position for 2:30 min. until Virtual Extra Light Body Fast Set has set. Any excess that has been displaced must be carefully removed with a suitable instrument, e.g. a scalpel.

Clinical try-in

Once the provisional restoration has been removed, the hybrid abutment or the hybrid abutment crown is screwed in manually using the dedicated screw. The geometry is checked in relation to the gingival margin (e.g. fit, gingival anaemia). If desired, the screw channel on the hybrid abutment can be sealed with a foam pellet. **Tip:** The inner aspect of the crown is isolated using glycerine gel (e.g. Try-in paste, Liquid Strip). The crown is placed on the hybrid abutment intraorally to check and, if necessary, adjust the proximal contacts.

Attention: Do not check the occlusal function at this stage. In order to check the occlusal function, the crown must be secured on the hybrid abutment with Virtual Extra Light Body Fast Set. Try-in paste must not be used for this purpose, as this material is not sufficiently resistant to compressive force. Virtual Extra Light Body Fast Set is applied to the inner aspect of the crown. The crown is pressed onto the hybrid abutment using the fingers until the final position is reached. The crown is held in the final position until Virtual Extra Light Body has set (approx. 2:30 min). Subsequently, excess material is removed.

Now the occlusion/articulation is checked. If necessary, adjustments are made using suitable grinding instruments (see the grinding instrument recommendations²). The crown is carefully removed from the hybrid abutment and the hybrid abutment/hybrid abutment crown from the implant. The implant site is rinsed with e.g. Cervitec Liquid (alcohol-free mouth rinse containing chlorhexidine) in order to clean it. Next, the temporary restoration is placed.

Completion of the ceramic structure

Depending on the desired processing technique and materials, the processing method to complete the ceramic structure is selected. Basically, there are three processing methods for completing the ceramic structure.

– Polishing technique (self glaze) on the blue restoration

The polishing technique is preferably used for the emergence profile of the hybrid abutment. For the hybrid abutment crown, a glaze firing is recommended.

For easier handling, screw the titanium bonding base onto a model analog. Secure the ceramic structure on the titanium bonding base with your fingers. **Attention:** Do not grind the titanium bonding base.

For polishing, observe the grinding instrument recommendations². Clean the restoration with ultrasound in a water bath or with the steam jet. Subsequently, attach the restoration to the crystallization pin as described under "Fixation of restorations on the IPS e.max CAD Crystallization Pin". Place the ceramic structure on the IPS e.max CAD Crystallization Tray and position the tray in the centre of the furnace. The firing program is selected based on the material (see "Crystallization and firing parameters").

– Staining technique on the blue restoration

The IPS e.max CAD Crystall./Glaze Spray is not recommended for glazing IPS e.max CAD Abutment Solutions, since it requires very targeted application. The glazing material must not touch the bonding surface to the titanium bonding base or the screw channel, as this may compromise the accuracy of fit.

Position the ceramic structure as described under "Fixation of restorations on the IPS e.max CAD Crystallization Pin". Subsequently, mix IPS e.max CAD Crystall./Glaze Paste/FLUO and apply it evenly on the areas to be glazed using a small brush. If the ready-to-use glaze needs to be diluted, it may be mixed with a small amount of IPS e.max CAD Crystall./Glaze Liquid. The glazing material must not touch the bonding surface to the titanium bonding base nor the screw channel, as this may compromise the accuracy of fit. Check the inner surface before firing

and carefully remove any contamination. On the hybrid abutment, do not apply any materials to the bonding surface to the crown, as this might compromise the fit of the crown. Do not apply the glaze too thickly. Avoid "pooling", especially on the occlusal surface of the hybrid abutment crown. Too thin a glaze layer may lead to an unsatisfactory gloss. If characterizations are desired, the ceramic structure can be customized using IPS e.max CAD Crystall./Shades and/or Stains before crystallization firing. Extrude the ready-to-use shades and stains from the syringe and mix them thoroughly. The shades and stains may be slightly thinned using IPS e.max CAD Crystall./Glaze Liquid. However, the consistency should still remain pasty. Apply the mixed shades and stains directly to the unfired glaze layer using a fine brush. More intensive shades are achieved by several staining procedures and repeated firing, not by applying thicker layers. Use the IPS e.max CAD Crystall./Shades Incisal to imitate the incisal area and create the effect of translucency on the hybrid abutment crown in the incisal and occlusal third. The cusps and fissures can be individualized using stains.

Then place the restoration in the centre of the IPS e.max CAD Crystallization Tray or place a maximum of 6 units on the tray and conduct the crystallization firing using the stipulated firing parameters (see Crystallization and firing parameters). Observe the directions under "How to proceed after firing".

Optional: Corrective firing

If additional characterizations or adjustments are required after crystallization, a corrective firing using IPS e.max CAD Crystall./Shades and Stains and Glaze can be conducted. Also use the IPS e.max CAD Crystallization Tray for the corrective firing cycle. Use the IPS e.max CAD Crystall./Add-On including the respective mixing liquid to carry out minor shape adjustments (e.g. proximal contact points).

– Staining technique on the tooth-coloured restoration

- Crystallization without the application of materials; separate Stains/Glaze firing either using IPS e.max CAD Crystall./ or IPS Ivocolor materials.
- Position the ceramic structure on the crystallization pin as described under "Fixation of ceramic structures on the IPS e.max CAD Crystallization Pin". For characterization and glazing, use either **IPS e.max CAD Crystall./Shades, Stains, Glaze** (see "Staining technique on the blue restoration", with a different firing program) or IPS Ivocolor.
- IPS e.max CAD Crystall./Shades, Stains, Glaze and IPS Ivocolor Shades, Essence, Glaze must not mixed with one another or applied one after the other.

When using IPS Ivocolor: To ensure better wetting, a small quantity of IPS Ivocolor Mixing Liquid may be slightly rubbed into the area that needs to be characterized. Mix the IPS Ivocolor Shades and Essences to the desired consistency using the respective IPS Ivocolor Liquids. More intensive shades are achieved by repeating the staining procedure and firing, not by applying thicker layers. Use the IPS Ivocolor Shades Incisal to imitate the incisal area and create the effect of translucency on the hybrid abutment crown in the incisal and occlusal third. The cusps and fissures can be individualized using Essences. On the hybrid abutment, only characterize the area of the emergence profile with IPS Ivocolor Shades and Essences. Under no circumstances must stains be applied to the bonding surface to the titanium bonding base and the screw channel, as this may compromise the accuracy of fit. Check the inner surface before firing and carefully remove any contamination. On the hybrid abutment, do not apply any materials to the bonding surface to the crown, as this may compromise the fit of the crown. Then secure the ceramic structure on the firing pin of the honey-comb tray with a little IPS Object Fix Putty or Flow and crystallize it using the stipulated firing parameters. Observe the directions under "How to proceed after firing".

Glaze firing is conducted with glaze paste or powder. On hybrid abutments, only the emergence profile is glazed. On hybrid abutment crowns, glaze is applied to the entire outer surface. For easier handling, the ceramic structure can be positioned on the titanium bonding base for glazing. For that purpose, secure the titanium bonding base on a model analog. Mix the glazing material (IPS Ivocolor Glaze Paste or Powder) with the IPS Ivocolor Mixing Liquid allround or longlife to achieve the desired consistency. Apply an even layer of glazing material covering all areas that are to be glazed.

The degree of gloss of the glazed surface is controlled via the consistency of the glazing material and the applied quantity, not by means of the firing temperature. For a higher degree of gloss, the glazing material must be applied in a correspondingly thicker layer. If required, the fluorescence may be increased by applying a fluorescing glazing material (Paste FLUO or Powder FLUO).

Under no circumstances must glazing material be applied to the bonding surface to the titanium bonding base or to the screw channel, as this may compromise the accuracy of fit. Check the inner surface before firing and carefully remove any contamination. On the hybrid abutment, do not apply any materials to the bonding surface to the crown, as this may compromise the fit of the crown.

Conduct the Stain/Glaze firing for IPS Ivocolor on a honey-comb tray using the stipulated firing parameters. Observe the directions under "How to proceed after firing".

Optional: Corrective firing

Mix IPS e.max Ceram Add-On Dentin or Incisal with IPS Build-Up Liquid soft or allround and apply on the corresponding areas. Fire with the stipulated parameters for the "Add-On after Glaze firing". Observe long-term cooling! If necessary, polish the adjusted areas to a high gloss after firing.

Completing the crown on the IPS e.max CAD hybrid abutment

The crown on the IPS e.max hybrid abutment can be completed using either the staining technique or the cut-back technique. For characterization and glazing, either the IPS e.max CAD Crystall./ materials or the IPS Ivocolor materials are used. The reduced areas (cut-back technique) are built up using IPS e.max Ceram layering materials. Use IPS Ivocolor materials for staining and glazing. Basically, the procedure for completing a crown is the same as that for a crown on a prepared tooth. For detailed information about the procedure, please refer to the IPS e.max CAD Instructions for Use.

Fixing the restoration on an IPS e.max CAD Crystallization Pin

1. Clean the ceramic structure with the steam jet to remove any contamination and grease residue. Avoid any contamination after cleaning.
2. Use the IPS e.max CAD Crystallization Pin XS for the crystallization of the ceramic structure.
3. Fill the inner surfaces of the ceramic structure with either IPS Object Fix Putty or Flow auxiliary firing paste. Immediately re-seal the IPS Object Fix Putty/Flow syringe after extruding the material. Once removed from the aluminium bag, the syringe is ideally stored in a resealable plastic bag or a container with a humid atmosphere.
4. Press the IPS e.max CAD Crystallization Pin XS only slightly into the IPS Object Fix Putty/Flow. **Important:** Do not press the pin in too deeply to make sure it does not touch the walls. This may lead to cracks in the ceramic structure.
5. Smooth out displaced auxiliary firing paste using a plastic spatula so that the pin is securely in place.
6. Prevent contamination of the outer surface/occlusal surface of the ceramic structure. Clean off any possible contamination with a brush dampened with water and dry.

Important: For crystallization, IPS e.max CAD restorations must not be placed directly on the IPS e.max CAD Crystallization Tray and Pins, i.e. without auxiliary firing paste.

How to proceed after firing

Remove the ceramic structure from the furnace once the firing cycle is complete (wait for the acoustic signal of the furnace) and allow it to cool to room temperature in a place protected from draft. Hot objects must not be touched with metal tongs. Remove the ceramic structure from the hardened IPS Object Fix Putty / Flow. Remove any residue with ultrasound in a water bath or with the steam jet. Residue must not be removed by blasting with Al₂O₃ or glass polishing beads. If the restoration needs to be adjusted by grinding², make sure that no overheating of the ceramic occurs. Finally, polish the ground areas to a high-gloss finish.

Permanently bonding the ceramic structure to the titanium bonding base

The contact surfaces must be meticulously prepared in order to ensure an optimum adhesive bond between the titanium bonding base and the ceramic structure.

	IPS e.max CAD ceramic structure (LS ₂)		Titanium bonding base
Blasting	–		Observe the manufacturer's directions.
Conditioning	Option 1	Option 2	
Etching	Etching the bonding surface to the titanium bonding base for 20 sec with IPS® Ceramic Etching Gel	Apply Monobond Etch & Prime® to the bonding surface to the titanium bonding base and allow to react for another 40 sec.	–
Silanizing	The bonding surface is conditioned with Monobond® Plus for 60 sec		The bonding surface is conditioned with Monobond® Plus for 60 sec
Adhesive cementation	Multilink® Hybrid Abutment		
Covering of the cement joint	Glycerine gel, e.g. Liquid Strip		
Curing	7 min auto-polymerization		
Polishing of the cement joint	Customary polishers for ceramic/resin materials		

– **Preparing the titanium bonding base**

- The titanium bonding base should be prepared according to the instructions of the manufacturer.
- Clean the titanium bonding base in an ultrasonic bath or with a steam cleaner and then dry it with blown air.
- Screw the titanium bonding base onto a model analog.
- Position the ceramic structure on the titanium bonding base and mark the relative position with a waterproof pen. This makes it easier to attain the correct relative end position when the parts are assembled at a later stage.
- After the bonding surface has been cleaned, it must not be contaminated under any circumstances as this would impair the bond.
- Apply Monobond Plus on the cleaned bonding surface and allow it to react for 60 seconds. After the reaction time, dry the remaining residue with water- and oil-free air.
- Seal the screw channel with a foam pellet or wax. Avoid any contamination of the bonding surface.

– **Preparing the ceramic structure**

- The ceramic structure must not be blasted.
- Clean the ceramic structure in an ultrasonic bath or with a steam cleaner and subsequently blow it dry with water- and oil-free air.
- Wax can be applied to protect the outer surfaces or the glazed areas.

For the conditioning of the ceramic structure, two options are available:

- **Option 1:** Conditioning of the bonding surfaces with IPS Ceramic Etching Gel and Monobond Plus
- **Option 2:** Conditioning of the bonding surfaces with Monobond Etch & Prime

Procedure when choosing option 1:

- Etch the bonding surface with 5% hydrofluoric acid gel (IPS Ceramic Etching Gel) for 20 seconds.
- Thoroughly rinse the bonding surface with running water and dry with water- and oil-free air.
- Apply Monobond Plus on the cleaned bonding surface and allow it to react for 60 seconds. After the reaction time, dry the remaining residue with water- and oil-free air.

Procedure when choosing option 2:

- Apply Monobond Etch & Prime to the bonding surface using a microbrush, scrub it in for 20 seconds and leave it to react for another 40 seconds.
- Then thoroughly rinse off Monobond Etch & Prime with water and dry the restorations with a strong stream of water- and oil-free air for about 10 seconds.

– **Cementing with Multilink Hybrid Abutment**

- Apply a thin layer of Multilink Hybrid Abutment directly from the mixing syringe to the bonding surface of the titanium bonding base and the bonding surface of the ceramic structure.
- Place the structure over the base in such a way that the position markings are aligned.
- Press the parts lightly and evenly together and check the correct relative position of the components (transition between the base and ceramic structure).
- Subsequently, firmly press the components together for 5 seconds.
- Carefully remove excess in the screw channel, e.g. with a microbrush or brush, using a rotating movement.
- **Note:** Do not remove circular excess resin cement before curing has started, i.e. 3 minutes after mixing. Use a suitable instrument for this purpose (e.g. Le Cron). Hold the components in place using light pressure.

- Apply glycerine gel (e.g. Liquid Strip) to the cement line to prevent the formation of an inhibition layer.
- Next, the composite cement is left to auto-polymerize completely within 7 min.
- **Important:** Do not move the components until Multilink Hybrid Abutment has completely cured. Hold them in place using e.g. diamond-coated tweezers to prevent movement.
- Cautiously polish the cement line with rubber polishers at low speed (< 5,000 rpm) to avoid overheating.
- Remove any remaining cement residue in the screw channel with suitable rotating instruments.

Insertion and aftercare

- Intraoral preparation

Remove the temporary restoration and clean the implant site. Then check the peri-implant tissue (emergence profile).

- Insertion

Do not use phenolic mouth rinses, as such products detrimentally affect the bond between the ceramic and the composite. Insert the hybrid abutment or the hybrid abutment crown intraorally into the implant. Manually screw in the matching implant screw and tighten it with a torque wrench (observe the manufacturer's directions).

- Hybrid abutment and separate crown

Insert a cotton or foam pellet into the screw channel of the hybrid abutment and seal the screw channel with a temporary composite (e.g. Telio® Inlay). This serves to ensure access to the screw at a later stage. Check the bonding surfaces for contamination/moisture and clean or dry with an air syringe, if necessary. Apply the luting material (e.g. SpeedCEM® Plus or Variolink® Esthetic) into the conditioned crown. Place the crown onto the hybrid abutment and secure in place in the final position.

Any cement excess is briefly light-cured using the "quarter technique" and can then be easily removed. Cover the cement joints with glycerine gel (e.g. Liquid Strip) and light-cure again with e.g. the Bluephase® curing light. Subsequently, rinse off the glycerine gel with water.

- Hybrid abutment crown

Insert a cotton or foam pellet into the screw channel of the hybrid abutment crown and seal the screw channel with a permanent composite (e.g. Tetric® Prime).

Check the occlusion and articulation. If the restoration is adjusted by grinding, the ground areas must subsequently be polished to a high gloss using silicone polishers (e.g. OpraGloss®). Also polish the restorations margins / cement joint. Finally, apply Cervitec® Plus (protective varnish) along the gingival margin.

Shade combination table

For the characterization and adjustment of the shade of IPS e.max CAD restorations, IPS e.max CAD Crystall./Shades, Stains or IPS Ivocolor Shades, Essences are used.

- IPS e.max CAD Crystall./Shades, Stains: For use on blue and tooth-coloured IPS e.max CAD restorations
- IPS Ivocolor Shades, Essences: For use on tooth-coloured IPS e.max CAD restorations

The combination table must be observed.

	BL1	BL2	BL3	BL4	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
IPS e.max CAD Crystall./ Shade	0				1					2				3				4		
IPS e.max CAD Crystall./ Shade Incisal	I1						I2			I1				I2						
IPS e.max CAD Crystall./ Stains	white, cream, sunset, copper, olive, khaki, mahogany																			

	BL1	BL2	BL3	BL4	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
IPS Ivocolor Shade	SD 0				SD 1	SD 2			SD 3	SD 4	SD 5			SD 6				SD 7	SD 6	
IPS Ivocolor Shade Incisal	SI 1						SI 2				SI 3									
IPS Ivocolor Essence	E 01	E 02	E 03	E 04	E 05	E 06	E 07	E 08	E 09	E 10										
	white	creme	lemon	sunset	copper	hazel	olive	khaki	terracotta	mahogany										
	E 11	E 12	E 13	E 14	E 15	E 16	E 17	E 18	E 19	E 20										
	cappuccino	espresso	terra	profundo	ocean	sapphire	anthracite	black	rose	coral										
			E 21		E 22			E 22												
			basic red		basic yellow			basic blue												

Crystallization and firing parameters

Ceramic furnaces without controlled (long-term) cooling function cannot be used. The ceramic furnace must be calibrated before the first crystallization and regularly every six months thereafter. Depending on the mode of operation, more frequent calibration may be required. Observe the manufacturer's directions.

Crystallization MO, Impulse, LT, MT, HT

with or without the application of IPS e.max CAD Crystall./ materials

Furnaces Programat	Stand-by temperature	Closing time	Heating rate	Firing temperature	Holding time	Heating rate	Firing temperature	Holding time	Vacuum 1	Vacuum 2	Long-term cooling	Cooling rate
	B [°C/°F]	S [min]	t1 [°C/°F/min]	T1 [°C/°F]	H1 [min]	t2 [°C/°F/min]	T2 [°C/°F]	H2 [min]	11 [°C/°F] 12 [°C/°F]	21 [°C/°F] 22 [°C/°F]	L [°C/°F]	tl [°C/°F/min]
P300 P500 P700	403/757	6:00	60/108	770/1418	0:10	30/54	850/1562	10:00	550/770 1022/1418	770/850 1418/1562	700/1292	0
P310 P510 P710	403/757	6:00	60/108	780/1436	0:10	30/54	860/1580	10:00	550/780 1022/1436	780/860 1436/1580	710/1310	0
CS/CS2/ CS3/CS4	Select the corresponding program											

Crystallization LT, MT, HT

with or without the application of IPS e.max CAD Crystall./ materials

Furnaces Programat	Stand-by temperature	Closing time	Heating rate	Firing temperature	Holding time	Heating rate	Firing temperature	Holding time	Vacuum 1	Vacuum 2	Long-term cooling	Cooling rate
	B [°C/°F]	S [min]	t1 [°C/°F/min]	T1 [°C/°F]	H1 [min]	t2 [°C/°F/min]	T2 [°C/°F]	H2 [min]	11 [°C/°F] 12 [°C/°F]	21 [°C/°F] 22 [°C/°F]	L [°C/°F]	tl [°C/°F/min]
P300 P500 P700	403/757	6:00	90/162	820/1508	0:10	30/54	840/1544	7:00	550/820 1022/1508	820/840 1508/1544	700/1292	0
P310 P510 P710	403/757	6:00	90/162	830/1526	0:10	30/54	850/1562	7:00	550/830 1022/1526	830/850 1526/1562	710/1310	0
CS/CS2/ CS3/CS4	Select the corresponding program											

Corrective firing/stain firing/glaze firing

with IPS e.max CAD Crystall./ materials

Furnaces Programat	Stand-by temperature	Closing time	Heating rate	Firing temperature	Holding time	Heating rate	Firing temperature	Holding time	Vacuum 1	Vacuum 2	Long-term cooling	Cooling rate
	B [°C/°F]	S [min]	t1 [°C/°F/min]	T1 [°C/°F]	H1 [min]	t2 [°C/°F/min]	T2 [°C/°F]	H2 [min]	11 [°C/°F] 12 [°C/°F]	21 [°C/°F] 22 [°C/°F]	L [°C/°F]	tl [°C/°F/min]
P300 P500 P700	403/757	6:00	90/162	820/1508	0:10	30/54	840/1544	3:00	550/820 1022/1508	820/840 1508/1544	700/1292	0
P310 P510 P710	403/757	6:00	90/162	830/1526	0:10	30/54	850/1562	3:00	550/830 1022/1526	830/850 1526/1562	710/1310	0
CS/CS2/ CS3/CS4	Select the corresponding program											

Firing parameters for the staining technique

with IPS Ivocolor Shade, Essence, Glaze

	Stand-by temperature	Closing time *	Heating rate	Firing temperature	Holding time	Vacuum 1	Vacuum 2	Long-term cooling **	Cooling rate
	B [°C/°F]	S [min]	t [°C/°F/min]	T [°C/°F]	H [min]	V1 [°C/°F]	V2 [°C/°F]	L [°C/°F]	tl [°C/°F/min]
Stain and Glaze firing	403/757	IRT/ 6:00	60/108	710/1310	1:00	450/842	709/1308	0	0

Firing parameters for the corrective firing (staining technique)

with IPS e.max Ceram Add-On

	Stand-by temperature	Closing time *	Heating rate	Firing temperature	Holding time	Vacuum 1	Vacuum 2	Long-term cooling **	Cooling rate
	B [°C/°F]	S [min]	t [°C/°F/min]	T [°C/°F]	H [min]	V1 [°C/°F]	V2 [°C/°F]	L [°C/°F]	tl [°C/°F/min]
Add-On after Glaze firing	403/757	IRT/ 6:00	50/90	700/1292	1:00	450/842	699/1290	0	0

* IRT standard mode

** Note: If the layer thicknesses exceed 2 mm, long-term cooling L to 500 °C/932 °F is required.

Note: Due to their geometry, the restorations may feature varying layer thicknesses. When the objects cool after the firing cycle, the different cooling speeds in the areas with different thicknesses may result in a build-up of internal tension. In the worst case, these internal tensions may result in fractures in the ceramic objects. By using slow cooling (long-term cooling L), these tensions can be minimized.

3 Safety information

Sterilization

The hybrid abutments or hybrid abutment crowns must be sterilized prior to insertion in the mouth. The local statutory provisions and hygiene standards that apply for dental practices must be observed.

Steam sterilization can be performed under 3x-fractionated pre-vacuum and observing the following parameters: Sterilization time 3 min; steam temperature 132 °C/270 °F; this correspond to a half-cycle exposure time of 2 min. The hybrid abutment or the hybrid abutment crown is to be used immediately. No storage after sterilization!

The user is responsible for the sterility of the hybrid abutment or hybrid abutment crown. It must be ensured that sterilization is only performed using suitable devices and materials as well as product-specific, validated methods. The devices used must be properly maintained and regularly serviced. It is the duty of the IPS e.max CAD Abutment Solutions users to inform their dentists that the ceramic structures require sterilization before being seated in the patient's oral cavity.

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclarvivadent.com, and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclarvivadent.com).
- The current Summary of Safety and Clinical Performance (SSPC) can be found in the download section of the Ivoclar Vivadent AG website (www.ivoclarvivadent.com).

Warnings

- IPS Ceramic Etching Gel contains hydrofluoric acid. Contact with skin, eyes and clothing must be prevented at all costs, since the material is extremely toxic and corrosive. The etching gel is intended for extraoral use only and must not be applied intraorally (inside the mouth).
- Monobond Etch & Prime is corrosive. Contact with skin and mucous membranes must be prevented. Monobond Etch & Prime is intended for extraoral use only and must not be applied intraorally (inside the mouth).
- Do not inhale ceramic dust during finishing. Use an extraction unit and wear a face mask.
- Risks associated with magnetic fields (e.g. MRI – Magnetic Resonance Imaging): Please note warnings, cautions and precautions of the abutment or implant manufacturer.
- Observe the Safety Data Sheet (SDS).

Disposal information

- Remaining stocks and extracted restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks. Some of these risks are listed below:

- Chipping / fracture
- Cement excess may lead to irritation of the soft tissue / gingiva.
- Decementation

4 Shelf life and storage

This product does not require any special storage conditions.

5 Additional information

Keep material out of the reach of children!

Not all products are available in all countries.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions.

¹ e.g. CEREC/inLab, PlanMill. The complete list is available at www.ivoclarvivadent.com.

CEREC/inLab and PlanMill are not registered trademarks of Ivoclar Vivadent AG.

² Ivoclar Vivadent flow chart "Recommended grinding instruments for extraoral and intraoral use".

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.