

Hybrid abutment**Hybrid abutment crown****Description**

Multilink® Hybrid Abutment is a self-curing luting composite designed for the permanent cementation of structures made of lithium disilicate glass-ceramic, zirconium oxide or PMMA on bonding bases made of titanium/titanium alloy for the fabrication of hybrid abutments or hybrid abutment crowns.

Shades

- HO 0 (High Opacity)
- MO 0 (Medium Opacity)

Working time

The working and curing times depend on the ambient temperature. Once Multilink Hybrid Abutment has been dispensed from the automix syringe, the following times apply:

	At room temperature 23 °C ± 1 °C (73 °F ± 2 °F)
Working time	approx. 2 minutes
Curing time (including working time)	approx. 7 minutes

The working and curing times become shorter at higher temperatures, while they become longer at lower temperatures.

Mixing ratio

Multilink Hybrid Abutment is extruded from the automix syringe in an optimum 1:1 mixing ratio.

Composition

The monomer matrix is composed of dimethacrylate and HEMA. The inorganic fillers are composed of barium glass, ytterbium trifluoride, spheroid <mixed oxide and titanium oxide. The particle size ranges between 0.15 and 3.0 µm. The average particle size is 0.9 µm. The total volume content of inorganic fillers is approx. 36 %.

Indication

Extraoral, permanent cementation of structures made of lithium disilicate glass ceramic, zirconium oxide or PMMA on bonding bases made of titanium/titanium alloy.

Multilink® Hybrid Abutment

EN Instructions for Use

- Self-curing dental luting composite for extraoral cementation

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Contraindication

The product is contraindicated

- for intraoral luting;
- if the stipulated working technique cannot be applied;
- if the patient is known to be allergic to components of Multilink Hybrid Abutment.

Side effects

Systemic side effects are not known. In rare cases allergic reactions to individual components have been reported.

Interactions

Phenolic substances (e.g. wintergreen oil, eugenol) inhibit polymerization. Therefore, materials containing these components must not be used. Oxidative disinfectants (e.g. hydrogen peroxide) may interact with the initiator system and impair the setting process. As a result, the automix syringe must not be disinfected with oxidative agents. The syringe can be disinfected by wiping it with medical alcohol for example.

Application

1. Preparation of the bonding base (e.g. Viteo® Base Ti)

- The bonding base should be prepared according to the instructions of the manufacturer.
- The bonding base is cleaned in an ultrasonic bath or with a steam cleaner and then dried with blown air.
- The bonding base is screwed to the model analog.
- The structure is placed on the bonding base and the relative position is marked with a waterproof pen. This makes it easier to attain the correct relative end position when the parts are assembled at a later stage.
- The emergence profile of the bonding base must not be blasted or modified in any way!
- **If the manufacturer recommends that the contact surfaces be sandblasted, the following procedure should be observed:**
 - Hard modelling wax is applied to protect the emergence profile, since this type of material is easy to remove later on.
 - The screw channel is also sealed with wax.
 - The bonding area surface is carefully blasted according to the instructions of the manufacturer.
 - An instrument and a steam jet are used for cleaning. The wax has to be removed meticulously.
- The bonding base is cleaned in an ultrasonic bath or with a steam cleaner and then dried with blown air.
- After cleaning, the bonding surface must not be contaminated under any circumstances, since this may negatively influence the bonding results.
- Monobond Plus is applied to the clean bonding surface and allowed to react for 60 s. After this reaction time, any residue is dried with air that is free of water and oil.
- The screw channel is sealed with a foam pellet or wax. The bonding surface must not be contaminated in the process.
- If using a Viteo Base Ti, insert the Viteo Screw Channel Pin into the screw channel of the titanium base. The Channel Pin can be shortened with a scalpel.

2. Preparation of the structures

2.1 made of lithium disilicate glass ceramics (e.g. IPS e.max® CAD/Press)

- The ceramic structure must not be blasted when it is prepared for cementation.
- The ceramic structure is cleaned in an ultrasonic bath and with a steam cleaner and then dried with blown air.
- After cleaning, the bonding surface must not be contaminated under any circumstances, since this may negatively influence the bonding results.
- Wax can be applied to protect the outer surfaces and the glazed areas.
- There are two options to prepare the bonding surface:
 - a) The bonding surface is etched with 5% hydrofluoric acid (IPS Ceramic Etching Gel) for 20 s. Subsequently, the bonding surface is thoroughly rinsed under running water and then dried with oil-free air. Monobond Plus is applied to the clean bonding surface and allowed to react for 60 s. After this reaction time, any residue is dried with air that is free of water and oil.
 - b) Apply Monobond Etch & Prime® with a microbrush to the bonding surface and agitate it into the surface for 20 s. Allow it to react for another 40 s. Then thoroughly rinse off Monobond Etch & Prime with water and dry the restoration with a strong steam of water- and oil free air for approximately 10 seconds.

2.2 made of zirconium oxide (ZrO₂)

- In order to facilitate its handling, the ceramic structure is attached to a holder (e.g. Microbrush) before it is blasted.
- The outer surfaces of the ceramic structure are protected with a wax coating.
- In order to better control the blasting procedure, the bonding surface can be coloured with a permanent marker.
- Carefully blast the bonding surface of the ceramic structure according to the instructions of the manufacturer.
- The ceramic structure is cleaned with a steam cleaner or in an ultrasonic bath.
- After cleaning, the bonding surface must not be contaminated under any circumstances, since this may negatively influence the bonding results.
- Monobond Plus is applied to the clean bonding surface and allowed to react for 60 s. After this reaction time, any residue is dried with air that is free of water and oil.

2.3 made of PMMA (e.g. Telio® CAD A16)

- The Telio CAD structure must not be blasted when it is prepared for cementation.
- The Telio CAD structure is cleaned in an ultrasonic bath or with a steam cleaner and then dried with blown air.
- After cleaning, the bonding surface must not be contaminated under any circumstances, since this may negatively influence the bonding results.
- Apply a thin layer of SR Connect onto the bonding surface using a disposable brush and agitate for 30 s. Subsequently polymerize for 40 s using a curing device (e.g. Bluephase Style).
- **The reaction time of 30 s for SR Connect must be observed. Do not use SR Connect to condition Telio CAD A16 for longer since this may compromise the fit.**

- Do not disperse SR Connect!
- 3. Cementation with Multilink Hybrid Abutment**
- The cleaned and conditioned components (structure, bonding base) are laid out ready for cementation.
- **The subsequent cementation procedure must be carried out quickly and without interruption. The working time of Multilink Hybrid Abutment is approx. 2 min. at 23 °C/73 °F (± 1°C/2 °F).**
- As a general rule, a new mixing tip is attached to the Multilink Hybrid Abutment syringe prior to each use.
- A thin layer of Multilink Hybrid Abutment is directly applied from the mixing syringe to the bonding surface of the base and **the bonding surface of the ceramic or PMMA structure.**
- The mixing tip is left on the Multilink Hybrid Abutment syringe until the next use. The remaining cement polymerizes in the tip and functions as a seal.
- The structure is placed on the base in such a way that the position markings are aligned.
- The parts are lightly and evenly pressed together and the correct relative position of the components is checked (transition between the base and ceramic structure).
- Subsequently, the parts are tightly pressed together for 5 s.
- Excess in the screw channel is carefully removed with, e.g. a Microbrush or brush using rotary movements.
- **Important: Excess must not be removed before curing has started, i.e. 3 minutes after mixing. For the purpose, a suitable dental lab instrument (e.g. Le Cron) is used and the components are held in place with light pressure in the process.**
- Glycerine gel is applied (e.g. Liquid Strip) on the cementation joint to prevent the formation of an inhibition layer. The glycerine gel must be applied cautiously to avoid blending it with or displacing the composite. The gel must be left on the cementation joint until polymerization is complete.
- Next, the luting composite is auto-polymerized within 7 min.
- **Important: The parts should not be moved until Multilink Hybrid Abutment has completely cured. They can be held immobile with e.g. diamond-coated tweezers.**
- After the completion of auto-polymerization, the glycerine gel is rinsed off with water
- **The cementation joint should be cautiously polished with rubber polishers at a low speed (< 5,000 rpm) to avoid overheating.**
- Any cement residue left in the screw channel is removed with suitable rotating instruments.
- The restoration is steam cleaned.

Additional notes

- Multilink Hybrid Abutment should be processed at room temperature. Low(er) temperatures may cause difficulty in dispensing and mixing the material and may prolong the working and curing time.
- Unpolymerized material may stain clothing.

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Notes on disinfection and sterilization

Hybrid abutments and hybrid abutment crowns cemented with Multilink Hybrid Abutment can be disinfected and sterilized according to the instructions of use for the individual components.

Warnings

Contact of uncured Multilink Hybrid Abutment with skin/mucous membrane and eyes must be prevented. In an unpolymerized state, Multilink Hybrid Abutment may cause slight irritation and lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Shelf life and storage

- Multilink Hybrid Abutment has to be stored at 2 – 28 °C (36 – 82 °F)
- Do not use Multilink Hybrid Abutment after the indicated expiration date.
- Expiry date: See note on syringe/packaging.
- In order to ensure a tight seal of Multilink Hybrid Abutment automix syringes, the mixing tip is left on the syringe after use.


Keep out of the reach of children!

For use in dentistry only.


The material has been developed solely for use in dentistry. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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