UTS 3D®

EN Cleaning and Sterilizing Instructions

Transferbow and 3D-Registration Joint

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Página 42

Cleaning and Sterilizing Instructions

List of Parts

1.0 Universal transferbow (complete)

- 1.1 Earpiece mounting
- 1.2 Earpiece rod screw
- 1.3 Side arms
- 1.4 Side arm support
- 1.5 Scale
- 1.6 Multifunctional mounting
- 1.7 Width setting screw
- 1.8 Registration joint connection
- 1.9 Connection for nosepiece
- 1.10 Nosepiece holder screw
- 1.11 Registration joint connection screw
- 1.12 Thread for screw 3.1

2.0 Individually adjustable nosepiece (complete)

- 2.1 Nosepiece holder
- 2.2 Nosepiece rod
- 2.3 Nosepad holder
- 2.4 Nosepiece rod screw

3.0 Reference indicator

- 3.1 Reference indicator screw

4.0 Axis pins

- 4.1 Thread
- 4.2 Adapter pin for articulator

5.0 Intercondylar distance (S,M,L)

6.0 Earpieces (complete)

- 6.1 Earpiece rod
- 6.2 Threaded openings for axis pins
- 6.3 Earpiece
- 6.4 Threaded pin for earpieces
- 6.5 Earpiece support

7.0 3D-Registration joint

- 7.1 Connector
- 7.2 Tommy screw
- 7.3 Bite fork connector
- 7.4 Bite fork connector screw
- 7.5 3D-bite fork connector

8.0 3D-Bite fork

- 8.1 Connector

9.0 Earpiece, black, 'large'

- 9.1 Threaded connector



10.0 Nosepads

- 10.1 Nosepad, green (medium)
- 10.2 Nosepad, blue (soft)

11.0 3D-Bite fork "S" (small)

- 11.1 3D-bite fork "S" connector

Components with contact to the patient

REPROCESSING (CLEANING, DISINFECTION, AND STERILIZATION) OF COMPONENTS (SEMI-CRITICAL CONTACT)

Fundamental points

All components with semi-critical contact to the patient are to be cleaned, disinfected, and sterilized prior to each application. This is also necessary for prior to initial use, after delivery of the unsterile components (cleaning and disinfection after removal of the protective packaging, sterilization after packaging). Effective cleaning and disinfection are indispensable prerequisites for effective sterilization of the components.

You are responsible for the sterility of the components with semi-critical contact to the patient. Therefore please ensure that cleaning, disinfection and sterilization measures are used, which are sufficiently validated for the specific devices and products, that the devices used (WD, sterilizer) are regularly serviced and inspected and that the validated parameters are applied to each cycle.

Additionally, please pay attention to the legal provisions valid for your country as well as to the hygienic instructions of the doctor's practice or of the hospital. This applies particularly to the different guidelines regarding the inactivation of prions (not relevant for USA).

For all further components (without contact to the patient) standard surface disinfection may be applied (not included in these instructions).

CLEANING AND DISINFECTION

Basics

If possible, an automated procedure (WD (Washer Disinfector)) should be used for cleaning and disinfection of the components. A manual procedure – even in case of application of an ultrasonic bath – should only be used if an automated procedure is not available; in this case, the significantly lower efficiency and reproducibility of a manual procedure must be taken into consideration.

The pre-treatment step is to be performed in both cases.

Pre-treatment

Please remove coarse impurities from the components directly after use (within a maximum of 2 h).

Procedure:

1. Disassemble the components as possible (see chapter "Specific Aspects").
2. Rinse the components at least 1 min under running water¹ (temperature < 35 °C/95 °F). If applicable (see chapter "Specific Aspects"): Rinse through all cavities within the components five times using a disposable syringe (minimum volume 5 ml).
3. Manually remove all visible impurities with a clean and soft brush (or a clean, soft and fluff-free cloth), which is only used for this purpose. Never use metal brushes or steel wool. If applicable (see chapter "Specific Aspects"): Brush holes clean with a conical interdental brush.
4. Rinse again at least 1 min under running water. If applicable (see chapter "Specific Aspects"): Rinse through all cavities within the components five times using a disposable syringe (minimum volume 5 ml).

Automated cleaning/disinfection (WD (Washer Disinfector))

Pay attention to following points when selecting the WD:

- fundamentally approved efficiency of the WD (for example CE marking according to EN ISO 15883 or DGHM or FDA approval/clearance/registration)
- possibility for an approved program for thermal disinfection (A0 value > 3000 or – in case of older devices – at least 5 min at 90 °C/194 °F; in case of chemical disinfection danger of remnants of the disinfectant on the components)
- fundamental suitability of the program for instruments as well as sufficient rinsing steps in the program
- post-rinsing only with sterile or low contaminated water (max. 10 germs/ml, max. 0.25 endotoxin units/ml), for example purified/highly purified water
- only use filtered air (oil-free, low contamination with microorganisms and particles) for drying
- regular maintenance and check/calibration of the WD

¹ In case a cleaning and disinfection detergent is used for this purpose (e.g. for reasons of personnel's safety) please consider that this should be aldehyde-free (otherwise fixation of blood impurities), possess a fundamentally approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking), be suitable for the disinfection of instruments made of metallic or plastic material, and be compatible with the components (see chapter "Material Resistance").

Please consider that a disinfectant used in the pre-treatment step ensures the personnel's safety only, but cannot replace the disinfection step later to be performed after cleaning.

² Still remaining remnants of registration material can be removed easily after cooling-down in a refrigerator.

Pay attention to following points when selecting the cleaning detergent:

- fundamental suitability for the cleaning of instruments made of metallic or plastic material
- additional application – in case of non-application of a thermal disinfection – of a suitable disinfectant with approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) compatible with the used cleaning detergent
- compatibility of the used detergents with the components (see chapter “Material Resistance”)

Pay attention to the instructions of the detergent manufacturers regarding concentration, temperature and soaking time as well as post-rinsing.

Procedure:

1. Transfer the disassembled components in the WD (pay attention that the components have no contact). If applicable (see chapter “Specific Aspects”): Use a basket for small parts.
2. Start the program.
3. Remove the components from the WD after end of the program.
4. Check and pack the components immediately after the removal (see chapters “Check”, “Maintenance”, and “Packaging”, if necessary after additional post-drying at a clean place).

Validation was carried out on an equivalent product by an independent accredited and recognized test laboratory by application of the WD G 7836 CD, Miele & Cie, Gütersloh, (thermal disinfection) and the cleaning detergent Neodisher medizym (Dr Weigert GmbH & Co. KG, Hamburg) according to the specified procedure.

Manual cleaning and disinfection

Pay attention to following points during selection of the cleaning and disinfection detergents:

- fundamental suitability for the cleaning and disinfection of instruments made of metallic or plastic material
- in case of application of an ultrasonic bath: suitability of the cleaning detergent for ultrasonic cleaning (no foam development)
- application of a disinfectant with approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) compatible with the used cleaning detergent
- compatibility of the used detergents with the components (see chapter “Material Resistance”)

Combined cleaning/disinfection detergents should not be used.

Only in case of extremely low contamination (no visible impurities) combined cleaning/disinfection could be used.

- Pay attention to the instructions of the detergent manufacturers regarding concentration, temperature and soaking time as well as post-rinsing. Please use only freshly prepared solutions as well as only sterile or low contaminated water (max. 10 germs/ml), low endotoxin contaminated water (max. 0.25 endotoxin units/ml), for example purified/highly purified water, and a soft, clean, and fluff-free cloth and/or filtered air for drying, respectively.

Procedure:

Cleaning

1. Soak the disassembled components for the given soaking time in the cleaning solution so that the components are sufficiently covered. Pay attention that there is no contact between the components and no remaining air in the cavities. Assist cleaning by careful brushing with a soft brush and ultrasonic treatment (duration at least once the required soaking time in addition to brushing). If applicable (see chapter “Specific Aspects”): Rinse all cavities of the components at least five times at the beginning and at the end of the soaking time using a disposable syringe (minimum volume 5 ml). Brush holes clean using a conical interdental brush.
2. Then, remove the components of the cleaning solution and post-rinse them at least three times intensively (at least 1 min) with water. If applicable (see chapter “Specific Aspects”): Rinse all cavities of the components at least five times at the beginning and at the end of the soaking time using a disposable syringe (minimum volume 5 ml).
3. Check the components (see chapters “Check” and “Maintenance”).

Disinfection

4. Soak the disassembled components for the given soaking time in the disinfectant solution so that the components are sufficiently covered. Pay attention that there is no contact between the components and no remaining air in the cavities. If applicable (see chapter “Specific Aspects”): Rinse all cavities of the components at least five times at the beginning and at the end of the soaking time using a disposable syringe (minimum volume 5 ml).
5. Then, remove the components from the disinfectant solution and post-rinse them at least five times intensively (at least 1 min) with water. If applicable (see chapter “Specific Aspects”): Rinse all cavities of the components at least five times at the beginning and at the end of the soaking time using a disposable syringe (minimum volume 5 ml).

6. Dry and pack the components immediately after the removal (see chapter "Packaging"), if necessary after additional post-drying at a clean place.

Validation was carried out on an equivalent product by an independent accredited and recognized test laboratory by application of the cleaning detergent Cidezyme/Enzol and the disinfectant Cidex OPA (Johnson & Johnson GmbH, Norderstedt) according to the specified procedure.

CHECK

Check all components after cleaning or cleaning/disinfection, respectively, for corrosion, damaged surfaces, and impurities. Do not further use damaged components (for limitation of the numbers of re-use cycles see chapter "Reusability"). Components which are still contaminated are to be cleaned and disinfected again.

MAINTENANCE

Instrument oils or grease must not be used.

Do not mount prior to packaging and sterilization (no sterilization of the other components without contact to the patient admitted).

PACKAGING

Please insert the cleaned and disinfected components into single-use sterilized packages (single packaging), which fulfill the following requirements (material/process):

- EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance)
- suitable for steam sterilization (temperature resistance up to at least 142 °C (288 °F), sufficient steam permeability)
- sufficient protection of the components as well as of the sterilized packages to mechanical damage

STERILIZATION

Please only perform listed sterilization procedures for sterilization; other sterilization procedures must not be applied.

Do not sterilize the other components without contact to the patient.

Steam sterilization

- fractionated vacuum/dynamic air removal procedure ³ or gravity displacement procedure ⁴ (with sufficient product drying ⁵)
- steam sterilizer according to EN 13060/EN 285 or ANSI AAMI ST79 (for USA: FDA clearance)
- validated according to EN ISO 17665 (valid IQ/OQ (commissioning) and product specific performance qualification (PQ))
- maximum sterilization temperature 138 °C (280 °F); plus tolerance according to EN ISO 17665
- sterilization time (exposure time at the sterilization temperature):

Area	fractionated vacuum/dynamic air removal	gravity displacement
USA	at least 4 min at 132 °C (270 °F), drying time at least 20 min ⁵	at least 30 min at 121 °C (250 °F), drying time at least 20 min ⁵
other countries	at least 3 min ⁶ at 132 °C (270 °F) / 134 °C (273 °F)	at least 20 min at 121 °C (250 °F)

³ at least three vacuum steps

⁴ The less effective gravity displacement procedure must not be used in case of availability of the fractionated vacuum procedure.

⁵ The required drying time depends directly on parameters for which the user has the sole responsibility (load configuration and density, sterilizer conditions ...), therefore, the drying time is to be determined by the user. Nevertheless, drying times shorter than 20 min must not be applied.

⁶ respectively 18 min (inactivation of prions, not relevant for USA)

Validation was carried out on an equivalent product by an independent accredited and recognized test laboratory by application of the steam sterilizer HST 6x6x6 (Zirbus technology GmbH, Bad Grund) and the fractionated vacuum/dynamic air removal procedure as well the gravity displacement procedure.

The flash/immediate use sterilization procedure must not be used.

Do not use dry heat sterilization, radiation sterilization, formaldehyde and ethylene oxide sterilization, as well as plasma sterilization.

STORAGE

Please store the components after sterilization in the sterilized packages at a dry and dust-free place.

MATERIAL RESISTANCE

Please take care that the listed substances are not ingredients of the cleaning or disinfection detergent:

- organic, mineral, and oxidizing acids (minimum admitted pH-value 6.5)
- lyes (maximum admitted pH-value 8.5, neutral/enzymatic cleaner recommended)
- organic solvents (for example: acetone, ether, alcohol, benzine)
- oxidizing agents (for example: peroxide)
- halogens (chlorine, iodine, bromine)
- aromatic, halogenated hydrocarbons

Please do not clean any components using metal brushes or steel wool.

Please do not use acidic neutralization agents or any rinse aids.

Please do not expose any components to temperatures higher than 142 °C (288 °F)!

REUSABILITY

The components can be reused – in case of adequate care and if they are undamaged and clean – 50 times. The user is responsible for each further use as well as for the use of damaged and dirty components (no liability in case of disregard).

Specific Aspects

Art. no.	Designation	Rinsing volume	Brush	Specific/additional procedure in case of:	
				Pretreatment	Manual cleaning/disinfection
6.3	Earpiece	5 ml (disposable syringe)	conical interdental brush (3–6.5 mm)	after removal from the support bracket, brushing outside and inside	five-times rinsing inside, brushing outside and inside
2.3	Nosepad holder	5 ml (disposable syringe with needle)	conical tooth brush (3–6.5 mm)	after removal from the nosepiece rod, brushing outside and inside	five-times rinsing, brushing outside and inside
10.0	Nosepads	n.a.	n.a.	not admitted	disposable part (for single use only)
8.0	3D-Bite fork	n.a.	tooth brush	after removal from the 3D registration joint, brushing outside and inside	five-times rinsing, brushing outside and inside
1.7	Width setting screw	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)
3.0	Reference indicator	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)
6.1	Earpiece rod	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)
4.0	Axis pins (two)	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)
1.9	Connection for nosepiece	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)
2.4	Nosepiece rod screw	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)
2.2	Nosepiece rod	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)
2.1	Nosepiece holder	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)
7.1	Connector	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)
1.11	Registration joint connection screw	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)
7.3	Bite fork connector	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)
7.4	Bite fork connector screw	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)
7.2	Tommy screw	n.a.	n.a.	not admitted	soaking not admitted (standard surface disinfection)

		Packaging	Sterilization	Maximum admitted cycle number	Recommended classification according to RKI/BfArM/KRINKO guideline (only Germany, with respect to intended use)
Automated cleaning/ disinfection	Maintenance				
small parts basket	lubrication not admitted	standard procedure	standard procedure	50	semi-critical B
small parts basket	lubrication not admitted	standard procedure	standard procedure	50	semi-critical B
not admitted	n.a.	n.a.	not admitted	n.a.	semi-critical B
small parts basket	lubrication not admitted	standard procedure	standard procedure	50	semi-critical B
not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient
not admitted	n.a.	n.a.	not admitted	n.a.	without contact to the patient

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.
Ivoclar Vivadent AG