

# Special feature

## Ivocerin® – a milestone in composite technology



[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

Ivoclar Vivadent AG

Bendererstr. 2 | FL 9494 Schaan | Liechtenstein | Tel.: +423 235 35 35 | Fax: +423 235 33 60

Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2 | D-73479 Ellwangen, Jagst | Deutschland | Tel. +49 7961 889 0 | Fax +49 7961 6326 | [www.ivoclarvivadent.de](http://www.ivoclarvivadent.de)

ivoclar  
vivadent®  
passion vision innovation



# Contents

page 4	from DZW 38/13, published on 18/09/2013
	<b>1 A new class of photoinitiators for the composite technology</b> Dr Thomas Hirt
page 5-6	from DZW 39/13, published on 25/09/2013
	<b>2 State of the art: Photopolymerization in dentistry</b> Prof. Dr Norbert Moszner
page 7-8	from DZW 40/13, published on 02/10/2013
	<b>3 The photoinitiators Ivocerin and camphorquinone in comparison</b> Dr Peter Burtscher
page 9-12	from DZW 41/13, published on 09/10/2013
	<b>4 Light initiator for a new filling material</b> Dipl.-Ing. Karin Vogel
page 13-15	from DZW 42/13, published on 16/10/2013
	<b>5 Materials science evaluation of a new bulk fill material</b> Joanna-Claire Todd
page 16-18	from DZW 43/13, published on 23/10/2013
	<b>6 Laboratory simulation tests with a new bulk fill material and their clinical significance</b> Dr Siegward Heintze
page 19-20	from DZW 44/13, published on 30/10/2013
	<b>7 A composite that results in a reliable cure even if polymerized incorrectly</b> Dr Siegward Heintze
page 21-23	from DZW 45/13, published on 06/11/2013
	<b>8 The new bulk fill material in clinical use</b> Dr med. dent. Arnd Peschke

# 1 A new class of photoinitiators for the composite technology

Dr Thomas Hirt talks about Ivocerin

The 2-mm increment technique has been the standard application method for direct dental composites for many years. Over time, the properties of resin composites have been consistently improved and the bond strength of dental adhesives optimized. In contrast, the light initiators used in composite filling materials and the light-emitting efficiency of the photo-polymerization process remained unchanged for many years. However, further developments have now also been achieved in this regard. In a series of scientific articles, we present the fundamentals of photo-polymerization and information on the new photoinitiator *Ivocerin* (Ivoclar Vivadent). We describe the development, working mechanism and application of a new class of photoinitiators for the visible light spectrum and light polymerization in general. *Ivocerin* has enabled us to develop a new bulk fill material (*Tetric EvoCeram Bulk Fill*, Ivoclar Vivadent) that can be applied and cured in 4-mm bulk increments.

The new *Ivocerin* photoinitiator is characterized by high quantum efficiency, high absorption capacity and very good bleaching properties.

## Further gradual development

After obtaining very promising fundamental research data, we proceeded to optimize the synthesis of *Ivocerin* in order to reproduce the molecule in large amounts and at a high level of purity. At the same time, the first successes regarding a 4-mm depth of cure were attained with experimental composites. Furthermore, in vitro studies showed excellent results in terms of marginal seal, retention and esthetics. The results were comparable to those of *Tetric EvoCeram* placed in 2-mm




Fig. 1: Schematic absorption spectra of the different photoinitiators

Source: Ivoclar Vivadent



Fig. 2: Clinical use of *Tetric EvoCeram Bulk Fill* featuring *Ivocerin*: The preoperative situation shows a large amalgam filling and mesial proximal caries.



Fig. 3: Large, undercutting defect after removal of the filling and caries excavation



Fig. 4: Restoration of the defect with *Tetric EvoCeram Bulk Fill* (IV AG)



Fig. 5: The completed restoration is very well integrated into the surrounding tooth structure despite the deep undercutting defect.

Figs 2 to 5 provided by Dr Arnd Peschke, Schaan (Liechtenstein)



The author of this article, Dr Thomas Hirt, is Chief Technology Officer at Ivoclar Vivadent AG, Schaan, Liechtenstein. Contact details: [thomas.hirt@ivoclarvivadent.com](mailto:thomas.hirt@ivoclarvivadent.com)

increments. Ultimately, *Tetric EvoCeram Bulk Fill* was produced on the basis of *Tetric EvoCeram*, with the addition of *Ivocerin* and other photoinitiators. The monomer-filler system was also optimized to enhance the esthetics of the material, to reduce shrinkage stress and increase depth of cure. Three technical developments patented by Ivoclar Vivadent have been incorporated into *Tetric EvoCeram Bulk Fill*: *Ivocerin* as a photoinitiator (polymerization booster), a light sensitivity filter and a filler technology (shrinkage stress reliever).

## Interesting alternative

*Ivocerin* represents a very interesting alternative to the established photoinitiators camphorquinone-amine, phenyl propandione and acyl phosphine oxide (Lucirin TPO). *Ivocerin* is capable of absorbing light at a higher wavelength range than acyl phosphine oxide, and can therefore be activated by all commercially available (halogen, LED) polymerization lights (Fig. 1). In contrast to the camphorquinone-amine initiators, *Ivocerin* can be used to produce amine-free composite resins, which are colour-stable under artificial sunlight and do not interact with acid monomers. For *Tetric EvoCeram Bulk Fill* however, the combination of *Ivocerin* plus camphorquinone-amine has proven to be most successful at ensuring an excellent cure in increments of 4 mm (Figs 2 to 5 show a clinical case).

**Dr. Thomas Hirt, Schaan (Liechtenstein)**

The full version of this article is available in the R&D Report, No 19, published by Ivoclar Vivadent AG in July 2013.



## 2 State of the art: Photopolymerization in dentistry

Prof. Dr Norbert Moczner talks about chemical relationships and developments

Light-curing filling materials have been on the market since the 1970s. These materials are composed of monomers and fillers and they additionally contain initiators, stabilizers and additives which are dissolved in the monomer. As the polymerization reaction can be initiated “on demand”, these materials give users sufficient time to apply them in the cavity.

Since light transmission through composite material is only limited, a cavity usually needs to be filled in several increments, each of which is cured separately. Normally, the thickness of each increment must not exceed 2 mm, as complete curing cannot be ensured with thicker layers. In order to increase the admissible maximum layer thickness, all the

factors that influence the depth of cure need to be taken into account, including translucency, shading, and the initiators and their concentration, as well as exposure time and light intensity. An important factor in this context is the photoinitiator used. In the case of *Tetric EvoCeram Bulk Fill*, for example, the photoinitiator *Ivocerin* developed by Ivoclar Vivadent in cooperation with the Vienna University of Technology has made a major contribution to increasing the admissible layer thickness without detrimentally affecting the optical properties of the composite, such as translucency and shade.

### The chemistry of Ivocerin

Our search for an alternative visible light photoinitiator was crowned with success: in cooperation with Prof. R. Liska of the Vienna University of Technology we succeeded in developing tailor-made visible light photoinitiators based on germanium compounds. Even though it was known from the literature that radicals are formed during light-induced cleavage of e.g. organic germanium compounds [1], the potential of these com-



Prof. Dr Norbert Moczner is Head of Department: Basic Research – Polymer Chemistry at Ivoclar Vivadent AG, Schaan (Liechtenstein).  
Contact details:  
norbert.moczner@ivoclarvivadent.com

pounds as photoinitiators for the visible light spectrum had not yet been fully recognized when we started our work in this field. We were able to show for the first time that germanium compounds such as benzoyltrimethylgermane (Ge-1) or dibenzoyldiethylgermane (Ge-2, Fig. 1) represent very efficient visible light photoinitiators for methacrylate resins [2, 3]. Compared to camphorquinone (CQ), Ge-1 and Ge-2 demonstrate much more intensive absorption in the visible region. The quantum efficiency of the light-induced cleavage for Ge-2 was determined to be 0.85, whereas the quantum efficiency for camphorquinone-amine photoinitiators is below 0.10 and thus significantly lower. Moreover, quick decolouration of Ge-1 and Ge-2 was observed upon exposure to light. An analysis of dental composites containing Ge-1 or Ge-2 as photoinitiators [3, 4] also showed that they offer considerable advantages compared to camphorquinone amine-based materials. Apart from demonstrating quicker curing and excellent bleaching behaviour, these materials require a much lower concentration of the photoinitiator to achieve comparable mechanical properties. Similar to camphorquinone amine-based composite resins, Ge-1 or Ge-2-based materials also show high storage stability.

In cooperation with Prof. G. Gerscheidt and Prof. R. Saf of the Graz University of Technology as well as Prof. R. Liska of the Vienna University of Technology, the mechanism of photolysis in diacylgermanes exemplified by Ge-2 was studied using state-of-the-art methods [5]. In the study, light-induced cleavage of Ge-2 was investigated under various conditions.




Fig. 1: Structure of the germanium photoinitiators Ge-1 and Ge-2




Fig. 2: Cleavage of the germanium compound Ge-2




Fig. 3: Synthesis of the germanium photoinitiator Ivocerin

The examinations conducted confirmed that dibenzoyldiethylgermane Ge-2 forms benzoyl (B·) and germyl radicals (G·) as direct cleavage products (Fig. 2). This also proved that

Compared to conventional camphorquinone amine-based photoinitiator systems, the novel germanium photoinitiator *Ivocerin* is characterized by intensive absorption of light in

of radical formation in conjunction with *Ivocerin* is many times higher than that in conjunction with camphorquinone.

Against this background, *Tetric EvoCeram Bulk Fill* takes an exceptional position among the currently available composite resins. Apart from providing a reliable four-millimetre cure when irradiated for ten seconds at a light intensity of more than 1,000 milliwatt per square centimetre, *Ivocerin* allows the opacity of the composite resin to be established at a slightly higher level compared with other bulk-fill materials. The translucency of the material is such that it ensures a seamless blend with the surrounding tooth structure (enamel in particular) due to the composite's optimized optical properties in combination with a favourable refractive index.

**Prof. Dr Norbert Moszner,  
Schaan (Liechtenstein)**

*The full version of this article is available in the R&D Report, No 19, published by Ivoclar Vivadent AG in July 2013.*



Fig. 4: *Tetric EvoCeram Bulk Fill* is the first material containing *Ivocerin*. Figures: Ivoclar Vivadent

when dibenzoyldiethylgermane Ge-2 is used as a photoinitiator, the benzoyl (B·) and germyl (G·) radicals produced by photolytic fission are the radicals that initiate the polymerization reaction.

Based on the results of these mechanistic investigations and the evaluation of different synthesis methods and structural variations of germanium compounds, bis-(4-methoxybenzoyl)diethylgermane Ge-3 was selected as the optimum photoinitiator and protected by a patent under the name of *Ivocerin* (Fig. 3). We were able to efficiently synthesize Ge-3 ( $\lambda_{\text{max}} = 408 \text{ nm}$ ) in two stages. It showed the strongest absorption in the visible region [6]. Ge-3 is a solid substance which is not soluble in water and has an intensive yellow colour. Its melting point is at approx. 50 °C. The synthesis of Ge-3 is started by metallating protected 4-methoxybenzaldehyde using n-butyl lithium, which is then followed by a coupling reaction with dichlorodiethylgermane. In a second stage, the protecting group is split off (Fig. 3). The resulting *Ivocerin* is of high purity (> 96% when determined by High Performance Liquid Chromatography - HPLC).

Apart from chemical characterization and application tests, toxicological investigations [7] were conducted with *Ivocerin* and comprehensive patent protection measures were initiated.

Moreover, two mutagenicity tests were carried out. Neither the Ames test (in vitro) [8] nor the Mouse Micronucleus Assay (in vivo) [9] revealed any mutagenic effects of *Ivocerin*.

We were able to obtain comprehensive patent protection for the germanium compounds. The US and EP patents were granted in 2009 [10].

the visible region and high photoreactivity. Furthermore, it imparts an excellent bleaching behaviour to composite resins. The radicals required to initiate the polymerization reaction are created by light-induced cleavage of *Ivocerin*. Additional co-initiators or accelerators are not required. The quantum efficiency

## L I T E R A T U R E

- [1] Taraban MB, Maryasova VI, Leshina TV, Rybin LI, Gendin DV, Vyazankin NS. About the mechanism of the photolysis of benzoyltriethylgermane, Et<sub>3</sub>GeCOPh. *J Organomet Chem* 326 (1987) 347-355
- [2] Ganster B, Fischer UK, Moszner N, Liska R. New photocleavable structures, 4. Acylgermane-based photoinitiator for visible light curing. *Macromol Rapid Commun* 29 (2008) 57-62
- [3] Ganster B, Fischer UK, Moszner N, Liska R. New photocleavable structures, Diacylgermane-based photo-initiators for visible light curing. *Macromolecules* 41 (2008) 2394-2400
- [4] Ganster B, Fischer UK, Moszner N, Liska R, Rheinberger V. Benzoyl germanium derivatives as novel visible light photoinitiators for dental materials. *Dent Mater* 24 (2008) 901-907
- [5] Neshchadin D, Rosspeintner A, Griesser M, Gorelik V, Liska R, Hametner C, Ganster B, Saf R, Moszner N, Gescheidt G. Phototrig-
- gered reactivity of acylgermanes: Short-lived intermediates and follow-up products. *J Amer Chem Soc* 135 (2013) 17314-17321
- [6] Moszner N, Zeuner F, Lamparth I, Fischer UK. Benzoylgermanium Derivatives as novel visible-light photoinitiators for dental composites. *Macromol Mater Eng* 294 (2009) 877-886
- [7] Heppenheimer A. Cytotoxicity assay in vitro: Evaluation of test items in the XTT-Test with K69, RCC-CCR Study, 1191102, July 2008
- [8] Sokolowski A. Salmonella typhimurium and Escherichia coli reverse mutation assay with K-69, RCC-CCR Study 1120104, Oktober 2007
- [9] Honarvar N. Micronucleus assay in bone marrow cells of the mouse with K69, RCC-CCR, 1191101, Sept 2008
- [10] N. Moszner, U. K. Fischer, U. Salz, R. Liska, P. Burtscher, F. Zeuner, B. Ganster, H. Gruber. Polymerisierbare Zusammensetzungen mit Acylgermaniumverbindungen als Initiatoren. US-Patent 7,605, 190 B2 (20.10.2009), EP-Patent 1,905,415 B1

### 3 The photoinitiators Ivocerin and camphorquinone in comparison

Dr Peter Burtscher talks about the particular properties of the photoinitiators used in dentistry



Dr Peter Burtscher is Director of Organic Chemistry at Ivoclar Vivadent AG, Schaan, Liechtenstein. Contact details: peter.burtscher@ivoclarvivadent.com

All customary composite filling materials are polymerized with blue light. Therefore, the light-absorbing initiators in the composites have an inherent yellow colour, which is the complimentary colour to blue light (Figs 1 to 3). When composites are cured, the colour of the initiator largely disappears. However, a slight yellowish tinge always remains. As a result, absolutely white restorations cannot be created with materials that cure with blue light. Nevertheless, the light tinge is accepted by dental

practitioners, since natural teeth also have a yellow hue.

Figures 1 to 3 show the initiators in their pure form. Camphorquinone (CQ) and *Ivocerin* in particular, exhibit a very strong yellow colour. Figure 4 compares the absorption spectra of conventional initiators with the spectrum of *Ivocerin*. Of particular note is the strong absorption of *Ivocerin* in comparison to that of camphorquinone despite the fact that its concentration is considerably lower. Acyl phosphine oxide, for example, Lucirin TPO, predominantly absorbs light in the UV range. Hence its inherent yellowish colour is very light, as shown in Figure 2.

#### Investigations with Ivocerin

Extensive studies were undertaken to examine the suitability of *Ivocerin* in dental composite resins. The aim of these studies was to establish at which concentration *Ivocerin* would achieve the same curing results as camphorquinone. Generally, a concentration of 0.3 wt% camphorquinone is added to the monomer. At lower concentrations, reactivity declines. At higher concentrations, the depth of cure decreases, because more blue light is absorbed by the inherent colour of this initiator.

Camphorquinone cannot be used on its own. It always has to be used together with a co-initiator, usually in the form of a tertiary aromatic amine. In a comparative investigation of the reactivity of *Ivocerin* and a standard monomer containing 0.3 wt% camphorquinone and 0.6 wt% amine, monomers

containing 0.1 to 0.4 wt% *Ivocerin* were mixed and tested in a Tetric-like formulation.


Table 1 shows that 0.2 percent *Ivocerin* in the monomer achieves comparable results to the conventional amount of camphorquinone and amine. Higher concentrations of *Ivocerin* increase the strength of the composite.

#### Influence of the concentration on the depth of cure

Similarly to the flexural strength findings, Table 2 shows that a concentration of 0.2 wt% *Ivocerin* also achieves good results in terms of the depth of cure, that is, the same as the standard concentration of camphorquinone and amine. A concentration of 0.4 wt% *Ivocerin* in the monomer increases the depth of cure significantly.

#### Combination with other light-activated initiators

It was interesting to find out whether or not *Ivocerin* offers improved depth of cure in



Figs 1 to 3: Photoinitiators in pure form




Fig. 4: Absorption spectra of acyl phosphine oxide (e.g. Lucirin TPO), Camphorquinone and Ivocerin

Source: R&D Ivoclar Vivadent, 2012



combination with other initiators. The ideal choice was to combine it with camphorquinone and amine. In this case, the amine concentration had to be twice as high as that of camphorquinone.

The determination of the depth of cure is a quick method for examining the reactivity of a composite. Vickers hardness tests, however, provide more meaningful results. The composite material is polymerized in the recom-

mended layer thickness and Vickers hardness values are established at the top and bottom of the sample. According to a study by David Watts (University of Manchester), the material is adequately cured when the hardness established at the bottom of the sample corresponds to at least 80% of that measured at the surface [Watts D, Amer O, Come E. Characteristics of visible light-activated composite systems. *Br Dent J* 156 (1984) 209-215]. Therefore, the influence of the initiator system on the depth of cure (DOC, irradiation with Bluephase G1 in HP mode for 10 seconds) and the Vickers hardness (same irradiation) was evaluated.

The results show that the depth of cure and the hardness of 4-mm increments significantly improve in the cases where *Ivocerin* is used together with camphorquinone. However, it is also evident that higher concentrations of *Ivocerin* do not lead to further increases in the curing performance (Table 3).

**Summary**

The different investigations on *Ivocerin* show that this initiator can be used successfully as an alternative to camphorquinone-amine. Very good product properties are achieved when *Ivocerin* is used in combination with camphorquinone-amine. Excellent depth of cure is attained with 4-mm increments. Therefore, this initiator combination was used for the first time in the innovative *Tetric EvoCeram Bulk Fill* composite.

The fundamental research on *Ivocerin* has been incorporated into the development of *Tetric EvoCeram Bulk Fill*. Given the various influences of *Ivocerin* on the physical properties of a composite resin, the concentration of this initiator was deliberately kept low in this product. As a consequence, the colour of the filling material before polymerization is only moderately influenced and the shrinkage stress is reduced to a minimum. The resulting composite resin is available in three popular basic shades to cover the esthetic demands of posterior teeth. Furthermore, the filling material can be placed in 4-mm thick increments, which are polymerized in 10 seconds with the Bluephase Style curing light (1,100 milliwatts per square centimetre – mW/cm<sup>2</sup>).

**Dr. Peter Burtscher,  
Schaan (Liechtenstein)**

*The full version of this article is available in the R&D Report, No 19, published by Ivoclar Vivadent AG in July 2013.*

Initiator in monomer	Flexural strength (MPa)	Modulus of elasticity (MPa)
0.10 % Ivocerin	124 ± 8	7200 ± 300
0.15 % Ivocerin	124 ± 14	8500 ± 200
0.20 % Ivocerin	135 ± 9	9400 ± 600
0.40 % Ivocerin	140 ± 15	10300 ± 300
0.30 % CQ and 0.60 % amine	138 ± 10	9200 ± 600

Table 1: Comparison of the flexural strength and modulus of elasticity in relation to the different concentrations of *Ivocerin* and the standard initiator system

Initiator in monomer	Depth of cure (mm)
0.2 % Ivocerin	4.7
0.4 % Ivocerin	5.1
0.3 % CQ / 0.6 % amine	4.7

Table 2: Comparison of the depth of cure in relation to the different concentrations of *Ivocerin* and the standard initiator system

Initiator concentration in monomer	DOC (mm)	Vickers hardness at surface (N/mm <sup>2</sup> )	Vickers hardness at bottom, 4 mm	
			Value (N/mm <sup>2</sup> )	% of the surface hardness
0.2 % Ivocerin	4.7	489	333	68
0.2 % / 0.1 % CQ	4.9	511	378	74
0.2 % Ivocerin / 0.2 % CQ	5.1	524	414	79
0.4 % Ivocerin	5.1	536	437	82
0.4 % Ivocerin / 0.3 % CQ	5.6	538	462	86
0.6 % Ivocerin / 0.3 % CQ	5.6	546	465	85

Table 3: Depth of cure (DOC) and Vickers hardness of composites containing various concentrations of the initiator



## 4 Light initiator for a new filling material

Dipl. Ing. Karin Vogel talks about the development steps from Tetric EvoCeram to Tetric EvoCeram Bulk Fill



Dipl. Ing. Karin Vogel is a Senior Research Associate: Restoratives and Prosthetics at Ivoclar Vivadent AG, Schaan. Contact details: [karin.vogel@ivoclarvivadent.com](mailto:karin.vogel@ivoclarvivadent.com)

Since its introduction almost 10 years ago, Tetric EvoCeram has been producing excellent clinical results. The composite filling material meets high standards with regard to its longevity and esthetics. Furthermore, it shows very good shrinkage behaviour and surface characteristics (polishing properties and wear resistance).

Since the new bulk-fill material is intended for use without a capping layer, its surface properties have to meet the same requirements as those of Tetric EvoCeram. Alongside having to satisfy these typical composite requirements, the developers of the product were also faced with the challenge of engineering a tooth-coloured material that could be placed in 4-mm bulk increments. The incorporation of the Ivocerin light initiator has enabled the development of a composite resin that has an enamel-like appearance and the ability to be cured in bulk increments of 4-mm.

**Optimization of the monomer blend:** In addition to the new Ivocerin initiator described in detail in the previous chapters, Tetric EvoCeram Bulk Fill also contains camphorquinone and 2,4,6 trimethyl benzoyl diphenyl phosphine oxide (acyl phosphine oxide – comparable to Lucirin® TPO). These components are used in well-balanced proportions to ensure optimized reactivity, working time, depth of cure and strength.

**Light sensitivity (working time in ambient light):** Due to the incorporation of the light initiators, camphorquinone, acyl phosphine oxide and Ivocerin, Tetric EvoCeram Bulk Fill can be applied and cured in increments of 4 mm. However, it is of utmost importance to prevent premature polymerization as a result of the heightened conversion rate of the monomer. That is, an adequate working time, for placing and sculpting the restoration should be ensured. As light-cured composites generally contain photoinitiators that react to blue light, the blue light contained in ambient light or operatory light can also trigger premature polymerization in these materials (Fig. 1).

Tetric EvoCeram Bulk Fill contains a patented light sensitivity inhibitor, which prevents premature polymerization and enables the filling material to be manipulated for three minutes (200 s) under defined light conditions (ISO 4049:2009) of 8000 lux. Conventional phenolic inhibitors (MeHQ, BHT) require a concentration of at least 1000 ppm relative to

the monomer in order to delay a reaction to ambient light.

Just 1/10 of this amount is necessary in the case of the light sensitivity inhibitor. This is advantageous, as the small amount of stabilizer/inhibitor delays the polymerization process at low-level blue light, without impairing the depth of cure or any of the other polymerization properties.

**Optimization of the filler composition:** The filler technology incorporated in Tetric EvoCeram Bulk Fill is based on that of the clinically proven Tetric EvoCeram. In order to fulfil the desired composite resin requirements, Tetric EvoCeram Bulk Fill contains a number of different fillers (Figs 2 and 3). The fine primary




Fig. 1: Sensitivity to ambient light/working time of various bulk-fill composites established according to ISO 4049 (R&D Ivoclar Vivadent AG, June 2011)





Fig. 2 and 3: Scanning electron micrographs (magnification: 1000x) of the filler composition and surface structure of Tetric EvoCeram Bulk Fill (R&D Ivoclar Vivadent AG, 2011)




Fig. 4: SonicFill/Kerr (1000x)(R&D Ivoclar Vivadent AG, 2011)




Fig. 5: Venus Bulk Fill/Heraeus Kulzer (1000x)(R&D Ivoclar Vivadent AG, 2011)




Fig. 6: SDR Dentsply (1000x)(R&D Ivoclar Vivadent AG, 2011)

particles of the fillers are responsible for the composite's wear resistance and its excellent polishing properties, which are manifested in a smooth surface texture and high lustre. The largest particles of the inorganic fillers of *Tetric EvoCeram Bulk Fill*

the filler content and reduces the polymerization shrinkage. However, when these materials are also used as the top-most layer, they have a distinct effect on the filling's surface texture, polishing properties, its wear behaviour and plaque resistance.

measure 3 µm. The composite fillers have a maximum size of 50 µm. In the polymerized state, they behave like the smaller inorganic primary particles. The large filler particles do not protrude from the surface. Consequently, the filling can be polished to a high-gloss finish. Other manufacturers of bulk-fill composites mainly use coarse fillers (Figs 4 to 6). This increases

**Polishing properties:** The polishing properties are determined by the composition of the different-sized fillers. *Tetric EvoCeram Bulk Fill* contains the same type of fillers as *Tetric EvoCeram*. As a result, its polishing properties are just as good as those of the conventional filling material. The two composites (SonicFill from Kerr and SDR Flow from Dentsply) that contain the coarse inorganic fillers show very little shine, even after 30 s of polishing (Fig. 7).

**Wear behaviour:** Apart from affecting the polishing properties, the type of filler used also has a considerable influence on the wear behaviour of a composite resin. The materials containing the visibly coarser fillers are less resistant to wear than the composites featuring smaller filler particles (Fig. 8).

**Shrinkage and shrinkage stress and their influencing factors:** The composite filler is responsible for reducing polymerization shrinkage and lowering shrinkage stress. *Tetric*

*EvoCeram Bulk Fill* contains a special composite filler that relieves shrinkage stress. The volumetric shrinkage of *Tetric EvoCeram Bulk Fill* and *SonicFill* after 1 h is below 2% and therefore comparable to that of conventional low-shrinkage composites (Fig. 9). The

shrinkage of flowable bulk-filling composites is significantly higher at 3.3%. The composite filler acts as a shrinkage stress reliever. It is a “gentle giant”, which is capable of absorbing the shrinkage stress due to its low modulus of elasticity of 10,000 MPa com-

pared with 70,000 MPa in conventional glass fillers. This is a particularly important factor in bulk-fill materials. The shrinkage stress of the composite rises, as the thickness of the increments increases from 0.8 mm to 2.0 mm. Nevertheless, the in-




Fig. 7: Average surface gloss of five different composite filling materials compared with *Tetric EvoCeram Bulk Fill* after polishing with *optraPol next Generation* – in relation to the polishing time (R&D Ivoclar Vivadent AG, 2011)




Fig. 8: Average vertical wear of *Tetric EvoCeram Bulk Fill* in comparison with *Tetric EvoCeram* and other bulk-fill materials (R&D Ivoclar Vivadent AG, 2011)




Fig. 9: Volumetric shrinkage of bulk-fill composites (\*) measured with a mercury dilatometer after 60 min, or (°) according to the Archimedes principle (K. Vogel, Rheinberger V. Shrinkage and contraction force of bulk filling and microhybrid composites. AADR Abstract, 858, Florida 2012)






Fig. 10: Shrinkage force of Tetric EvoCeram Bulk Fill in different layer thicknesses (R&D Ivoclar Vivadent AG, 2011)

**Shrinkage stress compared with that of other bulk-fill and conventional composites**

Increment thickness	Sculptable bulk-fill composites			Flowable bulk-fill composites		Conventional composites		
	Tetric EvoCeram Bulk Fill	SonicFill / Kerr	x-tra fil / Voco	SDR / Dentsply	Venus Bulk Fill / Heraeus Kulzer	Filtek Supreme XTE / 3M Espe	Filtek Z250 / 3M Espe	Herculite XRV Ultra / Kerr
0.8 mm	1.21 MPa	1.38 MPa	1.16 MPa	1.11 MPa	1.34 MPa	1.51 MPa	1.22 MPa	1.64 MPa
2.0 mm	1.51 MPa	1.77 MPa	1.8 MPa	not measurable		1.93 MPa	1.63 MPa	1.60 MPa
4.0 mm	1.40 MPa	1.77 MPa	1.83 MPa			not indicated		

Table 1: Shrinkage stress of bulk-fill and conventional composites in different increment thicknesses (R&D Ivoclar Vivadent AG, February 2013)

crease within the first few seconds of curing is lower in the thicker layers (Fig. 10). This slow build-up reduces the risk of marginal gap formation. The shrinkage stress in 4-mm increments is lower than in 2-mm layers, and the slope of the curve is much shallower than that of increments of 0.8 and 2.0 mm.

**Shrinkage stress in comparison:** The shrinkage stress rises in relation to the increasing thickness of the increment from 0.8

to 2.0 mm. At a layer thickness of 4 mm, the shrinkage stress ceases to increase. Furthermore, 4-mm increments of bulk-fill composites do not show higher shrinkage stress levels than conventional 2-mm layers. *Tetric EvoCeram Bulk Fill* shows the lowest shrinkage stress of all the sculptable bulk-filling materials in increments of 2 and 4 mm (Table 1). The new *Ivocerin* initiator has enabled the development of an esthetic bulk-fill material that can be applied in 4-mm increments and cured

within 10 seconds. What's more, the material is comparable to conventional direct filling composites in terms of its surface quality, working time and shade blend with natural teeth.

**Dipl.-Ing. Karin Vogel, Schaan (Liechtenstein)**

*The full version of this article is available in the R&D Report, No 19, published by Ivoclar Vivadent AG in July 2013.*



# 5 Materials science evaluation of a new bulk fill material

**Joanna-Claire Todd talks about materials science investigations underpinning Tetric EvoCeram Bulk Fill and the new photoinitiator it contains**



Joanna-Claire Todd (BA, MPH) is a Research Associate at the Research and Development Department of Ivoclar Vivadent AG, Schaan, Liechtenstein. Contact details: joanna-claire.todd@ivoclarvivadent.com

Before the introduction of bulk fill composites, standard dental teaching recommended a maximum layer thickness of 2 mm for composite fillings [1, 2]. This was in order to minimise shrinkage stress and to ensure adequate depth of cure. Assuming correct, adequate curing with a suitably functioning curing unit, translucency and shade have the most significant effect on the curing depth. The darker and more opaque a composite, the lower the depth of cure [3]. The characteristics of *Tetric EvoCeram Bulk Fill* allow an enamel-like translucency of 15%, which ensures adequate polymerisation at depth. The inclusion of *Ivocerin* means that *Tetric EvoCeram Bulk Fill* can be applied in increments of up to 4 mm, without compromising cure or esthetics. The following internal tests conducted by Ivoclar Vivadent's Research and Development (R&D) department, as well as external investigations confirm the effectiveness of the photoinitiator formulation with *Ivocerin*.

### Vickers hardness

Vicker hardness testing was carried out by R&D Ivoclar Vivadent AG: Samples of each of the three *Tetric EvoCeram Bulk Fill* shades were cut and the Vickers hardness measured at the top and at a depth of 4 mm after removal of the inhibition layer. The values measured at the top were set to 100% and the values measured at 4 mm are expressed as a percentage of this value. Professor David Watts of the

University of Manchester, UK, defined an acceptable curing depth as when the bottom hardness corresponds to at least 80% of the surface hardness [4]. Various light intensities were employed and the curing times were adjusted accordingly to ensure a similar light output in each case. For each of the shades, the 4 mm hardness value exceeded 80% of the surface hardness under all curing settings (Figs 1 to 3).

### Depth of cure in comparison to other composites



In July 2012, Dr A. Rzanny and M. Facht of Universitätsklinikum Jena, Germany, investigated the depth of cure achieved with Bluephase and Bluephase Style in *Tetric EvoCeram Bulk Fill* and compared it to that of other composites: Rzanny et al calculated the depth of cure of the composites: *Tetric EvoCeram Bulk Fill* (IVA), Venus Bulk Fill (Universal, Heraeus Kulzer) and *Tetric EvoCeram* (A3) by using a Penetrometer after curing for 10 seconds with Bluephase G2 (1200 mW/cm<sup>2</sup>) or Bluephase Style (1100 mW/cm<sup>2</sup>) and established their Vickers hardness values. Depth of cure measurement according to ISO 4049 yielded the following results: There was no significant difference between curing lamps for any of the composites. Both bulk fill composites *Tetric EvoCeram Bulk Fill* and Venus Bulk Fill far exceeded the manufacturer indicated allowable increment thickness (4 mm) in terms of depth of cure (approx. 5 mm). *Tetric EvoCer-*

*am* is not a bulk fill composite and is intended to be applied in 2 mm increments (Fig 4).

Vickers hardness in this investigation: The Vickers hardness results for *Tetric EvoCeram Bulk Fill* all exceeded the 80% ratio necessary. When cured with Bluephase G2 the ratio was 87.6% after 24 hours and 83.6% after 7 days. When cured with Bluephase Style it was 80.3% after 24 hours and 87.5% after 7 days. The authors conclude that both Bluephase G2 and Bluephase Style are equally suitable for polymerising the three composites investigated.

### Evaluation of the depth of cure and surface microhardness

C. Sabatini of the Dental Biomaterials Research Laboratory, State University of New York, Buffalo, USA, evaluated the depth of cure and surface microhardness of a new bulk fill composite system in October 2012: By means of Knoop hardness tests, Sabatini investigated the depth of cure and surface mi-







Fig. 1 to 3: Tetric EvoCeram Bulk Fill Shades IVA ①, IVB ② und IVW ③ – 4 mm depth hardness as percentage of surface hardness, measured with different light intensities Source: R&D, Ivoclar Vivadent

Fig. 4: Depth of cure for various composites when cured with Bluephase G2 and Bluephase Style for 10 seconds Source: Dr. A. Rzanny, M. Facht, Universitätsklinikum Jena, Germany




Fig. 5: Bottom/top Knoop hardness ratios at 0.5 mm increments for Tetric EvoCeram Bulk Fill

Source: Sabatini, October 2012

crohardness of Tetric EvoCeram Bulk Fill, two further bulk fill products: x-tra fil (Voco) and SonicFill (Kerr), and Tetric EvoCeram as a control. Two light curing units were employed: Bluephase G2 (1200 mW/cm<sup>2</sup>) and Bluephase 20i Turbo (2000mW/cm<sup>2</sup>) with exposure times of 10 and 5 seconds respectively.

Results: Two way analysis of variance (ANOVA) revealed no difference in the average bottom/top hardness values for the type of polymerisation unit used. However significant differences were found between certain restorative composites cured with the same light source ( $p < 0.001$ ). When polymerised with Bluephase G2 there were no significant differences between x-tra fil, Tetric EvoCeram

Bulk Fill and Tetric EvoCeram but all were significantly different to SonicFill. When polymerised with Bluephase 20i, there were also no significant differences between x-tra fil, Tetric EvoCeram Bulk Fill and Tetric EvoCeram however Sonic Fill was significantly lower than Tetric EvoCeram Bulk Fill and Tetric EvoCeram but not significantly lower than x-tra fil. Notably both Tetric EvoCeram Bulk Fill at 4 mm and Tetric EvoCeram at 2 mm fulfilled the Watts criterion with all figures exceeding 80%. In this investigation, x-tra fil and SonicFill did not achieve 80%.

Figure 5 shows the bottom/top hardness ratios at different depths for Tetric EvoCeram Bulk Fill when cured with both Bluephase lights.

	x-tra fil / Voco (4 mm)	SonicFill / Kerr (5 mm)	Tetric EvoCeram Bulk Fill (4 mm)	Tetric EvoCeram (2 mm)
Bluephase G2	70.6 %	47.1 %	85.7 %	85.1 %
Bluephase 20i	69.4 %	55.6 %	86.9 %	81.4 %

Table: Average bottom/top hardness ratios at the recommended increment thickness per material Source: Sabatini, October 2012

This investigation showed that Tetric EvoCeram Bulk Fill achieved in excess of the necessary 80% bottom/top hardness ratio at a depth of 4 mm, independent of the light source. Tetric EvoCeram also achieved this at a depth of 2 mm.




Fig. 6: Microhardness (HV) of various composites at various depths and curing times

Source: S. Zawawi, Boston University, USA




Fig. 7: Mean modulus of elasticity of various composites when cured for 10 vs. 40 seconds

Source: S. Zawawi, Boston University, USA

### The effect of curing duration on the depth of cure and modulus of elasticity

In vitro testing was carried out to evaluate the effect of curing duration on the depth of cure and modulus of elasticity in bulk fill composites [5]. S. Zawawi et al. from the Department of Restorative Sciences and Biomaterials of Boston University, Boston (MA), USA, used Vickers hardness measurements to evaluate the effect of curing duration on the depth of cure and modulus of elasticity in the bulk fill composites: *Tetric EvoCeram Bulk Fill*, SDR (Dentsply) and Venus Bulk Fill (Heraeus Kulzer). The Bluephase 16i (1600 mW/cm<sup>2</sup>) curing light was used to cure the composites for either 10 or 40 seconds. The Vickers hardness values for the composites at the surface and at a depth of 2 mm, 4 mm and 6 mm were determined.

In this study, *Tetric EvoCeram Bulk Fill* exhibited higher microhardness than SDR and Venus Bulk Fill at all depths and curing times (Fig. 6).

The mean modulus of elasticity was also measured for each bulk fill composite with both 10 seconds and 40 seconds of curing (Fig. 7). There was no significant difference between *Tetric EvoCeram Bulk Fill* samples when cured for 10 seconds or 40 seconds.

However there were significant differences in modulus between the different materials at both 10 and 40 seconds. Whereas curing duration had no discernible effect on the modulus of elasticity for *Tetric EvoCeram Bulk Fill*, this was not the case for SDR or Venus Bulk Fill – where there was a clear difference (increase) between 10 and 40 seconds polymerisation time.

Conclusion: *Tetric EvoCeram Bulk Fill* achieves higher mechanical properties than the other products and is almost indifferent to the length of cure (10s vs. 40s). Notably the bottom/top ratio of the Vickers hardness for *Tetric EvoCeram Bulk Fill* at 4 mm (10 seconds = 95.5%, 40 seconds = 98.5%) and even at 6 mm (10 seconds = 85.7%, 40 seconds = 96.2%) far exceeds the 80 % level stipulated by Watts. The studies involving Vickers hardness and Knoop hardness measurements confirmed the efficacy of the *Ivocerin* formulation, which means that *Tetric EvoCeram Bulk Fill* offers an adequate depth of cure of at least 4 millimetres.

**Joanna-Claire Todd (BA, MPH),  
Schaan (Liechtenstein)**

The full version of this article is available in the R&D Report, No 19, published by Ivoclar Vivadent AG in July 2013.

## LITERATURE

- [1] Pilo R, Oelgiesser D, Cardash H. A survey of output intensity and potential for depth of cure among light-curing units in clinical use. *J Dent* 27 (1999) 235-259.
- [2] Sakaguchi R, Douglas W, Peters M. Curing light performance and polymerization of composite restorative materials. *J Dent* 20 (1992) 183-8.
- [3] Kawaguchi M, Fukushima T, Miyazaki M. The relationship between cure depth and transmission coefficient of visible-light-activated resin composites. *J Dent Res* 73 (1994) 516-521.
- [4] Watts D, Amer O, Combe E. Characteristics of visible light activated composite systems. *Br Dent J* (1984) 209-215.
- [5] Zawawi S, Brulat N, Nathanson D. Curing duration vs. depth of cure and modulus of bulk fill composites. Brazil (2012) IADR Abstract 121.

## 6 Laboratory simulation tests with a new bulk fill material and their clinical significance

**Dr Siegward Heintze reports on in vitro tests carried out on Tetric EvoCeram Bulk Fill containing the new light initiator Ivocerin**

For decades, students of dentistry were taught that direct restorations made of composite resin should be applied into the cavity in increments. The main reason was the low depth of cure of conventional composites (one to two millimetres [mm]). Additional reasons included the belief that the incremental technique could partially offset the shrinkage, or shrinkage stress, which takes place as the composite resin polymerizes. It was said that polymerization shrinkage may result in marginal gaps, debonding from the cavity floor and/or, if a strong bond

to the tooth is present, may also lead to movement of the remaining tooth structure (cusp movement/deflection), with the latter being believed to entail the risk of enamel/dentin crack formation, infractions and cusp fractures.

Imagine a composite material that fully cures in thick increments of up to 4 mm when exposed to the light of a standard curing device and in addition is very tolerant of operator errors, e.g. oblique positioning of the light probe or polymerization from a substantial distance to the filling? Such a material is indeed available and on the following pages we will show which errors this material is able to tolerate and how it behaves in the cavities of extracted teeth.

### Tetric EvoCeram Bulk Fill and marginal gaps

The aim was to confirm or refute the assumption that thick composite layers result in poorer marginal quality than composites applied using the conventional incremental technique. This question was investigated with an experimental setup that used an ex-




*The author of this article, Dr Siegward Heintze, is Head of Preclinic Department at Ivoclar Vivadent AG in Schaan (Liechtenstein). Contact details: siegward.heintze@ivoclar-vivadent.com*

tracted mandibular molar in which proximal occlusal cavities were prepared on both the mesial and distal side. The cavity was designed with a proximal depth of 4 mm, an occlusal depth of 3 mm and a lingual-buccal width of 5 mm (Figs 1 to 3). On one side of the tooth, *Tetric EvoCeram Bulk Fill* was applied in a single increment; on the other side, *Tetric EvoCeram* was applied in three increments: a horizontal gingival, an oblique buccal and an oblique lingual increment (Fig. 2). Each layer was light-cured for 10 seconds with a Bluephase G2 curing light (1,200 milliwatts per square centimetre - mW/cm<sup>2</sup>). Two adhesive systems were employed: the single-component etch & rinse system ExciTE F and the self-etching two-component system AdheSE.

Eight fillings were placed for each test group and all restoration margins were confined to the enamel. After the fillings had been placed, excess removed and the restoration surfaces polished, the teeth were first immersed in water for 24 hours and subsequently subjected to thermocycling (10,000 cycles between 5 degrees Celsius and 55 degrees Celsius) for 10 days. After that, replicas were produced and the marginal quality was evaluated using scanning electron microscopy (SEM) (Figs 4 and 5).


Figure 6 shows the percentage of regular margin in the axio-proximal part. It is immediately evident that a) there is no difference between the bulk and incremental technique and b) that the adhesive system, rather than the incremental technique, is the determining factor for marginal quality: The etch & rinse system, which included enamel acid



*Fig. 1: Extracted lower molar with two-surface cavities*





*Fig. 2: After placement of filling and removal of excess. On the left: Tetric EvoCeram Bulk Fill in one increment; on the right: Tetric EvoCeram in 3 increments*



*Fig. 3: Proximal view: above Tetric EvoCeram Bulk Fill, below Tetric EvoCeram*





Figs 4 and 5: SEM images (x200) of the margin of a filling placed with Tetric EvoCeram Bulk Fill and the etch & rinse system Excite F in an extracted molar after 10,000 temperature cycles. Fig. 4 Axio-proximal enamel. Fig. 5 Cervical enamel




Fig. 6: Percentage of regular margin (and standard deviation) after SEM analysis (x200) in fillings on extracted lower molars after 10,000 temperature cycles in relation to filling material and filling technique




Fig. 7: Two-surface filling with Tetric EvoCeram Bulk Fill on a molar: horizontal section of a micro-CT image (7-µm resolution; courtesy of Scanco Medical, Switzerland). Good adaptation to the margin with only few air bubbles Source: Scanco Medical

etching with 36 per cent phosphoric acid, resulted in significantly better margins (fewer gaps and irregularities) than the self-etching system.

Figure 7 shows the cross section of a dental filling that was X-rayed using high-resolution micro-CT. The filling shows a good adaptation to the cavity margin and hardly any major air bubbles.

Given that the laboratory tests did not reveal any difference between the fillings placed with Tetric EvoCeram Bulk Fill (one increment) and


Tetric EvoCeram (three increments), would it not make sense to refer to the results of clinical studies on conventionally layered composite fillings placed together with Excite F (or any other etch & rinse system) and AdheSE? It could be expected that Tetric EvoCeram Bulk Fill would achieve similar results.

Such studies do exist. Clinical results of up to 6 years are available for Tetric EvoCeram and Excite. In one study, almost 90 per cent of all the restorations were still intact after six years [1]; the main reasons for replace-

ments were marginal caries and fractures, with marginal caries mainly occurring in patients with a high caries risk. In another study carried out by the same group of researchers, Excite and Tetric Ceram were used to place Class II fillings in 59 patients and this study concluded with similar results [2]: after seven years only about 10 per cent of the fillings showed marginal staining and only 2 fillings exhibited marginal caries. These results tally with the above mentioned in-vitro results. Other studies that used Excite and AdheSE confirmed the clinical superiority of Excite over AdheSE with regard

to the marginal quality of posterior restorations [3-5].

If we take into account that a) clinical data on Class II restorations placed with Tetric EvoCeram and the adhesive systems Excite F and AdheSE are available from the Internal Clinic and that b) data from other published studies on other composite resins placed in combination with AdheSE and Excite confirm the reliable performance of the restorations, we may assume that the restorations placed in bulk will provide



**Fig. 8: Percentage of regular margins (and standard deviation), examined with magnifying glasses and an explorer, in extracted lower molars after 10,000 temperature cycles in relation to the filling material and filling technique** Figs 1 to 6 and 8: Ivoclar Vivadent AG, Schaan

clinical results that are equivalent to the results achieved with *Tetric EvoCeram* and the incremental technique, bearing in mind, how-

ever, the limitations of laboratory tests and other physical characteristics such as flexural strength, expansion, shrinkage, etc.

How come the difference in marginal staining observed between the self-etch two-component system and the etch & rinse system was clinically smaller than the laboratory data would have led us to expect? The answer to this question is simple: In the laboratory, the margins of the fillings are evaluated far too minutely under a microscope; minor irregularities and gaps that are irrelevant for clinical performance, are discovered and recorded. If the fillings are evaluated in-vitro in the same way as the dentist does in the patient – i.e. using magnifying glasses and a sharp explorer – the difference is less pronounced (Fig. 8) [6].

Restoration margins exhibiting gaps or irregularities seldom lead to marginal caries but, in the course of time, may cause marginal staining, which may be mistaken for marginal caries and result in the filling being unnecessarily replaced. A meta-analysis of 59 clinical studies on posterior composite fillings revealed that, on average, marginal staining was diagnosed seven times more frequently than marginal caries after ten years (21 per cent of the fillings versus 3 per cent of the fillings) [7]. Furthermore, the occurrence of marginal caries was unrelated to the type of adhesive system used. Marginal caries is mostly found in patients with a high caries activity.

**Conclusion**

*Tetric EvoCeram Bulk Fill*, incorporating the innovative photoinitiator system *Ivocerin*, is ideally suited for the direct restorative technique because it allows cavities to be filled in large increments of up to 4 mm. Compared with the three-increment technique, 4-mm cavities filled with a single increment of *Tetric EvoCeram Bulk Fill* did not lead to more marginal gaps. Several studies have shown that the adhesive system, rather than the composite, is the determining factor for marginal quality. *Tetric EvoCeram Bulk Fill* allows for effective posterior restorations without compromises in quality.

### Conclusion

**Dr. Siegwald Heintze, Schaan (Liechtenstein)**

*The full version of this article is available in the R&D Report, No 19, published by Ivoclar Vivadent AG in July 2013.*

**L I T E R A T U R E**

- [1] van Dijken JW, Pallesen U. A six-year prospective randomized study of a nano-hybrid and a conventional hybrid resin composite in Class II restorations. *Dent Mater* 29 (2013) 191-198
- [2] van Dijken JW, Pallesen U. Clinical performance of a hybrid resin composite with and without an intermediate layer of flowable resin composite: a 7-year evaluation. *Dent Mater* 27 (2011) 150-156.
- [3] Boeckler A, Boeckler A, Eppendorf K, Schaller HG, Gernhardt CR. A Prospec-

- tive, randomized clinical trial of a two-step self-etching vs two-step etch-and-rinse adhesive and SEM margin analysis: Four-year results. *J Adhes Dent* 6 (2012) 582–592.
- [4] Palaniappan S, Elsen L, Lijnen I, Peumans M, Van Meerbeek B, Lambrechts P. Nanohybrid and micro-filled hybrid versus conventional hybrid composite restorations: 5-year clinical wear performance. *Clin Oral Invest* 16 (2012) 181-190.
- [5] Barabanti N, Gagliani M, Roulet JF, Testori T, Özcan M, Cerutti A.

- Marginal quality of posterior microhybrid resin composite restorations applied using two polymerization protocols: 5-year randomised split mouth trial. *J Dent* 41 (2013) 436–442.
- [6] Heintze SD. Clinical relevance of tests on bond strength, microleakage and marginal adaptation. *Dent Mater* 29 (2013) 59-84.
- [7] Heintze SD, Rousson V. Clinical effectiveness of direct class II restorations - A metaanalysis. *J Adhes Dent* 14 (2012) 407-431.



## 7 A composite that results in a reliable cure even if polymerized incorrectly

Dr Siegwald Heintze reports on in vitro investigations examining the polymerization of Tetric EvoCeram Bulk Fill

A bulk-fill composite such as Tetric EvoCeram Bulk Fill containing the innovative Ivocerin initiator system presents an ideal material for the direct restorative technique, because it allows cavities to be filled in large increments of up to 4 millimetres (mm). However, how does the material behave in the presence of polymerization-related errors, such as an oblique position of the light probe or curing from a large distance to the filling?

Lecturers and opinion leaders often refer to practitioner-errors with the suggestion that such mistakes lead to the insufficient poly-

merization of composite materials. They demonstrate that a shadow is created if the light guide is held at an angle to the composite




Fig. 1: Oblique position of the light probe and 5-mm distance to the tooth




Fig. 2: Vickers hardness indentations in the 4-mm composite filling

material to cure a proximal filling (Fig. 1). Un-toward effects are also to be expected if the distance between the light probe and filling is too large or if the light probe is not positioned precisely over the filling. Interestingly, these statements are not backed up by study results – and hardly any studies on the above-mentioned various deviations from the standard position of the light probe can be found. Only the increase in distance to the composite has been investigated in a few studies and it was found that the depth of cure decreased with an increase in distance. However, these results were to some extent affected by the curing devices used [1]. How does the *Tetric EvoCeram Bulk Fill* composite material behave in the presence of polymerization-related errors? To investigate this issue, we drew up a test plan to examine a variety of scenarios:

**Test group 1:** Light probe in vertical position to the composite

**Test group 2:** Light probe at a 40-degree angle to the composite

**Test group 3:** Light probe in vertical position, at a distance of 5 mm to the composite

**Test group 4:** Light probe at a 40-degree angle, at a distance of 5 mm to the composite

**Test group 5:** Light probe in vertical position over the centre of the molar, so that both proximal sides are only partially illuminated

**Test group 6:** Same as group 5 but with a distance of 5 mm

**Test group 7:** Light probe in vertical position to the composite. The entire light emission window is covered with a layer of composite material (thickness: 0.8 mm).

The same experimental setup was used as for the tests on marginal adaptation (see Part 6 of this series, Figs 1 to 3). Prior to the application of the composite, a steel matrix was placed and secured with a wooden wedge. The cavity was lightly isolated with glycerine gel to ensure that the composite could be removed in a single piece. *Tetric EvoCeram Bulk Fill* was applied in a single increment and polymerized for 10 seconds with Bluephase G2 (1,200 milliwatts per square centimetre - mW/cm<sup>2</sup>) – in the manner specified for the test groups above. Four fillings were placed for each group. After 24 hours of dry storage at 37 °C (incubator), the composite specimens were invested in epoxy resin in a dark-room, reduced to the centre by grinding and then polished. Subsequently, the hardness profile of the entire proximal length was determined, by making indentations every 0.5 mm using a Vickers Hardness tester (Fig. 2). The data were averaged and processed graphically. The results came as a surprise (Fig. 3). None of the test groups 1 to 6 showed cervical




Fig. 3: Relative Vickers hardness (in per cent) of the cervical portion of the Tetric EvoCeram Bulk Fill filling relative to the Vickers hardness of the coronal part – depending on the position of the light probe Figures: Ivoclar Vivadent

hardness values below 80 per cent of the surface hardness. The exception was the group with the composite-covered light emission window. This group achieved merely 71 % of the coronal hardness value.

It can therefore be said that *Tetric EvoCeram Bulk Fill* undergoes complete curing even if the light probe is positioned incorrectly. Obviously, operators should use a curing light capable of generating the required light output. Additionally, the light performance should be checked regularly with a measuring device. A field test in general dental practices in Germany revealed that half of the curing lights tested failed to provide the expected performance [2]. It is equally important to ensure that the light emission window is not contaminated or covered with composite material. The same field study found that 37 % of

the light probes were contaminated with composite or adhesive.

To sum up, *Tetric EvoCeram Bulk Fill*, incorporating the new *Ivocerin* initiator system, allows 4-mm deep cavities to be filled with a single increment. The material is tolerant of errors related to the polymerization process. Even careless positioning of the light probe

still results in the sufficient curing of 4 mm layers.

**Dr. Siegwald Heintze,  
Schaan (Liechtenstein)**

*The full version of this article is available in the R&D Report, No 19, published by Ivoclar Vivadent AG in July 2013.*

## L I T E R A T U R E

[1] Ernst CP, Meyer GR, Müller J, Stender E, Ahlers MO, Willershausen B. Depth of cure of LED vs QTH light-curing devices at a distance of 7 mm. *J Adhes Dent* 6 (2004) 141–150.

[2] Ernst CP, Busemann I, Kern T, Willershausen B. Feldtest zur Lichtemissionsleistung von Polymerisationsgeräten in zahnärztlichen Praxen. *Dtsch Zahnärztl Z* 61 (2006) 466–471



## 8 The new bulk fill material in clinical use

**Dr med. dent. Arnd Peschke discusses the clinical experiences with and the performance of modern bulk fill materials**

There is no question that dentists would like to be able to place posterior restorations in thicker layers and thus in fewer working steps. The layering technique [1–3], in which thin layers of composite are placed in defined orientation – in order to compensate for shrinkage stress, is complicated and its advantages in terms of clinical quality are debatable [4, 5]. Furthermore, the low curing depth of composites has so far limited the possibility for dentists to efficiently adapt their application technique to their patients’ individual anatomic situation. Bulk-fill composites should overcome these limitations.

Bulk-fill composites enable a more efficient and ergonomic working technique due to the fact that they can be placed in thicker increments (usually 4 mm). However, bulk-fill composites with a conventional initiator system generally comprise larger filler particles in order to counteract shrinkage and polymerization stress and they demonstrate a high level of translucency in order to allow increased layer thickness. Both material properties have a significant influence on the clinical characteristics such as surface quality and esthetics (Fig. 1).

Furthermore, such materials often require longer polymerization times to ensure sufficient curing (Fig. 2), or they react very sensitively to ambient light – limiting their process-

ing time under operatory light (see Fig. 1 regarding the light sensitivity in the article by Karin Vogel, page 9 of DZW 41/13).

In comparison to the bulk-fill materials of other manufacturers, *Tetric EvoCeram Bulk Fill* shows a beneficial esthetic integration (Fig. 1). Materials with a higher translucency level do not demonstrate such good shade adaptation – a notable problem in situations with discoloured dentin.

Due to the addition of the innovative photoinitiator *Ivocerin*, *Tetric EvoCeram Bulk Fill* can be polymerized in ten seconds (Fig. 2) at a light intensity of  $\geq 1,000 \text{ mW/cm}^2$ . *Tetric EvoCeram Bulk Fill* is an exceptional product, featuring well-balanced properties, optimized for



Dr med. dent. Arnd Peschke is Director R&D Clinic at Ivoclar Vivadent AG in Schaan, Liechtenstein. Contact details: [arnd.peschke@ivoclarvivadent.com](mailto:arnd.peschke@ivoclarvivadent.com)

the posterior region. These properties are demonstrated by the excellent results attained in laboratory and preclinical tests, however only clinical studies can provide truly reliable information as to the final clinical performance.

### In-house clinical study

In an in-house clinical investigation of *Tetric EvoCeram Bulk Fill*, 35 restorations (11 Class I and 24 Class II) were placed in combination with a one-bottle etch & rinse adhesive system. The restorations were assessed according to FDI criteria [6, 7], using a semi-quantitative clinical evaluation (SQUACE method), which



Fig. 1: Standardized Class II cavities in resin teeth with cavity floor characterized with IPS Empress Direct Color Grey and filled with various bulk-fill composites in different levels of transparency. (ICDE Schaan, December 2012)




Fig. 2: Polymerization times of different bulk-fill materials according to the manufacturer's instructions (status: March 2013)

enables the documentation of marginal deficiencies in relation to the overall filling margin [8]. In the following paragraphs, the results of the clinical aspects – postoperative sensitivities, esthetic integration, polishability, occurrence of enamel cracks and initial marginal quality – are briefly described.

**Postoperative sensitivity**

In the majority of cases (33), defective fillings represented the preoperative situation. Only two cases involved primary caries. In 51 per cent of the cases (18), a superficial to medium caries progression towards the pulp was observed. Deep lesions were discovered in 46 per cent (16) of the patients. Only in one case (3 per cent), was Caries Profunda treatment (application of a calcium hydroxide compound) necessary, due to proximity to the pulp.

Despite the fact that the lesions were partly very close to the pulp, no occurrence of postoperative discomfort, related to the ap-

plied adhesive or the filling material, was reported at baseline or in the meantime [9].

**Esthetic integration**

The optical properties have been optimized to meet the requirements of posterior restorations so that the shade system could be reduced to just three shades (IVA, IVB and IVW). The Ivocerin photoinitiator allows a slightly higher composite opacity compared to other bulk-fill materials due to its high quantum yield. Therefore, the translucency of the material could be adjusted in such a way that its optical properties, combined with the composite's favourable refraction index, blend in optimally with the tooth structure (particularly enamel) (Fig. 1). Thus, virtually invisible restorations can be achieved with average-sized fillings and in cases where the dentin is not discoloured.

The examiner rated 77 per cent of the restorations as "perfect" [7, 10] and 23 per cent as "good" (minor shade deviations). These results document the material's highly

developed "chameleon effect", especially given the fact that shade IVA was used in 80 per cent of the restorations.

Those cases in which the examiner identified slight shade deviations mostly involved either deep cavities or discoloured dentin. In such situations it is advisable to apply a thin layer of opaque material, e.g. Tetric EvoFlow Dentin, onto the dentin to obtain optimum esthetic integration (see also Figs 3 to 6).

**Polishability**

The clinical test also confirmed the material's excellent polishing properties. In 77 per cent of the restorations placed in this clinical in-house study (n=35), the surface gloss was rated to be "clinically perfect" according to FDI criteria after polishing with the composite polishing instrument OptraPol Next Generation (Ivoclar Vivadent). After a period of twelve months, the results were still the same for 69 per cent.

**Occurrence of enamel cracks**

Thick increments can cause cracks in the enamel due to their higher volume shrinkage and especially when good adhesion to the enamel is established [11, 12]. Therefore, the clinical evaluation focused on this aspect in particular.

The majority of the treated teeth (71 per cent) already showed enamel cracks before the placement of the restoration. The position of these cracks was noted to enable them to be distinguished from subsequently occurring cracks. Cold-light optical fibres (light probe from Lercher) were used for crack detection. However, no additional cracks were identified at baseline. This confirms that the tensions caused by Tetric EvoCeram Bulk Fill are not higher than those of other posterior compos-

FDI criteria /evaluation	Excellent	Good (excellent after correction)	Acceptable	Inadequate (repairable)	Unacceptable (new restoration necessary)
<b>Number (% of all restorations)</b>					
Postoperative sensitivity	35 (100 %)	0 (0 %)	0 (0 %)	0 (0 %)	0 (0 %)
Anatomic shape	34 (97 %)	1 (3 %)	0 (0 %)	0 (0 %)	0 (0 %)
Surface/gloss/pores	27 (77 %)	2 (6 %)	6 (17 %)	0 (0 %)	0 (0 %)
Esthetics	25 (71 %)	10 (29 %)	0 (0 %)	0 (0 %)	0 (0 %)
Surface discoloration	35 (100 %)	0 (0 %)	0 (0 %)	0 (0 %)	0 (0 %)
Patient satisfaction	34 (97 %)	1 (3 %)	0 (0 %)	0 (0 %)	0 (0 %)
Material fracture	35 (100 %)	0 (0 %)	0 (0 %)	0 (0 %)	0 (0 %)
Tooth integrity	34 (97 %)	1 (3 %)	0 (0 %)	0 (0 %)	0 (0 %)
Proximal contacts	33 (94 %)	1 (3 %)	1 (3 %)	0 (0 %)	0 (0 %)
<b>% of the entire restoration margin</b>					
Marginal discoloration	99.9 %	0.1 %	0 %	0 %	0 %
Marginal deficits	99.1 %	0.9 %	0 %	0 %	0 %
Submargination	99.9 %	0 %	0.1 %	0 %	0 %

Table: Results of the 12-month assessment according to FDI criteria (n=35)

ites which are applied in 2-mm layers. The results of this clinical study correspond with the preclinical data presented earlier. As it can be assumed that the water absorption of the composite compensates for the polymerization stress in a relatively short time [13], this baseline result can be considered reliable.

### Initial marginal quality

*Tetric EvoCeram Bulk Fill* restorations demonstrated excellent marginal quality at baseline as well as at the 12-month recall. With regard to marginal discolouration, submargination and marginal irregularities,

99 per cent of the examined marginal sections were rated as “clinically perfect”. A summary of the clinical results at the 12-month recall is displayed in the table above.

### Conclusions

The preliminary clinical data of *Tetric EvoCeram Bulk Fill* confirm the good results of the preclinical tests. Due to the photoinitiator *Ivocerin*, the material can be applied in 4-mm layers, with esthetic properties that meet the requirements of posterior restorations.

The increased curing depth of *Tetric EvoCeram Bulk Fill* gives the user more freedom

without loss of quality: dentists can restore cavities quickly with horizontal 4-mm increments or – perhaps even more importantly – they can adapt their layering technique to individual anatomic situations, their ergonomic preferences (e.g. the centripetal build-up technique for restorations) [14] or to esthetic aspects.

**Dr med. dent. Arnd Peschke,  
Schaan (Liechtenstein)**

*An extended version of this article has been published in: R&D Report, Ivoclar Vivadent AG, No. 19, July 2013.*



Fig. 3: Cavity preparation on tooth 16



Fig. 4: Tooth 16 after one layer of *Tetric EvoCeram Bulk Fill* has been used. The esthetic integration is optimal.



Fig. 5: Adjacent cavities in teeth 14 and 15



Fig. 6: Teeth 14 and 15 restored with one layer of *Tetric EvoCeram Bulk Fill IVA*

Figs 3 to 6: Arnd Peschke

## L I T E R A T U R E

- [1] Krejci I, Sparr D, Lutz F. Three-layer light hardening procedure with traditional composites for Black Class II restorations. *Quintessenz* 38 (1987) 1217-1229.
- [2] Lutz F, Krejci I, Luescher B, Oldenburg TR. Improved proximal margin adaptation of class II composite resin restorations by use of light-reflecting wedges. *Quintessence Int* (1986) 659-664.
- [3] Lutz F, Krejci I, Oldenburg TR. Elimination of polymerisation stresses at the margins of posterior composite resin restorations: a new restorative technique. *Quintessence Int* 17 (1986) 777-784.
- [4] Versluis A, Douglas WH, Cross M, Sakaguchi RL. Does an incremental filling technique reduce polymerization shrinkage stresses? *J Dent Res* 75 (1996) 871-878.
- [5] G. M. Lösche. Marginal adaptation of Class II composite fillings: guided polymerization vs reduced light intensity. *J Adhes Dent* 1 (1999) 31-39.
- [6] Hickel R, Peschke A, Tyas M, et al. FDI World Dental Federation: clinical criteria for the evaluation of direct and indirect restorations-update and clinical examples. *J Adhes Dent* 12 (2010) 259-272.
- [7] Hickel R, Roulet JF, Bayne S, et al. Recommendations for conducting controlled clinical studies of dental restorative materials. Science Committee Project 2/98--FDI World Dental Federation study design (Part I) and criteria for evaluation (Part II) of direct and indirect restorations including onlays and partial crowns. *J Adhes Dent* 9 Suppl 1 (2007) 121-147.
- [8] Peschke A. The semi-quantitative evaluation of restorations (SQUACE) – a tool to increase sensitivity als Teil des Symposiums: New Criteria for the Clinical Evaluation of Dental Restorations (Hickel, R., Sarrett, D., Peschke, A.). In: Hickel R, ed. IADR/AADR/CADR 87th General Session and Exhibition. Miami Beach, Florida, USA: 3. April 2009.
- [9] Ivoclar Vivadent Report, Nr. 19, Juli 2013. [10] Hickel R, Peschke A, Tyas M, et al. FDI World Dental Federation: clinical criteria for the evaluation of direct and indirect restorations-update and clinical examples. *Clin Oral Investig* 14 (2010) 349-366.
- [11] Batalha-Silva S, de Andrada MA, Maia HP, Magne P. Fatigue resistance and crack propensity of large MOD composite resin restorations: direct versus CAD/CAM inlays. *Dent Mater* 29 (2013) 324-331.
- [12] Clark DJ, Sheets CG, Paquette JM. Definitive diagnosis of early enamel and dentin cracks based on microscopic evaluation. *J Esthet Restor Dent* 15 (2003) 391-401; discussion 401.
- [13] Versluis A, Tantbirojn D, Lee MS, Tu LS, DeLong R. Can hygroscopic expansion compensate polymerization shrinkage? Part I. Deformation of restored teeth. *Dent Mater* 27 (2011) 126-133.
- [14] Bichacho N. The centripetal build-up for composite resin posterior restorations. *Pract Periodontics Aesthet Dent* 6 (1994) 17-23; quiz 24.



# The *future* of composite *technology*.

Available *now!*



## The *fast posterior composite*

- **Bulk placement** up to 4 mm due to Ivocerin<sup>®</sup>, the new light initiator
- **Low shrinkage** and low shrinkage stress for superior margins
- **Sculptable consistency**, extended working time under operatory light



## Tetric EvoCeram<sup>®</sup> Bulk Fill

The sculptable bulk-fill composite

Experience the *future of composites* featuring Ivocerin:  
[www.ivoclarvivadent.com/bulkfill\\_en](http://www.ivoclarvivadent.com/bulkfill_en)



[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

Ivoclar Vivadent AG  
Bendererstr. 2 | 9494 Schaan | Liechtenstein | Tel.: +423 / 235 35 35 | Fax: +423 / 235 33 60

ivoclar  
vivadent<sup>®</sup>  
passion vision innovation