# Cervitec® Plus

## EN Instructions for Use

– Protective varnish containing chlorhexidine

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ivoclar  
vivadent  
*clinical*

## Description

Cervitec® Plus is a varnish that protects sensitive dentin in the case of exposed cervicals. Chlorhexidine and thymol are the components in Cervitec Plus.

## Composition

Ethanol, water, acrylate copolymer, vinyl acetate copolymer and Chlorhexidine diacetate	1 %
Thymol	1 %

## Indication

- Protection of exposed root surfaces
- Treatment of hypersensitive cervicals

## Contraindication

If patients are known to be allergic to any of the ingredients in Cervitec Plus, the material should not be applied.

## Warning

In rare, isolated cases, chlorhexidine may cause allergic reactions.

## Side effects

In individual cases, contact with the mucous membrane may result in short-term, reversible irritation. Chlorhexidine may cause reversible discolouration of the tooth structure and restorations.

## Application



Cervitec Plus can be used for patients of all age groups and is professionally applied by dentists, dental hygienists or dental prophylaxis assistants. Generally, Cervitec Plus is applied every three months. If intensive treatment is required, however, the varnish may also be applied more frequently. One Free Stand Single Dose Cervitec Plus is enough for the treatment of one set of teeth. One tube is enough for about 20 to 30 sets of teeth.

## Step-by-step application procedure

1. Clean tooth surfaces thoroughly.
2. Dry with cotton rolls and air syringe.
3. a) **Single Dose:**  
Remove the foil from the Single Dose and apply the varnish directly from the Single Dose.  
b) **Tube:**  
Press out 3 drops into, for example, a dappen dish and reseat the tube.

Hold the tube in a vertical position with the cap facing upwards for opening and closing. For pressing out the contents, hold the opening vertically down.

4. Apply a thin coat of varnish by means of Vivadent® Applicator or suitable brush (diagram); in proximal areas, Cervitec Plus is applied with dental floss.
5. Let the varnish dry. The varnish can be dispersed with air.
6. Remove cotton rolls after 30 seconds.
7. Ask patients not to rinse their mouth.

## Tips for patients

After the application of Cervitec Plus, patients should not eat or drink for 1 hour.

## Special note

Operators can eliminate the risk of cross-contamination by using Cervitec Plus in the Single-Dose vessel for just one patient. The applicator is for single use only and cannot be sterilized.

## Storage

- Close the tube immediately after use
- Storage at temperature 2-8 °C / 36-46 °F
- Highly flammable. Keep away from sources of ignition
- Expiry date: see label / packaging

## Keep material out of the reach of children!

### For use in dentistry only!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes.

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