Bluephase® Style M8

EN Instructions for Use

Page 2

SV Bruksanvisning

Sidan 14

DA Brugsanvisning

Side 26

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Sivu 38

NO Bruksanvisning

Side 50

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Pagina 62

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CE 0123

For dental use only!

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist.

Made in Austria

Ivoclar Vivadent AG, FL-9494 Schaan / Liechtenstein


ivoclar
vivadent®
clinical

Table of Contents

1. Product overview	3
– List of parts	
– Indicators on the charging base	
– Indicators on the handpiece	
– Operating the light	
2. Safety	5
– Intended use	
– Indication	
– Signs & symbols	
– Safety notes	
– Contraindication	
3. Start-up	7
4. Operation	9
5. Maintenance and cleaning	10
6. What if ...	12
7. Procedure in case of repair	13
8. Product specifications	13

1. Product overview

List of parts


Indicators on the charging base

- Indicator is black = Battery is charged
- Indicator lights up in blue with different light intensity = battery is being charged

Indicators on the handpiece


With the handpiece switched on, the current charging status is shown on the handpiece as follows:

- **Green:**
Battery fully charged
(Curing capacity of approximately 20 minutes)
- **Orange:**
Battery weak
Time can still be set and a polymerization time of approximately 3 minutes is left. Place the light into the charging base as soon as possible.
- **Red:**
Battery completely discharged
The light can no longer be called up and the curing time can no longer be set.

Operating the light

Intuitive 2-button operation

2. Safety

Intended use

Bluephase Style M8 is an LED polymerization light that produces energy-rich blue light. It is used for the polymerization of light-curing dental materials immediately at the dental unit. The intended use also includes the observation of the notes and regulations in these Instructions for Use.

Indication

Bluephase Style M8 is suitable for polymerizing light-curing dental materials in the wavelength range of 430–490 nm. These materials include restoratives, bonding agents / adhesives, bases, liners, fissure sealants, temporary materials as well as luting materials for brackets and indirect restorations such as ceramic inlays.

Signs and symbols

Contraindication

Symbols on the curing light

Double insulation
(Device complies with safety class II)

Protection against electrical shock
(BF type apparatus)

Observe Instructions for Use

Observe Instructions for Use

Caution

The curing light must not be disposed with the normal domestic waste. Information regarding the disposal of the device can be found on the respective national Ivoclar Vivadent homepage.

Recyclable

AC voltage

DC voltage

Safety notes

Bluephase Style M8 is an electronic device and a medical product, which is subject to IEC 60601-1 (EN 60601-1) and EMC standards, as well as the 93/42/EEC Medical Device Directive. The appliance complies with the relevant EU regulations.

NRTL (Nationally Recognized Testing Laboratory)

CE 0123

The curing light has been shipped from the manufacturer in a safe and technically sound condition. In order to maintain this condition and to ensure risk-free operation, the notes and regulations in these Instructions for Use must be observed. To prevent damage to equipment and risks for patients, users and third parties, the following safety instructions must be observed:

Contraindications

Materials whose polymerization is activated outside the wavelength range of 430–490 nm. If you are not sure about certain products, please ask the manufacturer of the corresponding material.

Use without light probe.

This unit should not be used or charged near highly flammable or ignitable substances.

Portable and mobile high-frequency communication devices may interfere with medical equipment. The use of mobile phones is not allowed during operation.

Usage and liability

- Bluephase Style M8 must only be employed for the intended use. Any other uses are contraindicated. Do not touch defective, open devices. Liability cannot be accepted for damage resulting from misuse or failure to observe the Instructions for Use.
- The user is responsible for testing Bluephase Style M8 for its use and suitability for the intended purposes. This is particularly important if other equipment is used in the immediate vicinity of the curing light at the same time.
- Use only original spare parts and accessories from Ivoclar Vivadent (see Accessories). The manufacturer does not accept any liability for damage resulting from the use of other spare parts or accessories.
- The light probe is an applied part and may warm up to a maximum of 45 °C (113 °F) at the interface to the handpiece during operation.

Operating voltage

Before switching on, make sure that

- a) the voltage indicated on the rating plate complies with the local power supply;
- b) the unit has acquired the ambient temperature.

If the battery is handled separately, contact with patients or third parties must be prevented. Do not touch the exposed contacts of the connection plug (power pack).

Rating plate on the charging base

Assumption of impaired safety

If it has to be assumed that safe operation is no longer possible, the power must be disconnected and the battery removed to avoid accidental operation. This may be the case, for example, if the device is visibly damaged or no longer works correctly. A complete disconnection from the power supply is only ensured when the power cord is disconnected from the power source.

Eye protection

Direct or indirect exposure of the eyes must be prevented. Prolonged exposure to the light of the curing device is unpleasant to the eyes and may result in injury. Therefore, using the supplied anti-glare cones is recommended. Individuals who are generally sensitive to light, who take photosensitizing drugs, have undergone eye surgery, or people who work with the apparatus or in its vicinity for long periods of time should not be exposed to the light of this device and should wear protective goggles (orange) that absorb light below 515 nm.

Battery

Caution: Use only original spare parts, particularly Ivoclar Vivadent batteries and charging bases. Do not short out the battery and do not touch the battery contacts. Do not store at temperatures above 40 °C / 104 °F (or 60 °C / 140 °F for a short period). Always store batteries charged. The storage period must not exceed 6 months. May explode if disposed of in fire.

Please note that lithium-polymer batteries may react with explosion, fire, and smoke development if handled improperly or if they are mechanically damaged. Damaged lithium-polymer batteries must no longer be used.

The electrolytes and electrolyte fumes released during explosion, fire and smoke development are toxic and corrosive. In case of accidental contact with the eyes or skin, immediately wash with copious amounts of water. Avoid inhalation of vapour. In case of indisposition, see a physician immediately.

3. Start-up

Heat development

As it is the case with all high-performance lights, the high light intensity results in a certain heat development. Prolonged exposure of areas near the pulp and soft tissues may result in irreversible damage. Therefore, the recommended curing times have to be observed, especially in areas near the pulp (adhesives 10 seconds). Uninterrupted curing times of more than 40 seconds on the same tooth surface, as well as direct contact with the gingiva, oral mucous membrane or skin have to be prevented. Polymerize indirect restorations at intermittent intervals of 20 seconds each or use external cooling with an air stream.

Check the delivery for completeness and any possible transportation damage (see delivery form). If parts are damaged or missing, contact your Ivoclar Vivadent representative.

Charging base

Before you switch on the device, make sure that the voltage mentioned on the rating plate complies with your local power supply.

Connect the power cord with the power supply and the power pack. The Power On indicator on the power pack lights up in green and the illuminated ring on the charging base flashes briefly.

Handpiece

Remove the handpiece from its packaging and clean the light probe and the handpiece (see chapter Maintenance and cleaning). Slightly turn the light probe to remove and reposition it.

After that, mount the anti-glare cone on the light probe.

Battery

Fully charge the battery before the first use! If the battery is fully charged, it features a curing capacity of approximately 20 minutes. Slide the battery straight into the handpiece until you hear and feel it click into place.

Gently place the handpiece into the corresponding rest on the charging base. If a hygiene sleeve is used, please remove it before you charge the battery. If possible, use the light always with a fully charged battery. This will prolong the service life. It is therefore recommended to place the handpiece into the charging base after each patient. If the battery is fully discharged, the charging time is 2 hours.

Since the battery is a consumable, it has to be replaced after its typical life cycle has expired after approximately 2.5 years. See battery label for the age of the battery.

1106000944	dd/mm/yy
#637 692	30/01/13
+	-

Battery charging status

The respective charging status is indicated on the handpiece as described on Page 4.

4. Operation


Disinfect contaminated surfaces of the curing light as well as the light probe and anti-glare cone before each use. The light probe may also be autoclaved. Furthermore, make sure that the stipulated light intensity permits adequate polymerization. For this purpose, check the light probe for contamination and damage and check the light intensity at regular intervals.

Selecting the curing time

The curing time can be individually set using the time selection button. Users may choose between 10, 15, 20 and 30 seconds.

Observe the Instructions for Use of the material applied when selecting the curing time. The curing recommendations for composite materials apply to all shades and, if not mentioned otherwise in the Instructions for Use of the relevant material, to a maximum layer thickness of 2 mm. Generally, these recommendations apply to situa-

tions where the emission window of the light probe is placed directly over the material to be polymerized. Increasing the distance between the light source and the material will require the curing time to be extended accordingly. For instance, if the distance to the material is 9 mm, the effective light output is reduced by approx. 50%. In this case, the recommended curing time has to be doubled.

Curing times		Light intensity 800 mW/cm ² ± 10%
Restorative materials	Composites • 2 mm ¹⁾ IPS Empress Direct / IPS Empress Direct Flow Tetric EvoCeram / Tetric EvoFlow Tetric / Tetric Basic White / Tetric N-Ceram / Tetric N-Flow	15 s
	• 4 mm ²⁾ Tetric EvoCeram Bulk Fill / Tetric N-Ceram Bulk Fill / Tetric Basic White	
	Heliomolar / Heliomolar HB / Heliomolar Flow / all conventional composites ¹⁾	20 s
	Compomers ³⁾ Compoglass F / Compoglass Flow	40 s
Indirect restorations / cementation materials	Variolink II Base ⁴⁾ / Variolink Veneer / Variolink N Base ⁴⁾ / Variolink N Clear Veneer / Variolink Ultra ⁵⁾ / Dual Cement ⁵⁾ / Variolink II ⁵⁾	per mm ceramic: 10 s per surface
Adhesives	AdheSE / AdheSE One F Excite F / Excite F DSC Tetric N-Bond / Tetric N-Bond Self-Etch Te-Econom Bond Heliobond (Syntac)	10 s
Temporary materials	System.link / Telio CS Link	15 s per surface
	System.inlay/onlay / Telio CS Inlay/Onlay	10 s
	Fermit / Fermit N	10 s
	Telio Stains Telio Add-On Flow	15 s 20 s
Miscellaneous	Helioseal / Helioseal F / Helioseal Clear Heliosit Orthodontic Vivaglass Liner	20 s
	IPS Empress Direct Color IPS Empress Direct Opaque MultiCore Flow / Multicore HB	30 s 30 s 40 s

1) Applies to a maximum layer thickness of 2 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (which may be the case e.g. with dentin shades)

2) Applies to a maximum layer thickness of 4 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (which may be the case, e.g. with dentin shades)

3) Applies to a maximum layer thickness of 3 mm

4) Applies to light-curing (exclusive use of base paste)

5) Applies to dual-curing

5. Maintenance and cleaning

Cure Memory function

The latest time setting used is automatically saved.

Start

The light is switched on by means of the start button. During curing a signal beep sounds every 10 seconds. Once the selected curing time has elapsed, the curing program is automatically terminated. If desired, the light can be switched off before the set curing time has elapsed by pressing the start button again.

Acoustic signals

Acoustic signals can be heard for the following functions:

- Start (Stop)
- Every 10 seconds
- Curing time change
- Error message

Light intensity

The light intensity is maintained at a consistent level during operation. If the supplied 10 mm light probe is used, the light intensity has been calibrated to $800 \text{ mW/cm}^2 \pm 10\%$.

If another light probe as the one supplied is used, the indicated light intensity may be directly affected.

In a light probe with parallel walls (10 mm), the diameter of the light entry and the diameter of the light emission window are the same. When using focussing light probes ($10 > 8 \text{ mm}$) or Pin-Point light probes ($6 > 2 \text{ mm}$), the diameter of the light entry is larger than that of the light emission window. The incident blue light is thus bundled to a smaller area. As a result, the emitted light intensity is increased.

Pin-Point light probes are suitable for spot-on polymerization, e.g. to tack veneers in place before removing excess material. For complete polymerization, the light probe must be changed.

For reasons of hygiene, it is recommended to use a disposable protective sleeve for each patient. Make sure to fit the protective sleeve closely to the light probe. Disinfect contaminated areas of the device and anti-glare cones (FD 366/Dürr Dental, Incidin Liquid/Ecolab) and sterilize the light probe before each use if not using disposable protective sleeves. Make sure that no liquids or other foreign substances enter the handpiece, charging base and particularly the power pack during cleaning (risk of electrical shock). Disconnect the charging base from the power source when cleaning it.

Housing

Wipe the handpiece and handpiece holder with a customary aldehyde-free disinfecting solution. Do not clean with highly aggressive disinfecting solutions (e.g. solutions based on orange oil or solutions with an ethanol content of more than 40%), solvents (e.g. acetone), or pointed instruments, which may damage or scratch the plastic. Clean contaminated plastic parts with a soap solution.

Light probe

Before you clean and/or disinfect the light probe, it should be pretreated. This applies to both manual and machine cleaning and disinfection.

Pretreatment

- Remove substantial contamination immediately after use and no later than 2 hours after use. For this purpose, thoroughly rinse the light probe under running water (for at least 10 seconds). You may also use a suitable aldehyde-free disinfecting solution to prevent blood residue from becoming irremovable.
- To remove contamination manually, use a soft brush or a soft cloth. Partially polymerized composite can be removed with alcohol and a plastic spatula, if necessary. Do not use sharp or pointed instruments. They may scratch the surface.

Cleaning and disinfection

For cleaning, immerse the light probe in a cleaning solution in such a way that it is sufficiently covered with liquid (ultrasound or careful brushing with a soft brush can support the effect). A neutral-enzymatic cleaning agent is recommended.

When cleaning and disinfecting, please make sure that the agents used are free of

- organic, mineral, and oxidizing acids (the minimum admissible pH value is 5.5);
- bases (the maximum admissible pH value is 8.5);
- oxidizing agents (e.g. hydrogen peroxides).

Afterwards, remove the light probe from the solution and thoroughly rinse it under running water (for at least 10 seconds). Cleaning in a thermal disinfectant is an effective alternative.

Sterilization

Thorough cleaning and disinfecting is imperative to ensure that the subsequent sterilization is effective. Use only autoclave sterilization for this purpose. Sterilization time (exposure time at sterilization temperature) is at least 4 minutes at 134°C (273°F); pressure should be 2 bar (29 psi). Dry the sterilized light probe using either the special drying program of your steam autoclave or hot air. The light probe has been tested for up to 200 sterilization cycles.

After that, check the light probe for damage. Hold the light probe against light. If individual segments appear black, glass fibres are broken. If this is the case, replace the light probe with a new one.

Disposal

The device must not be disposed of in normal domestic waste. Dispose unserviceable batteries and polymerization lights according to the corresponding legal requirements in your country. Batteries must not be incinerated!

6. What if ...?

Indicator	Causes	Error rectification
<p>all LEDs orange</p> 	<p>The device is overheated.</p>	<p>Allow the device to cool down and try again after a certain time. If the error reappears, please contact your dealer or your local Service Center.</p>
<p>all LEDs red</p> 	<p>Electronic component of the handpiece defective</p>	<p>Remove and reinsert the battery. If the error reappears, please contact your dealer or your local Service Center.</p>
<p>The charging base is not illuminated during charging</p>	<ul style="list-style-type: none"> – Power pack not connected or defective – Battery fully charged 	<p>Check if the power pack is correctly positioned in the charging base or if the power pack is connected to the power supply by means of the power cord (display on the power pack lights up in green if it works correctly).</p>
<p>No indicator activity of the device with the battery in place</p>	<p>Battery empty</p>	<p>Place the device in the charging base and charge for at least 2 hours.</p>
	<p>Battery contacts contaminated</p>	<p>Remove battery and clean the battery contacts.</p>

7. Procedure in case of repair

The warranty period for Bluephase Style M8 is 2 years from the date of purchase (battery 1 year). Malfunctions resulting from faulty material or manufacturing errors are repaired free of charge during the warranty period. The warranty does not provide the right to recover any material or non-material damage other than the ones mentioned. The apparatus must only be used for the intended purposes. Any other uses are contraindicated. The manufacturer does not accept any liability resulting from misuse and warranty claims cannot be accepted in such cases. This is particularly valid for:

- damage resulting from improper handling, especially incorrectly stored batteries (see chapter Product specifications: Transportation and storage conditions);
- damage to components resulting from wear under standard operating conditions (e.g. battery);
- damage resulting from external influences, e.g. blows, drop to the floor;
- damage resulting from incorrect set-up or installation;
- damage resulting from connecting the unit to a power supply whose voltage and frequency do not comply with the ones stated on the rating plate;
- damage resulting from improper repairs or modifications that have not been carried out by certified centers.

In case of a claim under warranty, the complete apparatus (handpiece, charging base, power cord and power pack) must be returned, carriage paid, to the dealer or directly to Ivoclar Vivadent, together with the purchase document. Use the original packaging with the corresponding cardboard inserts for transportation. Repair work may only be carried out by a certified Ivoclar Vivadent Service Center. In case of a defect that cannot be rectified, please contact your dealer or your local Service Center (see addresses on the reverse side). A clear description of the defect or the conditions under which the defect occurred will facilitate locating the problem. Please enclose this description when returning the apparatus.

8. Product specifications

Light source	LED
Wavelength range	430–490 nm
Light intensity	800 mW/cm ² ± 10%
Operation	3 min. on / 7 min. off (intermittent)
Light probe	10 mm, black, disinfectable and autoclavable
Anti-glare cones	autoclavable
Signal transmitter	acoustic after 10 seconds and every time the start button / time selection button is pressed
Dimensions of the handpiece (without light probe)	L = 180 mm W = 30 mm H = 30 mm
Weight of the handpiece	120 g (including battery and light probe)
Operating voltage handpiece	3.7 VDC with battery
Operating voltage Charging base	5 VDC
Power pack	Input: 100–240 VAC, 50-60 Hz 310 mA Output: 5 VDC / 2 A Manufacturer Friwo Type: FW7401M/05
Weight of the power pack	165 g
Operating conditions	Temperature +10 °C to +35 °C (+50 °F to +95 °F) Relative humidity 30% to 75% Ambient pressure 700 hPa to 1060 hPa
Dimensions of the charging base	D = 125 mm, H = 70 mm
Weight of the charging base	195 g
charging time	approximately 2 hours (with the battery empty)
Power supply of the handpiece	Li-Po battery (approx. 20 min. with a new, fully charged battery)
Transportation and storage conditions	Temperature –20 °C to +60 °C (–4 °F to +140 °F) Relative humidity 10% to 75% Ambient pressure 500 hPa to 1060 hPa The curing light has to be stored in a closed, roofed room and must not be exposed to hard shocks. Battery: – Do not store at temperatures above 40 °C / 104 °F (or 60 °C / 140 °F for a short period). Recommended storage temperature 15–30 °C (59–86 °F). – Always store the battery charged and not for longer than 6 months.
Delivery form	1 Handpiece 1 Battery 1 Charging base with power cord and power pack 1 Light probe 10 mm, black 3 Anti-glare cones 1 Set of Instructions for Use

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Appendix

Guidance and manufacturer's declaration – electromagnetic emission

The Bluephase Style M8 is intended for use in the electromagnetic environment specified below. The customer or the user of the Bluephase Style M8 should assure that it is used in such an environment.

Emision test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Bluephase Style M8 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Bluephase Style M8 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	N/A
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	N/A

Guidance and manufacturer's declaration – electromagnetic immunity

The Bluephase Style M8 is intended for use in the electromagnetic environment specified below. The customer or the user of the Bluephase Style M8 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of typical commercial or dental environment.
Surge IEC 61000-4-5	± 1 kV line - line ± 2 kV line - earth	± 1 kV line - line no prot. earth	Mains power quality should be that of typical commercial or dental environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of typical commercial or dental environment. If the user of the Bluephase Style M8 requires continued operation during power mains interruptions, it is recommended that the Bluephase Style M8 be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic Field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or dental environment.

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity

The Bluephase Style M8 is intended for use in the electromagnetic environment specified below. The customer or the user of the Bluephase Style M8 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should not be used closer to any part of the Bluephase Style M8, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			Recommended separation distance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	10 V	$d = 0.35 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m	$d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz $d = 0.70 \sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environ-

ment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Bluephase Style M8 is used exceeds the applicable RF compliance level above, the Bluephase Style M8 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Bluephase Style M8.

^b Over the frequency range 150 kHz to 80 MHz, field strength should be less than 10 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the "Bluephase Style M8"

The Bluephase Style M8 is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Bluephase Style M8 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Bluephase Style M8 as recommended below, according to the maximum output power of the communication equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 0.35 \sqrt{P}$	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
0.01	0.04	0.04	0.07
0.1	0.13	0.13	0.22
1	0.40	0.40	0.70
10	1.3	1.3	2.2
100	4.0	4.0	7.0

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 4: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Ivoclar Vivadent – worldwide

Ivoclar Vivadent AG

Bendererstrasse 2
9494 Schaan
Liechtenstein
Tel. +423 235 35 35
Fax +423 235 33 60
www.ivoclarvivadent.com

Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive
P.O. Box 367
Noble Park, Vic. 3174
Australia
Tel. +61 3 9795 9599
Fax +61 3 9795 9645
www.ivoclarvivadent.com.au

Ivoclar Vivadent GmbH

Tech Gate Vienna
Donau-City-Strasse 1
1220 Wien
Austria
Tel. +43 1 263 191 10
Fax: +43 1 263 191 111
www.ivoclarvivadent.at

Ivoclar Vivadent Ltda.

Alameda Caiapós, 723
Centro Empresarial Tamboré
CEP 06460-110 Barueri – SP
Brazil
Tel. +55 11 2424 7400
Fax +55 11 3466 0840
www.ivoclarvivadent.com.br

Ivoclar Vivadent Inc.

1-6600 Dixie Road
Mississauga, Ontario
L5T 2Y2
Canada
Tel. +1 905 670 8499
Fax +1 905 670 3102
www.ivoclarvivadent.us

Ivoclar Vivadent Shanghai

Trading Co., Ltd.
2/F Building 1, 881 Wuding Road,
Jing An District
200040 Shanghai
China
Tel. +86 21 6032 1657
Fax +86 21 6176 0968
www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

Calle 134 No. 7-B-83, Of. 520
Bogotá
Colombia
Tel. +57 1 627 3399
Fax +57 1 633 1663
www.ivoclarvivadent.co

Ivoclar Vivadent SAS

B.P. 118
F-74410 Saint-Jorioz
France
Tel. +33 4 50 88 64 00
Fax +33 4 50 68 91 52
www.ivoclarvivadent.fr

Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2
D-73479 Ellwangen, Jagst
Germany
Tel. +49 7961 889 0
Fax +49 7961 6326
www.ivoclarvivadent.de

Wieland Dental + Technik

GmbH & Co. KG
Schwenninger Strasse 13
D-75179 Pforzheim
Germany
Tel. +49 7231 3705 0
Fax +49 7231 3579 59
www.wieland-dental.com

Ivoclar Vivadent Marketing (India)

Pvt. Ltd.
503/504 Raheja Plaza
15 B Shah Industrial Estate
Veera Desai Road, Andheri (West)
Mumbai, 400 053
India
Tel. +91 22 2673 0302
Fax +91 22 2673 0301
www.ivoclarvivadent.in

Ivoclar Vivadent s.r.l.

Via Isonzo 67/69
40033 Casalecchio di Reno (BO)
Italy
Tel. +39 051 6113555
Fax +39 051 6113565
www.ivoclarvivadent.it

Ivoclar Vivadent K.K.

1-28-24-4F Hongo
Bunkyo-ku
Tokyo 113-0033
Japan
Tel. +81 3 6903 3535
Fax +81 3 5844 3657
www.ivoclarvivadent.jp

Ivoclar Vivadent Ltd.

12F W-Tower, 1303-37
Seocho-dong, Seocho-gu,
Seoul 137-855
Republic of Korea
Tel. +82 2 536 0714
Fax +82 2 596 0155
www.ivoclarvivadent.co.kr

Ivoclar Vivadent S.A. de C.V.

Av. Insurgentes Sur No. 863,
Piso 14, Col. Napoles
03810 México, D.F.
México
Tel. +52 55 5062 1000
Fax +52 55 5062 1029
www.ivoclarvivadent.com.mx

Ivoclar Vivadent BV

De Fruittuin 32
2132 NZ Hoofddorp
Netherlands
Tel. +31 23 529 3791
Fax +31 23 555 4504
www.ivoclarvivadent.com

Ivoclar Vivadent Ltd.

12 Omega St, Rosedale
PO Box 303011 North Harbour
Auckland 0751
New Zealand
Tel. +64 9 914 9999
Fax +64 9 914 9990
www.ivoclarvivadent.co.nz

Ivoclar Vivadent Polska Sp. z o.o.

Al. Jana Pawla II 78
00-175 Warszawa
Poland
Tel. +48 22 635 5496
Fax +48 22 635 5469
www.ivoclarvivadent.pl

Ivoclar Vivadent Marketing Ltd.

Prospekt Andropova 18 korp. 6/
office 10-06
115432 Moscow
Russia
Tel. +7 499 418 0300
Fax +7 499 418 0310
www.ivoclarvivadent.ru

Ivoclar Vivadent Marketing Ltd.

Qlaya Main St.
Siricon Building No.14, 2nd Floor
Office No. 204
P.O. Box 300146
Riyadh 11372
Saudi Arabia
Tel. +966 11 293 8345
Fax +966 11 293 8344
www.ivoclarvivadent.com

Ivoclar Vivadent S.L.U.

C/ Ribera del Loira nº 46, 5ª planta
28042 Madrid
Spain
Tel. + 34 913 757 820
Fax + 34 913 757 838
www.ivoclarvivadent.es

Ivoclar Vivadent AB

Dalvägen 14
S-169 56 Solna
Sweden
Tel. +46 8 514 939 30
Fax +46 8 514 939 40
www.ivoclarvivadent.se

Ivoclar Vivadent Liaison Office

: Tesviyike Mahallesi
Sakayik Sokak
Nisantas' Plaza No:38/2
Kat:5 Daire:24
34021 Sisli – Istanbul
Turkey
Tel. +90 212 343 0802
Fax +90 212 343 0842
www.ivoclarvivadent.com

Ivoclar Vivadent Limited

Ground Floor Compass Building
Feldspar Close
Warrens Business Park
Enderby
Leicester LE19 4SE
United Kingdom
Tel. +44 116 284 7880
Fax +44 116 284 7881
www.ivoclarvivadent.co.uk

Ivoclar Vivadent, Inc.

175 Piney Drive
Amherst, N.Y. 14228
USA
Tel. +1 800 533 6825
Fax +1 716 691 2285
www.ivoclarvivadent.us