ProArt CAD Model beige

EN Instructions for Use

- Model disc for the CAD/CAM technology

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Product description

ProArt CAD Model beige are discs made from polyurethane for the fabrication of models using the CAD/CAM technique.

The discs can be processed using a suitable milling machine (please observe the manufacturer's instructions).

Indications

The ProArt CAD Model beige disc is used for the CAD/CAM-supported production of acrylic models. ProArt CAD Model beige is not suitable for intraoral use!

Contraindications

The ProArt CAD Model beige is not approved for use in the patient's mouth.

Processing with CAD/CAM

ProArt CAD Model beige can be processed in CAD/CAM milling systems with a standard holder of 98.5 mm. Only use milling tools recommended for this milling system when processing the material.

Processing after milling

Finishing

The milled model bases and model dies can be separated from the milling disc using cross-cut tungsten carbide burs suitable for use with acrylic or suitable cutting discs.

Storage

- Protect against direct sunlight.
- Store out of the reach of children.
- Protect against moisture.

Disposal

ProArt CAD Model beige consists of polyurethane. Residues can be disposed of in the general household waste.

This material has been developed solely for use in dentistry. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. These regulations also apply if the material is mixed or used in conjunction with products of other manufacturers.

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