Intended use**Intended purpose**

Protection of the tooth structure, bite adjustments, templates

Use

For dental use only.

Description

- Therapeutic appliances to correct TMJ problems and make occlusal adjustments
- Occlusal splints
- Drill templates

Technical data

		Specifications	Typical average value
Flexural strength ¹⁾	[MPa]	≥ 100	135
Flexural modulus ²⁾	[MPa]	≥ 2800	3106
Ball indentation hardness ³⁾	[MPa]	≥ 140	176
Water absorption ¹⁾	[µg/mm ³]	≤ 40	21
Solubility ¹⁾	[µg/mm ³]	≤ 7.5	0.0018

¹⁾ method based on EN ISO 10477:2004²⁾ determined by bending tests³⁾ internal method**Indication**

Dental dysfunctions and parafunctions, partial edentulism.

Contraindication

If a patient is known to be allergic to PMMA, ProArt CAD Splint should not be used.

Limitations of use

Not suitable for use as a drill template without a guiding sleeve.

Composition

Polymethyl methacrylate

ProArt CAD Splint

ProArt CAD Splint

EN Instructions for Use - PMMA disc for the CAD/CAM technology	SL Navodila za uporabo - Disk PMMA za tehnologijo CAD/CAM	SQ Udhëzime përdorimi - Disk akriliku për teknologji CAD/CAM
DE Gebrauchsinformation - PMMA-Scheibe für die CAD/CAM-Technik	HR Upute za uporabu - PMMA disk za CAD/CAM tehniku	RO Instrucţiuni de utilizare - Disc PMMA pentru tehnica CAD/CAM
FR Mode d'emploi - Disque PMMA pour la technologie CAD/CAM	CS Návod k použití - PMMA disk pro CAD/CAM techniku	UK Інструкція з використання - Диск PMMA для технології CAD/CAM
IT Istruzioni d'uso - Disco PMMA per la tecnica CAD/CAM	SK Návod na používanie - Disk PMMA pre technológiu CAD/CAM	ET Kasutamishüend - PMMA ketas CAD/CAM tehnoloogia jaoks
ES Instrucciones de uso - Discos de PMMA para la tecnología CAD/CAM	HU Használati utasítás - Polimetil-metakrilát (PMMA) korong CAD/CAM technológiához	LV Lietošanas instrukcija - PMMA diska CAD/CAM tehnoloģijai
PT Instruções de Uso - Disco de PMMA para tecnologia CAD/CAM	SR Упутство за употребу - PMMA диск за CAD/CAM технологију	LT Naudojimo instrukcija - CAD/CAM technologijai skirtas PMMA diskas
SV Bruksanvisning - PMMA-disk för CAD/CAM teknik	MK Упатство за употреба - PMMA диск за CAD/CAM техника	
	BG Инструкции за употреба - Диск от PMMA за CAD/CAM технология	
DA Brugsanvisning - PMMA skive til CAD/CAM teknik		
FI Käyttöohjeet - PMMA-kielko CAD/CAM-tekniikkaan		
NO Bruksanvisning - PMMA-skive til CAD/CAM-teknikk		
NL Gebruiksaanwijzing - PMMA disc voor de CAD/CAM-technologie		
EL Οδηγίες Χρήσεως - Δίσκος PMMA για τεχνική CAD/CAM		
TR Kullanma Talimatı - CAD/CAM teknolojisini için PMMA disk		
RU Инструкции по применению - Диск PMMA для технологии CAD/CAM		
PL Instrukcja stosowania - Dysk z krytlu do technologii CAD/CAM		

Application

Processing with CAD/CAM equipment

- ProArt CAD Splint discs can be processed in CAD/CAM milling systems with a standard holder of 98.5 mm. Only use milling tools recommended for the respective milling system when processing the material. Please note that the improper use of milling tools, polishing brushes, steam cleaners and water baths can lead to the overheating of the object, which may result in damage to the material.

Processing after the milling procedure

Finishing

- Once milling is completed, use tungsten carbide burs to separate the object from the disc.
- Any finishing work or adjustments to the shape must be carried out using cross-cut tungsten carbide burs. When doing so, pay particular attention to the margins, occlusal surfaces and inner surfaces. If a patient try-in is carried out, the object should subsequently be disinfected using an agent that is suitable for PMMA. Please note that strong alcohol-based disinfectants are not suitable for this material.

Polishing

- When polishing the objects, pay particular attention to the margins, occlusal surfaces and inner surfaces. For pre-polishing, use rubber polishers and silicone wheels with different abrasive values from coarse to fine. Polish to a high gloss with a goat-hair brush, cotton or leather polishing buff and universal polishing paste.

Individualization

- Splints made of ProArt CAD Splint can be individualized by adding light-curing or cold-curing acrylic materials.
- a) With SR Nexco® (light-curing)*
 - Prepare the area to be built up by grinding and/or abrasive blasting (Al_2O_3 , 80-100 μm grit at 1-2 bar), clean with a steam cleaner and then dry with oilfree compressed air. Then apply SR® Connect according to the relevant instructions for use, cure, and finally build up with SR Nexco materials (see the SR Nexco Instructions for Use).
- b) With cold-curing acrylic resin*
 - Prepare the area to be built up by grinding and/or abrasive blasting (Al_2O_3 , 80-100 μm grit at 1-2 bar), clean with a steam cleaner and then dry with oilfree compressed air. For the processing of cold-curing acrylic resin, see the respective instructions for use.

Important processing restrictions

Adhering to the following framework conditions is essential to ensure the successful processing of ProArt CAD Splint:

- The instructions of the implant manufacturer must be observed when fabricating drill templates.
- The discs must be milled using a compatible CAD/CAM system.
- For individualization only use materials which are approved or recommended.

Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclarvivadent.com, and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclarvivadent.com).

Warnings

- Failure to observe the stipulated limitations of use and processing instructions may lead to failure.
- Do not inhale grinding dust during finishing.
- Use suction equipment and a face mask.

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Shelf life and storage

Protect from sunlight.

Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions.

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Ivoclar Vivadent Marketing Ltd.
Ivoclar Vivadent AB
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.