# ProArt CAD Wax blue

## EN Instructions for Use

- Wax disc for the CAD/CAM technology

## DE Gebrauchsinformation

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### Product description

ProArt CAD Wax blue are wax discs for the CAD/CAM fabrication of restorations for the casting technique. The discs can be processed using a suitable milling machine (please observe the manufacturer's instructions).

### Indications

The production of burn-out restorations for the casting technique

### Contraindications

- If the patient is known to be allergic to any of the components
- Not approved for use directly in the patient's mouth

### Important processing restrictions

Not suitable for the press technique as the investment ring may develop cracks.

### Notes regarding the design

For the design, please observe the manufacturer's instructions for the relevant alloy.

### Processing with CAD/CAM

ProArt CAD Wax blue can be processed in CAD/CAM milling systems with a standard holder of 98.5 mm. Only use milling tools recommended for this milling system when processing the material. The ambient working temperature should not exceed 35 °C / 95 °F.

### Processing after milling

#### Finishing

After milling, separate the processed object from the disc using a suitable cross-cut tungsten carbide bur. Any finishing work or adjustments to the shape must

be carried out using cross-cut tungsten carbide burs. The material can be enhanced using conventional modelling wax.

### Attaching the sprues

The rules how to choose and position the sprues correctly depend on the alloy used. Detailed information concerning the alloys can be found in the corresponding manufacturer's instructions.

### Investing

Line the inside of the casting ring with ceramic liner. Mix the investment material and invest under vacuum. Please observe the manufacturer's instructions of the investment material and the alloy to be used.

### Storage

- Protect from direct sunlight.
- Store out of the reach of children!
- Store at 2 – 28 °C (36 – 82 °F).

### Disposal

ProArt CAD Wax blue consists of pure synthetic wax. Residues can be used for other customary wax applications. The wax can be disposed of in the general household waste.

This material has been developed solely for use in dentistry. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. These regulations also apply if the material is mixed or used in conjunction with products of other manufacturers.

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