


Removable prosthetics

Crossbite set-up ^{with}
Typ posterior teeth



Contents

▶ Introduction	2
▶ Definition of a crossbite	3
▶ Models in the articulator – Crossbite analysis	4
▶ Modification of the molars	6
▶ Set-up – Maxillary anterior teeth	7
▶ Set-up – Mandible	8
▶ Set-up – Maxillary posterior teeth	9
▶ In function – Comparison between a normal bite and crossbite set-up	10

▶ Introduction

Edentulous patients may have different types of occlusion. Often, due to skeletal conditions or atrophies, it is a crossbite.

This guide shows how a crossbite can be identified in an edentulous jaw and how Typ teeth can be modified for an appropriate set-up. This procedure applies to all Typ posterior teeth from Ivoclar Vivadent.




► Definition of a crossbite

The orthodontist Dr Edward H. Angle (1855-1930) divided misalignments of the dentulous jaw into three classes. The assessment is based on the occlusal relation of the first molars in the maxilla and mandible. The buccal cusps assume the function of the working cusps.

- Class I: Neutroclusion
- Class II: Distocclusion
- Class III: Mesioclusion

The crossbite belongs to Angle class III (mesioclusion). Rather than the palatal cusp, the buccal cusp of the upper molar occludes with the central fossa of the lower molar.



For the edentulous jaw, a further distinction is made within Angle class III:

1. One-sided crossbite due to atrophy of the maxilla and mandible
2. Crossbite on both sides with protrusion of the mandible (prognathism)

► Models in the articulator – Crossbite analysis

The working models are aligned with each other in the usual manner using the wax bite or the intraoral registration, and then articulated. Subsequently, the model analysis of both jaws is carried out. The midline of the jaw is transferred to the model edge for the subsequent working steps.

Comprehensive model analysis

Marking:

Centre of the incisive papilla

Relevance:

- Anatomical midline of the upper jaw
- Labial positioning of the central incisors

Marking:

Deepest point of the vestibule

Relevance:

Starting point for measuring the vertical dimension and the incisal height of the central incisors

Marking:

First large pair of rugae

Relevance:

Labial positioning of the canine teeth at the tip of the rugae

Marking:

Palatal suture (raphe palatina), anatomical midline

Relevance:

Reference point for the transversal symmetry of the anterior set-up

Marking:

Post dam

Relevance:

Posterior palatal limit of the denture base

Marking:

Crest of the alveolar ridge

Relevance:

Provides orientation for a static set-up

Marking:

Distal half of the retromolar pad (trigonum retromolare)

Relevance:

- Positioning of the setting up template on the dorsal aspect (corresponds to the height of the occlusal plane)
- Dorsal positioning of the lateral wings of the horizontal guide

Marking:

Lingual limit of the retromolar pads

Relevance:

Pound's Line, consideration for tongue space

Marking:

Model midline transferred from maxilla model, anatomical midline

Relevance:

- Bilateral orientation for the anterior set-up
- Position of the symphysis fork of the horizontal guide

Marking:

Crest of the alveolar ridge

Relevance:



Course of the central fissure of the posterior teeth (static positioning)

Marking:

Deepest point of the vestibule

Relevance:

Starting point for measuring the total vertical dimension



If the upper dental arch is narrower than the lower arch by more than one premolar width, a crossbite set-up in the posterior area is indicated. In the case shown, a crossbite situation is present on the left due to atrophy, while a normal bite is present on the right.



► Modification of the molars

In the crossbite situation, the buccal cusps of the maxillary molars slide into the central fossa of the mandibular molars and assume the working function.



The buccal cusps of the **maxillary molars** therefore require modification. Carefully, grind the edges off with fine rubber polishers to create abrasion facets and then polish the surfaces to a high gloss again.



The maxillary molar should lie flush in the central fossa of the mandibular molar.




Important!

The reduced areas must be polished to a high gloss again, as the modification would otherwise affect the abrasion behaviour of the tooth.

► Set-up – Upper anterior teeth

The front teeth are set-up according to the model analysis or the information from the wax registration.

In general, the maxillary incisors are aligned with the curvature of the dental arch and point toward the lower mucolabial fold. Their labial curvature harmoniously blends in with the vertical anterior tooth arch.




The central incisors are set-up perpendicular to the occlusal plane, the laterals are positioned with a slight mesial inclination and the canines are placed so that the long axis is almost vertical.



► Set-up – Mandible


The lower canines are set up in relation to tooth 13 and 23. The direct axis of the lower canine points between the upper lateral incisor and the canine.




The first premolar should be placed in such a way that there is a smooth transition to the canine. This means that the distal slope of the incisal edge is flush with the mesial rise of the first premolar. The template is adjusted to the cusp tips of the first mandibular premolar and the distal thirds of the retromolar pads.



The posterior teeth are placed with contact to the template on both sides.




They are oriented according to the centre of the alveolar ridge and Pound's line.





► Set-up – Maxillary posterior teeth

The upper posterior teeth are set up as in a normal bite situation with a one-to-two-tooth relationship. In a crossbite there is a reversed cusp-fossa relationship starting from the molar. The upper buccal cusps are in contact with the lower central fossa.

* For an optimum transition into the crossbite, the second premolar should be placed in an edge-to-edge bite situation. In this case, it is important not to reduce the buccal cusp tip of the upper premolar too much since otherwise the guidance contacts are affected.




To avoid the typical crossbite curve, minor changes are necessary:



- 1 Move the upper second premolars slightly toward the palatal area
- 2 Move the upper first premolars slightly more to the buccal
- 3 Move the lower first molar slightly more to the buccal
- 4 Move the lower second premolar slightly more toward the lingual area

Finally, the contact points of the set-up are checked again.



► In function – Comparison between a normal bite and crossbite set-up

Centric normal bite

The teeth are placed in a one-to-two-tooth relationship. Consequently, the primary contacts in the centric position are located in the central fossae of the mandible and on the marginal ridges.



Centric crossbite

Up to and including the second premolar, the occlusal relation is normal. From the first molar onwards, the teeth for the crossbite are situated in a statically correct, reverse relationship to the antagonist teeth.



Laterotrusion normal bite side

The mesio-buccal surfaces of the maxillary premolars slide over the distal-buccal edges of the mandibular premolars.



Laterotrusion crossbite side

The buccal cusps (working cusps) of the maxillary molars slide over the lingual surfaces of the mandibular molars.



Recommended literature

Removable Denture Prosthetics, Basics according to the BPS® Concept, Denture Manufacturing Protocol. Ivoclar Vivadent AG, Schaan (2014)

Removable Denture Prosthetics, Basics according to the BPS® Concept, Clinical Protocol. Ivoclar Vivadent AG, Schaan (2014)

Deep Overbite and Crossbite, Special Edition, Ivoclar Vivadent AG, Schaan (2006)

BPS Full Dentures, a Systematic Approach, Verlag Neuer Merkur GmbH, Munich (2003)

Ivoclar Vivadent – worldwide

Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan
Liechtenstein
Tel. +423 235 35 35
Fax +423 235 33 60
www.ivoclarvivadent.com

Ivoclar Vivadent Pty. Ltd.
1 – 5 Overseas Drive
P.O. Box 367
Noble Park, Vic. 3174
Australia
Tel. +61 3 9795 9599
Fax +61 3 9795 9645
www.ivoclarvivadent.com.au

Ivoclar Vivadent GmbH
Tech Gate Vienna
Donau-City-Strasse 1
1220 Wien
Austria
Tel. +43 1 263 191 10
Fax: +43 1 263 191 111
www.ivoclarvivadent.at

Ivoclar Vivadent Ltda.
Alameda Caiapós, 723
Centro Empresarial Tamboré
CEP 06460-110 Barueri – SP
Brazil
Tel. +55 11 2424 7400
Fax +55 11 3466 0840
www.ivoclarvivadent.com.br

Ivoclar Vivadent Inc.
1-6600 Dixie Road
Mississauga, Ontario
L5T 2Y2
Canada
Tel. +1 905 670 8499
Fax +1 905 670 3102
www.ivoclarvivadent.us

**Ivoclar Vivadent Shanghai
Trading Co., Ltd.**
2/F Building 1, 881 Wuding Road,
Jing An District
200040 Shanghai
China
Tel. +86 21 6032 1657
Fax +86 21 6176 0968
www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.
Calle 134 No. 7-B-83, Of. 520
Bogotá
Colombia
Tel. +57 1 627 3399
Fax +57 1 633 1663
www.ivoclarvivadent.co

Ivoclar Vivadent SAS
B.P. 118
F-74410 Saint-Jorioz
France
Tel. +33 4 50 88 64 00
Fax +33 4 50 68 91 52
www.ivoclarvivadent.fr

Ivoclar Vivadent GmbH
Dr. Adolf-Schneider-Str. 2
73479 Ellwangen, Jagst
Germany
Tel. +49 7961 889 0
Fax +49 7961 6326
www.ivoclarvivadent.de

**Ivoclar Vivadent Marketing (India)
Pvt. Ltd.**
503/504 Raheja Plaza
15 B Shah Industrial Estate
Veera Desai Road, Andheri (West)
Mumbai, 400 053
India
Tel. +91 22 2673 0302
Fax +91 22 2673 0301
www.ivoclarvivadent.in

Ivoclar Vivadent Marketing Ltd.
The Icon
Horizon Broadway BSD
Block M5 No. 1
Kecamatan Cisauk Kelurahan Sampora
15345 Tangerang Selatan – Banten
Indonesia
Tel. +62 21 3003 2932
Fax +62 21 3003 2934

Ivoclar Vivadent s.r.l.
Via Isonzo 67/69
40033 Casalecchio di Reno (BO)
Italy
Tel. +39 051 6113555
Fax +39 051 6113565
www.ivoclarvivadent.it

Ivoclar Vivadent K.K.
1-28-24-4F Hongo
Bunkyo-ku
Tokyo 113-0033
Japan
Tel. +81 3 6903 3535
Fax +81 3 5844 3657
www.ivoclarvivadent.jp

Ivoclar Vivadent Ltd.
12F W-Tower
54 Seocho-daero 77-gil, Seocho-gu
Seoul, 06611
Republic of Korea
Tel. +82 2 536 0714
Fax +82 2 596 0155
www.ivoclarvivadent.co.kr

Ivoclar Vivadent S.A. de C.V.
Calzada de Tlalpan 564,
Col Moderna, Del Benito Juárez
03810 México, D.F.
México
Tel. +52 (55) 50 62 10 00
Fax +52 (55) 50 62 10 29
www.ivoclarvivadent.com.mx

Ivoclar Vivadent BV
De Fruittuinen 32
2132 NZ Hoofddorp
Netherlands
Tel. +31 23 529 3791
Fax +31 23 555 4504
www.ivoclarvivadent.com

Ivoclar Vivadent Ltd.
12 Omega St, Rosedale
PO Box 303011 North Harbour
Auckland 0751
New Zealand
Tel. +64 9 914 9999
Fax +64 9 914 9990
www.ivoclarvivadent.co.nz

Ivoclar Vivadent Polska Sp. z o.o.
ul. Jana Pawla II 78
00-175 Warszawa
Poland
Tel. +48 22 635 5496
Fax +48 22 635 5469
www.ivoclarvivadent.pl

Ivoclar Vivadent LLC
Prospekt Andropova 18 korp. 6/
office 10-06
115432 Moscow
Russia
Tel. +7 499 418 0300
Fax +7 499 418 0310
www.ivoclarvivadent.ru

Ivoclar Vivadent Marketing Ltd.
Qlaya Main St.
Siricon Building No.14, 2nd Floor
Office No. 204
P.O. Box 300146
Riyadh 11372
Saudi Arabia
Tel. +966 11 293 8345
Fax +966 11 293 8344
www.ivoclarvivadent.com

Ivoclar Vivadent S.L.U.
Carretera de Fuencarral nº24
Portal 1 – Planta Baja
28108-Alcobendas (Madrid)
Spain
Tel. +34 91 375 78 20
Fax +34 91 375 78 38
www.ivoclarvivadent.es

Ivoclar Vivadent AB
Dalvägen 14
169 56 Solna
Sweden
Tel. +46 8 514 939 30
Fax +46 8 514 939 40
www.ivoclarvivadent.se

Ivoclar Vivadent Liaison Office
: Tesvikiye Mahallesi
Sakayik Sokak
Nisantas' Plaza No:38/2
Kat:5 Daire:24
34021 Sisli – Istanbul
Turkey
Tel. +90 212 343 0802
Fax +90 212 343 0842
www.ivoclarvivadent.com

Ivoclar Vivadent Limited
Compass Building
Feldspar Close
Warrens Business Park
Enderby
Leicester LE19 4SD
United Kingdom
Tel. +44 116 284 7880
Fax +44 116 284 7881
www.ivoclarvivadent.co.uk

Ivoclar Vivadent, Inc.
175 Pineview Drive
Amherst, N.Y. 14228
USA
Tel. +1 800 533 6825
Fax +1 716 691 2285
www.ivoclarvivadent.us