IPS Empress® Direct

IPS Empress® Direct Effect

EN Instructions for Use

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English

Description

IPS Empress® Direct is a light-curing, radiopaque composite, which is suitable for the fabrication of restorations that fulfill the highest demands. IPS Empress® Direct Effect is a flowable, light-curing, highly esthetic composite for the direct restorative therapy. IPS Empress Direct and IPS Empress Direct Effect cure with light in the wavelength range of 400-500 nm.

Shades

IPS Empress Direct is available in the following shades:

Enamel shades

A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, D2, D3, Bleach-L, Bleach-XL

Dentin shades

A1, A2, A3, A3.5, A4, IVA5, IVA6, B1, B2, C3, D2, Bleach-L, Bleach-XL

Translucent shades

Trans 20, Trans 30, Trans Opal*

*Not radiopaque due to its high transparency

IPS Empress Direct Effect is available in the following shades:

Enamel shades

Bleach XL

Highly translucent Incisal shades

Trans Opal, Trans 30

Compositions

IPS Empress Direct

The monomer matrix is composed of dimethacrylates (20-21.5 wt%, opalescent shade 17 wt%). The fillers contain barium glass, ytterbium trifluoride, mixed oxide, silicon dioxide and copolymer (77.5-79 wt%, opalescent shade 83 wt%). Additional contents: additives, initiators, stabilizers and pigments (<1.0 wt%). The total content of inorganic fillers is 75-79 wt% or 52-59 vol% (opalescent shade 60.5 wt% or 45 vol%). The particle size of the inorganic fillers is between 40 nm and 3 µm with a mean particle size of 550 nm.

IPS Empress Direct Effect

The monomer matrix is composed of dimethacrylates (42 wt% in the case of Effect Trans 30 and Effect Bleach XL, 33 wt% in the case of Effect Trans Opal). The fillers comprise barium glass, ytterbium trifluoride, highly dispersed silicon dioxide, mixed oxide and copolymer (58 wt% in the case of Effect Trans 30 and Effect Bleach XL, 67 wt% in the case of Effect Trans Opal). Additives, initiators, stabilizers and pigments are additional contents (< 1.0 wt%). The total content of inorganic fillers is 57 wt% or 36 vol% in the case of Effect Trans 30 and Effect Bleach XL, and 49wt% / 34 vol% in the case of Effect Trans Opal). The particle sizes of the inorganic fillers range between 40 nm and 3 µm. The mean particle size is 550 nm.

Indications

IPS Empress Direct

- Anterior restorations (Class III, IV)
- Restorations in the posterior region (Classes I and II)
- Tooth neck restorations (Class V, e.g. in the case of cervical caries, root erosion, wedge-shaped defects)
- Corrections of tooth positions and tooth shape (e.g. diastema closure, closure of interdental “black triangles”, lengthening of the incisal edge)
- Direct veneers

IPS Empress Direct Effect

- As an intermediate or covering layer in the fabrication of esthetically demanding composite restorations
- Repair of composite restorations (in particular filling of voids, levelling out of porosities and minor chipping)
- As initial layer under Class I and II restorations
- Small restorations of all types
- Extended fissure sealing
- Splinting of mobile teeth
- Blocking out of undercuts
- Repair of composite and ceramic veneers

Contraindication

Placement of IPS Empress Direct and IPS Empress Direct Effect restorations is contraindicated

- if a dry working field cannot be established, or if the stipulated technique cannot be applied.
- if a patient is known to be allergic to any of the ingredients in IPS Empress Direct and IPS Empress Direct Effect.

Side effects

In rare cases, components of IPS Empress Direct and IPS Empress Direct Effect may lead to sensitization. IPS Empress Direct or IPS Empress Direct Effect should not be used in such cases. To avoid possible irritation of the pulp, areas close to the pulp should be protected with a suitable pulp/dentin protector (selectively apply a calcium hydroxide-based preparation in areas close to the pulp and coat with a suitable liner).

Interactions

Phenolic substances, such as materials containing eugenol/clove oil, inhibit the polymerization of methacrylate-based materials. Consequently, the application of such materials together with IPS Empress Direct and IPS Empress Direct Effect must be avoided. Discoloration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.

Application

1. Shade determination

Clean the teeth prior to shade determination. The enamel shade is selected with the tooth still moist. The dentin shade can be taken either by placing the tab next to the prepared tooth or next to the cervical area of the canine.

2. Isolation

Adequate relative or absolute isolation using suitable auxiliaries, such as OptraGate® or OptraDam® Plus, is required.

3. Cavity preparation

Cavity preparation is carried out according to the requirements of the adhesive technique, i.e. protecting the tooth structure. Do not prepare any sharp internal edges nor additional undercuts in caries-free areas. The cavity geometry is mainly determined by the extend of caries and/or the old restoration. Enamel edges in the anterior region are bevelled. In the posterior region, only sharp enamel edges should be rounded (finishing diamonds, 25-40 µm). Subsequently, all residue in the cavity is removed with water spray. The cavity is dried with water- and oil-free air.

4. Pulp protection / Base

Do not apply a base material when using an enamel/dentin bonding agent. Only cover very deep areas close to the pulp with a calcium hydroxide material (e.g. ApexCal®) and subsequently use pressure-resistant cement (e.g. glass ionomer cement, such as Vivaglass® Liner). Do not cover other cavity walls, since they can be used to support the bond with an enamel/dentin adhesive.

5. Placement of a matrix / interdental wedge

Use a wrap-around matrix for cavities affecting the proximal area or a sectional matrix and wedge.

6. Conditioning / Application of the bonding agent

Condition the tooth structure and apply the bonding agent according to the Instructions for Use of the product in use. Ivoclar Vivadent recommends using the universal bonding agent Adhese® Universal, a single-component, light-cured adhesive suitable for direct and indirect restorations and all etching techniques.

7. Application of IPS Empress Direct and IPS Empress Direct Effect

- Since a large number of IPS Empress Direct Dentin and Enamel shades in different translucencies for the creation of anatomical details are available, any layering technique can be creatively used to achieve highly esthetic restorations. Apply IPS Empress Direct in layers of max. 2 mm (or 1.5 mm Dentin shades) thickness and, for an optimum result, adapt with e.g. OptraSculpt®/OptraSculpt® Pad. Apply IPS Empress Direct Effect in layers of max. 2 mm thickness.
- Sufficient exposure to the curing light prevents incomplete polymerization. Recommended curing times (exposure time) and light intensities are listed in Table 1. Hold the light emission window as closely as possible to the surface of the restorative material.
- When using a metal matrix, additionally polymerize the restoration from the oral and vestibular aspect after removing the matrix, if no Bluephase® curing light was used or if the light probe could not be optimally positioned (e.g. far away from the material or with divergent radiation).
- In many cases, a flowable composite is now used for the initial layer

in order to create an even cavity floor and to simplify the adaptation of the subsequently used restorative material. The use of a flowable composite (e.g. Tetric EvoFlow®) as an initial thin layer is optional. This layer must be light-cured in a separate step (the instructions of the respective Instructions for Use must be observed).

8. Finishing / Checking the occlusion / Polishing

Remove excess material with suitable finishers or fine diamonds after polymerization. Remove proximal excess with diamond finishers, tungsten carbide finishers, finishing strips, or flexible finishing discs. Check the occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restorations. Use silicone polishers (e.g. OptraPol®) as well as polishing discs and polishing strips to polish the restorations to a durable high gloss.

Additional information

1. IPS Empress Direct can be combined with IPS Empress Direct Effect.
2. If necessary, additional IPS Empress Direct or IPS Empress Direct Effect material can be directly applied to the polymerized material. If the IPS Empress Direct or IPS Empress Direct Effect restoration has already been polished, it must first be roughened and wetted with Adhese Universal before fresh IPS Empress Direct or IPS Empress Direct Effect is applied.
3. IPS Empress Direct and IPS Empress Direct Effect should be at room temperature when they are applied. Cool temperatures render the material difficult to extrude.
4.  For single use only. If IPS Empress Direct and IPS Empress Direct Effect are applied from the Cavifil or the application canula directly in the mouth of the patient, the Cavifil or the application canula must be discarded after use to prevent cross-contamination between patients.
5. Syringes or Cavifils should not be disinfected with oxidizing disinfection agents.

Warning

Unpolymerized IPS Empress Direct and IPS Empress Direct Effect should not come in contact with the skin, mucous membrane, or eyes. Unpolymerized IPS Empress Direct and IPS Empress Direct Effect may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage

- Storage temperature 2-28 °C / 36-82 °F.
- Close syringes/Cavifils immediately after use. Exposure to light causes premature polymerization.
- Do not use IPS Empress Direct and IPS Empress Direct Effect after the indicated date of expiration.
- Shelf life: see information on Cavifils, syringes and packages.

Keep material out of children's reach.

For use in dentistry only.

The materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Table 1

Program \ Unit	Bluephase® C8	Bluephase® Style	Bluephase®	Bluephase® Style 20i
Turbo	–	–	–	5 s
High Power	15 s	10 s	10 s	10 s
Soft Start	–	–	15 s	–
light intensity		exposure time		
≥ 500 mW/cm ²		20 s		
≥ 1000 mW/cm ²		10 s		

Ivoclar Vivadent AG

Bendererstrasse 2, 9494 Schaan, Liechtenstein
Tel. +423 235 35 35, Fax +423 235 33 60
www.ivoclarvivadent.com

Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive, P.O. Box 367
Noble Park, Vic. 3174, Australia
Tel. +61 3 9795 9599, Fax +61 3 9795 9645
www.ivoclarvivadent.com.au

Ivoclar Vivadent Ltda.

Alameda Caiapós, 723, Centro Empresarial Tamboré
CEP 06460-110 Barueri – SP, Brazil
Tel. +55 11 2424 7400, Fax +55 11 3466 0840
www.ivoclarvivadent.com.br

Ivoclar Vivadent Inc.

1-6600 Dixie Road, Mississauga, Ontario, L5T 2Y2, Canada
Tel. +1 905 670 8499, Fax +1 905 670 3102
www.ivoclarvivadent.us

Ivoclar Vivadent Shanghai Trading Co., Ltd.

2/F Building 1, 881 Wuding Road, Jing An District
200040 Shanghai, China
Tel. +86 21 6032 1657, Fax +86 21 6176 0968
www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

Calle 134 No. 7-B-83, Of. 520, Bogotá, Colombia
Tel. +57 1 627 3399, Fax +57 1 633 1663
www.ivoclarvivadent.co

Ivoclar Vivadent SAS

B.P. 118, F-74410 Saint-Jorioz, France
Tel. +33 4 50 88 64 00, Fax +33 4 50 68 91 52
www.ivoclarvivadent.fr

Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2
D-73479 Ellwangen, Jagst, Germany
Tel. +49 7961 889 0, Fax +49 7961 6326
www.ivoclarvivadent.de

Ivoclar Vivadent Marketing (India) Pvt. Ltd.

503/504 Raheja Plaza, 15 B Shah Industrial Estate
Veera Desai Road, Andheri (West), Mumbai, 400 053, India
Tel. +91 22 2673 0302, Fax +91 22 2673 0301
www.ivoclarvivadent.in

Ivoclar Vivadent s.r.l.

Via Isonzo 67/69, 40033 Casalecchio di Reno (BO), Italy
Tel. +39 051 6113555, Fax +39 051 6113565
www.ivoclarvivadent.it

Ivoclar Vivadent K.K.

1-28-24-4F Hongo, Bunkyo-ku, Tokyo 113-0033, Japan
Tel. +81 3 6903 3535, Fax +81 3 5844 3657
www.ivoclarvivadent.jp

Ivoclar Vivadent Ltd.

12F W-Tower, 1303-37
Seocho-dong, Seocho-gu, Seoul 137-855
Republic of Korea
Tel. +82 2 536 0714, Fax +82 2 596 0155
www.ivoclarvivadent.co.kr

Ivoclar Vivadent S.A. de C.V.

Av. Insurgentes Sur No. 863, Piso 14, Col. Napoles
03810 México, D.F., México
Tel. +52 55 5062 1000, Fax +52 55 5062 1029
www.ivoclarvivadent.com.mx

Ivoclar Vivadent BV

De Fruittuinen 32, 2132 NZ Hoofddorp, Netherlands
Tel. +31 23 529 3791, Fax +31 23 555 4504
www.ivoclarvivadent.com

Ivoclar Vivadent Ltd.

12 Omega St, Rosedale, PO Box 303011 North Harbour
Auckland 0751, New Zealand
Tel. +64 9 914 9999, Fax +64 9 914 9990
www.ivoclarvivadent.co.nz

Ivoclar Vivadent Polska Sp. z o.o.

Al. Jana Pawla II 78, 00-175 Warszawa, Poland
Tel. +48 22 635 5496, Fax +48 22 635 5469
www.ivoclarvivadent.pl

Ivoclar Vivadent Marketing Ltd.

Prospekt Andropova 18 korp. 6/
office 10-06, 115432 Moscow, Russia
Tel. +7 499 418 0300, Fax +7 499 418 0310
www.ivoclarvivadent.ru

Ivoclar Vivadent Marketing Ltd.

Qlaya Main St., Siricon Building No.14, 2nd Floor
Office No. 204, P.O. Box 300146, Riyadh 11372
Saudi Arabia
Tel. +966 11 293 8345, Fax +966 11 293 8344
www.ivoclarvivadent.com

Ivoclar Vivadent S.L.U.

C/ Ribera del Loira nº 46, 5ª planta
28042 Madrid, Spain
Tel. +34 913 757 820, Fax +34 913 757 838
www.ivoclarvivadent.es

Ivoclar Vivadent AB

Dalvägen 14, S-169 56 Solna, Sweden
Tel. +46 8 514 939 30, Fax +46 8 514 939 40
www.ivoclarvivadent.se

Ivoclar Vivadent Liaison Office

: Tesvikiye Mahallesi , Sakayik Sokak
Nisantas' Plaza No:38/2, Kat:5 Daire:24
34021 Sisli – Istanbul, Turkey
Tel. +90 212 343 0802, Fax +90 212 343 0842
www.ivoclarvivadent.com

Ivoclar Vivadent Limited

Ground Floor Compass Building, Feldspar Close
Warrens Business Park, Enderby
Leicester LE19 4SE, United Kingdom
Tel. +44 116 284 7880, Fax +44 116 284 7881
www.ivoclarvivadent.co.uk

Ivoclar Vivadent, Inc.

175 Pineview Drive, Amherst, N.Y. 14228, USA
Tel. +1 800 533 6825, Fax +1 716 691 2285
www.ivoclarvivadent.us