**EN Instructions for Use**

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**Product description**

Transparent, two-component addition silicone. Especially suitable for the fabrication of silicone keys (C&B technique) in dental laboratories.

**Composition**

Polysiloxane

**Properties**

- Light-curing composite resins can be cured through the transparent silicone key.
- Hard after setting (70 Shore A)
- After setting, the silicone can be finished with burs and/or cutting instruments.
- Precise results.

**Application**

- Fabrication of keys in the C&B technique in conjunction with the wax-up
- Especially suitable for use in the SR Nexco® Flask
- The cartridge tips of Transil F exactly fit into the openings of the SR Nexco Flask intended for this purpose.

**Limitations of use**

Latex gloves inhibit the setting of polysiloxanes. Therefore, Transil F must not come into direct or indirect contact with latex gloves. It is recommended that operators wear vinyl gloves.

**Side effects**

Polysiloxanes may provoke irritation or other allergic reactions.

**Processing**

Twist off the cap of the cartridge and mount the mixing tip. Subsequently, insert the cartridge into the dispenser. The base and catalyst are auto-mixed in a 1:1 ratio in the mixing tip. Extrude the material evenly until the key has been completed. After use, leave the used mixing tip on the cartridge. It serves as a sealing cap. Use a new mixing tip for each application. Transil F sets at room temperature (23 °C / 73 °F). The working time is 20–40 seconds and the setting time is 180–260 seconds. These time indications

apply from the start of mixing at 23 °C / 73 °F. Higher temperatures will reduce the working and setting times, lower temperatures will prolong them.

**Warning**

If uncured material comes into contact with the skin or eyes, rinse with copious amounts of water. If irritation persists, seek medical attention.

**Storage**

- Storage temperature 2 – 28°C (36 – 82°F)
- Expiration date: see label on packaging
- Do not use Transil F after the expiry date.
- Keep away from sources of heat.

**Keep material out of children's reach!****For use in dentistry only.**

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions.

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