Bluephase® PowerCure

The intelligent curing light

EN Instructions for Use

Page 2

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Pagina 56

Appendix

Page 110

DE Gebrauchsinformation

Seite 20

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Página 74

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Page 38

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Página 92

C US

CE 0123

For dental use only!
Made in Austria

Manufacturer
Ivoclar Vivadent AG
9494 Schaan/Liechtenstein
www.ivoclarvivadent.com

Rx ONLY

In USA distributed by
Ivoclar Vivadent Inc.,
175 Pineview Drive,
Amherst, NY 14228, USA

ivoclar
vivadent®
clinical

Dear Customer

Optimum polymerization is an important requirement for all light cured materials in order to consistently produce high quality restorations. The curing light selected also plays a decisive role in this respect. Therefore, we would like to thank you for having purchased Bluephase® PowerCure.

Bluephase® PowerCure is a high-quality medical device which has been designed according to the latest standard of science and technology in compliance with the relevant industry standards.

These Instructions for Use will help you safely start up the device, make full use of its capabilities and ensure a long service life.

Should you have any further questions, please do not hesitate to contact us (see addresses on the reverse page).

Your Ivoclar Vivadent Team

Table of Contents

1. Product Overview	4
1.1 List of parts	
1.2 Indicators on the charging base	
1.3 Indicators on the handpiece	
1.4 Operating the light	
2. Intended Purpose and Safety	7
2.1 Normal use	
2.2 Indications	
2.3 Signs & symbols	
2.4 Safety notes	
2.5 Contraindications	
3. Start-Up	10
4. Operation	13
5. Maintenance and Cleaning	16
6. What if ...?	18
7. Warranty /Procedure in Case of Repair	19
8. Product Specifications	19

1. Product Overview

1.1 List of parts


1.2 Indicators on the charging base

Charging base with integrated radiometer and Click & Cure function:

- Indicator is dark: Battery is charged
- Indicator shows battery symbol: Battery is charging
- Indicator shows light intensity: Measurement has been carried out

1.3 Indicators on the handpiece

1.4 Operating the light

Activating/deactivating polyvision

The polyvision function is activated by long pressing (>2 s) the program or time selection button (see 4. Operation). To confirm that polyvision is activated, the handpiece will beep and vibrate for a short time and the polyvision indicator will light up. Renewed long pressing of the program or time selection button results in the polyvision function being deactivated; the handpiece will not vibrate.

With the handpiece switched on, the current charging status is shown on the handpiece as follows:

- **No indicators lit up on the handpiece:**

- Battery sufficiently charged**

- Curing capacity of minimum 20 minutes in the High Power program, approx. 15 minutes in the Turbo program and approx. 8 minutes in the 3sCure program.

- **Battery symbol on the handpiece is flashing orange:**

- Battery weak**

- Time/intensity can still be set and a polymerization time of approx. 6 minutes in the High Power program, approx. 4 minutes in the Turbo program, or 30 cycles in the 3sCure program are left. Place the device into the charging base as soon as possible.

- **Battery symbol on the handpiece is flashing orange and a red "x" is shown:**

- Battery completely discharged**

- The light can no longer be activated and the curing time can no longer be set. However, the handpiece can be used in the Click & Cure corded mode.

2. Intended Purpose and Safety

2.1 Normal use

Bluephase PowerCure is an LED curing light that produces blue light. It is used for the polymerization of light-curing dental materials directly in the oral cavity of patients. The curing light is intended for use by the dentist or dental assistant at the dental practice, medical practice or hospital. Intended use also assumes the observation of the notes and regulations in these Instructions for Use. The curing light must only be operated by trained dental personnel.

The 3sCure program (3,000 mW/cm²) must only be used for direct restorations in the posterior region (Class I & II) in combination with Tetric PowerFill, Tetric PowerFlow or Adhese Universal. Do not use the 3sCure program in cases of deep cavities (caries profunda).

Normal use: light-curing of dental materials

2.2 Indications

With its broadband Polywave® LED technology, Bluephase PowerCure is suitable for the polymerization of all light-curing dental materials curing in the wavelength range of 385 – 515 nm. These materials include restoratives, bonding agents/adhesives, bases, liners, fissure sealants, temporaries, as well as luting materials for brackets and indirect restorations, such as ceramic inlays.

2.3 Signs and symbols

Contraindication

Symbols on the curing light

Double insulation
(device complies with safety class II)

Protection against electrical shock
(BF type apparatus)

Observe Instructions for Use

Observe Instructions for Use

The Instructions for Use are stored on the website of Ivoclar Vivadent AG (www.ivoclarvivadent.com).

Caution

The curing light must not be disposed of in the normal domestic waste. Information regarding disposal of the device can be found on the respective national Ivoclar Vivadent homepage.

Recyclable

AC voltage

DC voltage

2.4 Safety notes

Bluephase PowerCure is an electronic device and a medical product which is subject to IEC 60601-1 (EN 60601-1) and EMC directives IEC 60601-1-2 (EN60601-1-2), as well as the 93/42/EEC Medical Device Directive. The curing light complies with the relevant EU regulations.

CE 0123

The curing light has been shipped from the manufacturer in a safe and technically sound condition. In order to maintain this condition and to ensure risk-free operation, the notes and regulations in these Instructions for Use must be observed. To prevent damage to equipment and risks for patients, users and third parties, the following safety instructions must be observed.

Safety notes for the "3sCure program" (3,000 mW/cm²):
– Avoid direct exposure of the gingiva, oral mucous membrane or skin.

2.5 Contraindications

Materials whereby polymerization is activated outside the wavelength range of 385–515 nm (no such materials known to date). If you are not sure about certain products, please ask the manufacturer of the corresponding material.

Do not charge or use the appliance near flammable or combustible substances.

Never use without light guide.

The use of a light guide other than the one provided in the delivery form is not admissible.

Using the device stacked on or close to other equipment should be avoided because correct function can be disrupted. If such use is unavoidable, the device needs to be monitored and checked for correct function.

Portable and mobile high-frequency communication devices may interfere with medical equipment. The use of mobile phones during operation is not allowed.

Caution – The use of controls or adjustment devices or performing procedures other than those specified herein may result in hazardous radiation exposure.

Never use without eye protection for patients and users.

If the light emission window cannot be optimally placed in relation to the composite restoration, the restoration must be polymerized using a conventional method. If soft tissue exposure to the curing light cannot be avoided, the 3sCure program must not be used, as exposure may result in damage of the soft tissues.

Warning

This unit should not be used near flammable anaesthetics or mixtures of flammable anaesthetics with air, oxygen or nitric oxide.

In the case of serious adverse events in connection with the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan, www.ivoclarvivadent.com, and your local health authorities.

Usage and liability

- Bluephase PowerCure must only be employed for its intended use. Any other uses are contraindicated. Do not touch defective, open devices. Liability cannot be accepted for damage resulting from misuse or failure to observe the Instructions for Use.
- The user is responsible for testing Bluephase PowerCure for its use and suitability for the intended purposes. This is particularly important if other equipment is used in the immediate vicinity of the curing light at the same time.
- Use only original spare parts and accessories from Ivoclar Vivadent. The manufacturer does not accept any liability for damage resulting from the use of other spare parts or accessories.
- The light guide is an applied part and may warm up to a maximum of 45 °C (113 °F) at the interface to the handpiece during operation.
- Keep out of the reach of children!
- For use in dentistry only!

Operating voltage

Before switching on, make sure that

- a) the voltage indicated on the rating plate complies with the local power supply and
- b) the unit has acquired the ambient temperature.

Do not touch the exposed contacts of the connection plug (power pack). If the battery or power pack are used separately (i.e. during start-up), contact with patients or third parties must be prevented.

Assumption of impaired safety

If it has to be assumed that safe operation is no longer possible, the power must be disconnected and the battery removed to avoid accidental operation. This may be the case, for example, if the device is visibly damaged or no longer works correctly. A complete disconnection from the power supply is only ensured when the power cord is disconnected from the power source. Ensure that the device can be quickly and easily disconnected at any time.


Eye protection

Direct or indirect exposure of the eyes must be prevented. Prolonged exposure to the light is unpleasant for the eyes and may result in injury. To optimize user safety, the device has been equipped with intelligent anti-glare protection. For this, the "polyvision" function must be activated (see 4. Operation). With the polyvision function activated, the Bluephase PowerCure automatically senses if the handpiece is outside of the mouth and automatically switches off the light if it has been activated inadvertently. If a protective sleeve is used, this function is not available.

In addition, we also recommend using the anti-glare protective equipment supplied. Individuals who are sensitive to light, who take photosensitizing drugs, have undergone eye surgery, or work with the apparatus (or in its vicinity) for long periods of time, should not be exposed and should wear protective orange goggles that absorb light below a wavelength of 515 nm. The same applies to patients.

Battery

Caution: Use only original spare parts for Bluephase PowerCure, particularly Ivoclar Vivadent batteries and charging bases. Do not short circuit battery. Do not touch the battery contacts. Do not store at temperatures above 40 °C / 104 °F (or 60 °C / 140 °F for a short period). Always store batteries charged. The storage period must not exceed 6 months. May explode if disposed of in fire.

Please note that lithium-ion batteries may react with explosion, fire and smoke development if handled improperly or mechanically damaged. Damaged lithium-ion batteries must no longer be used.

The electrolytes and electrolyte fumes released during explosion, fire and smoke development are toxic and corrosive. Do not touch leaking batteries with bare hands. In case of contact with the eyes or skin, immediately wash with copious amounts of water. Avoid inhalation of fumes. In case of indisposition, see a physician immediately. Remove electrolyte residue from surfaces by washing/wiping with a moist cloth. Wash contaminated pieces of clothing immediately.

Heat development

As is the case with all high-performance lights, the high light intensity results in a certain heat development. Prolonged exposure of areas near the pulp and soft tissues may result in irreversible damage. Therefore, this high-performance curing light must only be operated by trained professionals.

Generally, the recommended curing times, particularly in areas near the pulp, must be observed. Using more than two 3sCure cycles on the same tooth within a period of 30 seconds and direct contact with the gingiva, oral mucous membrane or skin must be strictly avoided. In addition, the light emission window must remain positioned exactly on the material to be cured during the entire duration of the curing cycle (e.g. by holding it in place using a finger). If soft tissue exposure cannot be avoided, work with reduced light intensity. Polymerize indirect restorations at intermittent intervals of 10 seconds in the High Power program and 1 x 5 seconds each in the Turbo program or use external cooling with an air stream. The instructions regarding curing programs and curing times must be observed (see Selecting the curing program).

After several curing cycles on the same tooth, there is a risk that the pulp suffers damage caused by the increased temperature!

3. Start-Up

Check the delivery for completeness and any possible transportation damage (see List of parts). If parts are damaged or missing, contact your Ivoclar Vivadent representative.

Charging base

Before you switch on the device, make sure that the voltage mentioned on the rating plate complies with your local power supply.

Slide the connection plug of the power pack into the socket on the underside of the charging base. Tilt it slightly and apply slight pressure until you hear and feel it snap into place. Place the charging base on a suitable, flat table top.

Connect the power cord with the power supply and the power pack. Make sure that the power cord is easily accessible at all times and can easily be disconnected from the power supply. The charging base will briefly display "Bluephase PowerCure" as ticker text on the screen.

Handpiece

Unpack the handpiece from the packaging and detach the light guide by pulling it out. Then clean the handpiece and light guide (see Maintenance and Cleaning). After cleaning, reinsert the light guide.

For reasons of hygiene, we recommend using a disposable protective sleeve for each patient (see Maintenance and Cleaning). Make sure to fit the protective sleeve snugly to the light guide. Sleeves enclosed in the delivery form may be used or sleeves can be purchased per country specific regulations. Then, attach the anti-glare cone or anti-glare shield to the light guide.

Battery

The battery must be fully charged before it is used for the first time! When fully charged, the battery features a curing capacity of at least 20 minutes in the High Power program, approx. 15 minutes in the Turbo program and approx. 8 minutes in the 3-second curing program (3sCure).

Slide the battery straight into the handpiece until you hear and feel it click into place.

Gently place the handpiece in the corresponding rest in the charging base without using any force. If a hygiene sleeve is used, please remove it before you charge the battery. If possible, always use the light with a fully charged battery. This will prolong the service life. We therefore recommend placing the handpiece into the charging base after each patient. If the battery is fully discharged, the charging time is 2 hours.

The battery is an expendable part that typically needs to be replaced after approximately 2.5 years. See battery label for the age of the battery.


Battery: charging status

The current charging status is displayed on the handpiece as described on page 6.

Click & Cure corded operation

Bluephase PowerCure can be used in corded operation mode at any time and particularly when the battery is completely empty.

For corded operation, release the battery by pressing on the release button on the lower part of the handpiece and then pull the battery out of the handpiece completely.

Then remove the power pack from the underside of the charging base. Do not pull on the power cord.

Insert the connection plug straight into the handpiece until you hear and feel it click into place.

During corded operation, the charging base cannot charge the battery, since it is not connected to a power source.

A complete disconnection from the power supply is only ensured when the power cord is disconnected from the power socket.

Measuring the light intensity

The integrated radiometer allows the light intensity (mW/cm^2) to be measured easily and quickly while the charging base is connected.

To measure the light intensity, place the tip of the light guide, without protective cover, flush into the marked recess on the upper side of the charging base. Then activate the light and read the value displayed on the screen. Measurement accuracy is in the range of $\pm 10\%$. If the light intensity is below $400 \text{ mW}/\text{cm}^2$, the screen will display "LOW". Measure the light intensity of the curing light only with the accompanying charging base supplied in the delivery form.

4. Operation

Disinfect contaminated surfaces of the curing light as well as light guides and anti-glare cones before each use. Additionally, the light guide can be sterilized using the autoclaves intended for this purpose (see chapter Maintenance and Cleaning). Furthermore, make sure that the emitted light intensity permits adequate polymerization. For that purpose, inspect the light guide for contamination and damage and check the light intensity at regular intervals (see paragraph Measuring the light intensity).

Repeated consecutive curing cycles (more than 2 curing cycles in succession) on the same tooth entail the risk that the pulp may become damaged due to the increase in temperature.

The PreCure program must not be used for conventional light-curing!

Selecting the curing program and curing time

Bluephase PowerCure is equipped with 5 selectable curing times and 4 curing programs for different indications. Use the Time/Program selection button to set the desired curing time and/or the light intensity.

H* (High Power Program), 1,200 mW/cm²:

The following curing times can be selected in the High Power program: 10, 15 or 20 seconds.

T* (Turbo Program), 2,100 mW/cm²:


The curing time is set to 5 seconds in the Turbo program and cannot be changed.

3s* (3sCure program), 3,050 mW/cm²:

In the 3sCure program, the curing time is set to 3 seconds and cannot be changed. The 3sCure program can be activated no more than twice in direct succession. After the two curing cycles have been conducted, a 30-second break is necessary for safety reasons before the next curing cycle can be activated.

PRE (PreCure Program), 950 mW/cm²:

The PreCure program is used to tack cure light-curing, adhesive luting composites (Variolink Esthetic) in order to facilitate the removal of excess material. If other luting composites are used, the distance from the light guide to the composite must be increased or several curing cycles must be conducted. The curing time of the PreCure program is preset to 2 seconds and cannot be altered.

Observe the Instructions for Use of the material applied when selecting the curing time and intensity. The curing recommendations for composite materials apply to all shades and, if not mentioned otherwise in the Instructions for Use, to a maximum layer thickness of 2 mm. Generally, these recommendations apply to situations where the emission window of the light guide is placed directly over the material to be polymerized. Increasing the distance between the light source and the material will require the curing time to be extended

accordingly. For instance, if the distance to the material is approx. 8 mm, the effective light output is reduced by approx. 50%. In this case, the recommended curing time has to be doubled.

- 1) The information provided herein applies to the 10>9 mm light guide supplied in the delivery form.
- 2) The information regarding heat development and burn hazards must be taken into consideration (see Safety notes).

Curing times		HIGH POWER 1,200 mW/cm ²	TURBO 2,000 mW/cm ²	3sCure 3,000 mW/cm ²
3sCure System	Tetric PowerFill / Tetric PowerFlow 4 mm	10 seconds	5 seconds	3 seconds ¹⁾
	Adhese Universal	10 seconds	---	3 seconds ¹⁾
Restorative materials	Composites • 2 mm ²⁾ IPS Empress Direct / IPS Empress Direct Effect Tetric EvoCeram / Tetric EvoFlow	10 seconds	5 seconds	---
	• 4 mm ³⁾ Tetric EvoCeram Bulk Fill / Tetric EvoFlow Bulk Fill	10 seconds	5 seconds	---
Indirect restorations / cementation materials	Variolink Esthetic LC ⁴⁾ / Variolink Esthetic DC ⁵⁾	per mm ceramic: 10 seconds per surface	per mm ceramic: 5 seconds per surface	---
	Multilink Automix ⁵⁾ / SpeedCEM Plus ⁵⁾	per mm ceramic: 20 seconds per surface	per mm ceramic: 2 x 5 seconds per surface	---
Temporary materials	Telio CS Link Telio CS Inlay/Onlay	10 seconds per surface 10 seconds	10 seconds per surface 5 seconds	---
Miscellaneous	Helioseal F Plus	10 seconds	---	---
	IPS Empress Direct Color IPS Empress Direct Opaque	20 seconds 20 seconds	5 seconds 2 x 5 seconds	---
	MultiCore Flow / Multicore HB	20 seconds	2 x 5 seconds	---

1) The 3sCure program may only be used on the occlusal surfaces of Class I & II restorations in the posterior region and may not be used at all in cases of caries profunda, or deep cavities.

2) Applies to a maximum layer thickness of 2 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (which may be the case, e.g. with certain dentin shades)

3) Applies to a maximum layer thickness of 4 mm and provided that the Instructions for Use of the respective material do not state any other recommendation (which may be the case, e.g. with certain dentin shades)

4) Applies to light-curing

5) Applies to optional light-curing

Polyvision - Automatic assistant with intelligent anti-glare protection

Bluephase PowerCure features a fully automated "polyvision" assistant for safe polymerization.

This feature enables the curing light to automatically sense if the handpiece is inadvertently moved from its position during the curing procedure. To prevent any associated reduction in the amount of energy being transferred, the device will start to vibrate to alert the user to the improper use and will automatically extend the curing time by 10%. If the handpiece is significantly shifted away from its initial position (e.g. out of the oral cavity), the light will automatically switch off so that the curing process can be restarted and carried out correctly.

In addition, polyvision also acts as intelligent anti-glare protection. Polyvision ensures that the curing light cannot be used outside the mouth. The light can only be activated once the light guide is positioned directly over the material to be polymerized. This protects the operator and patient from being blinded by the light. If a protective sleeve is used, this function is not available. Contamination of the light guide may affect the functionalities of the polyvision assistant.

The automatic assistant is there to support the user. It does not obviate the need for monitoring by the user.

If you do not want to use the assistant, you can deactivate it at any time by long pressing (> 2 seconds) the Time or Program selection button. The symbol on the handpiece (see 1.3 Indicators on the handpiece) then disappears.

Measuring the light intensity

We recommend checking the light intensity emitted by the curing light at regular intervals using the integrated radiometer in order to be sure that the materials are adequately cured and the composite fillings are of a high and lasting quality (see 3. Start-Up).

Cure Memory function

The last settings used, together with the combination of curing program and curing time, are automatically saved.

Start/Stop

The light is switched on by means of the Start/Stop button. It is recommended that the emission window of the light guide is placed directly on the material to be polymerized. Once the selected curing time has elapsed, the curing program is automatically terminated. If desired, the light can be switched off before the set curing time has elapsed by pressing the Start/Stop button again.

Acoustic signals

Acoustic signals can be heard for the following functions:

- Start (Stop)
- Every 10 seconds
- Curing time and program change
- Inserting battery
- Error message (when anti-glare protection is activated or curing cycle aborted)
- Error message when attempting to activate the 3sCure program during the 30-second waiting time that follows after the program has been used twice in succession.

Light intensity

The light intensity is maintained at a consistent level during operation. If the supplied 10>9 mm light guide is used, the light intensity has been calibrated to 1,200 mW/cm² ± 10% in the High Power program. In the Turbo program, the light intensity is calibrated to 2,100 mW/cm² ± 10%. In the 3-second program (3sCure), the light intensity is calibrated to 3,050 mW/cm² ± 10%. If another light guide than the one supplied is used, it directly influences the emitted light intensity.

When focussing light guides (10>9 mm) or Pin-Point light guides (6>2 mm) are used, the diameter of the light entry is larger than the diameter of the light emission window. The incident blue light is therefore concentrated over a smaller area. As a result, the emitted light intensity is increased.

Pin-Point light guides are suitable for spot-on polymerization, e.g. to tack veneers in place before removing excess material. For complete polymerization, the light guide must be changed.

5. Maintenance and Cleaning

For reasons of hygiene, we recommend using a disposable protective sleeve for each patient. Make sure to fit the protective sleeve closely to the light guide. Use either the sleeves enclosed in the delivery form or other suitable and approved single-use protective sleeves. Disinfect contaminated surfaces of the device and anti-glare cones (FD 366/Dürr Dental, Incidin Liquid/Ecolab, CaviCide/Metrex Research). Sterilize the light guide before each use if disposable protective sleeves are not used. Make sure that no liquids or other foreign substances enter the handpiece, charging base and particularly the power pack during cleaning (risk of electrical shock). Disconnect the charging base from the power source when cleaning it.

Cleaning the housing

Wipe the handpiece and handpiece holder with a customary aldehyde-free disinfecting solution. Do not clean with highly aggressive disinfecting solutions (e.g. solutions based on orange oil or with an ethanol content of more than 40%), solvents (e.g. acetone), or pointed instruments, which may damage or scratch the plastic. Clean dirty plastic parts with a soapy solution.

Pre-treating the light guide

Before cleaning and/or disinfecting the light guide, pretreat it. This applies to both automated and manual cleaning and disinfection:

- Remove substantial contamination immediately after use or 2 hours after that at the latest. For this purpose, thoroughly rinse the light guide under running water (for at least 10 seconds). Alternatively, use a suitable aldehyde-free disinfecting solution to remove any adherent blood.
- To remove contamination manually, use a soft brush or soft cloth. Partially polymerized composite can be removed with alcohol and a plastic spatula, if necessary. Do not use sharp or pointed objects, as they may scratch the surface.

Manually cleaning and disinfecting the light guide

For manual cleaning, place the light guide in an instrument cleaning solution (e.g. ID 212 forte/ Dürr Dental) for the recommended reaction time (15 minutes). Make sure that the light guide is sufficiently immersed in cleaning solution (cleaning bath with sieve insert and lid). Observe the instructions for use of the disinfectant manufacturer when using cleaning and disinfectant solution.

When cleaning and disinfecting, please make sure that the agents used are free of:

- organic, mineral and oxidizing acids (the minimum admissible pH value is 5.5)
- alkaline solutions (the maximum admissible pH value is 11)
- oxidizing agent (e.g. hydrogen peroxide)

After the cleaning process, remove the light guide from the solution and thoroughly rinse it under running water (20 +/- 2 °C / 36 +/- 4 °F) for at least 10 seconds.

Mechanical cleaning and disinfection of the light guide (Disinfector/CDU (cleaning and disinfection unit))

As an alternative, cleaning and disinfection can be carried out mechanically: Neodisher® MediClean forte, Dr. Weigert, 0.5 %, according to cleaning program, e.g. cleaning 55 °C (131 °F) (+5 °C (41 °F)/-0 °C), 5–10 min, disinfection 90 °C (194 °F) (+5 °C (41 °F)/-0 °C), 5–10 min.

Sterilization of the light guide

Thorough cleaning and disinfecting is imperative to ensure that the subsequent sterilization is effective. Use only autoclave sterilization for this purpose: 3x pre-vacuum, sterilization time (exposure time at sterilization temperature) is 4 minutes at 134 °C (273 °F); pressure should be 2 bar (29 psi). Use nationally approved sterilization pouches. Dry the sterilized light guide (10 min) using either the special drying program of your steam autoclave or hot air. The light guide has been tested for up to 200 sterilization cycles.

Checking the light guide

After that, check the light guide for damage. Hold it against the light. If individual segments appear black, glass fibres are broken. If this is the case, replace the light guide with a new one. If you can still see signs of contamination on the light guide, the cleaning and disinfecting procedure must be repeated.

Disposal

The curing light must not be disposed of as normal household waste. Dispose unserviceable batteries and curing lights according to the corresponding legal requirements in your country.

Batteries must not be incinerated.

6. What if ...?

Indicator	Causes	Error rectification
<p>Red "x" lights up</p> 	<p>The device is overheated.</p> <p>Electronic component of the handpiece is defective.</p>	<p>Allow the device to cool down and try again after a certain time. If the error persists, please contact your dealer or your local Service Centre.</p> <p>Remove and reinsert the battery. If the error persists, please contact your dealer or your local Service Centre.</p>
<p>Red "x" and battery symbol light up</p> 	<p>Battery empty</p> <p>Battery contacts dirty</p>	<p>Place the device in the charging base and charge it.</p> <p>Remove battery and clean the battery contacts.</p>
<p>3s flashes</p> 	<p>In the 3sCure program, the curing time is set to 3 seconds and cannot be changed. The 3sCure program can be used no more than twice in direct succession.</p>	<p>After two consecutive curing cycles, an extended waiting time of 30 seconds is necessary for safety reasons before the next curing cycle can be started.</p> <p>If for compelling reasons an additional curing cycle is required before the prescribed interval of 30 seconds has elapsed, the 3sCure program can be activated again by selecting it anew with the program selection button.</p>
<p>The charging base is not illuminated during charging</p>	<ul style="list-style-type: none"> - Power pack not connected or defective - Battery fully charged 	<p>Check if the power pack is correctly positioned in the charging base or if the power pack is connected to the power supply by means of the power cord.</p>

7. Warranty / Procedure in Case of Repair

The warranty period for Bluephase PowerCure is 3 years from the date of purchase (battery: 1 year). Malfunctions resulting from faulty material or manufacturing errors are repaired free of charge during the warranty period. The warranty however does not confer rights to compensation for any material or non-material damage. The apparatus must only be used for the intended purposes. Any other uses are contraindicated. The manufacturer does not accept any liability resulting from misuse and warranty claims cannot be accepted in such cases. This is particularly valid for:

- Damage resulting from improper handling, especially incorrectly stored batteries (see Product Specifications: Transportation and storage conditions).
- Damage to components resulting from wear under standard operating conditions (e.g. battery).
- Damage resulting from external influences, e.g. blows, drops to the floor.
- Damage resulting from incorrect set-up or installation.
- Damage resulting from connecting the unit to a power supply, the voltage and frequency of which do not comply with the ones stated on the rating plate.
- Damage resulting from improper repairs or modifications that have not been carried out by certified Service Centres.

In case of a claim under warranty, the complete apparatus (handpiece, charging base, battery, power cord and power pack) must be returned, carriage paid, to the dealer or directly to Ivoclar Vivadent, together with the purchase document. Use the original packaging with the corresponding cardboard inserts for transportation. Repair work may only be carried out by a certified Ivoclar Vivadent Service Centre. In case of a defect that cannot be rectified, please contact your dealer or your local Service Centre (see addresses on the reverse side). A clear description of the defect or the conditions under which the defect occurred will facilitate locating the problem. Please enclose this description when returning the apparatus.

The product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

8. Product Specifications

Light source	Ivoclar Vivadent Polywave® LED
Wavelength range	385–515 nm
Light intensity	3sCure program: 3,050 mW/cm ² ± 10 % Turbo program: 2,100 mW/cm ² ± 10 % High Power program: 1,200 mW/cm ² ± 10 % PreCure program: 950 mW/cm ² ± 10 %
Operation	3 min on / 7 min off (intermittent)
Light guide	10>9 mm, autoclavable
Signal transmitter	Beeps every 10 seconds and every time the Start/Stop button or the Time/Program selection button is activated or the anti-glare protection is enabled or the curing process is aborted and 30-second waiting time after the 3sCure program has been used twice in succession
Dimensions of the handpiece	L = 170 mm, B = 30 mm, H = 30 mm
Weight of the handpiece	135 g (including battery and light guide)
Operating voltage of the handpiece	3.7 VDC with battery 5 VDC with power pack
Operating voltage of the charging base	5 VDC
Power supply	Input: 100–240 VAC, 50–60 Hz max 1 A Output: 5 VDC / 3 A Manufacturer: EDAC POWER ELEC. Type: EM1024B2
Operating conditions	Temperature +10 °C to +35 °C Relative humidity 30 % to 75 % Ambient pressure 700 hPa to 1060 hPa
Dimensions of the charging base	D = 110 mm, H = 55 mm
Weight of the charging base	145 g
Charging time	Approx. 2 hours (with the battery empty)
Power supply of the handpiece	Li-Ionen battery
Transportation and storage conditions	Temperature –20 °C to +60 °C Relative humidity 10 % to 75 % Ambient pressure 500 hPa to 1060 hPa The curing light has to be stored in closed, roofed rooms and must not be exposed to severe jarring. Battery: – Do not store at temperatures above 40 °C / 104 °F (or 60 °C / 140 °F for a short period). Recommended storage temperature 15–30 °C / 59–86 °F – Keep the battery charged and store no longer than 6 months.
Delivery form	1 Charging base with power cord and power pack 1 Handpiece 1 Handpiece support 1 Light guide 10>9 mm 1 Anti-glare shield 3 Anti-glare cones 1 Pack of sleeves (1x 50 pcs) 1 Instructions for Use

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Appendix

Bluephase PowerCure is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 3th (see the following tables) and IEC 60601-1-2:2014 4th Edition (according clause 7 and 8.9, tables 4 to 9). Bluephase PowerCure is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the attached EMC information.

Warning

Portable wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walki-talkie etc. can affect the Bluephase PowerCure and should be kept at least a distance of 30 cm (12 inches) to any part of Bluephase PowerCure.

Guidance and manufacturer's declaration – electromagnetic emission

The following tables are guidelines according to the medical standard IEC 60601-1-2.

Bluephase PowerCure is intended for use in the electromagnetic environment specified below. The customer or the user of Bluephase PowerCure should assure that it is used in such an environment.

Emision test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	Bluephase PowerCure uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	Bluephase PowerCure is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	N/A
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	N/A

Table: According to IEC 60601-1-2

Guidance and manufacturer's declaration – electromagnetic immunity

Bluephase PowerCure is intended for use in the electromagnetic environment specified below. The customer or the user of Bluephase PowerCure should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of typical commercial or dental environment.
Surge IEC 61000-4-5	± 1 kV line - line ± 2 kV line - earth	± 1 kV line - line no prot. earth	Mains power quality should be that of typical commercial or dental environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of typical commercial or dental environment. If the user of Bluephase PowerCure requires continued operation during power mains interruptions, it is recommended that Bluephase PowerCure be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or dental environment.

Table: According to IEC60601-1-2

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

Bluephase PowerCure is intended for use in the electromagnetic environment specified below. The customer or the user of Bluephase PowerCure should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should not be used closer to any part of Bluephase PowerCure, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			Recommended separation distance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	10 V	$d = 0.35 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m	$d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz $d = 0.70 \sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

Table: According to IEC60601-1-2

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic

environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Bluephase PowerCure is used exceeds the applicable RF compliance level above, Bluephase PowerCure should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating Bluephase PowerCure.

^b Over the frequency range 150 kHz to 80 MHz, field strength should be less than 10 V/m.

Recommended separation distances between portable and mobile RF communications equipment and "Bluephase Style PowerCure"

Bluephase PowerCure is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Bluephase PowerCure can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Bluephase PowerCure as recommended below, according to the maximum out-put power of the communication equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 0.35 \sqrt{P}$	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
0.01	0.035	0.12	0.23
0.1	0.11	0.38	0.73
1	0.35	1.2	2.3
10	1.1	3.8	7.3
100	3.5	12	23

Table: According to IEC60601-1-2

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 4: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Ivoclar Vivadent – worldwide

Ivoclar Vivadent AG
Benderstrasse 2
9494 Schaan
Liechtenstein
Tel. +423 235 35 35
Fax +423 235 33 60
www.ivoclarvivadent.com

Ivoclar Vivadent Pty. Ltd.
1 – 5 Overseas Drive
P.O. Box 367
Noble Park, Vic. 3174
Australia
Tel. +61 3 9795 9599
Fax +61 3 9795 9645
www.ivoclarvivadent.com.au

Ivoclar Vivadent GmbH
Tech Gate Vienna
Donau-City-Strasse 1
1220 Wien
Austria
Tel. +43 1 263 191 10
Fax +43 1 263 191 111
www.ivoclarvivadent.at

Ivoclar Vivadent Ltda.
Alameda Caiapós, 723
Centro Empresarial Tamboré
CEP 06460-110 Barueri – SP
Brazil
Tel. +55 11 2424 7400

Ivoclar Vivadent Inc.
1-6600 Dixie Road
Mississauga, Ontario
L5T 2Y2
Canada
Tel. +1 905 670 8499
Fax +1 905 670 3102
www.ivoclarvivadent.us

Ivoclar Vivadent Shanghai Trading Co., Ltd.
2/F Building 1, 881 Wuding Road,
Jing An District
200040 Shanghai
China
Tel. +86 21 6032 1657
Fax +86 21 6176 0968
www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.
Calle 134 No. 7-B-83, Of. 520
Bogotá
Colombia
Tel. +57 1 627 3399
Fax +57 1 633 1663
www.ivoclarvivadent.co

Ivoclar Vivadent SAS
B.P. 118
74410 Saint-Jorioz
France
Tel. +33 4 50 88 64 00
Fax +33 4 50 68 91 52
www.ivoclarvivadent.fr

Ivoclar Vivadent GmbH
Dr. Adolf-Schneider-Str. 2
73479 Ellwangen, Jagst
Germany
Tel. +49 7961 889 0
Fax +49 7961 6326
www.ivoclarvivadent.de

Ivoclar Vivadent Marketing (India) Pvt. Ltd.
503/504 Raheja Plaza
15 B Shah Industrial Estate
Veera Desai Road, Andheri (West)
Mumbai, 400 053
India
Tel. +91 22 2673 0302
Fax +91 22 2673 0301
www.ivoclarvivadent.in

Ivoclar Vivadent Marketing Ltd.
The Icon
Horizon Broadway BSD
Block M5 No. 1
Kecamatan Cisaug Kelurahan
Sampora
15345 Tangerang Selatan – Banten
Indonesia
Tel. +62 21 3003 2932
Fax +62 21 3003 2934
www.ivoclarvivadent.com

Ivoclar Vivadent s.r.l.
Via del Lavoro 47
40033 Casalecchio di Reno (BO)
Italy
Tel. +39 051 6113555
Fax +39 051 6113565
www.ivoclarvivadent.it

Ivoclar Vivadent K.K.
1-28-24-4F Hongo
Bunkyo-ku
Tokyo 113-0033
Japan
Tel. +81 3 6801 1301
Fax +81 3 5844 3657
www.ivoclarvivadent.jp

Ivoclar Vivadent Ltd.
4F TAMIYA Bldg.
215 Baumoe-ro, Seocho-gu
Seoul, 06740
Republic of Korea
Tel. +82 2 536 0714
Fax +82 2 6499 0744
www.ivoclarvivadent.co.kr

Ivoclar Vivadent S.A. de C.V.
Calzada de Tlalpan 564,
Col Moderna, Del Benito Juárez
03810 México, D.F.
México
Tel. +52 (55) 50 62 10 00
Fax +52 (55) 50 62 10 29
www.ivoclarvivadent.com.mx

Ivoclar Vivadent BV
De Fruittuinen 32
2132 NZ Hoofddorp
Netherlands
Tel. +31 23 529 3791
Fax +31 23 555 4504
www.ivoclarvivadent.com

Ivoclar Vivadent Ltd.
12 Omega St, Rosedale
PO Box 303011 North Harbour
Auckland 0751
New Zealand
Tel. +64 9 914 9999
Fax +64 9 914 9990
www.ivoclarvivadent.co.nz

Ivoclar Vivadent Polska Sp. z o.o.
Al. Jana Pawła II 78
00-175 Warszawa
Poland
Tel. +48 22 635 5496
Fax +48 22 635 5469
www.ivoclarvivadent.pl

Ivoclar Vivadent LLC
Prospekt Andropova 18 korp. 6/
office 10-06, 115432 Moscow
Russia
Tel. +7 499 418 0300
Fax +7 499 418 0310
www.ivoclarvivadent.ru

Ivoclar Vivadent Marketing Ltd.
Qlaya Main St.
Siricon Building No.14, 2nd Floor
Office No. 204, P.O. Box 300146
Riyadh 11372
Saudi Arabia
Tel. +966 11 293 8345
Fax +966 11 293 8344
www.ivoclarvivadent.com

Ivoclar Vivadent S.L.U.
Carretera de Fuencarral nº24
Portal 1 – Planta Baja
28108-Alcobendas (Madrid)
Spain
Tel. +34 91 375 78 20
Fax +34 91 375 78 38
www.ivoclarvivadent.es

Ivoclar Vivadent AB
Dalvägen 14
169 56 Solna
Sweden
Tel. +46 8 514 939 30
Fax +46 8 514 939 40
www.ivoclarvivadent.se

Ivoclar Vivadent Liaison Office
: Tesvikiyeh Mahallesi
Sakayik Sokak
Nisantas' Plaza No:38/2
Kat:5 Daire:24
34021 Sisli – Istanbul
Turkey
Tel. +90 212 343 0802
Fax +90 212 343 0842
www.ivoclarvivadent.com

Ivoclar Vivadent Limited
Compass Building
Feldspar Close
Warrens Business Park
Enderby
Leicester LE19 4SD
United Kingdom
Tel. +44 116 284 7880
Fax +44 116 284 7881
www.ivoclarvivadent.co.uk

Ivoclar Vivadent, Inc.
175 Pineview Drive
Amherst, N.Y. 14228
USA
Tel. +1 800 533 6825
Fax +1 716 691 2285
www.ivoclarvivadent.us

701090