Adhese® 100

EN Instructions for Use

Light-curing single-component dental adhesive

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2, 9494 Schaan/Liechtenstein
www.ivoclarvivadent.com

Rx ONLY

CE 0123

English

Description

Adhese® 100 is a light-curing single-component dental adhesive for enamel and dentin and is compatible with all etching techniques (self-etch, selective-enamel-etch and etch & rinse techniques).

Composition

Adhese 100 contains methacrylates, ethanol, water, highly dispersed silicon dioxide, initiators and stabilizers.

Indications

- a. Indications when conventionally light-cured (500 to 1,400 mW/cm²)**
- Direct-placed light-curing composite and compomer restorations
 - Direct-placed core build-ups with light-, self- and dual-curing composites
 - Repair of fractured composite and compomer restorations
 - Adhesive cementation of indirect restorations with light- and dual-curing luting composites
 - Sealing of prepared tooth surfaces before temporary / permanent cementation of indirect restorations
 - Desensitization of hypersensitive cervical areas
- b. Indications when light-cured using Bluephase® PowerCure in the 3sCure curing mode**
- Restorations in the posterior region of permanent dentition (Class I and II, including the replacement of individual cusps) when light-cured from the occlusal aspect

Safety notes

- Do not place light in direct contact with unprotected gingiva, mucous membrane or skin.
- The 3sCure curing mode must not be used in case of caries profunda and very deep cavities.

Contraindications

- Do not use Adhese 100 if the patient is known to be allergic to any of the materials' ingredients or if the stipulated working technique cannot be employed.
- Applications in which sufficient illumination cannot be ensured (e.g. cementation of root canal posts).
- Direct pulp cappings.

Side effects

In rare cases, components of Adhese 100 may lead to sensitization. In such cases, Adhese 100 should no longer be used.

Interactions

- Peroxide-based bleaching agents and disinfectants may inhibit the polymerization of adhesives. Such materials must not be used \leq 14 days prior to the application of Adhese 100.
- Phenolic substances (e.g. eugenol, salicylic derivatives) may inhibit the polymerization of Adhese 100.
- Before any treatment with Adhese 100, completely rinse off all iron- and aluminium-based styptic agents with water.
- Fluoridation agents must not be used prior to the application of Adhese 100.

Procedure for direct and indirect restorations

1. Isolation

Adequate relative or absolute isolation using suitable auxiliaries, such as OptraGate® or OptraDam® Plus is required.

2. Pulp protection / cavity liner

The 3sCure curing mode must not be used for restorations in areas close to the pulp. In very deep cavities, areas close to the pulp must be selectively coated with a calcium hydroxide liner (e.g. ApexCal®); and subsequently covered with pressure-resistant cement (e.g. glass ionomer cement such as Vivaglass® Liner).

3. Conditioning with phosphoric acid gel (optional)

The bond to enamel can be improved by selectively etching the enamel or by applying the "etch & rinse" technique. Unprepared enamel surfaces must be conditioned with phosphoric acid (e.g. Total Etch). Please observe the instructions for use for the phosphoric acid gel.

a. Selective-enamel-etch technique

Apply phosphoric acid gel (e.g. Total Etch) onto the enamel and allow it to react for 15–30 seconds. Then rinse thoroughly with a vigorous stream of water for at least 5 seconds and dry with oil- and water-free compressed air until the etched enamel surfaces appear chalky white.

b. Etch & rinse technique

Apply phosphoric acid gel (e.g. Total Etch) onto the prepared enamel first, and then onto the dentin. The etchant should be left to react on the enamel for 15–30 seconds and dentin for 10–15 seconds. Then rinse thoroughly with a vigorous stream of water for at least 5 seconds and dry with oil- and water-free compressed air until the etched enamel surfaces appear chalky white.

4. Handling bottle

- Dispense the desired amount of Adhese 100 into a mixing well and apply it using a disposable applicator (e.g. Vivadent® Applicator Brush).
- Protect the adhesive from light (e.g. VivaPad®).
- Carefully close the bottle immediately after each use.
- A new disposable applicator must be used for each application.

Notes

- Contamination on the bottle can be removed, e.g. by wiping with alcohol or other non-oxidizing disinfectants.
- Do not use the bottle intraorally. For hygienic reasons, suitable application aids should be used (e.g. VivaPad®).

5. Application of the adhesive

- Starting with the enamel, completely coat the tooth surfaces to be treated with Adhese 100.
- The adhesive must be scrubbed into the tooth surface for at least 20 seconds. This time must not be shortened. Applying the adhesive on the tooth surface without scrubbing is inadequate.
- Disperse Adhese 100 with oil- and moisture-free compressed air until a glossy, immobile film layer results.

Important: Avoid pooling, since this can compromise the fitting accuracy of the permanent restoration.

6. Light-curing the adhesive

- Adhese 100 can either be conventionally light-cured for 10 s (500 to 1,400 mW/cm²), or for 3 s at 2,700 to 3,300 mW/cm² in case of Class I and II restorations when light-cured from the occlusal aspect (e.g. using Bluephase PowerCure).

The instructions for use of the curing light must be observed.

light intensity (mW/cm ²)	exposure time
500–900	10 s
900–1400	10 s
1800–2200	—
2700–3300	3 s

7. Application of the restorative or luting composite

- Proceed according to the instructions for use of the restorative and/or luting composite used.

Procedure for intraoral repairs of fractured composite fillings/composite restorations with light-curing materials

- Roughen the surface to be repaired using diamond finishers and thoroughly clean with water spray and dry with oil- and water-free compressed air.
- Proceed according to the direct/indirect application for Adhese 100.

Procedure for the conditioning of Tetric® CAD restorations

- Proceed according to the instructions for use for Tetric CAD.

Warnings for Adhese 100

- Adhese 100 is irritant. Avoid contact with the skin, mucous membrane and eyes.
- In case of skin contact, immediately wash with soap and copious amounts of water.
- If the material should accidentally come in contact with the eyes, immediately rinse with copious amounts of water and consult a physician / ophthalmologist.
- In rare cases, contact with skin may lead to sensitization towards the ingredients.
- Commercial medical gloves do not provide protection against the sensitizing effects of methacrylate.
- Do not use Adhese 100 as a primer for ceramic restorative materials. A suitable ceramic primer must be used (e.g. Monobond Etch & Prime®).

Warnings for phosphoric etching gel

Phosphoric etching gel is corrosive. Avoid contact with the eyes, mucous membrane and skin (protective goggles are recommended for both the patients and practitioners). If the material should accidentally come in contact with the eyes, immediately rinse with copious amounts of water and consult a physician/ ophthalmologist.

Shelf life and storage

- Expiration date: see imprint on packaging.
- Storage temperature 2–28 °C / 36–82 °F.

Keep out of the reach of children!

For use in dentistry only!

The product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damage resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

The Material Safety Data Sheet is available online at www.ivoclarvivadent.com

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