# Tetric® 100 Bulk Fill

## EN Instructions for Use

Light-curing resin-based dental restorative material

Manufacturer:  
Ivoclar Vivadent AG  
Bendererstrasse 2, 9494 Schaan/Liechtenstein  
www.ivoclarvivadent.com

Rx ONLY

CE 0123

### English

#### 1 Intended use

##### Intended purpose

Direct restorations of posterior teeth

##### Patient target group

- Patients with permanent teeth

##### Intended users / Special training

- Dentists
- No special training required.

##### Use

For dental use only.

##### Description

Tetric® 100 Bulk Fill is a light-curing, radiopaque composite (300% Al) for the direct restorative treatment of posterior teeth (according to ISO 4049:2019 Type 1, Class 2, Group 1). Tetric 100 Bulk Fill is also suitable for restoring occlusal surfaces.

100% aluminum has a radiopacity equivalent to that of dentin and 200% aluminum is equivalent to enamel.

Tetric 100 Bulk Fill cures with light in the wavelength range of 400–500 nm and can be applied in layers of up to 4 mm.

When light-cured at a light intensity of  $\leq 1,200$  mW/cm<sup>2</sup>, Tetric 100 Bulk Fill can be used for:

- Reconstructive build-ups

When light-cured at a light intensity of  $\leq 2,000$  mW/cm<sup>2</sup>, Tetric 100 Bulk Fill can be used for:

- Restorations in the posterior region (Class I and II).

When light-cured with the Bluephase® PowerCure in the 3sCure curing mode (3,000 mW/cm<sup>2</sup>), Tetric 100 Bulk Fill can be used for:

- Restorations in the posterior region of permanent dentition (Classes I and II) when light-cured from the occlusal aspect

##### Indication

Missing tooth structure in posterior teeth (Classes I and II).

##### Contraindications

Do not use if a patient is known to be allergic to any of the ingredients of Tetric 100 Bulk Fill.

##### Limitations of use

- Due to esthetic reasons, Tetric 100 Bulk Fill is not suitable for Class III and IV restorations.
- Avoid direct exposure of the gingiva, mucous membrane or skin when using the 3sCure mode.
- The 3sCure mode must not be used in case of caries profunda and very deep cavities.
- A dry working field cannot be established.
- The stipulated working procedures cannot be applied.
- The product is not intended to be reprocessed or reused.

##### Side effects

In rare cases, components of Tetric 100 Bulk Fill may lead to sensitization. The product must not be used in such cases. In order to avoid irritation of the pulp, provide areas close to the pulp with suitable pulp/dentin protection. Selectively apply a calcium hydroxide-containing material in areas close to the pulp and cover with a suitable cavity liner.

##### Interactions

Phenolic substances such as eugenol/clove oil inhibit the polymerization of methacrylate-based materials. Consequently, the application of such materials together with Tetric 100 Bulk Fill must be avoided.

Discolouration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.

##### Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

##### Composition

Barium glass, copolymer, Si-Zr mixed oxide, Bis-GMA, ytterbium trifluoride, Bis-PMA, UDMA, Bis-EMA

Total content of inorganic fillers: 53–54 vol%

Particle size of inorganic fillers: between 0.11 µm and 15.46 µm.

#### 2 Application

##### 1. Shade determination

Clean the teeth prior to shade determination. The shade is selected with the tooth still moist using a shade guide (e.g. the Tetric 100 Bulk Fill and Tetric 100 Bulk Flow shade guide).

##### 2. Isolation

Adequate relative or absolute isolation is required.

The following auxiliary aids can be used for this purpose:

- OptraGate®
- OptraDam Plus

##### 3. Cavity preparation

The cavity is prepared according to the principles of the adhesive technique, i.e. by preserving as much of the tooth structure as possible. Do not prepare sharp internal edges and angles. Do not prepare additional undercuts in caries-free areas. The dimensions of the cavity are generally determined by the extent of the caries or the size of the old restoration. For the posterior region, only the sharp enamel edges should be rounded (finishing diamonds, 25–40 µm).

Caries-free cervical defects are not prepared, only cleaned with pumice or other suitable cleaning pastes with the help of rubber cups or rotary brushes. Remove all residue in the cavity with water spray. Dry the cavity with water- and oil-free air.

##### 4. Pulp protection / Base

The 3sCure curing mode must not be used for restorations in areas close to the pulp. Do not apply a base material when using an enamel/dentin bonding agent. For very deep areas close to the pulp – cover the area selectively with a calcium hydroxide-containing material followed by a pressure-resistant cement (e.g. a glass ionomer cement). Do not cover other cavity walls, since they can be used to support the bond with an enamel/dentin adhesive.

##### 5. Placement of matrix / interdental wedge

Use a wrap around matrix for cavities affecting the proximal area or a sectional matrix band and secure it with wedges.

## 6. Conditioning / Application of the bonding agent

Condition and apply the bonding agent according to the instructions for use of the product in use. Ivoclar Vivadent recommends using the universal bonding agent Adhese® 100, a single-component, light-cured adhesive for direct and indirect bonding procedures that features compatibility with all etching techniques.

## 7. Application of Tetric 100 Bulk Fill

- In order to achieve optimum results, Tetric 100 Bulk Fill should be applied in increments of max. 4 mm (e.g. Cavifil Injector) and adapted to the cavity walls with a suitable instrument (e.g. OptraSculpt®).
- Prevent incomplete polymerization of the restoration by ensuring sufficient exposure to the curing light.
- For the recommendations regarding exposure time per increment and light intensity see Table.

Light intensity	Exposure time	
500 – 900 mW/cm <sup>2</sup>	20 s	
900 – 1,400 mW/cm <sup>2</sup>	10 s	e.g. Bluephase® 100
1,800 – 2,200 mW/cm <sup>2</sup>	5 s	
2,700 – 3,300 mW/cm <sup>2</sup>	3 s	e.g. Bluephase® PowerCure

-  **The instructions for use of the curing light must be observed.**

- When using a metal matrix, additionally polymerize the composite material from the buccal and lingual/palatal aspect after removing the matrix, if no Bluephase® curing light is used.
- If the light guide could not be positioned ideally, e.g. at a distance to the composite or at a divergent light emission angle, light-cure the composite material again (not in 3sCure mode) observing the limitations of use mentioned above.
- Optionally, a flowable composite (e.g. Tetric 100 Flow or Tetric 100 Bulk Flow) can be used as an initial layer. Cure this layer separately according to the respective instructions for use.

## 8. Finishing / Checking the occlusion / Polishing

Remove excess material with tungsten carbide or diamond finishers after polymerization. Check the occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restoration. Use polishers (e.g. OptraGloss®), polishing discs and polishing strips to polish the restoration to a high gloss.

### Application notes

- Tetric 100 Bulk Fill can be used in combination with Tetric 100 Bulk Flow and Tetric 100 Flow. Use the 3sCure mode of Bluephase PowerCure for Tetric 100 Bulk Fill or Tetric 100 Bulk Flow materials only.
- In the case of repairs, additional Tetric 100 Bulk Fill can be directly applied to the polymerized material. If the Tetric 100 Bulk Fill restoration has already been polished, it must first be roughened and wetted with Adhese 100 before fresh Tetric 100 Bulk Fill is applied.
- Tetric 100 Bulk Fill should have ambient temperature when applied. Cool temperatures render the material difficult to extrude.
-  For single use only. If Tetric 100 Bulk Fill is directly applied from the Cavifil into the oral cavity of the patient, the Cavifil should only be used for one patient due to hygienic reasons (prevention of cross-contamination between patients).
- Syringes or Cavifils should not be disinfected with oxidizing disinfection agents.
- The recommended increment thickness is based on hardness profile measurements.

## 3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: [www.ivoclarvivadent.com](http://www.ivoclarvivadent.com), and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website ([www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)).
- The current Summary of Safety and Clinical Performance (SSCP) can be found in the download section of the Ivoclar Vivadent AG website ([www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)).

### Warnings

Unpolymerized Tetric 100 Bulk Fill may have a slight irritating effect and may lead to a sensitization against methacrylates. Unpolymerized Tetric 100 Bulk Fill should not come in contact with skin, mucous membrane and eyes. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

### Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

### Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks. Some of these risks are listed below:

- Failure of the adhesive bond (loss of the filling)
- Postoperative sensitivity
- Inclusion of air bubbles during the placement of the filling
- Wear of the filling
- Heat development during the curing procedure
- Chipping, fractures
- Ingestion of material

## 4 Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)
- Close syringes/Cavifils immediately after usage. Exposure to light leads to premature polymerization
- Do not use the product after the indicated date of expiration.
- Date of expiration: see note on Cavifils, syringes and packages.

## 5 Additional information

Keep out of the reach of children!

Not all products are available in all countries.

The material has been developed solely for use in dentistry and should be used strictly according to the Instructions for Use. Liability cannot be accepted for damage resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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