Tetric® 100 Bulk Flow

EN Instructions for Use

Light-curing, flowable resin-based dental restorative material

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2, 9494 Schaan/Liechtenstein
www.ivoclarvivadent.com

Rx ONLY

CE 0123

English

Description

Tetric® 100 Bulk Flow is a flowable, light-curing radiopaque composite for the direct restorative treatment in posterior teeth. As its opacity increases during polymerization, Tetric 100 Bulk Flow is also suitable for discoloured tooth structure. It is applied as initial layer in increments of up to 4 mm in Class I and II restorations. Tetric 100 Bulk Flow cures with light in the wavelength range of 400 – 500 nm.

Shades

Tetric 100 Bulk Flow is available in the Universal shade TMA.

Composition

The monomer matrix is composed of dimethacrylates (28 wt%). The fillers include barium glass, ytterbium trifluoride and copolymers (71 wt%). Additional contents: additives, initiators, stabilizers and pigments (<1.0 wt%). The total content of inorganic fillers is 68.2 wt.% or 46.4 vol.% respectively. The particle sizes of the inorganic fillers range between 0.1 µm and 30 µm.

Indications

Conventional application (Light intensity ≤ 2,000 mW/cm²)

- As initial layer / first increment in Class I and II composite restorations in permanent posterior teeth
- Restoration of deciduous teeth

Light-curing using the 3sCure mode of Bluephase® PowerCure (Light intensity 3,050 mW/cm²)

- As initial layer / first increment in Class I and II composite restorations in permanent posterior teeth when light-cured from the occlusal aspect

Contraindications

The placement of Tetric 100 Bulk Flow restorations is contraindicated

- If a dry working field cannot be established, or the stipulated working procedures cannot be applied.
- If a patient is known to be allergic to any of the ingredients of Tetric 100 Bulk Flow.
- Tetric 100 Bulk Flow must not be used as a luting composite.

Safety note

- Do not place light in direct contact with unprotected gingiva, mucous membrane or skin.
- The 3sCure curing mode must not be used in case of caries profunda and very deep cavities.

Side effects

In individual cases, components of Tetric 100 Bulk Flow may lead to sensitization. Tetric 100 Bulk Flow should not be used in such cases. To avoid possible irritations of the pulp, areas close to the pulp should be protected with a suitable pulp/dentin protector (apply a calcium hydroxide-based preparation to areas in close proximity of the pulp and cover it with an adequate lining).

Interactions

Materials containing eugenol/clove oil inhibit the polymerization of composite materials. Consequently, application of such materials together with Tetric 100 Bulk Flow must be avoided. Discolouration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.

Application

1. Shade determination

Clean the teeth prior to shade determination. The shade is selected with the tooth still moist. The shade of the shade tab corresponds with the shade of the composite after the change in opacity has taken place, e.g. after polymerization.

2. Isolation

Adequate relative or absolute isolation using suitable auxiliaries, such as OptraGate® or OptraDam® Plus, is required.

3. Cavity preparation

Cavity preparation is carried out according to the requirements of the adhesive technique, i.e. protecting the tooth structure. Do not prepare sharp, internal edges and angles or additional undercuts in caries-free areas. The dimensions of the cavity are generally determined by the extent of the caries or the size of the old restoration. Slightly bevel or round out the enamel margins using finishing diamonds (grit size 25 – 40 µm). Subsequently, remove all residue in the cavity with water spray and dry with water- and oil-free air.

4. Pulp protection / Base

The 3sCure mode must not be used for restorations in areas close to the pulp. Do not apply a base material when using an enamel/dentin bonding agent. Only cover very deep areas close to the pulp selectively with a calcium hydroxide material (e.g. ApexCal®) and subsequently use a pressure-resistant cement (e.g. a glass ionomer cement, such as Vivaglass® Liner). Do not cover other cavity walls, since they can be used to support the bond with an enamel/dentin adhesive.

5. Apply matrix / interdental wedge

Use a wrap around matrix for cavities affecting the proximal area or a sectional matrix and wedge it.

6. Conditioning / Application of the bonding agent

Condition and apply the bonding agent according to the Instructions for Use of the product in use. Ivoclar Vivadent recommends using the universal bonding agent Adhese® 100, a single-component, light-cured adhesives for direct and indirect bonding procedures that features compatibility with all etching techniques.

7. Application of Tetric 100 Bulk Flow

- In order to achieve optimum results, Tetric 100 Bulk Flow should be applied in increments of max. 4 mm.
- Tetric 100 Bulk Flow must be covered with a layer of methacrylate-based universal or posterior composite (e.g. Tetric 100 Universal, Tetric 100 Bulk Fill, IPS Empress® Direct). Processing and finishing should be conducted according to the respective Instructions for Use.
- Proximal contacts must be established by means of a matrix system. The filling material does not exert any pressure on the matrix band. The matrix band can be shaped with a suitable hand instrument (e.g. a ball-type condenser) before and/or during light-curing.
- Sufficient exposure to the curing light prevents incomplete polymerization.

For the recommendations regarding exposure time per increment and light intensity see table:

Light intensity	Exposure time	
500 – 900 mW/cm ²	20 s	e.g. Bluephase® Style M8
900 – 1,400 mW/cm ²	10 s	e.g. Bluephase® 100
1,800 – 2,200 mW/cm ²	5 s	e.g. Bluephase® Style 20i
2,700 – 3,300 mW/cm ²	3 s	e.g. Bluephase® PowerCure

The instructions for use of the curing light must be observed.

- When using a metal matrix, additionally polymerize the composite material from the buccal and lingual/palatal aspect after removing the matrix, if no Bluephase® curing light is used.
- If the light guide could not be positioned ideally, e.g. at a distance to the composite or at a diverging scattering angle, the composite material should also be light-cured again (no 3sCure mode), but observing the above mentioned safety notes.

8. Finishing / Checking the occlusion / Polishing

After polymerization, remove excess material with tungsten-carbide or diamond finishers. Check the occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restorations. Use silicone polishers (e.g. OptraPol®) as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Additional notes

1. Tetric 100 Bulk Flow can be used in combination with Tetric 100 Flow, Tetric 100 Universal and Tetric 100 Bulk Fill. Use the 3sCure mode of Bluephase PowerCure for Tetric 100 Bulk Fill or Tetric 100 Bulk Flow materials only.
2. Tetric 100 Bulk Flow should have ambient temperature when applied. Cool temperatures render the material difficult to extrude.
3.  For single use only. If Tetric 100 Bulk Flow is applied directly into the oral cavity of the patient, the application tips should only be used for one patient due to hygienic reasons (prevention of cross-contamination between patients).
4. Syringes should not be disinfected with oxidizing disinfection agents.
5. The use of any other application cannula may make the material hard to extrude.
6. The recommended increment thickness is based on hardness profile measurements.

Warning

Unpolymerized Tetric 100 Bulk Flow should not come in contact with skin, mucous membrane and eyes. Unpolymerized Tetric 100 Bulk Flow may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)
- Close syringes immediately after usage. Exposure to light leads to premature polymerization.
- Do not use Tetric 100 Bulk Flow after the indicated date of expiration.
- Date of expiration: see information on syringes and packages.

Keep out of the reach of children!

For use in dentistry only!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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