IPS e.max® CAD

EN Instructions for Use

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Intended use

Intended purpose

All-ceramic material for dental use

Use

For dental use only.

Description

IPS e.max CAD is a tried-and-tested lithium disilicate glass-ceramic block (LS₂) for the fabrication of fixed anterior and posterior restorations.

Types of restorations:

- Veneers
- Inlays
- Onlays (e.g. occlusal veneers, partial crowns)
- Crowns
- Three-unit bridges up to the second premolar as the terminal abutment

For hybrid abutments and hybrid abutment crowns, separate instructions for use are available.

IPS e.max CAD can be processed in an authorized CAD/CAM machine¹ in the intermediate crystalline state (≥ 130 MPa). After wet-processing of the block, the restoration is crystallized in a ceramic furnace.² Due to the resulting change in the microstructure, a biaxial flexural strength of ≥ 360 MPa and the respective optical properties are achieved. The mean biaxial flexural strength from 10 years of quality measurements is 530 MPa. (According to ISO 6872:2015)

Indication

- Missing tooth structure
- Partial edentulism

Contraindication

- Patients with substantially reduced residual dentition
- Bruxism
- If the patient is known to be allergic to any of the material's constituents

Limitations of use

- Inlay, cantilever and Maryland bridges
- Width of the pontic: anterior region > 11 mm, premolar region > 9 mm
- Temporary cementation of IPS e.max CAD restorations
- Complete veneering of molar crowns
- Very deep sub-gingival preparations
- Any other uses that are not included in the range of applications

Additional limitations of use for minimally invasive crowns:

- Layer thicknesses below 1 mm
- Preparations with sharp edges
- Preparations that are not anatomically supported and feature varying layer thicknesses
- Conventional and self-adhesive cementation
- Build-up materials other than composite resin
- Absence of canine guidance
- Crowns on implants

Processing restrictions

Failure to observe the following information will compromise the results achieved with IPS e.max CAD:

- Falling short of the required minimum thicknesses
- Milling of the blocks in a non-compatible CAD/CAM system
- Crystallization in a non-approved and non-calibrated ceramic furnace
- Mixing of IPS e.max CAD Crystall./Glaze, Shades and Stains with other dental ceramics (e.g. IPS Ivocolor® Glaze, Shades and Essences).
- Layering with a veneering ceramic other than IPS e.max Ceram

IPS e.max CAD Crystall./Glaze Spray:

- Intraoral application
- Application on IPS e.max CAD restorations, if they are to be veneered with IPS e.max Ceram.
- Restorations whose inner surfaces cannot be reliably and precisely covered with IPS Objekt Fix Putty/Flow (e.g. inlays)

System requirements

IPS e.max CAD must be processed with an authorized CAD/CAM system.¹

Side effects

There are no known side effects to date.

Composition

- Lithium disilicate glass-ceramic (LS₂) – Silicate-based glass-ceramic (LS₂), Type II / Class 3 according to ISO 6872:2015 (CTE 25–500 °C/77 – 932 °F: $10.1 \pm 0.5 \times 10^{-6}$ K)
- IPS e.max CAD Crystall./Glaze Spray – Glaze powder, propellant: isobutane

Application

Processing techniques and compatible materials

- **On the blue restoration**
 - *Polishing technique (self glaze):* e.g. with OptraFine®, followed by crystallization firing without individual characterization and glazing
 - *Staining technique*
 - a) Glazing with IPS e.max CAD Crystall./Glaze Spray, followed by speed crystallization (speed crystallization and glaze firing in one step)
 - b) Staining and glazing with IPS e.max CAD Crystall./Glaze Paste/Fluo, followed by crystallization and stains/glaze firing in one step
 - c) Staining and glazing with IPS e.max CAD Crystall./Glaze Spray, followed by crystallization and stains/glaze firing in one step
- **On the tooth-coloured restoration**
 - a) *Staining technique:* Characterization/glaze firing of the tooth-coloured restorations with either IPS e.max CAD Crystall./ or IPS Ivocolor materials
 - b) *Cut-back technique:* Incisal firing using IPS e.max Ceram materials. Stain/glaze firing using IPS Ivocolor materials
 - c) *Layering technique:* Dentin/incisal firing using IPS e.max Ceram materials. Stain/glaze firing using IPS Ivocolor materials

 The respective instructions for use must be observed.

Block concept

Degree of translucency	Block sizes/shades	Processing technique				Types of restorations											
		Polishing technique 	Staining technique 	Cut-back technique 	Layering technique 	Occlusal veneer ¹⁾	Thin veneer ¹⁾	Veneer	Inlay, onlay	Partial crown	Crown	3-unit bridge	Multi-unit bridge	Hybrid abutment	Hybrid abutment crown		
HT High Translucency	I12, C14, B40, B40 L available in 20 shades	✓	✓	✓		✓	✓ ⁴⁾	✓ ⁴⁾	✓ ⁴⁾	✓ ⁴⁾	✓ ⁴⁾						
MT Medium Translucency	C14, available in 7 shades	✓	✓	✓		✓	✓ ⁴⁾	✓ ⁴⁾	✓ ⁴⁾	✓ ⁴⁾	✓ ⁴⁾			✓	✓		
LT Low Translucency	I12, C14, A14, C16, A16, B32 in 20 shades ⁵⁾	✓	✓	✓				✓ ⁴⁾		✓ ⁴⁾	✓ ⁴⁾	✓ ²⁾		✓	✓		
MO Medium Opacity	C14, A14 in MO 0 – 4 ⁵⁾			✓							✓ ³⁾			✓			
I Impulse	C14 in O1 and O2	✓	✓	✓		✓	✓	✓									

- 1) The cut-back technique must not be used for thin and occlusal veneers.
- 2) Only up to the second premolar as the distal abutment
- 3) Up to the second premolar
- 4) Max. 2 units if the IPS e.max CAD Speed Crystallization Tray is used
- 5) The shade range may vary depending on the translucency level/block size or the CAD/CAM machine.

Compatible luting materials

	Adhesive cementation e.g. Variolink® Esthetic, Multilink® Automix	Self-adhesive cementation e.g. SpeedCEM® Plus	Conventional cementation e.g. Vivaglass® CEM
Preparation requirements	Non-retentive preparation	Retentive preparation (preparation angle of 4–8°, preparation height of at least 4 mm)	Retentive preparation (preparation angle of approx. 4–8°, preparation height of at least 4 mm)
Veneers	✓	–	–
Inlays, onlays (e.g. occlusal veneers, partial crowns)	✓	–	–
Minimally invasive crowns	✓	–	–
Crowns	✓	✓	✓
Three-unit bridges up to the second premolar as the terminal abutment	✓	✓	✓

For further information, please refer to the flow chart “Pretreatment and cementation”.

Shade selection

Clean the teeth prior to shade determination. Determine the shade based on the moist and unprepared tooth or the neighbouring teeth. Determine the shade of the prepared tooth based on the largest discoloured area of the preparation. The IPS e.max Shade Navigation App helps to select the most suitable block.

www.ipsemax.com/sna

Preparation

Prepare the tooth based on the guidelines for all-ceramic restorations and make sure to observe the minimum layer thicknesses:

- No angles or sharp edges
- Shoulder preparation with a rounded internal angle and/or pronounced chamfer
- The indicated dimensions reflect the minimum thickness for IPS e.max CAD restorations.
- The thickness of the edge of the preparation, particularly in anterior teeth, must be at least 1.0 mm in order to ensure optimum processing in the CAD/CAM machine.

Minimally invasive posterior crown (adhesive cementation mandatory)	Minimally invasive anterior crown (adhesive cementation mandatory)	Inlay	Onlay	Thin veneer
Veneer	Posterior crown/bridge abutment in the premolar region	Anterior crown/bridge abutment in the anterior region	Occlusal veneer (table top)	Partial crown

Dimensions in mm

Minimum layer thickness of restorations when using the staining technique

Cementation	Mandatory adhesive cementation				Optional adhesive, self-adhesive or conventional cementation			
	Thin veneer	Inlay	Onlay (e.g. occlusal veneer, partial crown)	Minimally invasive crown in the anterior and posterior region	Crown		Bridge	
					Anterior region	Posterior region	Anterior region	Posterior region
Minimum layer thickness IPS e.max CAD – Polishing technique								
Minimum layer thickness IPS e.max CAD – Staining technique								
Incisal/occlusal	0.5	1.0 Depth of the fissures	1.0	1.0	1.5	1.5	1.5	1.5
Circular	0.4	1.0 Isthmus width	1.0	1.0	1.2	1.5	1.2	1.5
Connector dimension	–	–	–	–	–	–	16 mm ² In general: height ≥ width	

The design of bridge connectors should be extended in a vertical direction rather than in a horizontal direction.

Cut-back and layering technique

When using the cut-back or layering technique, a reduced, tooth shape-supporting framework is built up to full contour using the IPS e.max Ceram layering materials.

In the fabrication of veneered or partially veneered restorations, the available space in large preparations must be filled by appropriately dimensioning the high-strength IPS e.max CAD component and not by adding additional IPS e.max Ceram layering material.

Cementation	Mandatory adhesive cementation				Optional adhesive, self-adhesive or conventional cementation				
	Veneer	Inlay	Onlay	Partial crown	Crown			Bridge	
					Anterior region	Premolar region	Molar region	Anterior region	Premolar region
Minimum layer thickness IPS e.max CAD – Cut-back technique									
Incisal/occlusal	0.4	–	–	1.3	0.8	1.0	1.3	0.8	1.0
Circular	0.6	–	–	1.5	1.2	1.5	1.5	1.2	1.5
Connector dimension	–	–	–	–	–	–	–	16 mm ² In general: height ≥ width	
Minimum layer thickness IPS e.max CAD – Layering technique									
Incisal/occlusal	–	–	–	–	0.8	1.0	–	–	–
Circular	–	–	–	–	0.8	0.8	–	–	–
Design type	–	–	–	–	Supporting the tooth shape		–	–	–

IPS e.max CAD is the high-strength component of the restoration and must, therefore, always make up at least 50% of the total layer thickness of the restoration. The total layer thickness of the restoration (depending on the type of restoration) consists of:

Total layer thickness of the restoration	0.8	1.0	1.2	1.5	1.8	2.0	2.5	3.0
Minimum thickness of the IPS e.max CAD framework	0.4	0.5	0.6	0.8	1.0	1.1	1.3	1.6
Maximum layer thickness of the IPS e.max Ceram veneer	0.4	0.5	0.6	0.7	0.8	0.9	1.2	1.4

Dimensions in mm

Finishing

Observe the grinding instrument recommendations⁴ and the minimum and maximum layer thicknesses when finishing and contouring the ceramic structures. Adjustments by grinding should be carried out while the restoration is still in the pre-crystallized (blue) state, at low speed and exerting only light pressure to prevent delamination and chipping at the margins. Overheating of the ceramic must be avoided. Smooth out the attachment point of the block, paying particular attention to the proximal contacts. If required, carry out individual shape adjustments. Finish the functional areas (occlusal contact surfaces) of the restoration with a fine diamond to smooth out the surface structure created by the CAD/CAM process. Do not "post-separate" bridge framework connectors. This may result in undesired predetermined breaking points, which will subsequently compromise the stability of the all-ceramic restoration. If desired, the restoration can be carefully tried in intraorally in the blue state in order to adjust the occlusion/articulation. Always clean the restoration either with ultrasound in a water bath or with the steam jet before crystallization. Make sure that the restoration has been thoroughly cleaned and any milling additive residue of the CAD/CAM machine has been removed before further processing it. Residue of the milling additive remaining on the surface may result in bonding problems and discolouration. The restoration must not be blasted with Al_2O_3 or glass polishing beads.

Completion

– Polishing technique (self glaze) on the blue restoration

For polishing, observe the grinding instrument recommendations⁴. Avoid overheating of the restoration. Use diamond rubber polishers⁵ for pre-polishing, and high-gloss rubber polishers⁵ for high-gloss polishing. Clean the restoration with ultrasound in a water bath or with the steam jet. Subsequently, attach the restoration to the crystallization pin as described under "Fixation of restorations on the IPS e.max CAD Crystallization Pin". Place the ceramic structure either on the IPS e.max CAD Speed Crystallization Tray or on the IPS e.max CAD Crystallization Tray and position the tray in the centre of the furnace. The firing program is selected based on the material and firing tray used (see "Crystallization and firing parameters").

– Staining technique on the blue restoration

a) **Glazing** with IPS e.max CAD Crystall./Glaze **Spray**, followed by speed crystallization (speed crystallization and glaze firing in one step).

Position the restoration as described under "Fixation of restorations on the IPS e.max CAD Crystallization Pin". Shake the IPS e.max CAD Crystall./Glaze Spray vigorously immediately before use until the mixing ball in the container is moving freely (approximately 20 seconds). If the spray is not sufficiently shaken, mainly the propellant is discharged with a spraying burst. As a result, the restoration is insufficiently coated with glazing powder. Observe a distance of 10 cm between the nozzle and the surface to be sprayed. Hold the spray can as upright as possible during spraying. Spray the restoration from all sides with short bursts while simultaneously rotating the restoration so that an even, covering layer is created. Shake the can again between individual sprays. Subsequently, spray the restoration a second time from all sides with short bursts while simultaneously rotating the restoration so that an even, covering layer is created. Shake the can again between individual sprays. Wait briefly until the glaze layer is dry and has assumed a whitish colour. Areas that do not show an even layer have to be sprayed again. Then position max. 2 restorations on the IPS e.max CAD Speed Crystallization Tray and conduct the speed crystallization firing using the stipulated firing parameters (speed crystallization). Observe the directions under "How to proceed after firing". For adjustments, please observe the information under "Corrective firing".

b) **Staining and glazing** using IPS e.max CAD Crystall./Glaze **Paste/Fluo**, followed by crystallization and stains/glaze firing in one step

Position the restoration as described under "Fixation of restorations on the IPS e.max CAD Crystallization Pin". Subsequently, apply IPS e.max CAD Crystall./Glaze Paste/Fluo evenly on the outer surfaces of the restoration using a brush. If the ready-to-use glaze needs to be diluted, it may be mixed with a small amount of IPS e.max CAD Crystall./Glaze Liquid. Do not apply the glaze too thickly. Avoid "pooling", especially on the occlusal surface. Too thin a glaze layer may lead to an unsatisfactory gloss. If characterizations are desired, the restoration can be individualized using IPS e.max CAD Crystall./Shades and/or IPS e.max CAD Crystall./Stains before crystallization firing. Extrude the ready-to-use shades and stains from the syringe and mix them thoroughly. The shades and stains may be slightly thinned using IPS e.max CAD Crystall./Glaze Liquid. However, the consistency should still remain pasty. Apply the mixed shades and stains directly to the unfired glaze layer using a fine brush.

Then place the restoration in the centre of the IPS e.max CAD Crystallization Tray or place a maximum of 6 units on the tray and conduct the crystallization firing using the stipulated firing parameters for the respective translucency. Observe the directions under "How to proceed after firing". For adjustments, please observe the information under "Corrective firing".

c) **Staining and glazing** using IPS e.max CAD Crystall./Glaze **Spray**, followed by crystallization and stain/glaze firing in one step.

Position the restoration as described under "Fixation of restorations on the IPS e.max CAD Crystallization Pin". Extrude the ready-to-use shades or stains from the syringe and mix them thoroughly. The shades and stains may be slightly thinned using IPS e.max CAD Crystall./Glaze Liquid. However, the consistency should still remain pasty. Apply the mixed shades and stains directly to the blue restoration using a fine brush. Spray the restoration with IPS e.max CAD Crystall./Glaze Spray from all sides to coat it with an even layer. Proceed according to the instructions under a).

Then place the restoration in the centre of the IPS e.max CAD Crystallization Tray or place a maximum of 6 units on the tray and conduct the crystallization firing using the stipulated firing parameters for the respective translucency. Observe the directions under "How to proceed after firing". For adjustments, please observe the information under "Corrective firing".

Corrective firing

If additional characterizations or adjustments are required after crystallization, a corrective firing using IPS e.max CAD Crystall./Shades and Stains and Glaze can be conducted. Also use the IPS e.max CAD Crystallization Tray for the corrective firing cycle. Use the IPS e.max CAD Crystall./Add-On including the respective mixing liquid to carry out minor shape adjustments (e.g. proximal contact points). The adjustments can be carried out in both the crystallization and corrective firing sequences.

– On the tooth-coloured restoration, crystallization without the application of materials

a) **Staining technique:** Stain/Glaze firing of the tooth-coloured restoration with either IPS e.max CAD Crystall./ or IPS Ivocolor materials.

Position the restoration as described under "Fixation of restorations on the IPS e.max CAD Crystallization Pin". For characterization and glazing, use either IPS e.max CAD Crystall./Shades, Stains, Glaze (see b under "Corrective/stains/glaze firing") or IPS Ivocolor.

When using IPS Ivocolor: To ensure better wetting, a small quantity of IPS Ivocolor Mixing Liquid may be slightly rubbed into the area that needs to be characterized. Mix the IPS Ivocolor Shades and Essences to the desired consistency using the respective IPS Ivocolor Liquids. More intensive shades are achieved by repeating the staining procedure and firing, not by applying thicker layers. Use the IPS Ivocolor Shades Incisal to imitate the incisal area and create the effect of translucency on the crown in the incisal and occlusal third. The cusps and fissures can be customized using the Essences. Then place the restoration in the centre of the IPS e.max CAD Crystallization Tray or place a maximum of 6 units on the tray and conduct the crystallization firing using the stipulated firing parameters for the respective translucency.

IPS e.max CAD Crystall./Shades, Stains, Glaze and IPS Ivocolor Shades, Essence, Glaze must not be mixed with one another or applied one after the other. Observe the directions under "How to proceed after firing".

b) and c) Cut-back and layering technique: using **IPS e.max Ceram materials**. Stain/glaze firing using IPS Ivocolor materials
When performing the cut-back, the minimum layer thicknesses of the framework must be observed. See points "Finishing" and "Minimum layer thicknesses" and as described under "Fixing the restoration on an IPS e.max CAD Crystallization Pin". Then place the restoration in the centre of the IPS e.max CAD Crystallization Tray or place a maximum of 6 units on the tray and conduct the crystallization firing using the stipulated firing parameters for the respective translucency. Observe the directions under "How to proceed after firing".

Option A: Wash firing using IPS e.max Ceram: If enough space is available, conduct the wash firing with the required IPS e.max Ceram Transpa Incisal and/or Impulse material. Use the IPS Build-Up Liquids allround or soft to mix the materials. Apply the wash in a thin coat on the entire framework.

Option B: Wash firing using IPS Ivocolor: If space is limited or to enhance the in-depth chroma effect, the wash firing can be conducted using IPS Ivocolor Shade, Essence and Glaze. Mix the paste or powder with the IPS Ivocolor Mixing Liquid allround or longlife to achieve the desired consistency and apply it in a thin coat on the entire framework.

The following applies for both options: The wash (foundation) must be fired before the actual layering procedure is started. Place the honey-comb firing tray in the furnace and conduct the wash (foundation) firing using the respective parameters.

1st/2nd Dentin / Incisal firing: With the IPS e.max Ceram layering materials, both the anatomical shape is completed and a customized esthetic appearance is achieved. The IPS e.max Ceram layering materials can be mixed with either IPS Build-Up Liquid allround or soft. If required, a second firing is conducted.

Stain/glaze firing using IPS Ivocolor materials

Preparing for stain and glaze firing: Finish the restoration using diamonds and give it a true-to-nature shape and surface texture, such as growth lines and convex/concave areas. Areas which should exhibit a higher gloss after glaze firing can be smoothed out and pre-polished using silicone discs. If gold and/or silver dust was used to visualize the surface texture, the restoration has to be thoroughly cleaned with steam. Make sure to remove all gold or silver dust in order to avoid any discolouration.

Stain firing is conducted using IPS Ivocolor Shades and/or Essences, and glaze firing is conducted using IPS Ivocolor Glaze Powder/Fluo or Paste/Fluo. Depending on the situation, the firing sequences may be conducted together or separately. The firing parameters are identical.

Clean the restoration with the steam jet and dry it with oil-free compressed air. Subsequently, mix the IPS Ivocolor Shades and Essences to the desired consistency using the respective IPS Ivocolor Liquids. In order to promote wetting of the stains and glaze materials, lightly moisten the surface with IPS Ivocolor Mixing Liquid. Subsequently, apply the glazing material in an even layer on the entire restoration. Customize the cusps and fissures using IPS Ivocolor Essence. Apply minor shade adjustments on the applied glazing material using IPS Ivocolor Shades. More intensive shades are achieved by repeating the staining procedure and firing, not by applying thicker layers. The degree of gloss of the glazed surface is controlled via the consistency of the IPS Ivocolor Glaze and the applied quantity, not by means of the firing temperature. In order to achieve a higher degree of gloss, avoid excessively diluting the glaze and/or apply more glazing material. Conduct the stain and glaze firing in a ceramic furnace on the corresponding honey-comb tray using the stipulated firing parameters for the cut-back and layering technique (stain firing using IPS Ivocolor).

Fixing the restoration on an IPS e.max CAD Crystallization Pin

1. Select the largest possible IPS e.max CAD Crystallization Pin (S, M, L) that best "fills" the inside of the restoration, but does not come into contact with the surrounding crown walls.
2. Subsequently, fill the inside of the restoration with IPS Object Fix Putty or Flow up to the restoration margin. Immediately re-seal the IPS Object Fix Putty/Flow syringe after extruding the material. Once removed from the aluminium bag, the syringe is ideally stored in a resealable plastic bag or a container with a humid atmosphere.
3. Press the selected IPS e.max CAD Crystallization Pin deeply into the IPS Object Fix Putty or Flow material so that it is adequately secured.
4. Smooth out displaced auxiliary firing paste using a plastic spatula, so that the pin is securely in place and the restoration margins are optimally supported. The IPS e.max CAD Crystallization Pin can be placed directly using a small amount of IPS Object Fix Flow.
5. Prevent contamination of the outer restoration surface. Clean off any possible contamination with a brush dampened with water and then dry.

Important: For crystallization, IPS e.max CAD restorations must not be placed directly on the IPS e.max CAD Crystallization Tray and Pins, i.e. without auxiliary firing paste.

How to proceed after firing

Remove the restoration from the furnace once the firing cycle is complete (wait for the acoustic signal of the furnace) and allow it to cool to room temperature in a place protected from draft. Hot objects must not be touched with metal tongs. Remove the restoration from the hardened IPS Object Fix Putty / Flow. Remove any residue with ultrasound in a water bath or with the steam jet. Residue must not be removed by blasting with Al₂O₃ or glass polishing beads. If the restoration needs to be adjusted by grinding, make sure that no overheating of the ceramic occurs. Finally, polish the ground areas to a high-gloss finish.

Shade combination table

For the characterization and adjustment of the shade of IPS e.max CAD restorations, IPS e.max CAD Crystall./Shades, Stains or IPS Ivocolor Shades, Essences are used.

- IPS e.max CAD Crystall./Shades, Stains: For use on blue and tooth-coloured IPS e.max CAD restorations
- IPS Ivocolor Shades, Essences: For use on tooth-coloured IPS e.max CAD restorations

The combination table must be observed.

	BL1	BL2	BL3	BL4	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
IPS e.max CAD Crystall./ Shade	0				1					2				3				4		
IPS e.max CAD Crystall./ Shade Incisal	I1						I2			I1				I2						
IPS e.max CAD Crystall./ Stains	white, cream, sunset, copper, olive, khaki, mahogany																			

	BL1	BL2	BL3	BL4	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
IPS Ivocolor Shade	SD 0				SD 1	SD 2			SD 3	SD 4	SD 5			SD 6				SD 7	SD 6	
IPS Ivocolor Shade Incisal	SI 1						SI 2				SI 3									
IPS Ivocolor Essence	E 01 white		E 02 creme		E 03 lemon		E 04 sunset		E 05 copper		E 06 hazel		E 07 olive		E 08 khaki		E 09 terracotta		E 10 mahogany	
	E 11 cappuccino		E 12 espresso		E 13 terra		E 14 profundo		E 15 ocean		E 16 sapphire		E 17 anthracite		E 18 black		E 19 rose		E 20 coral	
	E 21 basic red						E 22 basic yellow				E 22 basic blue									

Crystallization and firing parameters

Ceramic furnaces without controlled (long-term) cooling function cannot be used. The ceramic furnace must be calibrated before the first crystallization and regularly every six months thereafter. Depending on the mode of operation, more frequent calibration may be required. Observe the manufacturer's directions.

Crystallization MO, Impulse, LT, MT, HT

with or without the application of IPS e.max CAD Crystall./ materials

Furnaces Programat	Stand-by temperature B [°C/°F]	Closing time S [min]	Heating rate t1 [°C/°F/min]	Firing temperature T1 [°C/°F]	Holding time H1 [min]	Heating rate t2 [°C/°F/min]	Firing temperature T2 [°C/°F]	Holding time H2 [min]	Vacuum 1 11 [°C/°F] 12 [°C/°F]	Vacuum 2 21 [°C/°F] 22 [°C/°F]	Long-term cooling L [°C/°F]	Cooling rate tL [°C/°F/min]
P300 P500 P700	403/757	6:00	60/108	770/1418	0:10	30/54	850/1562	10:00	550/770 1022/1418	770/850 1418/1562	700/1292	0
P310 P510 P710	403/757	6:00	60/108	780/1436	0:10	30/54	860/1580	10:00	550/780 1022/1436	780/860 1436/1580	710/1310	0
CS/CS2/ CS3/CS4	Select the corresponding program											

Crystallization LT, MT, HT

with or without the application of IPS e.max CAD Crystall./ materials

Furnaces Programat	Stand-by temperature B [°C/°F]	Closing time S [min]	Heating rate t1 [°C/°F/min]	Firing temperature T1 [°C/°F]	Holding time H1 [min]	Heating rate t2 [°C/°F/min]	Firing temperature T2 [°C/°F]	Holding time H2 [min]	Vacuum 1 11 [°C/°F] 12 [°C/°F]	Vacuum 2 21 [°C/°F] 22 [°C/°F]	Long-term cooling L [°C/°F]	Cooling rate tL [°C/°F/min]
P300 P500 P700	403/757	6:00	90/162	820/1508	0:10	30/54	840/1544	7:00	550/820 1022/1508	820/840 1508/1544	700/1292	0
P310 P510 P710	403/757	6:00	90/162	830/1526	0:10	30/54	850/1562	7:00	550/830 1022/1526	830/850 1526/1562	710/1310	0
CS/CS2/ CS3/CS4	Select the corresponding program											

Speed crystallization (observe the block concept)

Max. 2 units *with or without* the application of IPS e.max CAD Crystall./Glaze Spray on an IPS e.max CAD Speed Crystallization Tray

Furnaces Programat	Stand-by temperature B [°C/°F]	Closing time S [min]	Heating rate t1 [°C/°F/min]	Firing temperature T1 [°C/°F]	Holding time H1 [min]	Heating rate t2 [°C/°F/min]	Firing temperature T2 [°C/°F]	Holding time H2 [min]	Vacuum 1 11 [°C/°F] 12 [°C/°F]	Vacuum 2 21 [°C/°F] 22 [°C/°F]	Long-term cooling L [°C/°F]	Cooling rate tl [°C/°F/min]
P300 P500 P700	403/757	1:00	110/198	800/1472	0:00	50/90	850/1562	3:00	690/800 1274/1472	800/850 1472/1562	700/1292	40/72
P310 P510 P710	403/757	0:30	120/216	850/1562	0:00	70/126	870/1598	3:30	690/850 1274/1562	850/870 1562/1598	705/1301	0
CS/CS2/ CS3/CS4	Select the corresponding program											

Corrective firing/stain firing/glaze firing

with IPS e.max CAD Crystall./ materials

Furnaces Programat	Stand-by temperature B [°C/°F]	Closing time S [min]	Heating rate t1 [°C/°F/min]	Firing temperature T1 [°C/°F]	Holding time H1 [min]	Heating rate t2 [°C/°F/min]	Firing temperature T2 [°C/°F]	Holding time H2 [min]	Vacuum 1 11 [°C/°F] 12 [°C/°F]	Vacuum 2 21 [°C/°F] 22 [°C/°F]	Long-term cooling L [°C/°F]	Cooling rate tl [°C/°F/min]
P300 P500 P700	403/757	6:00	90/162	820/1508	0:10	30/54	840/1544	3:00	550/820 1022/1508	820/840 1508/1544	700/1292	0
P310 P510 P710	403/757	6:00	90/162	830/1526	0:10	30/54	850/1562	3:00	550/830 1022/1526	830/850 1526/1562	710/1310	0
CS/CS2/ CS3/CS4	Select the corresponding program											

Firing parameters for the staining technique

with IPS Ivocolor Shade, Essence, Glaze

	Stand-by temperature B [°C/°F]	Closing time* S [min]	Heating rate t↗ [°C/°F/min]	Firing temperature T [°C/°F]	Holding time H [min]	Vacuum 1 V1 [°C/°F]	Vacuum 2 V2 [°C/°F]	Long-term cooling** L [°C/°F]	Cooling rate tl [°C/°F/min]
Stain and Glaze firing	403/757	IRT/ 6:00	60/108	710/1310	1:00	450/842	709/1308	0	0

* IRT normal mode

** Note: If the layer thicknesses exceed 2 mm, long-term cooling to 500 °C/932 °F is required.

Note: Due to their geometry, the restorations may feature varying layer thicknesses. When the objects cool after the firing cycle, the different cooling speeds in the areas with different thicknesses may result in a build-up of internal tension. In the worst case, these internal tensions may result in fractures in the ceramic objects. By using slow cooling (long-term cooling L), these tensions can be minimized. For monolithic restorations (staining technique) featuring layer thicknesses of more than 2 mm, long-term cooling L must be used.

Firing parameters for the cut-back and layering technique

with IPS e.max Ceram/IPS Ivocolor Shade, Essence, Glaze

	Stand-by temperature B [°C/°F]	Closing time* S [min]	Heating rate t↗ [°C/°F/min]	Firing temperature T1 [°C/°F]	Holding time H1 [min]	Heating rate t↗ [°C/°F/min]	Firing temperature T2 [°C/°F]	Holding time H2 [min]	Vacuum 1 11 [°C/°F] 12 [°C/°F]	Vacuum 2 21 [°C/°F] 22 [°C/°F]	Long-term cooling L [°C/°F]	Cooling rate tl [°C/°F/min]
Wash firing (foundation)	403/757	IRT/ 04:00	90/162	650/1202	00:00	20/36	730/1346	02:00	400/752 650/1202	650/1202 729/1344	0	0
1 st Dentin and Incisal firing	403/757	IRT/ 04:00	90/162	650/1202	00:00	20/36	730/1346	02:00	400/752 650/1202	650/1202 729/1344	0	0
2 nd Dentin and Incisal firing	403/757	IRT/ 04:00	90/162	650/1202	00:00	20/36	730/1346	02:00	400/752 650/1202	650/1202 729/1344	0	0
Stain firing using IPS Ivocolor	403/757	IRT/ 06:00	60/108	710/1310	01:00	-	-	-	450/842	709/1308	0	0
Glaze firing using IPS Ivocolor	403/757	IRT/ 06:00	60/108	710/1310	01:00	-	-	-	450/842	709/1308	0	0
Add-On with Glaze firing	403/757	IRT/ 06:00	60/108	710/1310	01:00	-	-	-	450/842	709/1308	0	0
Add-On after Glaze firing	403/757	IRT/ 06:00	50/90	700/1292	01:00	-	-	-	450/842	699/1290	0	0

* IRT normal mode

- ¹ e.g. PrograMill, CEREC/inLab, PlanMill. The complete list is available at www.ivoclarvivadent.com. CEREC/inLab, PlanMill, CEREC SpeedFire are not registered trademarks of Ivoclar Vivadent AG.
- ² e.g. Programat CS4, CS3, CS2, CS, CEREC SpeedFire
- ³ The availability of shades may vary depending on the degree of translucency/block size or CAD/CAM machine.
- ⁴ Ivoclar Vivadent flow chart "Recommended grinding instruments for extraoral and intraoral use".
- ⁵ OpraFine F
- ⁶ OpraFine P

Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclarvivadent.com, and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclarvivadent.com).

Warnings

- IPS Natural Die Material Separator contains hexane. Hexane is highly flammable and detrimental to health. Avoid contact of the material with the skin and eyes. Do not inhale vapours and keep away from sources of ignition.
- Do not inhale ceramic dust during finishing. Use an extraction unit and wear a face mask.
- IPS e.max CAD Crystall./Glaze Spray: Do not inhale the spray mist. Container is pressurized. Protect from direct sunlight and temperatures over 50°C/122 °F. Do not forcibly open or burn the can, even after it has been completely emptied. Observe the Safety Data Sheet (SDS).

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Shelf life and storage

This product does not require any special storage conditions.

Additional information

Keep out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions.

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.