Description

Telio® Inlay and Telio® Onlay are light-curing, single-component materials for temporary restorations, generally not requiring the application of a temporary cement. Both materials are available in a transparent (T) and universal (Uni) shade and are supplied in syringes and Cavifils.

Composition– **Telio Inlay**

The monomer matrix consists of monofunctional and difunctional methacrylates (43.4 wt%). The fillers are highly dispersed silicon dioxide and copolymers (55.8 wt%).

Fluoride (1500 ppm), initiators, stabilizers and pigments (0.3 wt%) are additional ingredients.

– **Telio Onlay**

The monomer matrix consists of monofunctional and difunctional methacrylates (36.3 wt%). The fillers are highly dispersed silicon dioxide and copolymers (62.3 wt%).

Fluoride (1500 ppm), initiators, stabilizers and pigments (0.3 wt%) are additional ingredients.

Indication– **Telio Inlay**

- Temporization of deep inlay preparations with parallel walls
- Sealing of implant screw access holes
- Relining of prefabricated, temporary crowns and bridges made of polycarbonate or methacrylates

– **Telio Onlay**

- Temporization of large, shallow cavities (onlay preparations). Please note that after curing, Telio Onlay is less elastic than Telio Inlay, which may complicate removal if undercuts are present.

Contraindication

Telio Inlay and Telio Onlay must not be used:

- for the fabrication of temporary crowns and bridges
- for the fabrication of long-term provisional restorations (worn longer than 6 weeks)
- if the patient is known to be allergic to any of the ingredients of Telio Inlay or Telio Onlay

Side effects

After longer periods of service, the material may demonstrate discoloration. However, this does not compromise the performance of the material. In rare cases, Telio Inlay and Telio Onlay may cause a sensitizing reaction in patients with a hypersensitivity to any of the materials' ingredients. In these cases, the materials should no longer be used.

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Interactions

Substances containing eugenol/clove oil may inhibit the polymerization of Telio Inlay and Telio Onlay. Therefore, Telio Inlay and Telio Onlay should not be applied in combination with zinc oxide-eugenol cements.

Usage

Application

1. Prepare, place base or liner and take impression as usual.
2. Clean cavity with water spray.
3. Isolate prepared teeth; use a matrix if necessary.
4. Application of Telio Desensitizer (optional):
 - a) The dentin surfaces should be clean and dry.
 - b) Apply Telio Desensitizer and brush the material into the dentin for 10 seconds using a suitable instrument (brush, application brush).
 - c) Cautiously disperse the desensitizer and dry it with oil-free air.
5. Place a sufficient quantity of Telio Inlay or Telio Onlay into the cavity using a spatula or another suitable instrument (e.g. OptraSculpt). If you use Cavifils, the material can be applied directly into the cavity by means of an application instrument (e.g. Cavifil® Injector). Spread out the material towards the preparation margin and remove any excess.
6. Increments of up to 4 mm can be cured in 10 seconds per surface with a standard curing light (at approx. 650 mW/cm², e.g. Bluephase using the LOW mode).
7. Remove matrix if required.
8. Check occlusion and adjust if required.

Removal

1. Insert a suitable instrument (e.g. probe/scaler) into the temporary restoration. Remove Telio Inlay or Telio Onlay from the cavity.
2. Subsequently, clean the cavity (with an oil- and fluoride-free cleaning paste, e.g. Proxyl® fluoride free).
3. Try in the final restoration and cement it as usual.

Notes

-  If Telio Inlay / Telio Onlay is applied from the Cavifil directly into the oral cavity of the patient, we recommend using the Cavifil only once due to hygienic reasons (prevention of cross-contamination between patients).
- Do not use oxidizing disinfectants to disinfect syringes and Cavifils.
- Telio Inlay / Telio Onlay remains well in place in standard inlay preparations. By placing a wooden wedge prior to the application of the material, tighter proximal contacts and thus improved retention can be achieved.
- Telio Inlay / Telio Onlay bonds with light-curing base/lining materials because of the similar chemical compositions. If such materials are utilized, isolate the base with glycerine gel (e.g. Liquid Strip). This prevents the base/lining from being accidentally removed together with Telio Inlay / Telio Onlay.
- To make adjustments by grinding or remove excess material, use suitable polishers (e.g. OptraPol®) or tungsten carbide finishers. Excess material may also be removed with a scalpel or a scaler.
- Telio Inlay / Telio Onlay can be contoured more easily by wetting the instrument with Telio Desensitizer or with an unfilled bonding agent (e.g. Heliobond).

Warning

Contact of uncured material with the skin, mucous membrane or eyes should be avoided, as this may lead to a slight irritation of the skin or, in very rare cases, to a sensitizing reaction. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage

- Do not use Telio Inlay / Telio Onlay after the indicated date of expiration.
- Shelf life: see note on Cavifils, syringes or packages.
- Close syringes / Cavifils immediately after usage. Exposure to light causes premature polymerization.
- Storage temperature: 2–28 °C / 36–82 °F

Keep the material out of the reach of children.

For use in dentistry only.

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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