### EN Instructions for Use

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## Instructions for Use

### Description

The OptraSculpt® system consists of two instruments and disposable elastic modelling tips in six different shapes. OptraSculpt allows easy application and manipulation of all composite materials due to the low adhesion of the modelling tips to composite.

### Shapes

#### Shapes of the handles

### Modelling tips

#### Recommended applications:

- 1 Cylinder**  
Initial increments of Class I & II restorations.
- 2 Chisel**  
Marginal ridges, cervical margins of Class V fillings.
- 3 Spatula**  
Class V fillings, all facial and buccal surfaces (or smooth surfaces in general). Also recommended for the removal of excess cement.
- 4 Sphere**  
Developing concavity of the lingual surface on anterior restorations.

Rx ONLY

Date information prepared:  
2015-12-14/Rev. 3  
590853/WE3

 **Manufacturer**  
Ivoclar Vivadent AG  
FL-9494 Schaan/Liechtenstein  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

CE 0123

  
ivoclar  
vivadent®  
clinical

#### 5 Point

Creating detailed anatomy such as pits and fissures, lingual contours, refining margins.

#### 6 Pyramid

Creating detailed anatomy such as pits, fissures, and other detailed morphology.

#### Indication

Sculpting and contouring anatomy in direct placement, composite resin restorations.

#### Use

**NOTE** – Brush NEW instruments with soapy water, rinse and dry. Then, carry out a complete hygiene cycle (see "Disinfection, cleaning and sterilization").

Choose and attach the appropriate modelling tip to the instrument.

If you do not dispense the material directly from the capsule, use the end with the metal spatula to place the material into the preparation. Use the elastic modelling tips to adapt, contour and sculpt the material.

The different shapes can be used as follows:

#### Cylinder

The cylinder-shaped modelling tips are used for placing and levelling the first increment in Class I and II cavities. They provide optimum adaptation to the tooth structure.

#### Chisel

The chisel-shaped modelling tips can be used for Class V restorations, labial and buccal surfaces as well as occlusal surfaces.

#### Spatula

Spatula-shaped modelling tips are particularly suitable for contouring the labial surfaces of anterior restorations (Class III and IV) and Class V restorations.

#### Sphere

Sphere-shaped modelling tips can be used to develop concavities on the lingual surface of anterior restorations (Class III and IV)

#### Point

The pointed modelling tips have been designed for contouring and sculpting the occlusal anatomy (Class I and II) as well as the lingual surfaces of anterior restorations (Class III and IV).

#### Pyramid

It is recommended to use the pyramid-shaped tips to accentuate fissures (Class I and II) and to create detailed morphology.

#### Disinfection, cleaning and sterilization

##### General notes:

Clean and disinfect every instrument before each use. It is strongly recommended that the instruments are additionally sterilized in a steam sterilizer. Do not use damaged instruments. Only suitably trained staff should carry out repair and maintenance work.

Do not place the instruments in NaCl solution (to avoid the risk of pitting and/or stress corrosion cracking).

Protective gloves and goggles complying with the requirements of directive 89/686/EEC must be worn to handle any used and contaminated instruments.

Disinfect and clean instruments immediately after use.

**Generally, processing in a washer-disinfector unit is preferable to manual cleaning.**

#### **Pre-cleaning:**

Carefully remove large unhardened residues from the instruments using a lint-free cellulose tissue. For preliminary disinfection, place the instruments in a disinfection bath (with a lid and sieve insert) filled with aldehyde-free alkaline disinfectant (e.g. immersion in ID 212 forte, Dürr Dental, for 15 min.). Make sure that the instruments are completely covered in disinfectant and that the disinfectant is free of bubbles.

Only use a soft brush for the manual removal of residues. Never use a metal brush or steel wool.

For automated processing, carefully rinse the instruments under running tap water after preliminary disinfection.

#### **Automated cleaning and disinfection:**

- Washer-disinfector unit (e.g. Miele thermal washer-disinfector G7882)
  - Cleaning and disinfection program (e.g. Vario TD)
  - Detergent and disinfectant (e.g. Neodisher FA, Dr. Weigert)
- If the instruments are only cleaned in a washer-disinfector, subsequent thermal disinfection in a steam sterilizer is mandatory.

#### **Manual cleaning and disinfection:**

- Soft brush (e.g. Miltex nylon cleaning brush for instruments)
- Detergent and disinfectant (e.g. immersion in ID 212 forte, Dürr Dental, for at least 15 min.)

After cleaning and disinfection, thoroughly rinse the instruments in the sieve insert under running tap water.

#### **Drying:**

- Compressed air or clean, lint-free cellulose tissue

#### **Sterilization:**

- Suitable sterile packaging (e.g. aluminium sterile container system)
- Autoclave (e.g. Systec HX-320)
- Program: steam sterilization with a fractionated pre-vacuum process at 134 °C (273 °F) for at least 4 min.

#### **Maintenance, checking and inspection:**

After cleaning and disinfecting the instruments, visually inspect them to see if they are clean, undamaged and function properly. If macroscopically visible residual debris is present, repeat the cleaning and disinfection process.

#### **Storage**

- Store instruments in dry conditions to avoid condensation build-up.
- Store OptraSculpt modelling tips in the original packaging at room temperature, protected from dust, moisture, compression and contamination.
- Expiry date: see printed inscription and label.

#### **Notes**

-  The OptraSculpt handle is autoclavable and designed for multiple use and at least 200 sterilization cycles in the autoclave. See “Disinfection, cleaning and sterilization” for the stipulated procedures.
- Before autoclaving the handle, remove the modelling tips.
-  The OptraSculpt modelling tips are intended for single use only. They must not be cleaned, disinfected, sterilized and reused.
- Since small parts generally involve a risk of ingestion and aspiration, it is advisable to use a rubber dam.
-  The Instructions for Use and detailed information on the reprocessing of handheld instruments are available on the website of Ivoclar Vivadent AG ([www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)).

#### **Keep out of the reach of children!**

#### **For use in dentistry only!**

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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**Ivoclar Vivadent Limited**
**Ivoclar Vivadent, Inc.**