Ivotion

EN Instructions for Use

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Date information prepared:
2019-06-18 / Rev.0
740209/WW

ivoclar
digital®

Please read these Instructions for Use carefully and familiarize yourself with the use of the material. The following materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the product for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

Product description

Ivotion is a PMMA-based tooth- and gingiva-shaded disc for the CAD/CAM technique. It is used for the fabrication of permanent restorations for removable denture prosthetics using a subtractive fabrication procedure.

Composition

Polymethyl methacrylate, pigments

Typical material properties

Denture base material	
Chem. description	PMMA-based denture resin
Flexural strength	≥ 65 MPa
Modulus of elasticity	≥ 2000 MPa
Residual MMA quantity	≤ 2.2 %
Water absorption	≤ 32 µg/mm ³
Water solubility	≤ 1.6 µg/mm ³
Tooth material	
Chem. description	Highly cross-linked PMMA-based material
Flexural strength	≥ 80 MPa
Water absorption	≤ 40 µg/mm ³
Water solubility	≤ 7.5 µg/mm ³

Purpose

Removable prosthetics

Indication

Removable dentures

Contraindication

Proven allergy to one of the components.

Side effects

Systemic side effects are not known to date. In individual cases local allergic reactions to methyl methacrylate materials have been reported.

CAD/CAM processing

For the fabrication of maxillary and mandibular dentures, there is one disc each with a different internal geometry. When designing the dentures, the minimum wall thickness of 2 mm in the area of the denture base and 2.5 mm in the tooth area must not be underscored. If the minimum values are not observed, this may lead to failure (e.g. fracture of the restoration) or negatively affect the tooth shade.

Ivotion discs may only be processed in CAD/CAM milling systems that allow the correct transfer of a predefined disc position within a specially adjusted workflow of CAD/CAM, milling machine and disc holder (e.g. 3Shape – CAM – PM7 and PrograMill Disc Holder in conjunction with Digital Denture holder rings). Use only the milling tools recommended for the respective milling system. Please note that the improper use of milling tools, polishing brushes or steam jets may expose the material to overheating, which in turn may result to damage of the material.

Processing after milling

Finishing

After milling in the CAD/CAM device, separate the restoration from the holding bar with a fine tungsten carbide bur or a diamond separating disc. Smooth out the attachment points with fine tungsten carbide burs to prevent the risk of injury by sharp edges. Make sure not to damage the anatomical surface structure in the process. Adjustments of shape may be made with cross-cut tungsten carbide burs or diamond burs. Make sure that the minimum wall thickness is maintained even after finishing.

Polishing

Manual polishing of the restoration is carried out with rotary instruments and polishing paste.

In order to achieve an optimum surface gloss, please observe the following procedure:

- Use the corresponding handpiece speed and little pressure for polishing to prevent heat development.
- Pre-polishing is done with the handpiece / goat hair brushes and fine pumice / Universal Polishing Paste.
- Interdental areas and occlusal surfaces should be given special attention.
- High-gloss polishing is achieved using a goat hair brush, cotton or leather buff, as well as Universal Polishing Paste.

Repair / extension / relining

Dentures made of Ivotion can be repaired, extended and relined with conventional self-curing denture base resin (e.g. ProBase® Cold). The practical procedure is the same as that for conventionally fabricated dentures. The corresponding details can be found in the respective Instructions for Use.

Cleaning and disinfection

Clean the dentures under running water using a brush. The dentures are disinfected with a suitable, customary disinfectant.

Important processing restrictions

The following stipulations must be followed for successful work with Ivotion:

- Observing the necessary minimum wall thickness (2 mm denture base, 2.5 mm tooth area)
- Milling of the discs using a compatible CAD/CAM system
- Staining / layering using materials that are approved or recommended by Ivoclar Vivadent

General notes

- When using Ivotion, make sure to select the correct disc for the respective jaw (see icon for the maxilla or mandibula on the disc/packaging).
- The tooth area of Ivotion follows the A-D shade system.
- Make sure that the positioning of the restoration in the CAD/CAM software is as precise as possible. If this is not the case, the pink and white esthetics may slightly vary.
- Prevent overheating of the material during finishing or with the alcohol torch.
- CNC processing without coolant liquid is recommended.
- Do not use any worn down or blunt cutters.
- Contact with solvents or monomer may lead to white discolouration.
- Do not inhale grinding dust during finishing – use exhaust air discharge and mouth protection.
- Varnishing with a glaze varnish is not recommended.
- After having been removed from the package, the discs must be protected from direct sunlight to prevent slight shade changes in the denture base material.
- Keep material out of children's reach.
- Failure to observe the stipulated limitations of use and processing instructions may lead to failure.

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.