# ProArt Print Model

## EN Instructions for Use

- Stereolithography model material

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## English

### Product description

ProArt Print Model is a light-curing resin for stereolithography printing. The material can be processed in PrograPrint printers (Please observe the corresponding Operating Instructions).

### Composition

> 95% (meth-)acrylate

< 5% initiators, stabilizers, colouring agents, pigments, fillers

### Indications

Production of jaw models

### Contraindications

The material should not be used if the user is known to be allergic to any of the ingredients of PrograPrint Model. The material is not approved for use in the patient's mouth.

### Safety notes

- Contains (meth-)acrylates.
- Irritating to eyes, respiratory organs and skin.
- Skin contact may cause sensitization.
- Avoid skin contact with unpolymerized resin and cleaning solutions. Commercial medical gloves do not provide protection against the sensitizing effect of (meth-)acrylates.
- Do not inhale grinding dust.

### Storage

- Store the material in a cool, dark and well-ventilated place.
- Storage temperature 2–28 °C (36–82 °F)
- Do not use the product after the indicated date of expiration.
- Store opened resin containers in the closed PrograPrint cartridge.
- To avoid premature curing, liquid resin should be exposed to daylight only briefly.

### Preparation for the printing procedure

ProArt Print Model is suitable for processing in PrograPrint printers. Before inserting a new resin container into the PrograPrint cartridge, shake it vigorously for two minutes. Subsequently, remove the screw cap of the resin container and snap on the dispensing valve included in the delivery. Make sure to correctly align the valve. The valve must be in the "Closed" position. The resin container can now be inserted into the PrograPrint cartridge. If the printer is not used for several weeks, components of the material may deposit in the PrograPrint Pool. If this is the case, the material must be manually mixed. Use the silicone wiper/scrapper included in the PrograPrint printer delivery for this purpose. Moreover, remove the resin container from the cartridge and shake it vigorously. During this procedure, the valve must remain on the cartridge and be in the "Closed" position. In order to prevent leakage during shaking, a towel can be placed over the valve. For details about the insertion of the PrograPrint cartridge into the PrograPrint printer and the operation of the device, please refer to the Operating Instructions of the PrograPrint printer. To ensure retraceability of the material batch, the batch number of each new resin container must be additionally documented.

### Processing after the printing procedure

Once the printing procedure is complete, the objects printed with ProArt Print Model must be cleaned and post-cured. For this purpose, the objects are left on the PrograPrint Stage. Cleaning and post-curing should be performed immediately after printing.

#### *Cleaning and post-curing*

- The PrograPrint Clean cleaning unit is suitable for cleaning the printed objects. For detailed information about the operation of the unit and the parameters to be used, please refer to the Operating Instructions of PrograPrint Clean.

- After the objects have been cleaned, cautiously dry them with a weak stream of oil-free compressed air. Should sticky areas or insufficiently cleaned spots still be recognizable after air-drying, the object must be re-cleaned. In individual cases (e.g. hollow model constructions positioned horizontally on the PrograPrint Stage) it may be recommendable to detach the object from the PrograPrint Stage to ensure optimum cleaning.
- Use the post-curing unit PrograPrint Cure to post-cure the printed objects. For detailed information about the operation of the unit and the parameters to be used, please refer to the Operating Instructions of PrograPrint Cure.
- Once the printed objects have been successfully cleaned and post-cured, they may be removed from the build platform. For this purpose, the tools included in the PrograPrint printer delivery can be used. Detach the objects including the support structure from the PrograPrint Stage with the help of these tools.

#### Finishing

- Use a side cutter to separate the supports structures from the object. Any remaining residues of the support structures should be removed using cross-cut tungsten carbide burs. These burs can also be used for shape adjustments.

#### Disposal

Liquid resin must not be disposed of with household refuse. Its disposal must take place in accordance with the respective local legal provisions. Prevent contamination of sewage systems, surface water and groundwater. Once the empty container has been thoroughly cleaned, it can be disposed of with household waste.

#### Keep material out of children's reach!

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

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**Ivoclar Vivadent Marketing Ltd.**
**Ivoclar Vivadent AB**
**Ivoclar Vivadent Limited**
**Ivoclar Vivadent, Inc.**