ProArt Print Wax

EN Instructions for Use

- Burn-out resin for stereolithography printing

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English

Product description

ProArt Print Wax is a light-curing resin for stereolithography printing. The material can be processed in PrograPrint printers (observe the corresponding Operating Instructions).

Composition

- > 60% (meth-)acrylate
- > 25% burning-out agent
- < 2% initiators, stabilizers, colouring agents, pigments

Indications

For the fabrication of burn-out restorations in the press technique.

Contraindications

In the case of a known allergy to any of the ingredients, the product must not be used. Not approved for use directly in the patient's mouth.

Safety notes

- Contains (meth-)acrylates
- Irritating to eyes, respiratory organs and skin
- Skin contact may cause sensitization
- Avoid skin contact with unpolymerized resin and cleaning solution. Commercial medical gloves do not provide protection against the sensitizing effect of (meth-)acrylates.
- Do not inhale grinding dust.
- Carry out burn-out procedure under extraction.

Storage

- Store the material in a cool, dark and well-ventilated place.
- Storage temperature 2–28 °C (36–82 °F)
- Do not use the products after the indicated date of expiration.
- Store opened resin containers in the closed PrograPrint cartridge.
- To avoid premature curing, liquid resin should be exposed to daylight only briefly.

Notes regarding the design

When designing the restorations, make sure to provide them with an adequate minimum wall thickness. For information on the minimum wall thickness for the press ceramic materials from Ivoclar Vivadent, please refer to the Instructions for Use of the ceramic being used.

Preparation for the printing procedure

ProArt Print Wax is suitable for processing in PrograPrint printers. Before inserting a new resin container into the PrograPrint cartridge, shake it vigorously for two minutes. Subsequently, remove the screw cap of the resin container and snap on the dispensing valve included in the delivery. Make sure to correctly align the valve. The valve must be in the "Closed" position. The resin container can now be inserted into the PrograPrint cartridge. If the printer is not used for several weeks, components of the material may deposit in the PrograPrint Pool. In this case, the resin must be manually mixed. Use the silicone wiper/scrapper included in the PrograPrint printer delivery for this purpose. Moreover, remove the resin container from the cartridge and shake it vigorously. During this procedure, the valve must be left on the cartridge and be in the "Closed" position. In order to prevent leakage during shaking, a towel can be placed over the valve. For details about the insertion of the PrograPrint cartridge into the PrograPrint printer and the operation of the device, please refer to the Operating Instructions of the PrograPrint printer. To ensure the traceability of the material batch, the batch number of each new resin container must be additionally documented.

Processing after the printing procedure

Once the printing process is complete, the ProArt Print Wax objects must be cleaned and post-cured. For this purpose, the objects are left on the PrograPrint Stage.

Cleaning and post-curing

- The PrograPrint Clean cleaning unit is suitable for cleaning the printed objects. The objects should be cleaned immediately after the printing process. For detailed information on the operation of the cleaning unit and the parameters to be used, please refer to the Operating Instructions of PrograPrint Clean.
- After the objects have been cleaned, cautiously dry them with a weak stream of oil-free compressed air. Should sticky areas or insufficiently cleaned spots still be recognizable after air-drying, the object must be re-cleaned.
- After that, the cleaned restoration should be stored in the air for approx. one hour. This drying procedure ensures a high accuracy of fit.
- Use the PrograPrint Cure unit to post-cure the printed objects. For detailed information on the operation of the unit and the parameters to be used, please refer to the Operating Instructions of PrograPrint Cure.
- Once the printed objects have been successfully cleaned and post-cured, they may be removed from the build platform. For this purpose, the tools included in the PrograPrint printer delivery can be used. Detach the objects including the support structures from the PrograPrint Stage with the help of these tools.

Finishing

- Use a side cutter to separate the support structures from the objects. Any remaining residues of the support structures should be removed using cross-cut tungsten carbide burs. These burs can also be used for shape adjustments. Conventional modelling wax can be used to add adjustments to the resin.

Note

Printed ProArt Print Wax objects should be further processed immediately after they have been post-cured to ensure a good accuracy of fit. If the printed objects cannot be further processed straight away, they must be stored in a sealed and preferably airtight plastic bag. The objects should undergo further processing within eight hours.

Printed objects must not be exposed to temperatures above 30 °C / 86 °F (e.g. cleaning with a steam cleaner) to safeguard their accuracy of fit and burn-out properties.

Spruing, investing

- The rules on how to select and attach the sprues and conduct the investment procedure depend on the press ceramic to be used. If using IPS e.max Press, please refer to the Instructions for Use of that material. IPS PressVest Premium investment material is especially suitable for the press process.

Disposal

Liquid resin must not be disposed of with the household refuse. Its disposal must take place in accordance with the respective local legal provisions. Prevent contamination of sewage systems, surface water and groundwater. Once the empty container has been thoroughly cleaned, it can be disposed of with the household waste.

Keep material out of the reach of children!

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

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Ivoclar Vivadent Marketing Ltd.
Ivoclar Vivadent AB
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.