# Cervitec® F

## EN Instructions for Use

- Protective varnish

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- Barniz protector

### Date information prepared:

2020-05-07/Rev.0

745294/3spr/USA

### Rx only

Made in Liechtenstein  
Ivoclar Vivadent AG  
Bendererstrasse 2  
9494 Schaan/Liechtenstein

  
ivoclar .  
vivadent®  
clinical

  
ivoclar .  
vivadent®  
clinical

## Description

Cervitec® F is a dental varnish for the mechanical protection of the tooth structure. The varnish system is characterized by its good moisture tolerance during application.

## Composition

Alcohol, Aqua, VA/Crotonates Copolymer, Cetylpyridinium chloride, Chlorhexidine diacetate, Ammonium fluoride, Aroma, Saccharin

Contains ammonium fluoride.

Fluoride content: 1400 ppm

## Indication

For treatment of dentinal hypersensitivity secondary to exposed dentin and root cervical surfaces.

## Side effects

Contact with the mucous membrane may result in short-term, reversible irritation.

If the varnish comes into contact with the eyes, immediately rinse with copious amounts of water and consult an ophthalmologist.

If users are known to be allergic to any of the ingredients in Cervitec F, the material should not be applied.

Chlorhexidine may cause reversible discolouration of the tooth structure and restorations.

In rare, isolated cases, chlorhexidine may cause allergic reactions.

## Application

Cervitec F is applied by the dentist, dental hygienists, dental prophylaxis assistants or professionally instructed personnel. Cervitec F is suitable for all age groups.

Generally, Cervitec F is applied every three months.

If required, it can be applied at shorter intervals.

Dosing tube 7 g: The contents are sufficient for about 30 full set of teeth.

Single Dose units: The contents are sufficient for 1 full set of teeth.

## Step-by-step description of the procedure

1. Clean the tooth surfaces.
2. Isolate the application area with cotton rolls and cotton swabs.  
The use of a saliva ejector or air syringe is optional.
3. a) Single Dose unit: Remove the foil from the Single Dose and apply the varnish directly from the Single Dose.  
b) Dosing tube: Dispense the required amount into e.g. a dappen dish and close the tube again.
4. Apply the varnish in a thin layer using e.g. Vivabrush G.
5. Allow the varnish to dry for 1 minute and then remove the cotton rolls.
6. No rinsing after application, only spitting.

Eating and drinking should be avoided for 1 hour after application of Cervitec F.

## Special note

Operators can eliminate the risk of cross-contamination by using Cervitec F in the Single-Dose vessel for just one patient. The applicator is for single use only and cannot be sterilized.

## Storage

- Carefully close the dispensing tube immediately after use.
- The product is highly flammable and should therefore be kept away from ignition sources.
- Do not use the product after the indicated date of expiration.
- Store the product at 2–28 °C (36–82 °F).

## Keep out of the reach of children!

### For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions of the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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## Composition

Alcool, Aqua, VA/Crotonates Copolymer, Cetylpyridinium chloride, Chlorhexidine diacetate, Ammonium fluoride, Aroma, Saccharin

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