Instructions for Use Mode d'emploi Instrucciones de Uso Instruções de Uso Kullanma Talimatı Инструкция

- Self-adhesive self-curing dental resin cement with light-curing option
- Composite de collage, autopolymérisant avec option photopolymérisation
- Cemento de resina dental autoadhesivo, autopolimerizable com opção de fotopolimerização
- Composto de cimentação auto-adesivo, autopolimerizável, com opção fotopolimerizável
- Diş hekimliğinde kullanılan, kendinden adeziv ve kendinden polimerize olan, ıskli polimerizasyon opsiyonlu kompozit siman
- Самоадгезивный, самотвердеющий стоматологический композитный цемент с дополнительной возможностью фотополимеризации

Description

Multilink Speed is a self-adhesive, self-curing dental resin cement with light-curing option for the cementation of indirect restorations made of high-strength all-ceramic, composite, metal-ceramic and metal. When Multilink Speed is used, additional dentin bonding agents are not required. The automix syringe enables a time-saving direct application of the luting material into the restoration.

Shades

Multilink Speed is available in the shade Transparent.

Working and setting times

The following times apply as soon as Multilink Speed has been dispensed from the automix syringe (the working and setting times depend on the ambient temperature):

	At room temperature 23 °C ± 1 °C / 73 °F ± 1 °F	Intra-orally approx. 37 °C ± 1 °C / 98 °F ± 1 °F
Working time	150 ± 30 sec	110 ± 30 sec
Setting time (without working time)	240 ± 60 sec	160 ± 40 sec

Attention

After Multilink Speed has been dispensed from the automix syringe, the setting process can be accelerated by intensive illumination (operating light, ambient light).

Composition

The monomer matrix is composed of dimethacrylates and acidic monomers. The inorganic fillers are barium glass, ytterbium trifluoride, co-polymer and highly dispersed silicon dioxide. Additional contents are initiators, stabilizers and colour pigments (< 1 %). The primary particle size of the inorganic fillers is between 0.1 µm and 7 µm. The mean particle size is 5 µm. The total content of inorganic fillers is approx. 40 vol%.

Indication

- 1) Permanent cementation on natural teeth in conjunction with restorations made of:
 - metal and metal-ceramic (inlays, onlays, crowns, bridges, endodontic posts)
 - high-strength all-ceramic: zirconium oxide, lithium disilicate and aluminium oxide ceramics (crowns, bridges, endodontic posts)
 - fibre-reinforced composite (endodontic posts)
- 2) Permanent cementation of crowns and bridges made of:
 - metal and metal-ceramic
 - high-strength all-ceramic: zirconium oxide, lithium disilicate and aluminium oxide ceramics

CE 0123

For dental use only.
Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein
Ivoclar Vivadent AG,
FL-9494 Schaan/Liechtenstein

ivoclar
vivadent®
clinical

- on implant abutments made of:
 - oxide ceramic (e.g. zirconium oxide)
 - metal (e.g. titanium)

Contraindication

Multilink Speed is contraindicated

- in situations where the preparation does not provide adequate retention (e.g. veneers);
- in situations where the shape of the implant abutment does not provide adequate retention;
- if a patient is known to be allergic to any of the ingredients of Multilink Speed;
- in general, Multilink Speed should not be applied to exposed pulp or to dentin which is close to the pulp.

Side effects

Systemic side effects are not known to date. In rare cases, allergic reactions to single components have been reported.

Interactions

Phenolic substances (e.g. eugenol, wintergreen oil) inhibit the curing process of composites. Consequently, the use of products, e.g. mouth rinses and temporary cements, which contain these components, should be avoided.

Disinfectants and bleaching agents (e.g. hydrogen peroxide) with an oxidative effect may interact with the initiator system, which in turn may impair the curing process. At least one week should be waited between the application of oxidative bleaching agents and disinfectants and adhesive cementation materials.

Alkaline jet media compromise the effect of Multilink Speed.

Application

Pre-treatment of the cavity or tooth preparation

1. Removal of the temporary restoration and thorough cleaning of the cavity

Remove possible residue of the temporary luting composite from the cavity or preparation with a polishing brush and an oil- and fluoride-free cleaning paste. Rinse with water spray. Then lightly dry with water- and oil-free air. Avoid overdrying.

Note:

Cleaning with alcohol can lead to overdrying of the dentin.

2. Try-in of the restoration and isolation

Next, the shade, accuracy of fit and occlusion of the restoration can be checked. Care should be taken when checking the occlusion of fragile and brittle ceramic objects before they are permanently seated, because there is a risk of fracture. If necessary, make corrections with fine diamonds at medium speed and with slight pressure and adequate water cooling. Polish ground surfaces. It is essential to ensure reliable isolation of the operative area with e.g. OptraGate®, cotton rolls, saliva ejector and absorbent pads. Dental hard tissue contaminated by saliva needs to be cleaned again as indicated under point 1.

Pre-treatment of implant abutment

1. Removal of temporary

Remove temporary restoration and, if present, temporary abutment. Rinse implant lumen and peri-implant gingiva.

2. Insert and check final abutment

According to the manufacturer's directions.

A retraction cord may be placed to be better able to check the accuracy of fit of the restoration and remove excess luting composite.

3. Try-in of the restoration and isolation

Next, the shade, accuracy of fit and occlusion of the restoration can be checked. Care should be taken when checking the occlusion of fragile and brittle ceramic objects before they are permanently seated, because there is a risk of fracture. If necessary, make corrections with fine diamonds at medium speed and with slight pressure. Polish ground surfaces. It is essential to ensure reliable isolation of the operative area – preferably with OptraGate, cotton rolls, saliva ejector and absorbent pads – when using adhesive cementation procedures with composites.

4. Pre-treatment of the abutment surface

- Seal screw channel of the abutment using, e.g. a temporary restorative material (System®.inlay). It is recommended to isolate the screw head from the temporary restorative material with cotton pellets.
- Rinse abutment with water spray.
- Dry abutment.

Pre-treatment and insertion of restoration

1. Cleaning and pre-treatment of the restoration

In general the manufacturer's directions must be followed.

Or else:

- Restorations made of metal and/or metal-based restorations
 - Sandblast the inner surfaces of the restoration (use the blasting parameters stated in the instructions for use of the relevant material) until a smooth surface is achieved.
 - If necessary, clean the restoration in an ultrasonic unit for approx. 1 minute.
 - Thoroughly rinse the restoration with water spray and dry with water-/oil-free air.
 - **IMPORTANT!** In order to achieve a reliable bond, do not clean the metal surfaces with phosphoric acid. In the case of contamination, the restoration surfaces may be cleaned with Ivoclean.
- Restorations made of zirconium oxide (e.g. IPS e.max® ZirCAD) or aluminium oxide ceramic
 - Sandblast the inner surfaces of the restorations (use the blasting parameters stated in the instructions for use of the relevant material, e.g. IPS e.max ZirCAD, max. 1 bar, Al₂O₃ 100 µm)
 - If necessary, clean the restoration in an ultrasonic unit for approx. 1 minute.
 - Thoroughly rinse the restoration with water spray and dry with water-/oil-free air.
 - **IMPORTANT!** In order to achieve a reliable bond, do not clean zirconium oxide surfaces with phosphoric acid. In the case of contamination, the restoration surfaces may be cleaned with Ivoclean.

- 1.3 Restorations made of lithium disilicate glass ceramic (e.g. IPS e.max® Press/CAD)
- Etch with 5% hydrofluoric acid (e.g. IPS® Ceramic Etching Gel) for 20 seconds or according to the instructions of the manufacturer.
 - Thoroughly rinse the restoration with water spray and dry with water-/oil-free air.
 - Apply Monobond® N or the silane Monobond-S to the pre-treated surfaces with a brush or microbrush and allow to react for 60 seconds. Subsequently, disperse the primer with a strong stream of air.

Note

- **IMPORTANT!** Contamination with saliva and blood must be avoided during and after the pre-treatment of the restoration (try-in).
- Clean contaminated restoration surfaces as follows, regardless of any previous conditioning:
 - After the try-in, rinse the restoration thoroughly with water spray and dry with oil-free air.
 - Apply Ivoclean with a microbrush or brush to the bonding surface of the restoration; make sure to achieve full coverage.
 - Allow Ivoclean to react for 20 seconds. Then, thoroughly rinse off with water spray and dry with oil-free air.
 - If the restoration is made of lithium disilicate glass-ceramic (e.g. IPS e.max Press/CAD), apply another coat of Monobond N or the silane Monobond-S to the cleaned surface using a brush or microbrush, allow to react for 60 seconds and then disperse the primer with a strong stream of air.

2. Application of Multilink Speed into the restoration

For each application, place a new mixing tip on the double syringe. Extrude Multilink Speed from the automix syringe and apply the desired quantity directly into the restoration. As the luting material will cure in the used mixing tip, it may serve as a seal for the contents of the syringe until needed once again (replace with a new tip just before the next use).

3. Placement of the restoration and removal of excess cement

a) solely self-curing

- Seat the restoration and retain it in place exerting uniform pressure.
- Remove excess material immediately with a microbrush/brush/foam pellet/dental floss or scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal or gingival margins).
- To remove excess material from implant abutments, special implant scalers/curettes should be used to avoid scratching of the implant surface.

b) self-curing with additional light-activation to tack-cure excess material or accelerate the curing process (quarter technique, indicated for cases with up to 2 bridge abutments = 3- to 4-unit bridge)

- Seat the restoration and fix/hold it in place exerting uniform pressure.
- Excess material is light-cured with the polymerization light (approx. 650 mW/cm², e.g. bluephase® LOW mode), for 1 second per quarter surface (mesio-oral, disto-oral, mesio-buccal, disto-buccal) at a distance of approx. 0-10mm. Subsequently, it can be easily removed with a scaler.

- Special implant scalers/curettes should be used to remove excess material from implant abutments to avoid scratching of the implant surface.
- Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal or gingival margins). After complete curing, the removal of excess can be difficult in these areas.
- Subsequently, light-cure all cement joints again for 20 seconds (approx. 1,200 mW/cm², e.g. bluephase HIGH mode).

Note

As with all composites, Multilink Speed is subject to oxygen inhibition. This means that the surface layer does not polymerize during curing, as it comes in contact with atmospheric oxygen. In order to avoid this, we recommend covering the restoration margins with glycerine gel/air block (e.g. Liquid Strip) immediately after the removal of excess. After complete polymerization, the glycerine gel/air block is rinsed off with water.

4. Finishing of the completed restoration

- Finish proximal areas with finishing and polishing strips.
- Check occlusion and function and adjust if necessary.
- Polish restoration margins with silicone polishers (OptraPol®, Astropol®, Politip®-F) or discs.

Special note

Multilink Speed should have room temperature when it is applied. Cool temperatures can render the material difficult to extrude and increase the processing and curing times. Do not disinfect the Multilink Speed automix syringe with agents that have an oxidative effect.

Warning

Unpolymerized Multilink Speed paste is an irritant. Avoid contact with the skin, mucous membrane and eyes. If Multilink Speed comes into contact with the eyes, immediately rinse with copious amounts of water and seek medical advice. If the material comes into contact with the skin, rinse with copious amounts of water. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage

- Do not use Multilink Speed after the indicated date of expiration.
- Cool storage (2–8 °C/ 36–46 °F) for Multilink Speed.
- Store Multilink Speed automix syringe with mixing tip attached after use.
- Shelf life: see expiration date.

Keep material out of the reach of children!

For use in dentistry only!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.