1 Intended use**Intended purpose**

Application of dental filling materials from mixing capsules with contents of ≤ 0.5 g

Patient target group

- Patients with permanent teeth
- Patients with deciduous teeth

Intended users / Special training

- Dentists (clinical procedure)
- Dental assistants (preparation of the clinical procedure)

No special training required.

Use

For dental use only.

Description

The Capsule Applicator is used to apply dental filling materials directly from the capsule into the tooth cavity to be restored.

Indications

None

Contraindications

None

Limitations of use

None

Side effects

There are no known side effects to date.

Interactions

None

Clinical benefit

- Direct intraoral application of materials from capsules

Composition

Stainless steel

2 Application

- Immediately after mixing, insert the capsule (filled with restoration material) into the applicator. If required, the capsule can be rotated 360°.
- Repeatedly activate the lever (pre-click) until the material becomes visible in the transparent applicator nozzle.
- Position the capsule in the mouth of the patient as desired and uniformly activate the lever to extrude the contents of the capsule.
- Avoid forcing the contents from the capsule too quickly or applying extremely high forces during application.

Capsule Applicator

| | |
|---|---|
| EN Instructions for Use Applicator for dental filling materials in capsules | FI Käyttöohjeet Kapselannosteilija hammastäyteaineille |
| DE Gebrauchsinformation Applikator für zahnärztliche Materialien in Kapseln | NO Bruksanvisning Applikator for odontologiske materialer i kapsler |
| FR Mode d'emploi Applicateur pour matériaux d'obturation dentaires en capsules | NL Gebruiksaanwijzing Applikator voor tandheelkundige vulmaterialen in capsules |
| IT Istruzioni d'uso Applicatore per materiali odontoiatrici in capsula | EL Οδηγίες Χρήσεως Πιστόλι εφαρμογής οδοντιατρικών εμφρακτικών υλικών σε κάψουλες |
| ES Instrucciones de uso Aplicador para obturaciones dentales en cápsulas | TR Kullanma Talimatı Kapsül halindeki dental dolgu materyalleri için aplikatör |
| PT Instruções de Uso Aplicador para materiais restauradores dentários em cápsulas | RU Инструкция по применению Аппликатор для стоматологических пломбировочных материалов в капсулах |
| SV Bruksanvisning Applikator för dentala fyllnadsmaterial i kapsel | PL Instrukcja stosowania Applikator do materiałów wypełniających w kapsułkach |

Observe the instructions for Use of the material being used to ensure correct processing.

- To remove the capsule from the applicator, press the release button and pull the capsule upwards. Subsequently, transfer the Capsule Applicator to the reprocessing procedure.

3 Reprocessing

Disinfection, cleaning and sterilization

General notes:

Please observe the country-specific regulations and guidelines regarding the hygiene standards and reprocessing of medical devices in dental practices.

New capsule applicators are not sterile and must be subjected to a complete reprocessing cycle before they are used for the first time. Subsequently, they must be cleaned and disinfected after every use. Oxidizing disinfectant solutions are not suitable.

Sterilization in the steam sterilizer is explicitly recommended!

Do not use damaged Capsule Applicators. Maintenance and repair work should only be carried out by specialist staff. Protective gloves and goggles complying with the requirements of directive 89/686/EEC must be worn when handling any used and contaminated applicators.

Automated reprocessing is not recommended, since the movable parts of the applicator could be damaged in the process.

Restrictions on reprocessing:

The life span of the product is mainly determined by wear and damage through use. The operator is responsible for checking the proper function of the applicator and lubricating the movable parts after each cleaning cycle using a dental instrument lubricant. Damaged applicators must no longer be used.

Preparation at the point of use:

Carefully remove large unhardened residues from the applicator immediately after use using a lint-free cellulose tissue. The Capsule Applicators may either be placed in a removal container (dry removal) or directly immersed in a combination of detergent and disinfectant (wet removal), see below, for transportation to the reprocessing area. It is advisable to transport the applicators to the reprocessing area as soon as possible. Waiting periods of several hours between use and reprocessing should be absolutely avoided.

Preparation for cleaning/disinfection:

Pretreatment of the applicators is expressly recommended to prevent infections and protect the staff. To prevent residues from drying and ensure preliminary disinfection, place the Capsule Applicator in a disinfection bath (with sieve insert and lid) filled with aldehyde-free alkaline disinfectant (e.g. All-in-One 4 Enzyme Detergent, brand: One Cleaner). Make sure that the applicator is completely immersed in disinfectant and that the disinfectant is free of bubbles. The disinfectant that you use should be suitable for the product in question and tested for effectiveness (e.g. DGHM or FDA certification or CE mark). Please adhere to the disinfectant manufacturer's instructions (e.g. concentration rate, temperature and contact time) to ensure that the disinfectant is used correctly.

Manual cleaning and disinfection:

| Step 1 | | |
|--|-----------------------|--|
| <u>Solution</u> | <u>Time (minutes)</u> | <u>Temperature</u> |
| Water | 2 | Room temperature |
| <u>Instructions</u> Carefully rinse the soiled applicator under running tap water. Use a brush with soft bristles if necessary. | | |
| Step 2 | | |
| <u>Solution</u> | <u>Time (minutes)</u> | <u>Temperature</u> |
| Enzymatic detergent e.g. All-in-One 4 Enzyme Detergent, brand: One Cleaner | 20 | Room temperature |
| <u>Instructions</u> Immerse the product in the cleaning solution and allow it to soak for the time indicated. Thoroughly clean the applicator using a soft brush until you can no longer see any contamination on the surface. It is very important to make sure that all the joints and handles are effectively cleaned. If necessary, use a brush with a small diameter. Subsequently check all the surfaces for visible signs of contamination. | | |
| Step 3 | | |
| <u>Solution</u> | <u>Time (minutes)</u> | <u>Temperature</u> |
| Water | 2 | Corresponding with hot water supply |
| <u>Instructions</u> Rinse the applicator under warm running tap water as long as required. | | |
| Step 4 | | |
| <u>Solution</u> | <u>Time (minutes)</u> | <u>Temperature</u> |
| Air | If required | Room temperature |
| <u>Instructions</u> Allow the applicator to air-dry in a clean place. Blow dry the cavities and hollow spaces. Use a clean source of air with a filter. | | |

Maintenance, checking and inspection:

After cleaning and disinfection, check the intactness, cleanliness and function of the Capsule Applicator. If macroscopically visible residual debris is present on the applicator, repeat the cleaning and disinfection process.

Sterilization:

Requirements to carry out the validated procedure:

- suitable sterile packaging (autoclave bag)
- steam sterilizer (AMSCO STEAM STERILIZER) Model MGS-GDS-1 EAGLE 2013)

Sterilization program:

| | Cycle #1 | Cycle #2 |
|-------------------------|---------------------------------|---------------------------------|
| Cycle type | Pre-vacuum, in autoclave bag | Pre-vacuum, in autoclave bag |
| Temp. (C°) | 121 | 132 |
| Impulse | 3 | 3 |
| Exposure time (minutes) | 30 | 4 |
| Drying time (minutes) | 20 | 20 |

It is strongly recommended to sterilize semi-critical medical devices of Classes A & B!

Packaging:

Autoclave the Capsule Applicator in a suitable sterile packaging (e.g. autoclave bag or aluminium sterilization container). Sterilize the items in a steam sterilizer (the unit should be compliant with EN 285 and validated according to EN ISO 17665) using a fractionated pre-vacuum at 132 °C (273 °F) for at least 20 minutes. Please also observe the directions of the device manufacturer. If you sterilize several products at the same time in the steam sterilizer, make sure that you observe the device manufacturer's stipulations on the maximum permissible load to avoid overloading. The stipulated sterilization procedure has been validated by the manufacturer as appropriate for the purpose in question according to EN ISO 17664. The manufacturer does not assume responsibility for other sterilization procedures. If other sterilization processes are used, the operator is responsible for validating them to ensure that they are effective and fit for the purpose.

4 Safety information

- The Capsule Applicator is optimally coordinated with the requirements of Ivoclar Vivadent capsules. However, it can also be used with other suitable capsule systems.
-  The Capsule Applicator is autoclavable and designed for multiple use and should be able to undergo at least 200 autoclaving cycles. See "Disinfection, cleaning and sterilization" for the stipulated procedures
- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclarvivadent.com, and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclarvivadent.com).

Warnings

None

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

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ivoclar
vivadent®

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks. Some of these risks are listed below:

- Cross-contamination if the Capsule Applicator has not been sufficiently cleaned and/or sterilized.

5 Shelf life and storage

- Store the Capsule Applicator in a dry place to prevent condensation build-up.

6 Additional information

Keep out of the reach of children!

Not all products are available in all countries.

The product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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Ivoclar Vivadent Marketing Ltd.
Ivoclar Vivadent AB
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.