PCT® FlexVest

Instructions for Use

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Mode d'emploi

Istruzioni d'uso

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Bruksanvisning

Brugsanvisning

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Bruksanvisning

Productinformatie

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English

Instructions for Use

PCT FlexVest is a carbon-free, phosphate-bonded investment material for precision alloy castings. PCT FlexVest can be used with a conventional burnout technique or a rapid burnout technique.

Indication

- Precious alloys
- Non-precious alloys
- Investing with a metal investment ring
- Ringless investment

Contraindication

- Pressed ceramics (e.g. IPS Empress)

Working procedures

Wax pattern

Fabricate wax patterns for crown and bridge frameworks.

Resin pattern

Coat patterns made of modelling resin with a thin layer of wax.

Spruing

Attach sprues of an appropriate dimension to the bridge framework, using a direct or indirect method. (See Instructions for Use)

Investing with a metal investment ring

Apply **one layer** of PCT Flex Vest Liner if metal investment rings of sizes 1 and 3 are used and **two layers** if ring sizes 6 and 9 or a non-precious alloy is used.

Slightly wet PCT FlexVest Liner. Absorb excess water by pressing a piece of dry liner onto the wet material. Line the metal investment ring with damp liner. Allow the liner endings to overlap by approx. 0.5 cm. Next, press the liner into place. Make sure that PCT FlexVest Liner is flush with the **upper** rim of the metal investment ring. Carefully position the metal investment ring lined with the **PCT FlexVest Liner** on the base former.

Note: Any wax surface tension releasing agent needs to be completely dried off. Contact of the investment material with surfactant may lead to incomplete setting and, consequently, to uneven investment surfaces.

Processing procedure

Only use distilled water to modify the concentration of PCT FlexVest Liquid/Liquid H.E. It is advisable to use a dosing syringe to take the desired amount of liquid out of the container.

Note: No not use mixing liquid if it is crystallized.

The higher the concentration of PCT FlexVest Liquid/Liquid H.E. is, the higher is the expansion rate of the investment material.

Concentration of mixing liquids

Alloys	Ceramic Alloys		C&B/BioUniversal Alloys	
	PCT FlexVest Liquid	Distilled water	PCT FlexVest Liquid	Distilled water
High Gold and reduced gold Alloys	60–80%	40–20%	50–60%	50–40%
Pd Base Alloys	70–80%	30–20%	70–80%	30–20%
Non precious Alloys (Co Cr / Ni Cr)	80–100%	20–0%	80–100%	20–0%

Mixing ratio

60 g powder	13 ml mixture PCT FlexVest Liquid/Liquid H.E. // distilled water
90 g powder	20 ml mixture PCT FlexVest Liquid/Liquid H.E. // distilled water
160 g powder	35 ml mixture PCT FlexVest Liquid/Liquid H.E. // distilled water

Note: The mixing ratios described above represent recommendations. Adjustments may be required for various alloys and casting techniques.

Mixing

Prior to blending the mixture, rinse a clean mixing bowl with water. Place the premeasured PCT FlexVest Liquid/Liquid H.E. in the mixing bowl. Add the preweighed powder to the liquid. Stir the two components with a spatula to form a homogeneous mixture. Place the mixing bowl in a vacuum device. Set the vacuum value. After the set vacuum has been reached, mix for 90 seconds. After completing the mixing procedure, maintain the vacuum for an additional 15 seconds and then flood the mixing bowl manually.

Note: Mixing bowls that are not completely clean or severely scratched can absorb up to 2 ml of mixing liquid. Extended mixing times will result in smoother casting surfaces, similar to gypsum-bonded investment materials for C&B casting alloys. Investing Pour PCT FlexVest into the investment ring at a slow and steady rate, using slight vibration. Stop vibrating the investment ring as soon as it has been completely filled. Do not touch the investment.

Note: It is advisable to pour the most intricate and least accessible parts of the pattern first to prevent air from being trapped.

Working time

The working time is 5 minutes at a room temperature of 23 °C (73 °F).

Note: The powder and liquid should always be stored in a cool, dry place. Reseal the container of the mixing liquid immediately after use. The powder and mixing liquid should have the same temperature when in use (ideally between 17 and 22 °C (63 to 72 °F)).

Note: If the working temperature is > 24 °C (> 75 °F), the setting expansion will increase considerably, while low working temperatures of < 16 °C (< 61 °F) result in a reduction in the setting expansion.

Setting time

At a room temperature of 23 °C (73 °F), the setting time is approx. 20 to 30 minutes, depending on the size of the investment ring used (speed burnout). After the thermal reaction has ended, roughen the surface of the investment. This will aid in the removal of gases during the first stages of burnout.

Note: Allow the investment to set completely before proceeding with burnout!

Rapid burnout

Quick heating can be effected up to 8 hours after setting.

Place the investment ring in an upright position in a furnace preheated to max. 850 °C (1562 °F) with the sprue hole down and elevated so that the investment ring does not touch the base plate (spacer), or alternatively, place the investment ring horizontally. If a higher temperature is required, place the investment ring in the preheated furnace at 850 °C (1562 °F) and then raise to the final temperature. The maximum temperature for the investment material is 1050 °C (1922 °F). If alloys with lower temperatures are used, place the investment ring in the furnace at 850 °C (1562 °F) and then decrease the temperature according to the manufacturer's instructions.

Note: Please observe the instructions of the corresponding alloy manufacturer regarding the preheating temperatures.

Holding time

Investment ring size 1	30 – 45 min
Investment ring size 3	45 – 60 min
Investment ring size 6	60 – 75 min
Investment ring size 9	75 – 90 min

Note: If modelling resins which burn out without leaving residue are employed, the holding time should be increased by 10 min. Observe the manufacturer's instructions! If several investment rings are placed in the preheating furnace at the same time, the holding time has to be increased accordingly.

Conventional burnout / holding times

Allow the investment material to set for at least 1 hour. Place the investment ring into the furnace at room temperature. Increase the temperature as indicated on the chart below.

Temperature	Temperature increase	Holding time 1x	Holding time 3x	Holding time 6x	Holding time 9x
23–270°C 73–518°F	7–10°C/min 13–18°F/min	30–60 min	30–60 min	30–60 min	30–60 min
270–560°C 518–1040°F	7–10°C/min 13–18°F/min	20 min	20 min	30 min	60 min
from 560°/ 1004°F to final temperature*	7–10°C/min 13–18°F/min	30 min	45 min	60 min	60 min

* (max. final temperature is 1050 °C/1920°F)

Note: If modelling resins which burn out without leaving residue are employed, the holding time should be increased by 10 minutes. Observe the manufacturer's instructions! If several investment rings are placed in the preheating furnace at the same time, the holding time has to be increased accordingly.

Casting

Cast alloy according to the alloys manufacturer's instructions.

Cooling / Divesting

After casting the alloy, allow the investment ring to cool to room temperature. Do not quench in water. Prior to divesting, wet the investment ring in water to prevent dust from forming. Use plaster pliers. Do not use a hammer for divesting. Remove residual investment from the casting by means of an ultrasonic cleaning solution or careful application of jet polishing medium/AL₂O₃.

Storage

PCT FlexVest Powder and PCT FlexVest Liquid/Liquid H.E. should always be stored in a cool, dry place and the containers should remain tightly sealed when not in use. The storage and shipping temperature of PCT FlexVest Liquid/ Liquid H.E. and PCT FlexVest Powder is 12–28 °C (54–82 °F). The materials are ideally used at temperatures between 17 and 22 °C (63 – 72 °F).

PCT FlexVest is sensitive to temperatures below freezing point. Do not use PCT FlexVest Liquid/Liquid H.E. if it exhibits a jelly-like consistency or crystals. The shelf life of PCT FlexVest Powder and PCT FlexVest Liquid/ Liquid H.E. is two years if they are stored according to the storage conditions stipulated.

Warning

Investments contain quartz. **Do not breathe dust!** Dust may cause delayed lung injury (silicosis, lung cancer). Use scissors to cut the pouches open. Prevent dust from forming when pouring the pouches' contents into a mixing bowl. Rinse empty pouches with water prior to discarding them. Always remove dust from your working space by means of a damp cleaning method. Immerse the investment ring in water for a short while after it has completely cooled to prevent dust from forming during divesting. Use suction equipment with a fine dust filter when sandblasting the castings.

Material properties (ISO 9694)

	PCT FlexVest Liquid 100 %	PCT FlexVest Liquid H.E. 100 %
Mixing ratio (100 g : 22 ml)	100 g : 22 ml	100 g : 22 ml
Flow path	approx. 14 cm	approx. 14 cm
Start of setting time	approx. 10 min	approx. 10 min
Compressive strength	approx. 5 MPa	approx. 5 MPa
Setting expansion	1.5 %	1.5 %
Thermal expansion	1.0 %	1.3 %

Delivery forms*

PCT FlexVest Trial Kit (90 g x 4 / 100ml)

PCT FlexVest 4 kg (60 g x 67)

PCT FlexVest 4 kg (90 g x 45)

PCT FlexVest 4 kg (160 g x 25)

PCT FlexVest 5 kg

PCT FlexVest Liquid 1000 ml

PCT FlexVest Liquid H.E. 1000 ml

PCT FlexVest Dosing bottle

PCT FlexVest Liner 20 m

* The delivery forms can vary in the different countries.

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damage resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent AB
Ivoclar Vivadent UK Limited
Ivoclar Vivadent, Inc.
Dentigenix Inc.