Die Heliomolar[®] Familie:

Heliomolar[®]

Heliomolar[®] Flow

Heliomolar[®] HB

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6. Literatur

Arend J, van der Zee Y

Fluoride uptake in bovin enamel and dentin from fluoride-releasing composite resin
 Quintessence International 21 (1990) 541-544

Arends J, Ruben J

Fluoride release from a composite
 Quintessence International 19 (1988) 513-514

Clinical Research Associates

Comparative performance of 21 class 2 materials at 3 years
 Clinical Research Associates Newsletter, Volume 18, Issue 5 May 1994

Clinical Research Associates

Restorative resins, new products, one year clinical performance
 Clinical Research Associates Newsletter, Volume 18, Issue 11 November 1994

Feilzer AJ, De Gee AJ, Davidson CL

Curing contraction of composites and glass-ionomer cements
 J Prosthetic Dent 59 (1988) 297-300

Knibbs PJ, Smart ER

The clinical performance of a posterior composite resin restorative material, HELIOMOLAR RO: 3-year report
 J Oral Rehabil 19(1992) 231--237

Leinfelder KF, Mirshaahidi M, Cury C, O'Neal W

An in vitro wear device for determining wear of posterior composites
 J Dent Res 70 (1991) 345

Ludin SA, Andersson B, Koch G, Rasmusson CG
Class II composite resin restorations: A three year clinical study of six different posterior composites
 Swed Dent J 14 (1990) 104-115

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Restorative resins, new products, status report #1 -clinical characteristics
 CRA Newsletter 17 (1993) 0

Alternatives for class 3 restorations (results of CRA clinical trial of 21 materials at 3 years)
 CRA Newsletter 18 (1994) 1-4

Alhadainy HA, Abdalla AI

2-year clinical evaluation of dentin bonding systems
 Am J Dent 9 (1996) 77-79

Arends J, Ruben J, Dijkman AG

The effect of fluoride release from a fluoride-containing composite resin on secondary caries: an in vitro study
 Quintessence Int 21 (1990) 671-674

Arends J, Van der Zee Y

Fluoride uptake in bovine enamel and dentin from a fluoride-releasing composite resin
 Quintessence Int 21 (1990) 541-544

Ashe MJ, Tripp GA, Eichmiller FC, George LA, Meiers JC

Surface roughness of glass-ceramic insert - composite restorations: assessing several polishing techniques
 J Am Dent Assoc 127 (1996) 1495-1500

Bayne SC, Heymann HO, Edward J, Swift JR

Update on dental composite restorations
 J Am Dent Assoc 125 (1994) 687-701

Bollen CML, Lambrechts P, Quiryen M

Comparison of surface roughness of oral hard materials to the threshold surface roughness for bacterial plaque retention: A review of the literature
 Dent Mater 13 (1997) 258-269

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- Bouschlicher MR, Cobb DS, Boyer DB
Radiopacity of compomers, flowable and conventional resin composites for posterior restorations
Oper Dent 24 (1999) 20-25
- Bryant RW, Hodge KLV
A clinical evaluation of posterior composite resin restorations
Aust Dent J 39 (1994) 77-81
- Bryant RW, Marzbani N, Hodge KLV
Occlusal margin defects around different types of composite resin restorations in posterior teeth
Oper Dent 17 (1992) 215-221
- Castline JM, TU SJ
Measured light intensities at simulated dentin/composite bonding interfaces
J Dent Res 79 (2000) 280
- Christensen GJ
Restorative Dentistry: An update for Practitioners, Educators, Examining Boards
J Am Dent Assoc 126 (1995) 1165-1168
- Christensen GJ, Christensen RP
A new technique for resoration of worn anterior teeth - 1995
J Am Dent Assoc 126 (1995) 1543-1546
- Chung K
Effects of finishing and polishing procedures on surface texture of resin composites
Dent Mater 10 (1994) 325-330
- Condon JR, Ferracane JL
Assessing the effect of composite formulation on polymerization stress
J Am Dent Assoc 131 (2000) 497-503
- Condon JR, Ferracane JL
Polymerization contraction stress of commercial composites
J Dent Res 77 (1998) 639
- Condon JR, Ferracane JL
Evaluation of composite wear with a new multi-mode oral wear simulator
Dent Mater 12 (1996) 218-226
- Davidson CL, Abdalla AI
Effect of thermal and mechanical load cycling on the marginal integrity of Class II resin composite restorations
Am J Dent 6 (1993) 39-42
- De Backer J, Dermaut L
Visible light sources and posterior visible light cured resins: a practical mixture
Quintessence Int 17 (1986) 635-641
- De Goes MF, Garcia-Godoy F, Cardenas L
Linear and volumetric polymerization shrinkage fo the restorative materials
J Dent Res 78 (1999) 484
- Dietschi D, Ciucchi B, Holz J
A clinical trial of four light curing posterior composite resins: 9-month report
Quintessence Int 20 (1989) 641-652
- Dietschi D, Holz J
Klinische Untersuchung von vier lichthärtenden Kompositmaterialien für den Seitenzahnbereich - Ergebnisse nach zwei Jahren (1)
Quintessenz 42 (1991) 743-751
- Dietschi D, Magne P, Holz J
Recent trends in esthetic restorations for posterior teeth
Quintessence Int 25 (1994) 659-677
- Dijkman GEHM, Arends J
Secondary Caries in situ around Fluoride-Releasing Light-Curing Composites: A Quantitative Model Investigation on Four Materials
Caries Res 26 (1992) 351-357
- Dijkman GEHM, De Vries J, Lodding A, Arends J
Long-term fluoride release of visible light-activated composites in vitro: a correlation with in situ demineralisation data
Caries Res 27 (1993) 117-123
- Donly K, Gomez C
In vitro demineralization-remineralization of enamel caries at restoration margins utilizing fluoride-releasing composite resin
Quintessence Int 25 (1994) 355-358
- Fallo GJ, Wakefield CW
Effects of uncontrolled outdoor storage on the polymerization, manipulation, and appearance of visible light-cured composite resin and resin-modified glass ionomer materials
Military Medicine 161 (1996) 290-293
- Feilzer AJ, De Gee AJ, Davidson CL
Relaxation of Polymerization Contraction Shear Stress by Hygroscopic Expansion
J Dent Res 69 (1990) 36-39
- Feilzer AJ, De Gee AJ, Davidson CL
Curing contraction of composites and glass-ionomer cements
J Prosthet Dent 59 (1988) 297-300
- Ferracane JL
Current Trends in Dental Composites
Crit Rev Oral Biol Med 6 (1995) 302-318
- Ferracane JL, Condon JR
Post-cure heat treatments for composites: properties and fractography
Dent Mater 8 (1992) 290-295
- Ferracane JL, Mitchem JC, Condon JR, Todd R
Wear and marginal breakdown of composites with various degrees of cure
J Dent Res 76 (1997) 1508-1516
- Francci C, Perdigo J, Cardoso PEC, Meira JBC, Nunes MF

The effect of composite resin, adhesive system, and low-viscosity liner on microleakage

J Dent Res 78 (1999) 388

Freilich MA, Goldberg AJ, Gilpatrick RO, Simonsen RJ

Three-year occlusal wear of posterior composite restorations

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Dtsch Zahnärztl Z 49 (1994) 729-732

Fujimitsu T, Kato Hltoh K, Wakumoto S

The adaptation of composite to the dentin cavity wall

Dent Mater J 8 (1989) 141-146

Fukushima M, Setcos JC, Phillips RW

Marginal fracture of posterior composite resins

J Am Dent Assoc 117 (1988) 577-583

Fuqua S, Burgess JO, Chan DCN, Nunez A

Recurrent caries in OF restorations in a rodent model

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Dtsch Zahnärztl Z 46 (1991) 270-273

Herrgott AM, Ziemięcki TL, Dennison JB

An evaluation of different composite resin systems finished with various abrasives

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Quintessence Int 17 (1986) 25-27

Kanca III J

Visible light-activated composite resins for posterior use - a comparison of surface hardness and uniformity of cure. Update

Quintessence Int 16 (1985) 687-690

Kawai K, Leinfelder K

In vitro evaluation of OCA wear resistance of posterior composites

Dent Mater 11 (1995) 246-251

Knibbs PJ, Smart ER

The clinical performance of a posterior composite resin restorative material, Heliomolar R.O.: 3-year report

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Dtsch Zahnärztl Z 46 (1991) 400-406

Lambert D

Closing diastemas, replacing small amalgams using Heliomolar enforced microfill restoratives

Dental Products Report 0 (2000) 92-93

Leinfelder K

Posterior composite resins: The materials and their clinical performance

J Am Dent Assoc 126 (1995) 663-676

Leinfelder K

Using composite resin as a posterior restorative material

J Am Dent Assoc 122 (1991) 65-70

Leinfelder K, Beaudreau RW, Mazer R

An in vitro device for predicting clinical wear

Quintessence Int 20 (1989) 755-761

Liebenberg WH

Direct access to equivocal approximal carious lesions

Quintessence Int 27 (1996) 607-617

Lundin SA, Andersson B, Koch G, Rasmusson CG

Class II composite resin restorations: a three-year clinical study of six different posterior composites

Swed Dent J 14 (1990) 105-114

Mair LH

The silver sorption layer in dental composites: three year results

Dent Mater J 15 (1999) 408-412

Mair LH

An investigation into the permeability of composite materials using silver nitrate

Dent Mater 5 (1987) 109-114

Mair LH

Subsurface compression fatigue in seven dental composites

Dent Mater 10 (1994) 111-115

Mair LH, Vowles RW

The effect of thermal cycling on the fracture toughness of seven composite restorative materials

Dent Mater 5 (1989) 23-26

Margolis HC, Moreno EC, Murphy BJ

Effect of low levels of fluoride in solution on enamel demineralization in vitro

J Dent Res 65 (1986) 23-29

Mazer R, Leinfelder K

Evaluating a microfill posterior composite resin - a five year study

J Am Dent Assoc 123 (1992) 33-38

Mazer R, Leinfelder K, Russell CM

Degradation of microfilled posterior composite

Dent Mater 8 (1992) 185-189

Milosevic A

The influence of surface finish and in-vitro pellicle on contact-angle measurement and surface morphology of three commercially available composite restoratives

J Oral Rehab 19 (1992) 85-97

Momoi Y, McCabe JF

Hygroscopic expansion of resin based composites during 6 months of water storage

Br Dent J 176 (1994) 91-96

Mount GJ

The tensile strength of the union between various glass ionomer cements and various composite resins

Aust Dent J 34 (1989) 136-146

Munksgaard EC, Nolte J, Kristensen K

Adherence of chewing gum to dental restorative materials

Am J Dent 8 (1995) 137-139

NaBadalung DP, Nicholls JI, Brudvik JS

Frictional resistance of removable partial dentures with retrofitted resin composite guide planes

Int J Prosthodont 10 (1997) 116-122

Oden A, Ruyter IE, Oysaed H

Creep and recovery of composites for use in posterior teeth during static and dynamic compression

Dent Mater 4 (1988) 147-150

Palenik CJ, Setcos JC

Antimicrobial abilities of various dentine bonding agents and restorative materials

J Dent 24 (1996) 289-295

Papagiannoulis L, Tzoutzas J, Eliades G

Effect of topical fluoride agents on the morphologic characteristics and composition of resin composite restorative materials

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Comparison of two- and three-body wear of glass-ionomers and composites

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Dtsch Zahnärztl Z 53 (1998) 61-64

Pilo R, Brosh T, Lugassy H, Baharav H, Helft M

The effect of irradiation time on the shear strength of composites

Dent Mater 10 (1994) 338-342

Powers JM, Bakus ER, Goldberg AJ

In vitro color changes of posterior composites

Dent Mater 4 (1988) 151-154

Powers JM, McDowell GC, Lang BR

In vivo wear. Part II: Wear and abrasion of composite restorative materials

J Prosthodont 60 (1988) 242-249

Prati C

Early marginal microleakage in class II resin composite restorations

Dent Mater 5 (1989) 392-398

Rasmusson CG, Lundin SA

Class II restorations in six different posterior composite resins: Five-year results.

Swed Dent J 19 (1995) 173-182

Rees JS, Jacobsen PH

The polymerization shrinkage of composite resins

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Dtsch Zahnärztl Z 47 (1992) 176-178

Retief DH, Mandras RS, Russell CM, Denys FR

Evaluation of the Syntac bonding system

Am J Dent 6 (1993) 17-21

Retief H, Mandras R, Russell CM

Shear bond strength required to prevent microleakage at the dentin/restoration interface

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Quintessenz Zahntech 19 (1993) 1467-1473

Sakaguchi RL, Peters MC, Nelson SR, Douglas WH, Poort HW

Effects of polymerization contraction in composite restorations

J Dent 20 (1992) 178-182

Saunders WP

The marginal microleakage of resin-retained bridges in association with existing composite and amalgam restorations

Dent Mater 6 (1990) 20-23

Scheufele P

Suitability of different antagonist for pre-clinical in-vitro wear tests

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Phillip J 12 (1995) 93-95

Shinkai K, Suzuki S, Katoh Y

Effect of an adhesive bonding system on wear resistance of resin composite restorations

Quintessence Int 28 (1997) 687-693

Shinkai K, Suzuki S, Leinfelder K, Katoh Y

How heat treatment and thermal cycling affect wear of composite resin inlays

J Am Dent Assoc 125 (1994) 1467-1472

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A 3-year follow-up study of preformed beta-quartz glass-ceramic insert restorations

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ZWR 100 (1991) 648-653

Steinmetz MJ, Pruhs RJ, Brooks JC, Dhuru VB, Post AC

Rechargeability of fluoride releasing pit and fissure sealants and restorative composites

Am J Dent 10 (1997) 36-40

Suliman AA, Boyer DB, Lakes RS

Cups movement in premolars resulting from composite polymerization shrinkage

Dent Mater 9 (1993) 6-10

Suliman AA, Boyer DB, Lakes RS

Interferometric measurements of cusp deformation of teeth restored with composites

J Dent Res 72 (1993) 1532-1536

Tani Y, Goto H, Ida K

Wear of posterior composite resins

Dent Mater J 6 (1987) 165-174

Taylor DF, Bayne SC, Leinfelder KF, Davis S, Koch GG

Pooling of long term clinical data for posterior composites

Am J Dent 7 (1994) 167-174

Tjan AHL, Chan CA

The polishability of posterior composites

J Prosthet Dent 61 (1989) 138-146

Turner CW, Meiers JC

Repair of an aged, contaminated indirect composite resin with a direct, visible-light-cured composite resin

Oper Dent 18 (1993) 187-194

Versluis A, Douglas WH, Sakaguchi RL

Thermal expansion coefficient of dental composites measured with strain gauges

Dent Mater 12 (1996) 290-294

Von Beetzen M, Li J, Nicander I, Sundstrom F

Microhardness and porosity of Class 2 light-cured composite restorations cured with a transparent cone attached to the light-curing wand

Oper Dent 18 (1993) 103-109

Wassel RW, McCabe JF, Walls AWG

Wear rates of regular and tempered composites

J Dent 25 (1997) 49-52

Wassel RW, McCabe JF, Walls AWG

A two-body friction wear test

J Dent Res 73 (1994) 1546-1553

Wassel RW, McCabe JF, Walls AWG

Subsurface deformation associated with hardness measurements of composites

Dent Mater 8 (1992) 218-223

Wassel RW, McCabe JF, Walls AW

Wear characteristics in a two-body wear test

Dent Mater 10 (1994) 269-274

Watts DC, Cash AJ

Analysis of optical transmission by 400-500 nm visible light into aesthetic dental biomaterials

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A classification of dental composites according to their morphological and mechanical characteristics

Dent Mater 8 (1992) 310-319

Willems G, Lambrechts P, Braem MJA, Vuylsteke-Wauters M, Vanherle G

The surface roughness of enamel-to enamel areas compared with the intrinsic roughness of dental composites

J Dent Res 70 (1991) 1299-1305

Xu X, Burgess Jo, Al-Owyyed, Weathersby

Artificial caries with fluoride-releasing and non-fluoride-releasing restorative materials

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Lowe E
Heavy duty tips and tricks for using Heliomolar HB
Dental Town Magazine 0 (2001) 14-15

Jackson RD
Replacing failed amalgam fillings using Heliomolar HB direct adhesive resin
Dental Products Report 3 (2001) 70

Neme AL, Pink FE, Aksu MN, Maxson BB, Kashani N
In-vitro microleakage in class II packable resin composite restorations
J Dent Res 80 (2001) 105

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Ivoclar Vivadent AG