Heliomolar® Flow

| | |
|----------------------|-------------------|
| Instructions for Use | Bruksanvisning |
| Gebrauchsinformation | Bruksanvisning |
| Mode d'emploi | Käyttöohjeet |
| Istruzioni d'uso | Bruksanvisning |
| Instrucciones de uso | Productinformatie |
| Instruções de Uso | Οδηγίες Χρήσεως |

- Light-curing, resin-based dental restorative material
- Lichthärtendes zahnärztliches Füllungscomposite
- Matériau de restauration photopolymérisable
- Materiale dentale da otturazione fotoindurente
- Material de restauración dental fotopolimerizable
- Material fotopolimerizável de restauração dental
- Ljushärdande dentalt fyllnadsmaterial
- Lyshärdende dentalt plastfyldningsmateriale
- Valokovetteinen yhdistelmämuovi-täyttemateriaali
- Lichtuithardend tandheelkundig vulmateriaal
- Φωτοπολυμεριζόμενο οδοντιατρικό υλικό αποκαταστάσεων

For dental use only!

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein

Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

Complies with / entspricht:
ISO 4049 / EN ISO 4049

CE 0123

ivoclar
vivadent®
clinical

English

Instructions for Use

Description

Heliomolar Flow is a radiopaque microfilled composite that cures with light in the wavelength range of 400–500 nm (blue light).

Shades

Heliomolar Flow is available in the following 7 shades:

| Chromascope | A–D shades |
|-------------|------------|
| 105 | – |
| 110 | A1 |
| 140 | A2 |
| 210 | A3 |
| 410 | D3 |
| 110T | A1T |
| 420T | – |

Composition

The monomer matrix consists of Bis-GMA, urethane dimethacrylate and triethylene glycol methacrylate (40.5 wt%). The fillers are highly dispersed silicon dioxide, ytterbiumtrifluoride and copolymer (59 wt%). Additional contents are catalysts, stabilizers and pigments (0.5 wt%).

The total content of inorganic fillers is 30 vol% or 51 wt%. The particle size ranges between 0.04 and 0.2 µm.

Indication

- Class V restorations (cervical caries, root erosion, wedge-shaped defects)
- Anterior restorations (Class III)
- Restorative therapy of mini-cavities of all types
- Preventive resin restorations in molars and premolars
- Repair of composite/ceramic veneers
- Blocking out of undercuts
- Adhesive cementation of ceramic and indirect composite restorations

- Splinting of mobile teeth
- As the first layer of Class I or II restorations
- As the last layer applied on a viscous composite material

Contraindication

- Heliomolar Flow restorations are contraindicated
- if a dry operating field cannot be established, or if the stipulated application technique cannot be used;
 - if the patient is known to be allergic to any of the ingredients in Heliomolar Flow.

Side effects

In rare cases, components of Heliomolar Flow may cause a sensitizing reaction in patients with a hypersensitivity to any of its ingredients. In these cases, Heliomolar Flow should no longer be used.

Interaction

- Eugenol may hamper the polymerization of Heliomolar Flow. Therefore, cements containing zinc oxide eugenol must not be used in combination with Heliomolar Flow.
- Cationic mouthwashes, plaque disclosing agents and chlorhexidine may cause discolouration.

Application procedure

1. Shade selection

Clean the teeth prior to shade selection. The teeth should be slightly moist for shade selection.

2. Isolation

Provide adequate isolation, best with a rubber dam (e.g. OptraDam® Plus).

3. Cavity preparation

Prepare the cavity according to the principles of the adhesive technique, i.e. preserving the tooth structure. Do not prepare sharp internal edges or additional undercuts in areas free of caries. The dimensions of the cavity depend on the extent of the caries and/or the old restoration. Where indicated, the enamel margins can be bevelled. Caries-free cervical defects do not require preparation.

Merely clean with pumice and a suitable cleaning paste, using a rubber cup or rotary brush. Next, remove all residue with a water spray.

4. Pulp protection

If an enamel-dentin bonding agent is used, a base/liner is generally not required. In very deep cavities, cover the areas close to the pulp with a calcium hydroxide material.

5. Matrix (e.g. OptraMatrix®)/interdental wedge

6. Conditioning / Application of the bonding agent

Condition and apply the bonding agent according to the Instructions for Use of the product in use. We recommend using Syntac® or ExcITE® F (with phosphoric acid etching) or AdheSE® or AdheSE® One F self etching adhesive.

7. Application of Heliomolar Flow

Apply Heliomolar Flow in layers of max. 2 mm and adapt with a suitable instrument. Light-cure each layer with a polymerization light (e.g. bluephase®, bluephase® 20i) for 40 seconds. Hold the light emission window of the curing light as close to the restorative as possible.

Important: Heliomolar Flow should have room temperature for application. This ensures a smooth consistency. Do not expose Heliomolar Flow to intensive light during application, since exposure to light reduces the working time.

8. Finishing/Polishing

Remove excess material with suitable finishers (e.g. Astropol F) or fine diamonds after polymerization. Remove proximal excess with diamond or tungsten carbide finishers or finishing strips. Check occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restorations. Use silicone polishers (e.g. Astropol® P/Astropol HP, Astrobrush®), polishing discs, and polishing strips to polish the restorations to a durable high gloss.

Additional notes

Heliomolar Flow may also be applied on polymerized material. In order to attain appropriate adhesion between the interfaces, the polymerized material should be roughened and cleaned prior to the application of fresh Heliomolar Flow.

Warning

- Avoid contact of Heliomolar Flow with skin/mucous membrane and eyes. Uncured Heliomolar Flow may cause slight irritation and may result in sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage

- Close syringes immediately after use. Exposure to light causes premature polymerization.
- Do not use Heliomolar Flow after the indicated date of expiration.
- Storage temperature: 2–28 °C / 36–82 °F.
- Shelf life: see date of expiration on label/package.

Keep out of the reach of children.

For use in dentistry only.

Date information prepared

12/2009, Rev. 1

Manufacturer

Ivoclar Vivadent AG
FL-9494 Schaan / Liechtenstein

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent UK Limited
Ivoclar Vivadent, Inc.