Dual Cement

English

Instructions for Use	Bruksanvisning
Gebrauchsinformation	Bruksanvisning
Mode d'emploi	Käyttöohjeet
Istruzioni d'uso	Bruksanvisning
Instrucciones de uso	Productinformatie
Instruções de Uso	Οδηγίες Χρήσεως

- Dual-curing resin-based dental luting material
- Dualhärtendes Befestigungscomposite
- Composite de collage photo- et chemo-polymérisable "dual"
- Cemento composito ad indurimento duale (foto-/autoindurimento)
- Cemento de fijación dual de curado fotopolimerizable
- Cimento de polimerização dual
- Dual- och ljushärdande kompositcement
- Dualhærdende plastretentionsciment
- Kaksoikovetteinen yhdistelmämuovisementti
- Dualherdende sementeringskomposit
- Duaaluithardend bevestigingscomposit
- Διπλά πολυμεριζόμενη ρητινώδης κοιλία

Instructions for Use

Description

Dual Cement is a microfilled, light- and self-curing composite cement for the adhesive cementation of composite, ceramic, or Ceromer/FRC restorations. Dual Cement is distinguished for its excellent physical properties, high abrasion resistance, continuous fluoride release, as well as good radiopacity. Dual Cement is available in one highly translucent shade that optimally adapts itself to the restorations and the surrounding tooth structure.

Composition

The monomer matrix consists of urethane dimethacrylate and decandiol dimethacrylate. The inorganic fillers are silicon dioxide and ytterbium trifluoride. Additional contents are catalysts, stabilizers, and pigments. The total contents of inorganic fillers are 61 wt. % or 38 vol. %.

The particle size is between 0.04 and 0.3 µm, with a mean grain size of 0.2 µm.

Indication

Dual Cement has been designed for the adhesive cementation of lab-fabricated restorations made of composite, ceramics, or Ceromer/FRC.

Contraindication

- The use of Dual Cement is contraindicated
- if a dry operating field cannot be established or if the stipulated working technique cannot be applied, or
 - if patients are known to be allergic to any of the ingredients of Dual Cement.

For dental use only!

Caution: U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein

Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

CE 0123

ivoclar
vivadent®
clinical

558544/0503/3/WE3/G

Side effects

Systemic side effects are not known to date. In individual cases, allergic reactions were reported. Irritation of the pulp may occur without the use of a suitable pulp/dentin protection or without base/liner in deep or medium-sized cavities.

Interactions

Phenolic substances (e.g. eugenol) inhibit the polymerization. Therefore, materials containing such substances should not be used in conjunction with Dual Cement.

Application

Step-by-step

For detailed information, please also refer to the Instructions for Use of the products mentioned.

1. Removing the temporary

Remove any residue of the temporary cement from the cavity or the preparation using a polishing brush and a cleaning paste that is free of oil and fluoride. Subsequently, rinse with water spray.

Note: Use only temporary cements that are free of eugenol (e.g. Provilink®), since eugenol may inhibit the polymerization of Dual Cement. If Fermit®/Fermit®N is used for inlays, a temporary cement is usually not required.

2. Trying-in the restoration

Check exact fit and proximal contacts of the restoration. The occlusion should only be checked after cementation (risk of fracture).

3. Isolating

The use of a rubber dam is recommended.

4. Conditioning of the restoration

According to the instructions of the manufacturer.

5. Conditioning of the cavity/preparation

- 5.1 Etch the enamel margins, e.g. with Total Etch , for 30–60 seconds, rinse, and dry with water- and oil-free air. If necessary, etch dentin according to the instructions of the manufacturer of the dentin adhesive used.
- 5.2 Apply a suitable **dentin bonding** agent (e.g. Syntac®, Excite® – observe the corresponding instructions for Use).

6. **Mix Dual Cement in a 1:1 ratio.** Dual Cement must not be exposed to direct sunlight or other intensive light sources during this process. The working time of Dual Cement is 4–5 minutes (at 23 °C / 73 °F).

- 6.1 Cool storage (2–8 °C / 36–46 °F) prolongs the working time of Dual Cement.
 - 6.2 Excessive mixing may result in air being trapped in the material, which may hamper the polymerization of the material and subsequently weaken the mechanical properties.
7. Apply Dual Cement on the restoration and in the cavity (avoid entrapment of air).

8. **Place** the restoration, e.g. inlay.

9. **Remove excess** cement. The inlay can be 'secured' by removing occlusal excess and quickly curing only the marginal areas by means of a polymerization light. After that, remove the remaining excess.

10. Subsequently, **cure** Dual Cement with light (Astralis®) as thoroughly as possible for 40–60 seconds.

Note on oxygen inhibition

As all composites, Dual Cement is affected by oxygen inhibition, i.e. the surface layer (approx. 50 µm), which comes in contact with oxygen during polymerization, does not cure. This problem can be solved in two ways:

- 1) Leave minimum Dual Cement excess. In this way, only the excess material is affected by oxygen inhibition. This excess is easily removed during finishing of the margins.
- 2) Cover the restoration margins with Liquid Strip after having removed the excess material, but before polymerization. This will prevent oxygen inhibition.

11. Finishing and polishing

Remove excess material and finish the restoration with finishing diamonds, flexible disks, and suitable silicone rubber polishers (Politip P, Politip F). Check the occlusion.

Warning

Contact of Dual Cement with skin/mucous membrane or eyes should be avoided. Uncured Dual Cement may have a slight irritant effect and may lead to a sensitization against methacrylates.

Storage and shelf life

- Do not use Dual Cement after the indicated date of expiration.
- Storage in the refrigerator is required for the catalyst paste (2–8 °C / 36–46 °F).
- **Note:** Close syringes immediately after use. Exposure to light causes premature polymerization.
- Shelf life: see date of expiration

**Store out of the reach of children.
For use in dentistry only.**

Date information prepared:
05/2003

Manufacturer:
Ivoclar Vivadent AG
FL-9494 Schaan/Liechtenstein

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent Ltd
Ivoclar Vivadent S.A.
Ivoclar Vivadent AB
Ivoclar Vivadent, Inc.