

# Programat® P700 & P700/G2

## Software Revision Log

Version	Changes from previous version	Status
<b>V8.1</b> 25.08.2017	<b>Software optimized</b> <ul style="list-style-type: none"> <li>• New program structure</li> <li>• New programs (e.g. IPS e.max Ceram)</li> </ul>	<b>Currently delivered version</b>
<b>V8.0</b> 01.04.2016	<b>Software optimized</b> <ul style="list-style-type: none"> <li>• New programs for IPS Style and IPS Ivocolor</li> </ul>	<b>No longer available!</b>
<b>V7.1</b> 01.10.2014	<b>Software optimized</b> <ul style="list-style-type: none"> <li>• New programs for IPS InLine powder opaquer</li> <li>• Various software improvements</li> </ul>	<b>No longer available!</b>
<b>V7.0</b> 05.05.2014	<b>Software optimized</b> <ul style="list-style-type: none"> <li>• New service interval indicator</li> <li>• Various software improvements</li> <li>• New, special speed crystallization program for IPS e.max CAD (P85): The new program replaces the existing IPS e.max CAD Speed Crystallization. The program may only be used in conjunction with the IPS e.max CAD Speed Crystallization Tray (article code 647633AN). Please observe the processing information for the IPS e.max CAD Speed Crystallization.</li> </ul>	<b>No longer available!</b>
<b>V6.0</b> 04.04.2011	<b>Software optimized</b> <ul style="list-style-type: none"> <li>• New Ivoclar Vivadent Programs for IPS InLine One (Programs 121 – 126)</li> <li>• New Ivoclar Vivadent Programs for IPS e.max CAD-on (Programs 83 – 84)</li> </ul>	<b>No longer available!</b>
<b>V5.0</b> 01.05.2010	<b>Software optimized and new functions</b> <ul style="list-style-type: none"> <li>• Ivoclar Vivadent Firing parameters were changed</li> <li>• Coloured firing curve display</li> </ul>	<b>No longer available!</b>
<b>V4.1</b> 01.12.2009	<b>Software optimized and new functions</b> <ul style="list-style-type: none"> <li>• Improved voice output and optimized remaining time indicator</li> </ul>	<b>No longer available!</b>

# Programat® P700 & P700/G2

<b>V4.0</b> 25.08.2009	<b>New Software P700/G2: Software optimized and new features</b> <ul style="list-style-type: none"> <li>• Power Saving Technology</li> <li>• CSP: Cooling Shock Protection</li> <li>• Acoustic voice output in 12 languages</li> <li>• New icons and optimized user interface</li> <li>• and many more!</li> </ul>	<b>No longer available!</b>
<b>V3.0</b> 15.11.2007	<b>Software optimized and new functions</b> <ul style="list-style-type: none"> <li>• Updated Ivoclar standard firing programs P1-P120</li> <li>• Newly added languages: Chinese simplified, Chinese traditional, Swedish, Russian, Polish, Turkish</li> </ul>	<b>No longer available!</b>
<b>V1.04</b> 03.09.2007	<b>Software optimized and new functions</b> <ul style="list-style-type: none"> <li>• Extended USB memory stick recognition</li> </ul>	<b>No longer available!</b>
<b>V1.03</b> 04.06.2007	<b>Software optimized</b>	<b>No longer available!</b>
<b>V1.02</b> 01.06.2007	<b>New Function</b> <ul style="list-style-type: none"> <li>• Improved opening and closing movement with smooth start function</li> </ul>	<b>No longer available!</b>
<b>V1.01</b> 22.05.2007	<b>New Function</b> <ul style="list-style-type: none"> <li>• Protocols of firing programs include SW information</li> </ul>	<b>No longer available!</b>
<b>V1.00</b> 14.05.2007	<b>Basic Version</b>	<b>No longer available!</b>

Version 25.08.2017