

# Programat® EP 5000 & EP 5000/G2

## Software Revision Log

Version	Changes from previous version	Status
<b>V8.1</b> 25.08.2017	<b>Software optimized</b> <ul style="list-style-type: none"> <li>New program structure</li> <li>New programs (e.g. IPS e.max Ceram)</li> </ul>	<b>Currently delivered version</b>
<b>V8.0</b> 01.04.2016	<b>Software optimized</b> <ul style="list-style-type: none"> <li>New programs for IPS Style, IPS Ivocolor and IPS e.max Press MT</li> <li>Various software optimizations</li> </ul>	<b>No longer available!</b>
<b>V7.1</b> 30.01.2015	<b>Software optimized</b> <ul style="list-style-type: none"> <li>Various software improvements</li> </ul>	<b>No longer available!</b>
<b>V7.0</b> 01.10.2014	<b>Software optimized</b> <ul style="list-style-type: none"> <li>New programs for IPS InLine powder opaquer</li> <li>New programs for IPS e.max Press Multi</li> <li>Various software improvements</li> </ul>	<b>No longer available!</b>
<b>V6.9</b> 05.05.2014	<b>Software optimized</b> <ul style="list-style-type: none"> <li>New service interval indicator</li> <li>Software improvements for service technician</li> <li>Program structure prepared for speed crystallization program for IPS e.max CAD</li> <li>Various software improvements</li> </ul>	<b>No longer available!</b>
<b>V6.1</b> 06.07.2011	<b>Software optimized</b> <ul style="list-style-type: none"> <li>Optimized program selection for individual pressing programs</li> </ul>	<b>No longer available!</b>
<b>V6.0</b> 04.04.2011	<b>New Software EP 5000/G2: Software optimized and new features</b>  Compared to the predecessor model Programat EP 5000 following functions got integrated: <ul style="list-style-type: none"> <li>Power Saving Technology</li> <li>TSP: Thermo Shock Protection</li> <li>CSP: Cooling Shock Protection</li> <li>New program-structure (130 Ivoclar Vivadent programs, 500 individual programs)</li> <li>New Ivoclar Vivadent Program for IPS e.max Press Multi ingots</li> <li>New Ivoclar Vivadent Programs for IPS InLine One (Programs 121 – 126)</li> <li>New Ivoclar Vivadent Programs for IPS e.max CAD-on (Programs 83 – 84)</li> </ul>	<b>No longer available!</b>

# Programat® EP 5000 & EP 5000/G2

<b>V5.0</b> 01.05.2010	<b>Software optimized</b> <ul style="list-style-type: none"> <li>Firing parameters were changed</li> </ul>	<b>No longer available!</b>
<b>V3.3</b> 19.05.2009	<b>Software optimized</b> <ul style="list-style-type: none"> <li>Pressing parameters for IPS e.max Press HT for 100 g and 200 g IPS Investment ring system</li> <li>IPF (Intelligent Press Function) parameters for IPS e.max ZirPress (for IPS Investment ring system 300 g) and for IPS Inline PoM (for IPS investment ring system 300 g)</li> </ul>	<b>No longer available!</b>
<b>V3.2</b> 12.12.2008	<b>Software optimized</b> <ul style="list-style-type: none"> <li>Optimized setting of the Crack Detection System (CDS) for 300 g investment ring in program "IPS InLine PoM – 300 g"</li> </ul>	<b>No longer available!</b>
<b>V3.1</b> 14.10.2008	<b>New functions</b> <ul style="list-style-type: none"> <li>Selection option for 300 g investment ring is integrated</li> <li>Parameter for 300 g investment ring is integrated</li> </ul>	<b>No longer available!</b>
<b>V3.0</b> 15.11.2007	<b>Software optimized and new functions</b> <ul style="list-style-type: none"> <li>New STD press program "InLine PoM"</li> <li>Aesthetic press programs new with stop speed of 180 rather than 250</li> <li>Updated Ivoclar standard firing parameters P1-P120</li> <li>Newly added languages: Chinese simplified, Chinese traditional, Swedish, Russian, Polish, Turkish</li> </ul>	<b>No longer available!</b>
<b>V1.04</b> 28.08.2007	<b>Software optimized and new functions</b> <ul style="list-style-type: none"> <li>Extended USB memory stick recognition</li> </ul>	<b>No longer available!</b>
<b>V1.03</b> 04.06.2007	<b>Software optimized</b>	<b>No longer available!</b>
<b>V1.02</b> 01.06.2007	<b>New function</b> <ul style="list-style-type: none"> <li>Improved opening and closing movement with smooth start function</li> </ul>	<b>No longer available!</b>
<b>V1.01</b> 22.05.2007	<b>New function</b> <ul style="list-style-type: none"> <li>Protocols of firing programs include SW information</li> </ul>	<b>No longer available!</b>
<b>V1.00</b> 14.05.2007	<b>Basic version</b>	<b>No longer available!</b>

Version 25.08.2017