

Programat® EP 3000 & EP 3000/G2

Software Revision Log

Version	Changes from previous version	Status
V8.1 25.08.2017	Software optimized <ul style="list-style-type: none"> • New program structure • New programs (e.g. IPS e.max Ceram) 	Currently delivered version
V8.0 01.04.2016	Software optimized <ul style="list-style-type: none"> • New programs for IPS Style, IPS Ivocolor and IPS e.max Press MT 	No longer available!
V7.1 30.01.2015	Software optimized <ul style="list-style-type: none"> • Various software improvements 	No longer available!
V7.0 01.10.2014	Software optimized <ul style="list-style-type: none"> • New programs for IPS InLine powder opaquer • New programs for IPS e.max Press Multi • Various software improvements 	No longer available!
V6.9 05.05.2014	Software optimized <ul style="list-style-type: none"> • New service interval indicator • Software improvements for service technician • Program structure prepared for speed crystallization program for IPS e.max CAD • Various software improvements 	No longer available!
V6.3 07.02.2012	Software optimized <ul style="list-style-type: none"> • Optimized temperature calibration in conjunction with the ATK2 temperature checking set 	No longer available!
V6.1 06.07.2011	Software optimized <ul style="list-style-type: none"> • Optimized program selection for individual pressing programs 	No longer available!

Programat® EP 3000 & EP 3000/G2

V6.0 04.04.2011	<p>New Software EP 3000/G2: Software optimized and new features</p> <p>Compared to the predecessor model Programat EP 3000 following functions got integrated:</p> <ul style="list-style-type: none"> • Power Saving Technology • TSP: Thermo Shock Protection • New program-structure (130 Ivoclar Vivadent programs, 300 individual programs) • New Ivoclar Vivadent Program for IPS e.max Press Multi ingots • New Ivoclar Vivadent Programs for IPS InLine One (Programs 121 – 126) • New Ivoclar Vivadent Programs for IPS e.max CAD-on (Programs 83 – 84) 	No longer available!
V5.0 01.05.2010	<p>Software optimized</p> <ul style="list-style-type: none"> • Ivoclar Vivadent firing parameters were changed 	No longer available!
V1.3 19.05.2009	<p>Software optimized</p> <ul style="list-style-type: none"> • Pressing parameters for IPS e.max Press HT for 100g and 200g IPS Investment ring system 	No longer available!
V1.2 12.12.2008	<p>Software optimized</p> <ul style="list-style-type: none"> • Optimized setting of the Crack Detection System (CDS) for 300g investment ring in program "IPS InLine PoM - 300g" 	No longer available!
V1.1 14.10.2008	<p>New functions</p> <ul style="list-style-type: none"> • Selection option for 300g investment ring is integrated • Parameter for 300g investment ring is integrated 	No longer available!
V1.0 01.07.2008	<p>Basic version</p>	No longer available!

Version 25.08.2017