

IvoBase® Injector

Software Revision Log

Version	Changes from previous version	Status
V2.0 18.02.2014	Software optimized • Error 1936: Optimized “capsule filling level detection” and new error description	Currently delivered version
V1.1 06.03.2013	Software optimized • Error 1914: Optimized error detection and display of error messages • Error 1935: Optimized error detection and display of error messages • Transport protection: If an IvoBase Injector is shipped or transported, a transport protection function must now be activated in the settings menu	No longer available!
V1.0 07.12.2009	Basic Version of the IvoBase Injector	No longer available!

Stand: 18.02.2014