


**1 Identification**

- **Product identifier**
- **Trade name:** SR Adoro Incisal, Dentin, Deep Dentin, Neck, Transpa Incisal, Special Incisal, Inter Incisal, Occlusal Dentin, Mamelon, Effect, Opal Effect, Incisal Edge, Transpa, Transpa Clear, Add-On
- **Relevant identified uses of the substance or mixture and uses advised against**  
No further relevant information available.
- **Application of the substance / the mixture** Crowns- and bridges material
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
Ivoclar Vivadent Inc.  
175 Pineview Drive, Amherst, N.Y. 14228  
USA  
Tel. +1 800 533 6825  
Fax +1 716 691 2285
- 
- Ivoclar Vivadent Inc.  
1-6600 Dixie Road  
Mississauga, Ontario  
L5T 2Y2  
Canada  
Phone: +1 905 670 8499  
Fax: +1 905 670 3102
- **Information department:** Quality Assurance / Regulatory Affairs
- **Emergency telephone number:**  
24 Hour Emergency Assistance:  
Emergency-Call USA - Infotrac: 1-800-535-5053  
Emergency-Call Canada - Canutec: 1-613-996-6666
- 
- General SDS Assistance:  
US: 1-800-533-6825  
Canada: 1-800-263-8182

**2 Hazard(s) identification**

- **Classification of the substance or mixture**  
Eye Irrit. 2A H319 Causes serious eye irritation.  
Skin Sens. 1 H317 May cause an allergic skin reaction.

---

- **Label elements**
- **GHS label elements**  
The product is classified and labeled according to the Globally Harmonized System (GHS).
- **Hazard pictograms**
- 
- GHS07
- **Signal word** Warning
- **Hazard-determining components of labeling:**  
urethane dimethacrylate  
1,10-decandiol dimethacrylate

## Safety Data Sheet

acc. to OSHA HCS

Printing date 04/29/2015

Version number 4

Reviewed on 04/16/2015

**Trade name: SR Adoro Incisal, Dentin, Deep Dentin, Neck, Transpa Incisal, Special Incisal, Inter Incisal, Occlusal Dentin, Mamelon, Effect, Opal Effect, Incisal Edge, Transpa, Transpa Clear, Add-On**

(Contd. of page 1)

- **Hazard statements**

Causes serious eye irritation.  
May cause an allergic skin reaction.

- **Precautionary statements**

Wear protective gloves/protective clothing/eye protection/face protection.  
Wash contaminated clothing before reuse.  
If skin irritation or rash occurs: Get medical advice/attention.  
If eye irritation persists: Get medical advice/attention.

- **Classification system:**

- **NFPA ratings (scale 0 - 4)**



Health = 1  
Fire = 0  
Reactivity = 0

- **HMIS-ratings (scale 0 - 4)**



Health = 1  
Fire = 0  
Reactivity = 0

- **Other hazards**

- **Results of PBT and vPvB assessment**

· **PBT:** Not applicable.  
· **vPvB:** Not applicable.

### \* 3 Composition/information on ingredients

- **Chemical characterization: Mixtures**

· **Description:** Mixture of the substances listed below with nonhazardous additions.

- **Dangerous components:**

72869-86-4	urethane dimethacrylate	10-<20%
6701-13-9	1,10-decandiol dimethacrylate	2.5-<10%
43048-08-4	tricyclodocane dimethanol dimethacrylate	1-<2.5%

### 4 First-aid measures

- **Description of first aid measures**

- **After inhalation:**

Supply fresh air and to be sure call for a doctor.  
In case of unconsciousness place patient stably in side position for transportation.

- **After skin contact:** Rinse with water.

- **After eye contact:**

Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.

- **After swallowing:**

Rinse out mouth and then drink plenty of water.  
If symptoms persist consult doctor.

- **Information for doctor:**

· **Most important symptoms and effects, both acute and delayed** No further relevant information available.

(Contd. on page 3)

## Safety Data Sheet

acc. to OSHA HCS

Printing date 04/29/2015

Version number 4

Reviewed on 04/16/2015

**Trade name: SR Adoro Incisal, Dentin, Deep Dentin, Neck, Transpa Incisal, Special Incisal, Inter Incisal, Occlusal Dentin, Mamelon, Effect, Opal Effect, Incisal Edge, Transpa, Transpa Clear, Add-On**

- **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

(Contd. of page 2)

### 5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**  
CO<sub>2</sub>, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:** No special measures required.

### 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:** Pick up mechanically.
- **Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

### 7 Handling and storage

- **Handling:**
- **Precautions for safe handling**  
Only adequately trained personnel should handle this product.  
For use in dentistry only.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** Store only in the original receptacle.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:**  
Protect from exposure to the light.  
Protect from heat and direct sunlight.  
Keep receptacle tightly sealed.
- **Specific end use(s)** No further relevant information available.

### 8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**  
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.

(Contd. on page 4)

US

## Safety Data Sheet

acc. to OSHA HCS

Printing date 04/29/2015

Version number 4

Reviewed on 04/16/2015

**Trade name: SR Adoro Incisal, Dentin, Deep Dentin, Neck, Transpa Incisal, Special Incisal, Inter Incisal, Occlusal Dentin, Mamelon, Effect, Opal Effect, Incisal Edge, Transpa, Transpa Clear, Add-On**

(Contd. of page 3)

- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**  
Usual hygienic measures for dental practice.  
Wash hands before breaks and at the end of work.  
Keep away from foodstuffs, beverages and feed.  
Avoid contact with the skin.
- **Breathing equipment:** Not required.
- **Protection of hands:**



Protective gloves

After use of gloves apply skin-cleaning agents and skin cosmetics.

- **Material of gloves**  
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
- **Penetration time of glove material**  
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Not suitable are gloves made of the following materials:**  
Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.
- **Eye protection:** Safety glasses

### 9 Physical and chemical properties

- |  |   |
|--|---|
| · <b>Information on basic physical and chemical properties</b> |   |
| · <b>General Information</b>                                   |   |
| · <b>Appearance:</b>   |   |
| Form:  | Pasty   |
| Color:   | According to product specification            |
| Odor:  | Characteristic                                |
| Odour threshold:   | Not determined.                               |
| pH-value:  | Not applicable.                               |
| · <b>Change in condition</b>                                   |   |
| Melting point/Melting range:                                   | Undetermined.                                 |
| Boiling point/Boiling range:                                   | Undetermined.                                 |
| Flash point:   | Not applicable.                               |
| Flammability (solid, gaseous):                                 | Not applicable.                               |
| Auto igniting:   | Product is not selfigniting.                  |
| Danger of explosion:   | Product does not present an explosion hazard. |
| · <b>Explosion limits:</b>                                     |   |
| Lower:   | Not determined.                               |
| Upper:   | Not determined.                               |
| Vapor pressure:  | Not determined.                               |

(Contd. on page 5)

## Safety Data Sheet

acc. to OSHA HCS

Printing date 04/29/2015

Version number 4

Reviewed on 04/16/2015

**Trade name: SR Adoro Incisal, Dentin, Deep Dentin, Neck, Transpa Incisal, Special Incisal, Inter Incisal, Occlusal Dentin, Mamelon, Effect, Opal Effect, Incisal Edge, Transpa, Transpa Clear, Add-On**

(Contd. of page 4)

· <b>Density:</b>	Not determined.
· <b>Relative density</b>	Not determined.
· <b>Vapour density</b>	Not determined.
· <b>Evaporation rate</b>	Not determined.
· <b>Solubility in / Miscibility with Water:</b>	Nearly insoluble.
· <b>Partition coefficient (n-octanol/water):</b>	Not determined.
· <b>Viscosity:</b>	
· <b>Dynamic:</b>	Not determined.
· <b>Kinematic:</b>	Not determined.
· <b>Other information</b>	No further relevant information available.

### 10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability**  
Stable under normal handling and storage conditions.  
Photoreactive (Danger of polymerization).
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid**  
Protect from exposure to the light.  
Keep away from heat and direct sunlight.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** None under normal conditions of storage and use.

### 11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**
- **Primary irritant effect:**
- **on the skin:** No irritant effect.
- **on the eye:** Irritating effect.
- **Sensitization:** Sensitization possible through skin contact.
- **Additional toxicological information:**
- **Carcinogenic categories**

#### · IARC (International Agency for Research on Cancer)

9011-14-7 polymethylmethacrylate

3

#### · NTP (National Toxicology Program)

None of the ingredients is listed.

#### · OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

(Contd. on page 6)

## Safety Data Sheet

acc. to OSHA HCS

Printing date 04/29/2015

Version number 4

Reviewed on 04/16/2015

**Trade name: SR Adoro Incisal, Dentin, Deep Dentin, Neck, Transpa Incisal, Special Incisal, Inter Incisal, Occlusal Dentin, Mamelon, Effect, Opal Effect, Incisal Edge, Transpa, Transpa Clear, Add-On**

(Contd. of page 5)

### 12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Behavior in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:**  
Water hazard class 2 (Self-assessment): hazardous for water  
Do not allow product to reach ground water, water course or sewage system.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

### 13 Disposal considerations

- **Waste treatment methods**
- **Recommendation:**  
Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to official regulations.

### 14 Transport information

- |   |  |
|---|--|
| <ul style="list-style-type: none"> <li>· <b>UN-Number</b></li> <li>· <b>DOT, ADR, ADN, IMDG, IATA</b></li> </ul>  | Void   |
| <ul style="list-style-type: none"> <li>· <b>UN proper shipping name</b></li> <li>· <b>DOT, ADR, ADN, IMDG, IATA</b></li> </ul>                            | Void   |
| <ul style="list-style-type: none"> <li>· <b>Transport hazard class(es)</b></li> <li>· <b>DOT, ADR, ADN, IMDG, IATA</b></li> <li>· <b>Class</b></li> </ul> | Void   |
| <ul style="list-style-type: none"> <li>· <b>Packing group</b></li> <li>· <b>DOT, ADR, IMDG, IATA</b></li> </ul>   | Void   |
| <ul style="list-style-type: none"> <li>· <b>Environmental hazards:</b></li> <li>· <b>Marine pollutant:</b></li> </ul>                                     | No   |
| <ul style="list-style-type: none"> <li>· <b>Special precautions for user</b></li> </ul>   | Not applicable.  |
| <ul style="list-style-type: none"> <li>· <b>Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code</b></li> </ul>                        | Not applicable.  |
| <ul style="list-style-type: none"> <li>· <b>Transport/Additional information:</b></li> </ul>  | Product is not classified as a dangerous good for transport (ADR, IMDG, IATA). |

(Contd. on page 7)

US

## Safety Data Sheet

acc. to OSHA HCS

Printing date 04/29/2015

Version number 4

Reviewed on 04/16/2015

**Trade name: SR Adoro Incisal, Dentin, Deep Dentin, Neck, Transpa Incisal, Special Incisal, Inter Incisal, Occlusal Dentin, Mamelon, Effect, Opal Effect, Incisal Edge, Transpa, Transpa Clear, Add-On**

(Contd. of page 6)

· UN "Model Regulation": -

### 15 Regulatory information

· Safety, health and environmental regulations/legislation specific for the substance or mixture

· Sara

· Section 355 (extremely hazardous substances):

None of the ingredients is listed.

· Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

· TSCA (Toxic Substances Control Act):

All ingredients are listed.

· Proposition 65

· Chemicals known to cause cancer:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

· Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

· Carcinogenic categories

· EPA (Environmental Protection Agency)

None of the ingredients is listed.

· TLV (Threshold Limit Value established by ACGIH)

None of the ingredients is listed.

· NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

· GHS label elements

The product is classified and labeled according to the Globally Harmonized System (GHS).

· Hazard pictograms



GHS07

· Signal word Warning

· Hazard-determining components of labeling:

urethane dimethacrylate

1,10-decandiol dimethacrylate

· Hazard statements

Causes serious eye irritation.

(Contd. on page 8)

US

## Safety Data Sheet acc. to OSHA HCS

Printing date 04/29/2015

Version number 4

Reviewed on 04/16/2015

**Trade name: SR Adoro Incisal, Dentin, Deep Dentin, Neck, Transpa Incisal, Special Incisal, Inter Incisal, Occlusal Dentin, Mamelon, Effect, Opal Effect, Incisal Edge, Transpa, Transpa Clear, Add-On**

(Contd. of page 7)

*May cause an allergic skin reaction.*

· **Precautionary statements**

*Wear protective gloves/protective clothing/eye protection/face protection.*

*Wash contaminated clothing before reuse.*

*If skin irritation or rash occurs: Get medical advice/attention.*

*If eye irritation persists: Get medical advice/attention.*

· **National regulations:**

· **Other regulations, limitations and prohibitive regulations**

*The product is a medical device according to the Directive 93/42/EEC.*

*This product is classified as a medical device under US and Canadian regulations and has been reviewed by the US Food and Drug Administration and Health Canada.*

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

### 16 Other information

*This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.*

· **Date of preparation / last revision 04/29/2015 / 3**

· **Abbreviations and acronyms:**

*ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)*

*IMDG: International Maritime Code for Dangerous Goods*

*DOT: US Department of Transportation*

*IATA: International Air Transport Association*

*ACGIH: American Conference of Governmental Industrial Hygienists*

*EINECS: European Inventory of Existing Commercial Chemical Substances*

*ELINCS: European List of Notified Chemical Substances*

*CAS: Chemical Abstracts Service (division of the American Chemical Society)*

*NFPA: National Fire Protection Association (USA)*

*HMIS: Hazardous Materials Identification System (USA)*

*Eye Irrit. 2A: Serious eye damage/eye irritation, Hazard Category 2A*

*Skin Sens. 1: Sensitisation - Skin, Hazard Category I*

· **\* Data compared to the previous version altered.**