

# Evetric®

State-of-the-art universal nanohybrid composite

Discover the  
secret of  
invisible fillings.



ivoclar  
vivadent®  
passion vision innovation



# Evetric

State-of-the-art universal  
nanohybrid composite

- High-end universal nanohybrid composite resin
- Different degrees of translucency of the dentin and enamel shades for very natural-looking restorations of both posterior and anterior teeth
- Available in 10 shades

# Invisible fillings for highly esthetic results

- Integrated nano pigments for an exceptional chameleon effect
- Homogeneously distributed pigments for shades with an intensive effect

## Benefits for you

Improved shade blend with the surrounding natural dentition

## Benefits for your patients


Fillings that look like the natural tooth structure and are therefore almost invisible



Prof. Dr Vlatko Pandurić, University of Zagreb, Croatia

“ Evetric provides  
excellent shade  
adaptation ”





### Visible if clinically necessary

- Use of a functional filler for an outstanding radiopacity of 400% Al


#### Benefits for you

The composite resin is clearly distinguishable on X-rays to help you in your diagnosis.

#### Benefits for your patients

Unnecessary drilling can be avoided by reliably identifying restorations in future examinations.

Radiopacity [Al %]



\* These brands are not registered trademarks of Ivoclar Vivadent AG.  
Measurement acc. to ISO 4049, R&D Ivoclar Vivadent AG, Report April 2014

## Advanced technology makes work easier

- Integrated nano-modifier for excellent modelling and non-stick properties
- Long working time due to the composite's reduced sensitivity to light


### Benefits for you

Modelling without time pressure enables you to produce high-precision esthetic outcomes.

### Benefits for your patients

Excellent results and a natural-looking smile



## Sensitivity to ambient light [s]



\* These brands are not registered trademarks of Ivoclar Vivadent AG.  
Measurement acc. to ISO 4049, R&D Ivoclar Vivadent AG, Report April 2014



# Delivery Forms

Description	Quantity	Shade	Art. No.
 <b>Evetric Assortment</b> Shade Guide	8 x 3.5 g	A1, A2, A3, A3.5, B2, C2, T, A3.5 Dentin	640102AN
 <b>Evetric System Kit</b> Evetric Bond 1x 6 g	4 x 3.5 g	A2, A3	668409AN
<b>Evetric Refill</b>	1 x 3.5 g	A1	639668AN
<b>Evetric Refill</b>	1 x 3.5 g	A2	639304AN
<b>Evetric Refill</b>	1 x 3.5 g	A3	639305AN
<b>Evetric Refill</b>	1 x 3.5 g	A3.5	639306AN
<b>Evetric Refill</b>	1 x 3.5 g	A3.5 Dentin	639327AN
<b>Evetric Refill</b>	1 x 3.5 g	B2	639669AN
<b>Evetric Refill</b>	1 x 3.5 g	B2 Dentin	639673AN
<b>Evetric Refill</b>	1 x 3.5 g	B3	639670AN
<b>Evetric Refill</b>	1 x 3.5 g	C2	639671AN
<b>Evetric Refill</b>	1 x 3.5 g	T	639672AN
<b>Evetric Shade Guide</b>			640013



## Direct Restoratives

Evetric® forms a part of the “Direct Restoratives” product category. The products of this category cover the procedure involved in the direct restoration of teeth – from preparation to restoration care. The products are optimally coordinated with each other and enable successful processing and application.



### THESE ARE FURTHER PRODUCTS OF THIS CATEGORY:

#### Bluephase® Style

The curing light



The smallest LED for every use

- Every hand – comfortable to hold for men and women
- Every material – universal use due to Polywave® LED with broadband spectrum
- Every time – optional corded operation due to Click & Cure

#### Evetric® Bond

State-of-the-art self-etch adhesive



Discover the secret of convenient bonding

- Fast and esthetic due to two bonding options
- Optimal handling properties due to homogeneous film formation
- Reliable composite restorations of all classes

Would you like to know more about the products of the “Direct Restoratives” category? Simply get in touch with your contact person at Ivoclar Vivadent or visit [www.ivoclarvivadent.com](http://www.ivoclarvivadent.com) for more information.