



IPS d.SIGN® LIEFERFORMEN

Chromascop	556458	Basic Kit	Opaquer pastes, intensive opaquer pastes, dentin, incisal, transparent, add-on, and glazing materials, as well as all the necessary liquids. The Kit also contains a materials shade guide for the most popular Chromascop shades.
	556534	Deep Dentin Kit	Deep dentin materials in the 10 most popular Chromascop shades. They are used to intensify the chroma in layers of limited thickness.
	556547	Margin Kit	Margin materials in the 10 most popular Chromascop shades and 4 intensive materials for special effects. 1 margin add-on material is used for final adjustment of the ceramic shoulder.
	557242	Trial Kit	Starter Kit to try out the system.

A-D	558195	Basic Kit	Opaquer pastes, intensive opaquer pastes, dentin, incisal, transparent, add-on, and glazing materials, as well as all the necessary liquids. The Kit also contains a materials shade guide for the most popular A-D shades.
	558256	Deep Dentin Kit	Deep dentin materials in the 7 most popular A-D shades. They are used to intensify the chroma in layers of limited thickness.
	558266	Margin Kit	Margin materials in the 8 most popular A-D shades and 4 intensive materials for special effects. 1 margin add-on material is used for final adjustment of the ceramic shoulder.
	558196	Trial Kit	Starter Kit to try out the system.

Independent of any shade system	602171	Bleach Kit BL	4 accurately coordinated shades meet the highest demands related to contemporary "white aesthetics".
	556567	Impulse Kit 1	15 ready-mixed individual fluorapatite leucite ceramic materials ranging from Mamelon to Opal materials facilitate the application of true-to-nature effects.
	556568	Impulse Kit 2	14 ready-mixed Impulse materials based on the inspiration and philosophies of renowned experts, such as Enrico Steger, Donald F. Cornell, and Dr. Robert R. Winter.
	556613	Gingiva Kit	5 true-to-nature gingival materials that can be individually used for each patient situation in combination with the 4 gingival modifiers.
	556624	Essence Kit	4 ceramic stains in powder form. They may be applied alone or mixed with other materials.
	556630	Stains Kit	14 intensive stains in paste form to reproduce natural characteristics on the ceramic surface.
	561571	Shade Kit	7 dentin stains in paste form for subsequent shade adjustments of IPS d.SIGN restorations.

Art.No.	Description	Quantity	Shade
---------	-------------	----------	-------

IPS d.SIGN Basic Kit

The IPS d.SIGN Basic Kits have been reduced to the most popular shades (8 x A–D / 10 x Chromascop). Therefore, the most common basic tooth shades can be easily reproduced.

558195



IPS d.SIGN Basic Kit A–D

Contents:

- 8x 3 g IPS d.SIGN Opaquer; Shades: A2, A3, A3.5, B2, B3, C3, D2, D3
- 4x 3 g IPS d.SIGN Intensive Opaquer; Shades: white, violet, brown, incisal
- 8x 20 g IPS d.SIGN Dentin; Shades: A2, A3, A3.5, B2, B3, C3, D2, D3
- 1x 20 g IPS d.SIGN Cervical Dentin D2/D3
- 3x 20 g IPS d.SIGN Incisal; Gradations: TS1, TS2, TS3
- 1x 20 g IPS d.SIGN Transpa neutral
- 1x 3 g IPS d.SIGN Glaze
- 1x 20 g IPS d.SIGN Add-on
- 1x 15 ml IPS d.SIGN Opaquer Liquid
- 2x 60 ml IPS d.SIGN Build-up Liquid, Medium and Optimum 2
- 1x 15 ml IPS d.SIGN Glaze and Stain Liquid
- 1x 50 ml Model Sealer
- 1x 15 ml IPS Ceramic Separating Liquid
- 1x Opaquer A–D materials shade guide
- 1x Dentin A–D materials shade guide
- 1x Incisal A–D materials shade guide

556458



IPS d.SIGN Basic Kit Chromascop

Contents:

- 10x 3 g IPS d.SIGN Opaquer; Shades: 130, 140, 210, 220, 230, 310, 410, 420, 430, 510
- 4x 3 g IPS d.SIGN Intensive Opaquer; Shades: white, violet, brown, incisal
- 10x 20 g IPS d.SIGN Dentin; Shades: 130, 140, 210, 220, 230, 310, 410, 420, 430, 510
- 3x 20 g IPS d.SIGN Incisal; Gradations: S1, S2, S3
- 1x 20 g IPS d. SIGN Transpa neutral
- 1x 3 g IPS d.SIGN Glaze
- 1x 20 g IPS d.SIGN Add-on
- 1x 15 ml IPS d.SIGN Opaquer Liquid
- 2x 60 ml IPS d.SIGN Build-up Liquid, Medium and Optimum 2
- 1x 15 ml IPS d.SIGN Glaze and Stain Liquid
- 1x 50 ml IPS Model Sealer
- 1x 15 ml IPS Ceramic Separating Liquid
- 1x Opaquer materials shade guide
- 1x Dentin materials shade guide
- 1x Incisal materials shade guide
- 1x Chromascop shade guide

IPS d.SIGN Trial Kit

The starter kit will give you a first impression. All the important material components are contained in this assortment.

558196

IPS d.SIGN Trial Kit A–D

Contents:

- 1x 3 g IPS d.SIGN Opaquer; Shade: A3
- 1x 3 g IPS d.SIGN Dentin; Shade: A3
- 1x 3 g IPS d.SIGN Incisal; Gradaation: TS2
- 1x 3 g IPS d. SIGN Deep Dentin; Shade: A3
- 1x 3 g IPS d.SIGN Margin; Shade: A3
- 1x 3 g IPS d.SIGN Transpa; Shade: neutral
- 1x 3 g IPS d.SIGN Occlusal Dentin; Shade: orange
- 2x 3 g IPS d.SIGN Effect; Shades: 1, 3
- 1x 3 g IPS d.SIGN Mamelon; Shade: light
- 1x 3 g IPS d.SIGN Glaze
- 1x 5 ml IPS d.SIGN Build-up Liquid, Medium
- 1x 5 ml IPS d.SIGN Glaze and Stain Liquid
- 1x 5 ml IPS d.SIGN Margin Build-up Liquid
- 1x materials shade guide

Art.No.	Description	Quantity	Shade
---------	-------------	----------	-------

557242



IPS d.SIGN Trial Kit Chromascop

Contents:

- 1x 3 g IPS d.SIGN Opaquer; Shade: 210
- 1x 3 g IPS d.SIGN Dentin; Shade: 210
- 1x 3 g IPS d.SIGN Incisal; Gradation: S2
- 1x 3 g IPS d.SIGN Deep Dentin; Shade: 210
- 1x 3 g IPS d.SIGN Margin; Shade: 210
- 1x 3 g IPS d.SIGN Transpa; Shade: neutral
- 1x 3 g IPS d.SIGN Occlusal Dentin; Shade: orange
- 2x 3 g IPS d.SIGN Effect; Shades: 1, 3
- 1x 3 g IPS d.SIGN Mamelon; Shade: light
- 1x 3 g IPS d.SIGN Glaze
- 1x 5 ml IPS d.SIGN Build-up Liquid, Medium
- 1x 5 ml IPS d.SIGN Glaze and Stain Liquid
- 1x 5 ml IPS d.SIGN Margin Build-up Liquid
- 1x materials shade guide

IPS d.SIGN Deep Dentin Kit

With the IPS d.SIGN Deep Dentin materials, restorations that demonstrate true-to-nature shade effects can be achieved even with limited layer thicknesses. Particularly if space is limited, IPS d.SIGN Deep Dentin materials can be used to intensify the chroma. The kit contains the most popular shades.

558256



IPS d.SIGN Deep Dentin Kit A-D

Contents:

- 7x 20 g IPS d.SIGN Deep Dentin; Shades: A2, A3, A3.5, B2, B3, C3, D2/D3
- 1x Deep Dentin materials shade guide

556534



IPS d.SIGN Deep Dentin Kit Chromascop

Contents:

- 10x 20 g IPS d.SIGN Deep Dentin; Shades: 130, 140, 210, 220, 230, 310, 410, 420, 430, 510
- 1x Deep Dentin materials shade guide

IPS d.SIGN Margin Kit

The IPS d.SIGN Margin material permits the design of a ceramic shoulder. The kit contains the most popular shades (8 x A-D / 10 x Chromascop), the corresponding Intensive shades, as well as the required liquids.

558266



IPS d.SIGN Margin Kit A-D

Contents:

- 8x 20 g IPS d.SIGN Margin; Shades: A2, A3, A3.5, B2, B3, C3, D2/D3
- 4x 20 g IPS d.SIGN Intensive Margin; Shades: yellow, orange, orange-pink, opaque
- 1x 20 g IPS d.SIGN Margin add-on
- 1x 60 ml IPS d.SIGN Margin Build-up Liquid
- 1x 20 ml IPS Margin Sealer
- 1x 15 ml IPS Ceramic Separating Liquid
- 1x Margin materials shade guide

556548



IPS d.SIGN Margin Kit Chromascop

Contents:

- 10x 20 g IPS d.SIGN Margin; Shades: 130, 140, 210, 220, 230, 310, 410, 420, 430, 510
- 4x 20 g IPS d.SIGN Intensive Margin; Shades: yellow, orange, orange-pink, opaque
- 1x 20 g IPS d.SIGN Margin add-on
- 1x 60 ml IPS d.SIGN Margin Build-up Liquid
- 1x 20 ml IPS Margin Sealer
- 1x 15 ml IPS Ceramic Separating Liquid
- 1x Margin materials shade guide

Art.No.	Description	Quantity	Shade
---------	-------------	----------	-------

602171



IPS d.SIGN Bleach Kit BL

The accurately coordinated shades meet the highest demands related to contemporary “white aesthetics”.

Contents:

- 2x 3 g IPS d.SIGN Opaquer; Shades: BL1/BL2, BL3/BL4
- 4x 20 g IPS d.SIGN Dentin; Shades: BL1, BL2, BL3, BL4
- 1x 20 g IPS d.SIGN Incisal; Shade: BL
- 2x 20 g IPS d.SIGN Margin; Shades: BL1, BL4*
- 2x 20 g IPS d.SIGN Deep Dentin; Shades: BL1, BL4*
- 1x 20 g IPS d.SIGN Add-On; Shade: BL
- 1x IPS d.SIGN/IPS InLine materials shade guide BL
- 1x Shade guide Bleach BL

* The Margin and Deep Dentin materials are only available in shades BL1 and BL4. The shades BL2 and BL3 are achieved with the following mixing ratios:

- BL2: 2/3 BL1 : 1/3 BL4
- BL3: 1/3 BL1 : 2/3 BL4

556567



IPS d.SIGN Impulse 1 Kit

15 ready-mixed individual fluorapatite leucite ceramic materials ranging from Mamelon to Opal materials facilitate the application of true-to-nature effects.

Contents:

- 2x 20 g IPS d.SIGN Occlusal Dentin; Shades: orange, brown
- 3x 20 g IPS d.SIGN Mamelon; Shades: light, yellow-orange, salmon
- 6x 20 g IPS d.SIGN Effect; Gradations: Effect 1, Effect 2, Effect 3, Effect 4, Effect 5, Effect violet
- 3x 20 g IPS d.SIGN Transpa; Shades: blue, brown, grey, orange-grey
- 1x 20 g IPS d.SIGN Incisal Edge
- 1x Impulse 1 materials shade guide

556568



IPS d.SIGN Impulse 2 Kit

14 ready-mixed Impulse materials based on the inspiration and philosophies of renowned experts, such as Enrico Steger, Donald F. Cornell, and Dr. Robert R. Winter.

Contents:

- 3x 20 g IPS d.SIGN Brilliant Dentin; Shades: yellow, orange, white
- 2x 20 g IPS d.SIGN Palatal Dentin; Shades: yellow, red
- 2x 20 g IPS d.SIGN Insert; Shades: orange, grey
- 2x 20 g IPS d.SIGN Special Incisal; Shades: yellow, grey
- 1x 20 g IPS d.SIGN Inter Incisal; Shades: white-blue
- 4x 20 g IPS d.SIGN Cervical Incisal; Shades: yellow, orange-pink, khaki, orange
- 1x Impulse 2 materials shade guide

556613



IPS d.SIGN Gingiva Kit

These five lifelike Gingiva materials are used in combination with the four Gingiva Modifiers depending on the individual clinical situation. The shade gradations range from orange to reddish to bluish.

Contents:

- 1x 3 g IPS d.SIGN Gingiva Opaquer
- 5x 20 g IPS d.SIGN Gingiva; Gradations: G1, G2, G3, G4, G5
- 4x 20 g IPS d.SIGN Gingiva Modifier; Gradations: GM1, GM2, GM3, GM4
- 1x Gingiva materials shade guide

556624



IPS d.SIGN Essence Kit

IPS d.SIGN Essence materials are ceramic stains in powder form that may be mixed with IPS d.SIGN layering materials in order to intensify the shade. Furthermore, IPS d.SIGN Essence can be used for surface characterizations to reproduce the special characteristics of the natural tooth.

Contents:

- 4x 1 g IPS d.SIGN Essence; Shades: orange, yellow, grey, brown
- 1x 15 ml IPS d.SIGN Glaze and Stain Liquid 15 ml

Art.No.	Description	Quantity	Shade
---------	-------------	----------	-------



556630

IPS d.SIGN Stains Kit

These intensive stains in paste form are available in 14 different shades. They are used to reproduce natural characteristics on the ceramic surface. With the three primary colours basic blue, basic red, and basic yellow, custom shades may be developed.

Contents:

- 14x 1 g IPS d.SIGN Stains; Shades: white, mahogany, khaki, orange, grey, vanille, crackliner, olive, yellow, black, maroon, basic blue, basic red, basic yellow
- 1x 3 g IPS d.SIGN Glaze
- 1x 15 ml IPS d.SIGN Glaze and Stain Liquid



561571

IPS d.SIGN Shade Kit

These 7 dentin stains in paste form make shade adjustments to IPS d.SIGN restorations possible. They are coordinated with the 20 Chromascop shades and the 16 A–D shades.

Contents:

- 7x 3 g IPS d.SIGN Shade; Shades: 1, 2, 3, 4, 5, 6, 7

Refills

IPS d.SIGN Opaquer A–D Refill 3 g

558198	IPS d.SIGN Opaquer A-D	3 g	A1
558199	IPS d.SIGN Opaquer A-D	3 g	A2
558200	IPS d.SIGN Opaquer A-D	3 g	A3
558201	IPS d.SIGN Opaquer A-D	3 g	A3,5
558202	IPS d.SIGN Opaquer A-D	3 g	A4
558203	IPS d.SIGN Opaquer A-D	3 g	B1
558204	IPS d.SIGN Opaquer A-D	3 g	B2
558205	IPS d.SIGN Opaquer A-D	3 g	B3
558206	IPS d.SIGN Opaquer A-D	3 g	B4
558207	IPS d.SIGN Opaquer A-D	3 g	C1
558208	IPS d.SIGN Opaquer A-D	3 g	C2
558209	IPS d.SIGN Opaquer A-D	3 g	C3
558210	IPS d.SIGN Opaquer A-D	3 g	C4
558211	IPS d.SIGN Opaquer A-D	3 g	D2
558212	IPS d.SIGN Opaquer A-D	3 g	D3
558213	IPS d.SIGN Opaquer A-D	3 g	D4

IPS d.SIGN Opaquer Chromascop Refill 3 g

556460	IPS d.SIGN Opaquer	3 g	110
556461	IPS d.SIGN Opaquer	3 g	120
556462	IPS d.SIGN Opaquer	3 g	130
556463	IPS d.SIGN Opaquer	3 g	140
556464	IPS d.SIGN Opaquer	3 g	210
556465	IPS d.SIGN Opaquer	3 g	220
556466	IPS d.SIGN Opaquer	3 g	230
556467	IPS d.SIGN Opaquer	3 g	240
556468	IPS d.SIGN Opaquer	3 g	310
556469	IPS d.SIGN Opaquer	3 g	320
556470	IPS d.SIGN Opaquer	3 g	330
556471	IPS d.SIGN Opaquer	3 g	340
556472	IPS d.SIGN Opaquer	3 g	410
556473	IPS d.SIGN Opaquer	3 g	420
556474	IPS d.SIGN Opaquer	3 g	430
556475	IPS d.SIGN Opaquer	3 g	440
556476	IPS d.SIGN Opaquer	3 g	510
556477	IPS d.SIGN Opaquer	3 g	520
556478	IPS d.SIGN Opaquer	3 g	530
556479	IPS d.SIGN Opaquer	3 g	540
556480	IPS d.SIGN Opaquer	3 g	white
556481	IPS d.SIGN Opaquer	3 g	violet
556482	IPS d.SIGN Opaquer	3 g	brown
556483	IPS d.SIGN Opaquer	3 g	incisal





Art.No.	Description	Quantity	Shade
599632	IPS d.SIGN Opaquer F Refill 1 g IPS d.SIGN Opaquer	1 g	F
	IPS d.SIGN Dentin A–D Refill 20 g		
558215	IPS d.SIGN Dentin A-D	20 g	A1
558216	IPS d.SIGN Dentin A-D	20 g	A2
558217	IPS d.SIGN Dentin A-D	20 g	A3
558218	IPS d.SIGN Dentin A-D	20 g	A3,5
558219	IPS d.SIGN Dentin A-D	20 g	A4
558220	IPS d.SIGN Dentin A-D	20 g	B1
558221	IPS d.SIGN Dentin A-D	20 g	B2
558222	IPS d.SIGN Dentin A-D	20 g	B3
558223	IPS d.SIGN Dentin A-D	20 g	B4
558224	IPS d.SIGN Dentin A-D	20 g	C1
558225	IPS d.SIGN Dentin A-D	20 g	C2
558226	IPS d.SIGN Dentin A-D	20 g	C3
558227	IPS d.SIGN Dentin A-D	20 g	C4
558228	IPS d.SIGN Dentin A-D	20 g	D2
558229	IPS d.SIGN Dentin A-D	20 g	D3
558230	IPS d.SIGN Dentin A-D	20 g	D4
558684	IPS d.SIGN Cervical Dentin A-D	20 g	D2/D3
	IPS d.SIGN Dentin A–D Refill 100 g		
558232	IPS d.SIGN Dentin A-D	100 g	A1
558233	IPS d.SIGN Dentin A-D	100 g	A2
558234	IPS d.SIGN Dentin A-D	100 g	A3
558235	IPS d.SIGN Dentin A-D	100 g	A3,5
558236	IPS d.SIGN Dentin A-D	100 g	A4
558237	IPS d.SIGN Dentin A-D	100 g	B1
558238	IPS d.SIGN Dentin A-D	100 g	B2
558239	IPS d.SIGN Dentin A-D	100 g	B3
558240	IPS d.SIGN Dentin A-D	100 g	B4
558241	IPS d.SIGN Dentin A-D	100 g	C1
558242	IPS d.SIGN Dentin A-D	100 g	C2
558243	IPS d.SIGN Dentin A-D	100 g	C3
558244	IPS d.SIGN Dentin A-D	100 g	C4
558245	IPS d.SIGN Dentin A-D	100 g	D2
558246	IPS d.SIGN Dentin A-D	100 g	D3
558247	IPS d.SIGN Dentin A-D	100 g	D4
558685	IPS d.SIGN Cervical Dentin AD	100 g	D2/D3
	IPS d.SIGN Dentin A–D Refill 250 g		
563507	IPS d.SIGN Dentin A-D	250 g	A1
563508	IPS d.SIGN Dentin A-D	250 g	A2
563509	IPS d.SIGN Dentin A-D	250 g	A3
563510	IPS d.SIGN Dentin A-D	250 g	A3,5
563511	IPS d.SIGN Dentin A-D	250 g	A4
563512	IPS d.SIGN Dentin A-D	250 g	B1
563513	IPS d.SIGN Dentin A-D	250 g	B2
563514	IPS d.SIGN Dentin A-D	250 g	B3
563515	IPS d.SIGN Dentin A-D	250 g	B4
563516	IPS d.SIGN Dentin A-D	250 g	C1
563517	IPS d.SIGN Dentin A-D	250 g	C2
563518	IPS d.SIGN Dentin A-D	250 g	C3
563519	IPS d.SIGN Dentin A-D	250 g	C4
563520	IPS d.SIGN Dentin A-D	250 g	D2
563521	IPS d.SIGN Dentin A-D	250 g	D3
563522	IPS d.SIGN Dentin A-D	250 g	D4
563523	IPS d.SIGN Cervical Dentin A-D	250 g	D2/D3

Art.No.	Description	Quantity	Shade
---------	-------------	----------	-------



IPS d.SIGN Dentin Chromascop Refill 20 g			
558188	IPS d.SIGN Bleach Dentin	20 g	010
558189	IPS d.SIGN Bleach Dentin	20 g	020
558190	IPS d.SIGN Bleach Dentin	20 g	030
558191	IPS d.SIGN Bleach Dentin	20 g	040
556485	IPS d.SIGN Dentin	20 g	110
556486	IPS d.SIGN Dentin	20 g	120
556487	IPS d.SIGN Dentin	20 g	130
556488	IPS d.SIGN Dentin	20 g	140
556489	IPS d.SIGN Dentin	20 g	210
556490	IPS d.SIGN Dentin	20 g	220
556491	IPS d.SIGN Dentin	20 g	230
556492	IPS d.SIGN Dentin	20 g	240
556493	IPS d.SIGN Dentin	20 g	310
556494	IPS d.SIGN Dentin	20 g	320
556495	IPS d.SIGN Dentin	20 g	330
556496	IPS d.SIGN Dentin	20 g	340
556497	IPS d.SIGN Dentin	20 g	410
556498	IPS d.SIGN Dentin	20 g	420
556499	IPS d.SIGN Dentin	20 g	430
556500	IPS d.SIGN Dentin	20 g	440
556501	IPS d.SIGN Dentin	20 g	510
556502	IPS d.SIGN Dentin	20 g	520
556503	IPS d.SIGN Dentin	20 g	530
556504	IPS d.SIGN Dentin	20 g	540

IPS d.SIGN Dentin Chromascop Refill 100 g			
556506	IPS d.SIGN Dentin	100 g	110
556507	IPS d.SIGN Dentin	100 g	120
556508	IPS d.SIGN Dentin	100 g	130
556509	IPS d.SIGN Dentin	100 g	140
556510	IPS d.SIGN Dentin	100 g	210
556511	IPS d.SIGN Dentin	100 g	220
556512	IPS d.SIGN Dentin	100 g	230
556513	IPS d.SIGN Dentin	100 g	240
556514	IPS d.SIGN Dentin	100 g	310
556515	IPS d.SIGN Dentin	100 g	320
556516	IPS d.SIGN Dentin	100 g	330
556517	IPS d.SIGN Dentin	100 g	340
556518	IPS d.SIGN Dentin	100 g	410
556519	IPS d.SIGN Dentin	100 g	420
556520	IPS d.SIGN Dentin	100 g	430
556521	IPS d.SIGN Dentin	100 g	440
556522	IPS d.SIGN Dentin	100 g	510
556523	IPS d.SIGN Dentin	100 g	520
556524	IPS d.SIGN Dentin	100 g	530
556525	IPS d.SIGN Dentin	100 g	540

IPS d.SIGN Dentin Chromascop Refill 250 g			
563484	IPS d.SIGN Dentin	250 g	110
563485	IPS d.SIGN Dentin	250 g	120
563486	IPS d.SIGN Dentin	250 g	130
563487	IPS d.SIGN Dentin	250 g	140
563488	IPS d.SIGN Dentin	250 g	210
563489	IPS d.SIGN Dentin	250 g	220
563490	IPS d.SIGN Dentin	250 g	230
563491	IPS d.SIGN Dentin	250 g	240
563492	IPS d.SIGN Dentin	250 g	310
563493	IPS d.SIGN Dentin	250 g	320
563494	IPS d.SIGN Dentin	250 g	330
563495	IPS d.SIGN Dentin	250 g	340
563496	IPS d.SIGN Dentin	250 g	410
563497	IPS d.SIGN Dentin	250 g	420
563498	IPS d.SIGN Dentin	250 g	430
563499	IPS d.SIGN Dentin	250 g	440
563500	IPS d.SIGN Dentin	250 g	510
563501	IPS d.SIGN Dentin	250 g	520
563502	IPS d.SIGN Dentin	250 g	530
563503	IPS d.SIGN Dentin	250 g	540


Art.No.	Description	Quantity	Shade
IPS d.SIGN T-Incisal Refill 20 g			
558249	IPS d.SIGN T-Incisal	20 g	TS-1
558250	IPS d.SIGN T-Incisal	20 g	TS-2
558251	IPS d.SIGN T-Incisal	20 g	TS-3
IPS d.SIGN T-Incisal Refill 100 g			
558253	IPS d.SIGN T-Incisal	100 g	TS-1
558254	IPS d.SIGN T-Incisal	100 g	TS-2
558255	IPS d.SIGN T-Incisal	100 g	TS-3
IPS d.SIGN T-Incisal Refill 250 g			
563524	IPS d.SIGN T-Incisal	250 g	TS-1
563525	IPS d.SIGN T-Incisal	250 g	TS-2
563526	IPS d.SIGN T-Incisal	250 g	TS-3
IPS d.SIGN Incisal Refill 20 g			
556527	IPS d.SIGN Incisal	20 g	1
556528	IPS d.SIGN Incisal	20 g	2
556529	IPS d.SIGN Incisal	20 g	3
IPS d.SIGN Incisal Refill 100 g			
556531	IPS d.SIGN Incisal	100 g	1
556532	IPS d.SIGN Incisal	100 g	2
556533	IPS d.SIGN Incisal	100 g	3
IPS d.SIGN Incisal Refill 250 g			
563504	IPS d.SIGN Incisal	250 g	1
563505	IPS d.SIGN Incisal	250 g	2
563506	IPS d.SIGN Incisal	250 g	3
IPS d.SIGN Deep Dentin A-D Refill 20 g			
560512	IPS d.SIGN Deep Dentin A-D	20 g	A1
558258	IPS d.SIGN Deep Dentin A-D	20 g	A2
558259	IPS d.SIGN Deep Dentin A-D	20 g	A3
558260	IPS d.SIGN Deep Dentin A-D	20 g	A3,5
560513	IPS d.SIGN Deep Dentin A-D	20 g	A4
560514	IPS d.SIGN Deep Dentin A-D	20 g	B1
558262	IPS d.SIGN Deep Dentin A-D	20 g	B2
558263	IPS d.SIGN Deep Dentin A-D	20 g	B3
560515	IPS d.SIGN Deep Dentin A-D	20 g	B4
560516	IPS d.SIGN Deep Dentin A-D	20 g	C1
560517	IPS d.SIGN Deep Dentin A-D	20 g	C2
558264	IPS d.SIGN Deep Dentin A-D	20 g	C3
560518	IPS d.SIGN Deep Dentin A-D	20 g	C4
558261	IPS d.SIGN Deep Dentin A-D	20 g	D2/D3
560519	IPS d.SIGN Deep Dentin A-D	20 g	D4
IPS d.SIGN Deep Dentin Chromascoop Refill 20 g			
559682	IPS d.SIGN Bleach Deep Dentin	20 g	010
560502	IPS d.SIGN Deep Dentin	20 g	110
560503	IPS d.SIGN Deep Dentin	20 g	120
556536	IPS d.SIGN Deep Dentin	20 g	130
556537	IPS d.SIGN Deep Dentin	20 g	140
556538	IPS d.SIGN Deep Dentin	20 g	210
556539	IPS d.SIGN Deep Dentin	20 g	220
556540	IPS d.SIGN Deep Dentin	20 g	230
560504	IPS d.SIGN Deep Dentin	20 g	240
556541	IPS d.SIGN Deep Dentin	20 g	310
560505	IPS d.SIGN Deep Dentin	20 g	320
560506	IPS d.SIGN Deep Dentin	20 g	330
560507	IPS d.SIGN Deep Dentin	20 g	340
556542	IPS d.SIGN Deep Dentin	20 g	410
556543	IPS d.SIGN Deep Dentin	20 g	420
556544	IPS d.SIGN Deep Dentin	20 g	430
560508	IPS d.SIGN Deep Dentin	20 g	440
556545	IPS d.SIGN Deep Dentin	20 g	510
560509	IPS d.SIGN Deep Dentin	20 g	520
560510	IPS d.SIGN Deep Dentin	20 g	530
560511	IPS d.SIGN Deep Dentin	20 g	540



Art.No.	Description	Quantity	Shade
---------	-------------	----------	-------

IPS d.SIGN Margin A-D Refill 20 g

560530	IPS d.SIGN Margin A-D	20 g	A1
558268	IPS d.SIGN Margin A-D	20 g	A2
558269	IPS d.SIGN Margin A-D	20 g	A3
558270	IPS d.SIGN Margin A-D	20 g	A3,5
560531	IPS d.SIGN Margin A-D	20 g	A4
560532	IPS d.SIGN Margin A-D	20 g	B1
558272	IPS d.SIGN Margin A-D	20 g	B2
558273	IPS d.SIGN Margin A-D	20 g	B3
560533	IPS d.SIGN Margin A-D	20 g	B4
560534	IPS d.SIGN Margin A-D	20 g	C1
560535	IPS d.SIGN Margin A-D	20 g	C2
558274	IPS d.SIGN Margin A-D	20 g	C3
560536	IPS d.SIGN Margin A-D	20 g	C4
558271	IPS d.SIGN Margin A-D	20 g	D2/D3
560537	IPS d.SIGN Margin A-D	20 g	D4



IPS d.SIGN Margin Chromascop Refill 20 g

558192	IPS d.SIGN Bleach Margin	20 g	010
558193	IPS d.SIGN Bleach Margin	20 g	040
560520	IPS d.SIGN Margin	20 g	110
560521	IPS d.SIGN Margin	20 g	120
556551	IPS d.SIGN Margin	20 g	130
556553	IPS d.SIGN Margin	20 g	140
556554	IPS d.SIGN Margin	20 g	210
556555	IPS d.SIGN Margin	20 g	220
556556	IPS d.SIGN Margin	20 g	230
560522	IPS d.SIGN Margin	20 g	240
556557	IPS d.SIGN Margin	20 g	310
560523	IPS d.SIGN Margin	20 g	320
560524	IPS d.SIGN Margin	20 g	330
560525	IPS d.SIGN Margin	20 g	340
556558	IPS d.SIGN Margin	20 g	410
556559	IPS d.SIGN Margin	20 g	420
556560	IPS d.SIGN Margin	20 g	430
560526	IPS d.SIGN Margin	20 g	440
556561	IPS d.SIGN Margin	20 g	510
560527	IPS d.SIGN Margin	20 g	520
560528	IPS d.SIGN Margin	20 g	530
560529	IPS d.SIGN Margin	20 g	540
556562	IPS d.SIGN Margin	20 g	gelb
556563	IPS d.SIGN Margin	20 g	orange
556564	IPS d.SIGN Margin	20 g	orange-pink
556565	IPS d.SIGN Margin	20 g	Add-On
556566	IPS d.SIGN Margin	20 g	Opaque










IPS d.SIGN Bleach BL Refill 20 g

602924	IPS d.SIGN Opaquer Bleach BL	20 g	BL1/BL2
602925	IPS d.SIGN Opaquer Bleach BL	20 g	BL3/BL4
602953	IPS d.SIGN Dentin Bleach BL	20 g	BL1
602954	IPS d.SIGN Dentin Bleach BL	20 g	BL2
602955	IPS d.SIGN Dentin Bleach BL	20 g	BL3
602956	IPS d.SIGN Dentin Bleach BL	20 g	BL4
602957	IPS d.SIGN Incisal Bleach BL	20 g	BL
602947	IPS d.SIGN Margin Bleach BL	20 g	BL1
602948	IPS d.SIGN Margin Bleach BL	20 g	BL4
602950	IPS d.SIGN Deep Dentin Bleach BL	20 g	BL1
602951	IPS d.SIGN Deep Dentin Bleach BL	20 g	BL4
602958	IPS d.SIGN Add-On Bleach BL	20 g	BL



Art.No.	Description	Quantity	Shade
IPS d.SIGN Brilliant Dentin Refill 20 g			
556570	IPS d.SIGN Brilliant Dentin	20 g	yellow
556571	IPS d.SIGN Brilliant Dentin	20 g	orange
556572	IPS d.SIGN Brilliant Dentin	20 g	white
IPS d.SIGN Palatinal Dentin Refill 20 g			
556573	IPS d.SIGN Palatinal Dentin	20 g	yellow
556574	IPS d.SIGN Palatinal Dentin	20 g	red
IPS d.SIGN Mamelon Refill 20 g			
556575	IPS d.SIGN Mamelonmasse	20 g	light
556576	IPS d.SIGN Mamelonmasse	20 g	yellow-orange
556577	IPS d.SIGN Mamelonmasse	20 g	salmon
IPS d.SIGN Insert Refill 20 g			
556578	IPS d.SIGN Insert	20 g	orange
556579	IPS d.SIGN Insert	20 g	grey
IPS d.SIGN Inter Incisal Refill 20 g			
556580	IPS d.SIGN Inter Incisal	20 g	white-blue
IPS d.SIGN Special Incisal Refill 20 g			
556581	IPS d.SIGN Special Incisal	20 g	yellow
556582	IPS d.SIGN Special Incisal	20 g	grey
IPS d.SIGN Transpa Refill 20 g			
556583	IPS d.SIGN Transpa	20 g	blue
556584	IPS d.SIGN Transpa	20 g	brown-grey
556585	IPS d.SIGN Transpa	20 g	orange-grey
556586	IPS d.SIGN Transpa	20 g	neutral
IPS d.SIGN Transpa Refill 100 g			
563527	IPS d.SIGN Transpa	100 g	neutral
IPS d.SIGN Transpa Refill 250 g			
563528	IPS d.SIGN Transpa	250 g	neutral
IPS d.SIGN Occlusal Dentin Refill 20 g			
556587	IPS d.SIGN Occlusal Dentin	20 g	orange
556588	IPS d.SIGN Occlusal Dentin	20 g	brown
IPS d.SIGN Incisal Edge Refill 20 g			
556589	IPS d.SIGN Incisal Edge	20 g	
IPS d.SIGN Cervical Incisal Refill 20 g			
556590	IPS d.SIGN Cervical Incisal	20 g	yellow
556591	IPS d.SIGN Cervical Incisal	20g	orange-pink
556592	IPS d.SIGN Cervical Incisal	20 g	khaki
556593	IPS d.SIGN Cervical Incisal	20 g	orange
IPS d.SIGN Effect Refill 20 g			
556594	IPS d.SIGN Effect	20 g	1
556595	IPS d.SIGN Effect	20 g	2
556596	IPS d.SIGN Effect	20 g	3
556597	IPS d.SIGN Effect	20 g	4
574564	IPS d.SIGN Effect	20 g	5
556598	IPS d.SIGN Effect	20 g	violet

	Art.No.	Description	Quantity	Shade
		IPS d.SIGN Gingiva Opaquer Refill 3 g		
	556623	IPS d.SIGN Gingiva Opaquer	3 g	
		IPS d.SIGN Gingiva Refill 20 g		
	556615	IPS d.SIGN Gingiva	20 g	1
	556616	IPS d.SIGN Gingiva	20 g	2
	556617	IPS d.SIGN Gingiva	20 g	3
	556618	IPS d.SIGN Gingiva	20 g	4
	556802	IPS d.SIGN Gingiva	20 g	5
		IPS d.SIGN Gingiva Modifier Refill 20 g		
	556619	IPS d.SIGN Gingiva Modifier	20 g	1
	556620	IPS d.SIGN Gingiva Modifier	20 g	2
	556621	IPS d.SIGN Gingiva Modifier	20 g	3
	556622	IPS d.SIGN Gingiva Modifier	20 g	4
		IPS d.SIGN Essence Refill 1 g		
	556626	IPS d.SIGN Essence	1 g	orange
	556627	IPS d.SIGN Essence	1 g	yellow
	556628	IPS d.SIGN Essence	1 g	grey
	556629	IPS d.SIGN Essence	1 g	brown
		IPS d.SIGN Shade Refill 3 g		
	561572	IPS d.SIGN Shade	3 g	1
	561573	IPS d.SIGN Shade	3 g	2
	561574	IPS d.SIGN Shade	3 g	3
	561575	IPS d.SIGN Shade	3 g	4
	561576	IPS d.SIGN Shade	3 g	5
	561577	IPS d.SIGN Shade	3 g	6
	561578	IPS d.SIGN Shade	3 g	7
		IPS d.SIGN Stains Refill 1 g		
	556632	IPS d.SIGN Stains	1 g	basic blue
	556633	IPS d.SIGN Stains	1 g	white
	556634	IPS d.SIGN Stains	1 g	mahogany
	556635	IPS d.SIGN Stains	1 g	khaki
	556636	IPS d.SIGN Stains	1 g	orange
	556637	IPS d.SIGN Stains	1 g	grey
	556638	IPS d.SIGN Stains	1 g	vanille
	556639	IPS d.SIGN Stains	1 g	crackliner
	556640	IPS d.SIGN Stains	1 g	olive
	556641	IPS d.SIGN Stains	1 g	yellow
	561102	IPS d.SIGN Stains	1 g	basic yellow
	561103	IPS d.SIGN Stains	1 g	basic red
	561104	IPS d.SIGN Stains	1 g	black
	561105	IPS d.SIGN Stains	1 g	maroon
		IPS d.SIGN Glaze Paste Refill 3 g		
	556643	IPS d.SIGN Glaze Paste	3 g	
		IPS d.SIGN Add-On Refill 20 g		
	556642	IPS d.SIGN Add On	20 g	1
		IPS d.SIGN Opaquer Liquid Refill 15 ml		
	556645	IPS d.SIGN Opaquer Liquid	15 ml	
		IPS d.SIGN Margin Build-up Liquid Refill 60 ml		
	556648	IPS d.SIGN Margin Build-up Liquid	60 ml	

Art.No.	Description	Quantity	Shade
---------	-------------	----------	-------



IPS d.SIGN Build-Up Liquid Medium Refill 60 / 250 ml

This liquid is suitable for mixing IPS d.SIGN materials that will be used for smaller restorations or when smaller amounts of material are needed for 2nd dentin and incisal firings (corrective firing). The Build-Up Liquid Light is suitable for those who prefer a drier working consistency, thus avoiding the need for constant blotting since the liquid evaporates much quicker. In this way, IPS d.SIGN can be easily manipulated with an instrument.

Components: Water, butylene glycol, and additives

556644	IPS d.SIGN Build-Up Liquid Medium	60 ml	
558985	IPS d.SIGN Build-Up Liquid Medium	250 ml	



IPS d.SIGN Build-Up Liquid Light Refill 60 / 250 ml

This liquid is suitable for mixing IPS d.SIGN materials that will be used for smaller restorations or when smaller amounts of material are needed for 2nd dentin and incisal firings (corrective firing). The Build-Up Liquid Light is suitable for those who prefer a drier working consistency, thus avoiding the need for constant blotting since the liquid evaporates much quicker. In this way, IPS d.SIGN can be easily manipulated with an instrument.

Components: Water, butylene glycol, and additives

564045	IPS d.SIGN Build-Up Liquid Light	60 ml	
564043	IPS d.SIGN Build-Up Liquid Light	250 ml	



IPS d.SIGN Build-Up Liquid Premium Refill 60 / 250 ml

This liquid is suitable for mixing IPS d.SIGN materials that will be used for fabricating larger restorations where a longer working time and smooth consistency is desired. The Build-Up Liquid Premium is suitable for technicians who prefer a long working time and moist material coupled with high stability. In this way, IPS d.SIGN can be processed without the need for constant moistening. Since the material stays moist for a longer period of time controlled blotting is required as well as a longer closing time (pre-drying time) for the firing procedure.

Components: Water, butylene glycol-chloride solution

564046	IPS d.SIGN Build-Up Liquid Premium	60 ml	
564044	IPS d.SIGN Build-Up Liquid Premium	250 ml	



IPS d.SIGN Build-Up Liquid Optimum 1 Refill 60 / 250 ml

Excellently suitable for mixing IPS d.SIGN, with its outstanding firing behaviour and lowest possible shrinkage. The materials are easily suctioned off and condensed and are particularly suitable for the manual reduction of materials with the help of instruments.

Components: water, polyethylene glycol (PEG) and additives

575541	IPS d.SIGN Build-Up Liquid Optimum 1	60 ml	
575542	IPS d.SIGN Build-Up Liquid Optimum 1	250 ml	



IPS d.SIGN Build-Up Liquid Optimum 2 Refill 60 / 250 ml

Excellently suitable for mixing IPS d.SIGN, with its properties that support the build-up, as well as the optimized firing properties with low shrinkage. This liquid is particularly suitable for targeted material application without continuous suction.

Components: water, polyethylene glycol (PEG) and additives

575543	IPS d.SIGN Build-Up Liquid Optimum 2	60 ml	
575545	IPS d.SIGN Build-Up Liquid Optimum 2	250 ml	



IPS d.SIGN Glaze and Stain Liquid Refill 15 ml

IPS d.SIGN Glaze and Stain Liquid 15 ml

556646	IPS d.SIGN Glaze and Stain Liquid	15 ml	
--------	-----------------------------------	-------	--

Art.No.	Description	Quantity	Shade
---------	-------------	----------	-------



IPS Model Sealer
The IPS Model Sealer is used to isolate the stone model. It seals the stone pores and thus prevents the moisture of the ceramic from being absorbed by the stone.

IPS Model Sealer Refill 50 ml
593401 IPS Model Sealer 50 ml



IPS Margin Sealer
The IPS Margin Sealer is used to isolate the stone die during the fabrication of a ceramic shoulder using Margin materials.

IPS Margin Sealer Refill 20 ml
593402 IPS Margin Sealer Liquid 20 ml



IPS Ceramic Separating Liquid with Brush
The liquid is used to coat those areas that have previously been impregnated with IPS Model Sealer or IPS Margin Sealer. A clean separation of the ceramic material and the stone model or stone die is thus possible.

IPS Ceramic Separating Liquid with Brush Refill 15 ml
593403 IPS Ceramic Separating Liquid with Brush 15 ml

- Materials Shade Guide**
- 605191 Materials Shade Guide Opaquer
 - 556652 Materials Shade Guide Dentin
 - 556653 Materials Shade Guide Deep Dentin
 - 556654 Materials Shade Guide Margin
 - 556655 Materials Shade Guide Incisal
 - 556656 Materials Shade Guide Impulse 1
 - 556657 Materials Shade Guide Impulse 2
 - 605192 Materials Shade Guide Gingiva
 - 602984 Materials Shade Guide Bleach BL
 - 605190 Materials Shade Guide Opaquer A-D
 - 558278 Materials Shade Guide Dentin A-D
 - 558279 Materials Shade Guide Deep Dentin A-D
 - 558280 Materials Shade Guide Margin A-D
 - 559046 Materials Shade Guide T-Incisal
 - 561317 Materials Shade Guide Deep Dentin 2
 - 561318 Materials Shade Guide Margin 2
 - 561319 Materials Shade Guide Deep Dentin A-D 2
 - 561320 Materials Shade Guide Margin A-D 2



IPS d.SIGN Material box
557684 IPS d.SIGN Material box empty