
IPS Classic®

Metal-ceramic system



A sound concept is the key to success



IPS Classic®

Quality is the secret to success

A classic in dental lab applications: The well-proven IPS Classic metal-ceramic

Distinguished by a balanced distribution of particle sizes, the layering materials of the IPS Classic® metal-ceramic range offer excellent modelling properties and high stability, even after several firing cycles.

Exceptional shade matching results can be achieved due to the ideal level of opacity and coordinated colour saturation of the IPS Classic Paste Opaquer and Powder Opaquer.





Final restoration made of IPS Classic (Oliver Brix, Bad Homburg/Germany)

Overview of benefits

- Flexible handling due to coordinated Powder Opaquer
- Efficient and economical to process
- Highly esthetic results
- Long-term clinical experience

Consistent basis – flexible handling

Even if applied in thin layers, the new IPS Classic Powder Opaquer provides an ideal basis for achieving the selected shade in highly esthetic restorations. The powder opaquer can be applied conventionally with a brush or instrument or by means of the spray-on technique and therefore gives users more flexibility.

Individualized characterization

A comprehensive range of special materials including opalescent Incisal and Effect materials allows a high degree of individuality to create lifelike restorations.

An impressive all-round package

Alloys, Programat® ceramic furnaces and luting materials complement the coordinated range of Ivoclar Vivadent products – for the clinical long-term success of lab fabricated restorations.



This is a product from the "Fixed Prosthetics" product category. The products from this category are optimally coordinated with each other.



