

Politip

Exceptionally smooth composite and amalgam restorations

Popular
general purpose
polishing system



Politip

Exceptionally smooth composite and amalgam restorations

Standard solution suitable for a variety of materials

Politip is a proven two-step polishing system that has been on the market for several years. The popular polishers are suited for polishing composite and amalgam restorations.

First-class results in two steps


- Grey polishers: coarse grit size for pre-polishing
- Green polishers: fine grit size for high gloss polishing

Cost-effective

- Autoclavable and reusable
- Excellent price/performance ratio

Proven

- Long-term clinical experience
- 3 shapes to suit different treatment areas



Delivery forms:

Assortment/Professional Set

Polishing tips, coarse (F):
6 x shape B (Flame)
3 x shape C (Large Cup)
3 x shape D (Small Cup)

Assortment F

Polishing tips, coarse (F):
2 x shape B (Flame)
2 x shape C (Large Cup)
2 x shape D (Small Cup)

Polishing tips, fine grit (P):

6 x shape B (Flame)
3 x shape C (Large Cup)
3 x shape D (Small Cup)

Assortment P

Polishing tips, fine grit (P):
2 x shape B (Flame)
2 x shape C (Large Cup)
2 x shape D (Small Cup)

Refills

Each shape is available separately in packages of 6.



This product forms part of the "Direct Restoratives" product category. Products from this category are optimally coordinated with each other.

Ivoclar Vivadent AG
Bendererstr. 2
9494 Schaan
Liechtenstein
Tel. +423 235 35 35
Fax +423 235 33 60
www.ivoclarvivadent.com

Descriptions and data constitute no warranty of attributes. Printed in Liechtenstein
© Ivoclar Vivadent AG, Schaan Liechtenstein

538115/2015-06-18/EN

ivoclar
vivadent[®]
passion vision innovation