


SpeedCEM® Plus


Placement of crowns and bridges (oxide ceramic, metal-ceramic, lithium-disilicate)
 Befestigung von Kronen und Brücken (Oxidkeramik, Metallkeramik, Lithium-Disilikat)

Zirconium oxide, aluminium oxide, metal-ceramic, metal.
 Zirkoniumoxid, Aluminiumoxid, Metallkeramik, Metall.

Lithium-disilicate.
 Lithium-Disilikat.




Sandblast
Sandstrahlen



H₂O




AIR!
LUFT!




Try-in
Einprobe



Ivoclean




20 s



H₂O




AIR!
LUFT!



Try-in
Einprobe




Monobond
Etch & Prime




H₂O



AIR!
LUFT!




20 s





H₂O

40 s




AIR!
LUFT!


10 s



Self-curing



Quarter Technique



CE 0123



See Instructions for Use
Bitte Gebrauchsinformation beachten