

SpeedCEM® Plus


Cementation of root posts (metal, ceramic, FRC)
Befestigung von Wurzelstiften (Metall, Keramik, FRC)

Root Post / Wurzelstift


Metal / FRC / ceramic: Preparation according to the manufacturer's instructions

Metall / FRC / Keramik: Vorbehandlung gemäss Herstellerangaben

Mono-
bond
Plus



Preparation / Präparation



more information
www.cementation-navigation.com

CE 0123



See Instructions for Use
Bitte Gebrauchsinformation beachten

SpeedCEM® Plus

Cementation of root posts (metal, ceramic, FRC)


Befestigung von Wurzelstiften (Metall, Keramik, FRC)

Cementation / Befestigung

Adhesive: according to manufacturer's instructions

Adhäsiv: gemäss Herstellerangaben

SpeedCEM Plus



Core build-up / Stumpfaufbau

MultiCore Flow

