









Viteo® Base Ti

COMBINATION TABLE

Implant Manufacturer*	Implant System*	Implant Diameter	Torque Diameter (clinical insertion)	Viteo®			Press technique		CAD/CAM technique		
				Viteo® Base Ti	Viteo® Screw	Viteo® Base Trimmer	Viteo® Base Ti Press Sleeve	IPS e.max® Press	Elos Accurate Intra Oral Scanabutments*	Telio® CAD**	IPS e.max® CAD**
											
Dentsply Implants	Astra Tech OsseoSpeed™	3.5/4.0 mm	20Ncm	AT-OS 3.5/4.0 (MD)	AT-OS 3.5/4.0	MD	MD	✓	IO 3 A-B SA	–	MD
		4.5/5.0 mm	25Ncm	AT-OS 4.5/5.0 (MD)	AT-OS 4.5/5.0	MD	MD	✓	IO 3 A-C SA	–	MD
Biomet 3i™	OSSEOTITE® Certain®	3.4 mm	20Ncm	BI-OC 3.4 (SD)	BI-OC 3.4/4.1/5.0	SD	SD	✓	IO 7 B-A	SD	SD
		4.1 mm		SD		SD	✓	IO 7 B-B	SD	SD	
		5.0 mm		MD		MD	✓	IO 7 B-B	MD	MD	
Nobel Biocare®	NobelActive™	3.5 mm	35Ncm	NB-NA 3.5 (MD)	NB-NA 3.5	MD	MD	✓	IO 2 B-A SA	MD	MD
		4.3 mm		NB-NA 4.3 (MD)		MD	MD	✓	IO 2 B-B SA	MD	MD
		5.0 mm		NB-NA 5.0 (MD)		MD	MD	✓	IO 2 B-B SA	MD	MD
Straumann	BoneLevel®	3.3 mm	35Ncm	ST-BL 3.3 (SD)	ST-BL 3.3/4.1/4.8	SD	SD	✓	IO 4A-A SA	SD	SD
		4.1/4.8 mm		SD		SD	✓	IO 4A-B SA	SD	SD	
		4.1/4.8 mm		MD		MD	✓		MD	MD	
Straumann	Standard Implant	4.8 mm RN	35Ncm	ST-SI 4.8 RN (MD)	ST-SI 4.8/6.5	MD	MD	✓	IO 4B-A SA	MD	MD
Nobel Biocare®	Nobel Replace®	3.5 mm NP	35Ncm	NB-NR 3.5 (MD)	NB-NR 3.5	MD	MD	✓	IO 2A-A	MD	MD
		4.3 mm RP		NB-NR 4.3 (MD)		MD	MD	✓	IO 2A-B	MD	MD
		5.0 mm WP		NB-NR 5.0 (MD)		MD	MD	✓	IO 2A-C	MD	MD
Nobel Biocare®	Brannemark System	3.5 mm NP	35Ncm	NB-BS 3.5 (MD)	NB-BS 3.5	MD	MD	✓	IO 6A-A	MD	MD
		4.1 mm RP		NB-BS 4.1 (MD)		MD	MD	✓	IO 6A-B	MD	MD
Biomet 3i™	OSSEOTITE External Hex	3.4 mm	32Ncm	BI-OT 3.4 (MD)	BI-OT 3.4/4.1/5.0	MD	MD	✓	IO 7A-A	MD	MD
		4.1 mm		MD		MD	✓	IO 7A-B	MD	MD	
		5.0 mm		MD		MD	✓		MD	MD	
Dentsply Implants	Xive/Frialit	3.4 mm	24Ncm	DI-XI/FR 3.4 (SD)	DI-XI/FR 3.4/3.8/4.5/5.5	SD	SD	✓	IO 8A-B	SD	SD
		3.8 mm		SD		SD	✓	IO 8A-C	SD	SD	
		4.5 mm		SD		SD	✓	IO 8A-D	SD	SD	
Zimmer®	Screw Vent®	3.5 mm	30Ncm	ZI-SV 3.5 (SD)	ZI-SV 3.5/4.5/5.7	SD	SD	✓	IO 5A-A	SD	SD
		4.5 mm		SD		SD	✓	IO 5A-B	SD	SD	
		5.7 mm		MD		MD	✓	IO 5A-C	MD	MD	

* Not a registered trademark of Ivoclar Vivadent

** Blocks with a prefabricated interface for Viteo® Base Ti are available for PlanMill® 405*; other systems in preparation.

CE 0123

Rx ONLY
For dental use only!



Manufacturer:
Ivoclar Vivadent AG, 9494 Schaan/Liechtenstein
www.ivoclarvivadent.com











see instructions

Date information prepared: 2019-04, Rev. 4
687083/D/E

ivoclar
digital®

Viteo® Base Ti

KOMBINATIONSTABELLE

Implantat-Hersteller*	Implantat-System*	Implantat-Durchmesser	Empfohlener Drehmoment (klinische Eingliederung)	Viteo®			Presstechnik		CAD/CAM-Technik		
				Viteo® Base Ti	Viteo® Screw	Viteo® Base Ti Trimmer	Viteo® Base Ti Press Sleeve	IPS e.max® Press	Elos Accurate Intra Oral Scanabutments *	Telio® CAD **	IPS e.max® CAD **
											
Dentsply Implants	Astra Tech OsseoSpeed™	3.5/4.0 mm	20Ncm	AT-OS 3.5/4.0 (MD)	AT-OS 3.5/4.0	MD	MD	✓	IO 3 A-B SA	–	MD
		4.5/5.0 mm	25Ncm	AT-OS 4.5/5.0 (MD)	AT-OS 4.5/5.0	MD	MD	✓	IO 3 A-C SA	–	MD
Biomet 3i™	OSSEOTITE® Certain®	3.4 mm	20Ncm	BI-OC 3.4 (SD)	BI-OC 3.4/4.1/5.0	SD	SD	✓	IO 7 B-A	SD	SD
		4.1 mm		BI-OC 4.1 (SD)		SD	SD	✓	IO 7 B-B	SD	SD
		5.0 mm		BI-OC 5.0 (MD)		MD	MD	✓	IO 7 B-B	MD	MD
Nobel Biocare®	NobelActive™	3.5 mm	35Ncm	NB-NA 3.5 (MD)	NB-NA 4.3/5.0	MD	MD	✓	IO 2 B-A SA	MD	MD
		4.3 mm		NB-NA 4.3 (MD)		MD	MD	✓	IO 2 B-B SA	MD	MD
		5.0 mm		NB-NA 5.0 (MD)		MD	MD	✓	IO 2 B-B SA	MD	MD
Straumann	BoneLevel®	3.3 mm	35Ncm	ST-BL 3.3 (SD)	ST-BL 3.3/4.1/4.8	SD	SD	✓	IO 4A-A SA	SD	SD
		4.1/4.8 mm		ST-BL 4.1/4.8 (SD)		SD	SD	✓	IO 4A-B SA	SD	SD
		4.1/4.8 mm		ST-BL 4.1/4.8 (MD)		MD	MD	✓		MD	MD
Straumann	Standard Implant	4.8 mm RN	35Ncm	ST-SI 4.8 RN (MD)	ST-SI 4.8/6.5	MD	MD	✓	IO 4B-A SA	MD	MD
Nobel Biocare®	Nobel Replace®	3.5 mm NP	35Ncm	NB-NR 3.5 (MD)	NB-NR 4.3/5.0/6.0	MD	MD	✓	IO 2A-A	MD	MD
		4.3 mm RP		NB-NR 4.3 (MD)		MD	MD	✓	IO 2A-B	MD	MD
		5.0 mm WP		NB-NR 5.0 (MD)		MD	MD	✓	IO 2A-C	MD	MD
Nobel Biocare®	Brannemark System	3.5 mm NP	35Ncm	NB-BS 3.5 (MD)	NB-BS 4.1	MD	MD	✓	IO 6A-A	MD	MD
		4.1 mm RP		NB-BS 4.1 (MD)		MD	MD	✓	IO 6A-B	MD	MD
Biomet 3i™	OSSEOTITE External Hex	3.4 mm	32Ncm	BI-OT 3.4 (MD)	BI-OT 3.4/4.1/5.0	MD	MD	✓	IO 7A-A	MD	MD
		4.1 mm		BI-OT 4.1 (MD)		MD	MD	✓	IO 7A-B	MD	MD
		5.0 mm		BI-OT 5.0 (MD)		MD	MD	✓		MD	MD
Dentsply Implants	Xive/Frialit	3.4 mm	24Ncm	DI-XI/FR 3.4 (SD)	DI-XI/FR 3.4/3.8/4.5/5.5	SD	SD	✓	IO 8A-B	SD	SD
		3.8 mm		DI-XI/FR 3.8 (SD)		SD	SD	✓	IO 8A-C	SD	SD
		4.5 mm		DI-XI/FR 4.5 (SD)		SD	SD	✓	IO 8A-D	SD	SD
Zimmer®	Screw Vent®	3.5 mm	30Ncm	ZI-SV 3.5 (SD)	ZI-SV 3.5/4.5/5.7	SD	SD	✓	IO 5A-A	SD	SD
		4.5 mm		ZI-SV 4.5 (SD)		SD	SD	✓	IO 5A-B	SD	SD
		5.7 mm		ZI-SV 5.7 (MD)		MD	MD	✓	IO 5A-C	MD	MD

* Kein eingetragenes Warenzeichen der Ivoclar Vivadent

** Blöcke mit vorgefertigter Schnittstelle für Viteo® Base Ti sind verfügbar für PlanMill® 405*; weitere Systeme in Vorbereitung.