

SERVICE+ | ÜBERSICHT IMPLANTATSCHNITTSTELLEN

Hersteller*	Implantat-System*	Implantat-Plattform (ø mm)					Abutmentlevel	Einteilige Abutments	Implantatgetragene Suprakonstruktionen
Biomet 3i	Osseotite®	3,4		4,1	5,0	6,0		✓	✓
	Osseotite® Certain®	3,4		4,1	5,0	6,0		✓	✓
CAMLOG®	Camlog®	3,3	3,8	4,3	5,0	6,0		✓	✓
	Conelog®	3,3	3,8	4,3	5,0			✓	✓
	iSy®		3,8	4,3	5,0			✓	✓
Dentsply Implants	Astra Tech Osseo Speed™	3,5/4,0		4,5/5,0				✓	✓
	Frialit® + Xive®	3,4	3,8	4,5/5,5				✓	✓
	Ankylos®						Balance Basisaufbau C / schmal		✓
medentis medical	ICX	3,45	3,75	4,1	4,8			✓	✓
Nobel Biocare	Brånemark System®	NP 3,5		RP 4,1	WP 5,1			✓	✓
	NobelActive™	NP 3,5		RP 4,3/5,0				✓	✓
	NobelReplace™	NP 3,5		RP 4,3	WP 5,0	6,0		✓	✓
							Multi-unit		✓
Straumann	Bone Level	NC 3,3		RC 4,1	RC 4,8			✓	✓
	synOcta®				RN 4,8	WN 6,5		✓	✓
Zimmer	Screw-Vent®	3,5		4,5	5,7			✓	✓

*Die aufgeführten Implantat-Systeme und Marken sind keine eingetragenen Warenzeichen der Ivoclar Vivadent AG.