

# Programat® S1 1600

## Software Revision Log

Version	Changes from previous version	Status
<b>V6.0</b> 11.03.2019	<b>Software optimized</b> <ul style="list-style-type: none"> <li>• New program for IPS e.max ZirCAD Prime</li> <li>• Program order modified               <ul style="list-style-type: none"> <li>○ <b>PLEASE NOTE: Individual programs on program number 23 and 24 will be deleted in case of an update!</b></li> </ul> </li> <li>• Various software improvements</li> </ul>	<b>Currently</b> delivered version
<b>V5.0</b> 29.05.2017	<b>Software optimized – new programs</b> <ul style="list-style-type: none"> <li>• New programs for IPS e.max ZirCAD MT Multi/MT/LT</li> <li>• Extended function with Programat Dosto-Tray</li> <li>• 50 programs available (22 pre-programmed)               <ul style="list-style-type: none"> <li>○ <b>PLEASE NOTE: All individual programs will be deleted in case of an update!</b></li> </ul> </li> <li>• Various software improvements</li> </ul>	<b>No</b> longer available!
<b>V4.0</b> 15.09.2015	<b>Basic Version of the Programat S1 1600</b>	<b>No</b> longer available!

Version 11.03.2019