

Programat® CS2 / Programat® CS3

Program structure

Program	Description
1	Speed Crystallization program for IPS e.max CAD HT/MT/LT
2	Crystallization/Glaze program for IPS e.max CAD HT/MT/LT
3	Crystallization program for IPS e.max CAD Impulse/MO/HT/MT/LT
4	Corrective/Stain/Glaze program for IPS e.max CAD with IPS e.max CAD Crystall./materials
5	Stain/Glaze program for IPS e.max ZirCAD with IPS e.max CAD Crystall./materials
6	Stain/Glaze program für IPS Empress CAD with IPS e.max CAD Crystall./materials
7	Stain/Glaze program for IPS Empress CAD with IPS Empress Universal
8	Fusion/Crystallization program for the IPS e.max CAD-on technique
9	Characterization/Glaze program for the IPS e.max CAD-on technique

Valid from software:
Programat CS2 V4.0
Programat CS3 V3.0

Version: 1
Date of issue: 2017-06
en



Programat® CS2 / Programat® CS3

Programmstruktur

Programm	Beschreibung
1	Speed Kristallisationsprogramm für IPS e.max CAD HT/MT/LT
2	Kristallisations-/Glanzprogramm für IPS e.max CAD HT/MT/LT
3	Kristallisationsprogramm für IPS e.max CAD Impulse/MO/HT/MT/LT
4	Korrektur-/Malfarben-/Glanzprogramm für IPS e.max CAD mit IPS e.max CAD Crystall./Materialien
5	Malfarben-/Glanzprogramm für IPS e.max ZirCAD mit IPS e.max CAD Crystall./Materialien
6	Malfarben-/Glanzprogramm für IPS Empress CAD mit IPS e.max CAD Crystall./Materialien
7	Malfarben-/Glanzprogramm für IPS Empress CAD mit IPS Empress Universal
8	Füge-/Kristallisationsprogramm für IPS e.max CAD-on Technik
9	Charakterisierung-/Glanzprogramm für IPS e.max CAD-on Technik

Gültig ab Software:
Programat CS2 V4.0
Programat CS3 V3.0

Version: 1
Ausgabedatum: 2017-06
de



ivoclar
vivadent®
passion vision innovation