


## 1 Identification

- **Product identifier**
- **Trade name:** ***Tetric CAD***
- **Application of the substance / the mixture** *Composite block for the CAD/CAM technology*
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
*Ivoclar Vivadent Inc.*  
*175 Pineview Drive, Amherst, N.Y. 14228*  
*USA*  
*Tel. +1 800 533 6825*  
*Fax +1 716 691 2285*
  
- Ivoclar Vivadent Inc.*  
*1-6600 Dixie Road*  
*Mississauga, Ontario*  
*L5T 2Y2*  
*Canada*  
*Phone: +1 905 670 8499*  
*Fax: +1 905 670 3102*
- **Information department:** *Quality Assurance / Regulatory Affairs*
- **Emergency telephone number:**  
*24 Hour Emergency Assistance:*  
*Emergency-Call USA - Infotrac: 1-800-535-5053*  
*Emergency-Call Canada - Canutec: 1-613-996-6666*
  
- General SDS Assistance:*  
*US: 1-800-533-6825*  
*Canada: 1-800-263-8182*

## 2 Hazard(s) identification

- **Classification of the substance or mixture**  
*The product is not classified according to the Globally Harmonized System (GHS).*

---

- **Label elements**
- **GHS label elements** *Void*
- **Hazard pictograms** *Void*
- **Signal word** *Void*
- **Hazard statements** *Void*
- **Classification system:**
- **NFPA ratings (scale 0 - 4)**  
 *Health = 0*  
*Fire = 0*  
*Reactivity = 0*
- **HMIS-ratings (scale 0 - 4)**  

HEALTH	0
FIRE	0
REACTIVITY	0

*Health = 0*  
*Fire = 0*  
*Reactivity = 0*
- **Other hazards** *Avoid breathing sanding dust /swarf.*
- **Results of PBT and vPvB assessment**
- **PBT:** *Not applicable.*

## Safety Data Sheet acc. to OSHA HCS

Printing date 11/23/2017

Version number 1

Reviewed on 11/23/2017

**Trade name: Tetric CAD**

· **vPvB:** Not applicable.

(Contd. of page 1)

### 3 Composition/information on ingredients

- **Chemical characterization: Mixtures**
- **Description:** Mixture of the substances listed below with nonhazardous additions.
- **Dangerous components:** Void

### 4 First-aid measures

- **Description of first aid measures**
- **General information:** No special measures required.
- **After inhalation:** Grinding dust: Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:**  
Grinding dust: Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing:** Grinding dust: Rinse out mouth and then drink plenty of water.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

### 5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**  
CO<sub>2</sub>, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:** No special measures required.

### 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** No special measures required.
- **Methods and material for containment and cleaning up:** Pick up mechanically.
- **Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.
- **Protective Action Criteria for Chemicals**

· **PAC-1:**

CAS: 1317-61-9	iron oxid black (E172)	21 mg/m <sup>3</sup>
CAS: 1309-37-1	diiron trioxide (red)	15 mg/m <sup>3</sup>
CAS: 20344-49-4	iron hydroxide oxide	24 mg/m <sup>3</sup>
CAS: 13463-67-7	titanium dioxide	30 mg/m <sup>3</sup>

· **PAC-2:**

CAS: 1317-61-9	iron oxid black (E172)	230 mg/m <sup>3</sup>
----------------	------------------------	-----------------------

(Contd. on page 3)

## Safety Data Sheet

acc. to OSHA HCS

Printing date 11/23/2017

Version number 1

Reviewed on 11/23/2017

**Trade name: Tetric CAD**

(Contd. of page 2)

CAS: 1309-37-1	diiron trioxide (red)	360 mg/m <sup>3</sup>
CAS: 20344-49-4	iron hydroxide oxide	260 mg/m <sup>3</sup>
CAS: 13463-67-7	titanium dioxide	330 mg/m <sup>3</sup>
<b>· PAC-3:</b>		
CAS: 1317-61-9	iron oxid black (E172)	1,400 mg/m <sup>3</sup>
CAS: 1309-37-1	diiron trioxide (red)	2,200 mg/m <sup>3</sup>
CAS: 20344-49-4	iron hydroxide oxide	1,600 mg/m <sup>3</sup>
CAS: 13463-67-7	titanium dioxide	2,000 mg/m <sup>3</sup>

### 7 Handling and storage

- **Handling:**
- **Precautions for safe handling**  
Only adequately trained personnel should handle this product.  
For use in dentistry only.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** Store only in the original receptacle.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Protect from heat and direct sunlight.
- **Specific end use(s)** No further relevant information available.

### 8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**  
The following constituents are the only constituents of the product which have a PEL, TLV or other recommended exposure limit.  
At this time, the other constituents have no known exposure limits.
- **Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**  
Usual hygienic measures for dental practice and dental laboratories.  
Keep away from foodstuffs, beverages and feed.  
Wash hands before breaks and at the end of work.
- **Breathing equipment:**  
Mechanical processing:  
Use respiratory protective device against the effects of dust.
- **Recommended filter device for short term use:** Filter P1
- **Protection of hands:**  
Protective gloves should always be worn during mechanical and thermal processing.  
After use of gloves apply skin-cleaning agents and skin cosmetics.
- **Material of gloves**  
Mechanical processing:  
Leather gloves  
Strong gloves

(Contd. on page 4)

US

## Safety Data Sheet acc. to OSHA HCS

Printing date 11/23/2017

Version number 1

Reviewed on 11/23/2017

Trade name: **Tetric CAD**

(Contd. of page 3)

- **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

- **Eye protection:**

Always wear safety goggles during mechanical processing (grinding, sawing /cutting, drilling, milling).

### 9 Physical and chemical properties

- **Information on basic physical and chemical properties**

- **General Information**

- **Appearance:**

· <b>Form:</b>	Solid
· <b>Color:</b>	According to product specification
· <b>Odor:</b>	Odorless
· <b>Odor threshold:</b>	Not determined.

- **pH-value:** Not applicable.

- **Change in condition**

· <b>Melting point/Melting range:</b>	Undetermined.
· <b>Boiling point/Boiling range:</b>	Undetermined.

- **Flash point:** Not applicable.

- **Flammability (solid, gaseous):** Not determined.

- **Auto igniting:** Product is not selfigniting.

- **Danger of explosion:** Product does not present an explosion hazard.

- **Vapor pressure:** Not applicable.

- **Density at 20 °C (68 °F):** 1.96 g/cm<sup>3</sup> (16.356 lbs/gal)

- **Relative density** Not determined.

- **Vapor density** Not applicable.

- **Evaporation rate** Not applicable.

- **Solubility in / Miscibility with**

· **Water:** Insoluble.

- **Partition coefficient (n-octanol/water):** Not determined.

- **Viscosity:**

· **Dynamic:** Not applicable.

· **Kinematic:** Not applicable.

- **Other information** No further relevant information available.

### 10 Stability and reactivity

- **Reactivity** No further relevant information available.

- **Chemical stability** Stable under normal handling and storage conditions.

- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.

- **Possibility of hazardous reactions** No dangerous reactions known.

- **Conditions to avoid** Keep away from heat and direct sunlight.

- **Incompatible materials:** No further relevant information available.

- **Hazardous decomposition products:** None under normal conditions of storage and use.

US

(Contd. on page 5)

# Safety Data Sheet

## acc. to OSHA HCS

Printing date 11/23/2017

Version number 1

Reviewed on 11/23/2017

Trade name: **Tetric CAD**

(Contd. of page 4)

### 11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**
- **on the skin:**  
Irritant to skin and mucous membranes.  
Based on available data, the classification criteria are not met.
- **on the eye:** Irritating effect.
- **Sensitization:** Sensitization possible through skin contact.
- **Additional toxicological information:** No further relevant information available.
- **Carcinogenic categories**

- **NTP (National Toxicology Program)**

None of the ingredients is listed.

- **OSHA-Ca (Occupational Safety & Health Administration)**

None of the ingredients is listed.

### 12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Behavior in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:** Generally not hazardous for water
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

### 13 Disposal considerations

- **Waste treatment methods**
- **Recommendation:**  
Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to official regulations.

### 14 Transport information

- **UN-Number**
- **DOT, ADR/RID/ADN, ADN, IMDG, IATA** Void

- **UN proper shipping name**
- **DOT, ADR/RID/ADN, ADN, IMDG, IATA** Void

(Contd. on page 6)

US

## Safety Data Sheet

acc. to OSHA HCS

Printing date 11/23/2017

Version number 1

Reviewed on 11/23/2017

**Trade name: Tetric CAD**

(Contd. of page 5)

· <b>Transport hazard class(es)</b>	
· <b>DOT, ADR/RID/ADN, ADN, IMDG, IATA</b>	
· <b>Class</b>	Void
· <b>Packing group</b>	
· <b>DOT, ADR/RID/ADN, IMDG, IATA</b>	Void
· <b>Environmental hazards:</b>	
· <b>Marine pollutant:</b>	No
· <b>Special precautions for user</b>	Not applicable.
· <b>Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code</b>	Not applicable.
· <b>Transport/Additional information:</b>	Product is not classified as a dangerous good for transport (ADR, IMDG, IATA).
· <b>UN "Model Regulation":</b>	Void

### 15 Regulatory information

· **Safety, health and environmental regulations/legislation specific for the substance or mixture**

· **Sara**

· **Section 355 (extremely hazardous substances):**

None of the ingredients is listed.

· **Section 313 (Specific toxic chemical listings):**

None of the ingredients is listed.

· **TSCA (Toxic Substances Control Act):**

All ingredients are listed.

· **Proposition 65**

· **Chemicals known to cause cancer:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

· **Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

· **Carcinogenic categories**

· **EPA (Environmental Protection Agency)**

None of the ingredients is listed.

· **TLV (Threshold Limit Value established by ACGIH)**

None of the ingredients is listed.

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

· **GHS label elements** Void

· **Hazard pictograms** Void

· **Signal word** Void

(Contd. on page 7)

## Safety Data Sheet acc. to OSHA HCS

Printing date 11/23/2017

Version number 1

Reviewed on 11/23/2017

**Trade name: Tetric CAD**

(Contd. of page 6)

· **Hazard statements** *Void*

· **National regulations:**

· **Other regulations, limitations and prohibitive regulations**

*The product is a medical device according to the Directive 93/42/EEC.*

*This product is classified as a medical device under US and Canadian regulations and has been reviewed by the US Food and Drug Administration and Health Canada.*

· **Chemical safety assessment:** *A Chemical Safety Assessment has not been carried out.*

### 16 Other information

*This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.*

· **Date of preparation / last revision** 11/23/2017 / -

· **Abbreviations and acronyms:**

*ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)*

*IMDG: International Maritime Code for Dangerous Goods*

*DOT: US Department of Transportation*

*IATA: International Air Transport Association*

*ACGIH: American Conference of Governmental Industrial Hygienists*

*EINECS: European Inventory of Existing Commercial Chemical Substances*

*ELINCS: European List of Notified Chemical Substances*

*CAS: Chemical Abstracts Service (division of the American Chemical Society)*

*NFPA: National Fire Protection Association (USA)*

*HMIS: Hazardous Materials Identification System (USA)*

*PBT: Persistent, Bioaccumulative and Toxic*

*vPvB: very Persistent and very Bioaccumulative*

*NIOSH: National Institute for Occupational Safety*

*OSHA: Occupational Safety & Health*

*TLV: Threshold Limit Value*

*PEL: Permissible Exposure Limit*

*REL: Recommended Exposure Limit*